

## SecurityLINK Fingertight Fittings

### Tips for Care and Use





Thank you for purchasing SecurityLINK Fittings. Below are recommended instructions for the care and use of this product.

#### General Information

Each SecurityLINK Fitting is developed, produced and distributed according to a Quality Management System that is certified for compliance with ISO 9001:2015. This product has passed all Phenomenex test requirements and thus meets quality control standards.

#### Inspection

Upon receipt of SecurityLINK, please verify the product you received is the one that was ordered (i.e. length, ID, Tubing material). Additionally, remove the red fitting tips and inspect both fittings and tubing for any physical damage potentially caused during shipment.

Specifications	Tubing Options		
	PEEKsil™	PEEK-Lined Stainless Steel	Stainless Steel
			
<b>Description</b>	<ul style="list-style-type: none"> <li>• Polymer-sheathed fused silica tubing</li> <li>• Excellent chemical compatibility</li> <li>• Low adsorption</li> <li>• Fused silica core provides a consistent pathway with very tight tolerances</li> <li>• Flexible, making connections easier</li> <li>• Compatible with most organic solvents</li> <li>• Resistant to strong acids (except hydrofluoric acid)</li> <li>• Lower carry over between samples improved reproducibility</li> <li>• Recommended for hydrophobic samples</li> </ul>	<ul style="list-style-type: none"> <li>• Chemically inert to most commonly used solvents</li> <li>• Smooth ID surface finish</li> </ul>	<ul style="list-style-type: none"> <li>• Seamless 316 stainless steel</li> <li>• Capable of tight bend radii to facilitate routing in the instrument</li> </ul>
<b>Direct Replacement for</b>	PEEK-Lined Stainless Steel	Stainless Steel	
<b>pH Range</b>	0–10	1–14	1–14
<b>Autoclavable</b>	Yes	Yes	Yes
<b>Fitting Size</b>	1/16"	1/16"	1/16"
<b>ID Range</b>	25 µm to 100 µm	25 µm to 100 µm	100 µm to 254 µm
<b>Pressure Rating</b>	19,000 psi	19,000 psi	19,000 psi
<b>Temperature Limit</b>	120 °C	120 °C	120 °C
<b>Biocompatible</b>	Yes	Yes	No
<b>Tubing Bend Radius Limit</b>	1.00" (~25 mm)	0.25" (~6.3 mm)	0.25" (~6.3 mm)
<b>Solvent Compatibility</b>	There is a small piece of PEEK at the end of each SecurityLINK fittings. Therefore solvents known to impact PEEK may impact the overall lifetime. See page 2 for more details.		

## Solvent Compatibility

All three SecurityLINK tubing types contain a small piece of PEEK at the end of each fitting. PEEK exhibits excellent resistance to a wide range of organic and inorganic chemicals. PEEK is compatible with almost any of the solvents used in HPLC. The only solvent which will attack PEEK are concen-

trated nitric acid and sulphuric acids. Examples of commonly used solvents known to impact PEEK are: DMSO, THF, Methylenechloride and solutions containing above 20-30% nitric acid. Therefore the use of these solvents may impact the overall lifetime of SecurityLINK fittings.

### PEEK Chemical Resistance Chart

Chemical	Rating	Chemical	Rating
Acetaldehyde	A	Butanol	A
Acetic acid (20%)	A	Calcium hydroxide	A
Acetic acid (80%)	A	Carbon tetrachloride	A
Acetic acid (glacial)	A	Chlorine (gas)	A
Acetone	A	Chlorine (liquid)	C
Acetonitrile	A	Chloroacetic acid	A
Acrylic acid	A	Chlorobenzene	A
Ammonia, anhydrous	A	Chloroform	A
Ammonia (10%)	A	Cyclohexane	A
Ammonia (Liquid)	B	Cyclohexanone	A
Ammonium hydroxide	A	Diethylamine	A
Aqua regia	C	Diethylether	A
Aromatic hydrocarbons	A	Diethylformamide	A
Benzene	A	Dioxane	A
Benzoic acid	A	Ethanol	A
Benzaldehyde	A	Ether	A
Bromine/dibromoethane	C	Ethyl acetate	A
Bromine (dry)	C	Ethylene glycol	A
Bromine (wet)	C	Heptane	A
Boric acid	A	Hexane	A

A = Suitable, B = Marginal-dependent on application, C = Not recommended

Chemical	Rating	Chemical	Rating
Hydrobromic acid (100%)	C	Phosphoric acid (100%)	A
Hydrobromic acid (20%)	C	Phosphoric acid (40%)	A
Hydrochloric acid (100%)	A	Potassium hydroxide (dilute)	A
Hydrochloric acid (20%)	A	Potassium hydroxide (70%)	A
Hydrofluoric acid (100%)	C	Propanol	A
Hydrogen peroxide (100%)	A	Pyridine	A
Hydrogen peroxide (50%)	A	Sodium hydroxide (80%)	A
Hydrogen peroxide (10%)	A	Sodium hydroxide (20%)	A
Iso-octane	A	Sulphuric acid (100%)	C
Isopropanol	A	Sulphuric acid (75%)	C
Ketones	A	Sulphuric acid (40%)	A
Methanol	A	Tetrahydrofuran	B
Methyl ethylketone	A	Toluene	A
Methylene chloride	B	Trichloroacetic acid	A
Nitric acid (100%)	C	Trichloroethane	A
Nitric acid (20%)	A	Trichloroethylene	A
Pentane	A	Water	A
Perchloric acid	A	Water (distilled)	A
Phenol (dilute)	A	Xylene	A
Phenol (conc.)	C		

### The chemical resistance of PEEK at high temperature

Chemical	Temperature	Rating
Acetic acid	200 °C (420 °F)	B
Ethylene glycol	200 °C (420 °F)	B
Ethylene glycol (50%)	140 °C (284 °F)	A
Methylethylketone	200 °C (420 °F)	C
Nitrobenzene	200 °C (420 °F)	C
Phosphoric acid (50%)	200 °C (420 °F)	A
Sulphuric acid (50%)	200 °C (420 °F)	B

Chemical	Temperature	Rating
Sodium hydroxide solution	200 °C (420 °F)	A
Liquid ammonia	200 °C (420 °F)	B
Sulphur dioxide gas	200 °C (420 °F)	A
Hydrogen sulphide gas	200 °C (420 °F)	A
Carbon monoxide gas	200 °C (420 °F)	A
Ammonia gas	200 °C (420 °F)	A
Methane	200 °C (420 °F)	A

## Fitting Installation Instructions

- 1 Insert SecurityLINK into column port or system.



- 2 Fingertighten until first "CLICK" feedback is received.



## Installation Considerations

SecurityLINK is guaranteed up to 100 leak-free connections with a single click. **Just one click is needed to achieve a leak-free connection. Remember to stop tightening after the initial click.** A connection is defined as a single click of SecurityLINK to a column, system port, or union. It is suggested to keep count of how many connections the SecurityLINK achieves on a specific system and method when initially installed. This will give a benchmark on when to proactively replace a SecurityLINK before any tubing failure is observed. *The use of pliers to achieve the 'click' may damage the torque mechanism and the tubing underneath.*

**Do not Bend or remove SecurityLINK while under pressure** as it can cause damage to the tubing material and/or fitting material of your SecurityLINK, leading to leaking.

The bend radius limit of SecurityLINK differs by the type of tubing. **SecurityLINK PEEKsil has a minimum bend radius of 1.00" (25 mm) while the SecurityLINK Stainless Steel and PeekLined Stainless Steel have a minimum bend radius of 0.25" (6.33 mm).** One method to approximate the 1.00" bend radius of the SecurityLINK PEEKsil is to create a circle with your thumb and index fingers which is ~1.00". This bend should be safe for the tubing and reduce the likelihood of an accidental weak spot or break in the tubing where leaking can occur.

The SecurityLINK Fitting cannot be removed from the tubing. If the fitting head is too large for the system port we recommend the single-sided SecurityLINK.

## Fitting Cleaning

If contamination of the SecurityLINK is suspected, the tubing can be cleaned by removing the end-fitting in the down direction of the flow. If this includes the column in the path, first remove the column and put an appropriate union in place as to not harm the column. Place the open fitting tip in a clean beaker and flush with HPLC grade water followed by Methanol

## Fitting Compatibility

### Unions and Ports

SecurityLINK is compatible with unions that have a thru hole  $\leq 0.020$ " (~0.4 mm). Zero-Dead Volume Unions are NOT compatible with SecurityLINK.

Phenomenex unions AQ0-8507 and AQ0-1674 are recommended for SecurityLINK. Please note, that AQ0-1674 is only pressure rated to 5K psi.

### Columns and Systems

SecurityLINK incorporates 10-32 Fittings and designed for  $\frac{1}{16}$  in. ports with a thru hole  $\leq 0.020$ " (~0.4 mm).

## Product Warranties

Phenomenex products are warranted to be free of defects in material and workmanship. If you are unsatisfied for any reason, please give your Phenomenex Technical Representative a call. We'll do our best to solve the problem to your satisfaction. Should it become necessary to return the product, a Return Authorization Number must be obtained from Phenomenex first.



Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

For any additional questions visit:

[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

or submit an inquiry to:

[support.phxtechnical@zendesk.com](mailto:support.phxtechnical@zendesk.com)

### Terms and Conditions

Subject to Phenomenex standard terms and conditions which can be viewed at [www.phenomenex.com/TermsandConditions](http://www.phenomenex.com/TermsandConditions)

### Trademarks

SecurityLINK and BE-HAPPY are trademarks of Phenomenex. PEEKsil is a trademark of Trajan Scientific Australia Pty Ltd.

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.

**QUALITY  
MANAGEMENT SYSTEM  
CERTIFIED BY DNV GL  
= ISO 9001:2015 =**