



TN-1333

One Minute Screening of a Multiclass Drug Panel with Liquid Chromatography Coupled to High Resolution Mass Spectrometry

Lawrence Andrade, Ana Grenier, and Bryan Tackett, PhD

¹Dominion Diagnostics LLC, 211 Circuit Drive, North Kingstown, RI 02852 USA

²Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Introduction

Due to regulatory changes in clinical urine drug testing, preliminary screening is now required before definitive testing can be performed. In many instances there are gaps in commercially available screening technologies such as EIA or ELISA. HRMS has emerged as an alternative technology to fill in the gaps. The purpose of this study was to validate a bioanalytical (LC-MS) method for preliminary screening of 9 drugs in urine using a Kinetex 2.6 μm EVO C18 LC column. This ultra-fast screening method (cycle time ~ 1 min) includes analytes where traditional screening methods, such as EIA and ELISA, could not be used because commercial kits are not available, or are cost prohibitive. Complete hydrolysis was possible with a 5-minute incubation utilizing a novel, purified hydrolysis enzyme. Linearity, LLOQ, Precision and Accuracy, Specificity, Selectivity, and Matrix Effects were evaluated and met all specifications.

HRMS Conditions

Polarity:	Positive
Full Scan:	100-600 m/z
Sheath Gas Flow Rate:	65
Aux Gas Flow Rate:	20
Sweep Gas Flow Rate:	2
Spray Voltage:	3.5 kV
Collision Gas:	9 psi
Spray Current:	0 μA
Capillary Temperature:	300 $^{\circ}\text{C}$
S-lens RF Level:	55
Aux Gas Heater Temperature:	420 $^{\circ}\text{C}$
Inclusion:	10 ppm
Resolution:	140 K

Sample Preparation

In a 96-well plate, 75 μL of samples, calibrators, and controls were diluted with internal standard solution (40 μL), and buffered enzyme solution (300 μL). The plate was capped and incubated for 5 minutes at room temperature followed by centrifugation at 4000 rpm for 7 minutes. 7.5 μL was injected. Positive identity was confirmed by accurate mass, retention, time and isotope pattern.

LC Conditions

Column:	Kinetex™ 2.6 μm EVO C18	
Dimensions:	20 x 2.1 mm	
Part No.:	OOM-4725-AN	
Mobile Phase:	A: 0.1 % Formic Acid in Water B: 0.1 % Formic Acid in Acetonitrile	
Gradient:	Time (min)	%B
	0	5
	0.5	60
	0.6	99
	0.7	99
	0.8	5
Flow Rate:	0.6 mL/min	
Injection Volume:	7.5 μL	
Temperature:	30 $^{\circ}\text{C}$	
LC System:	Waters® ACQUITY® I Class UPLC	
Detection:	HRMS	
Detector:	Q Exactive™ Orbitrap™	

Table 1. Full Scan: Exact Mass.

Compound	Chemical Formula	Extracted Mass	Internal Standard
Carisoprodol	$\text{C}_{12}\text{H}_{24}\text{N}_2\text{O}_4$	261.1809	Carisoprodol-D ₇
Fentanyl	$\text{C}_{22}\text{H}_{28}\text{N}_2\text{O}$	337.2274	Fentanyl-D ₅
Gabapentin	$\text{C}_9\text{H}_{17}\text{NO}_2$	172.1332	Gabapentin-D ₁₀
Meprobamate	$\text{C}_9\text{H}_{18}\text{N}_2\text{O}_4$	219.1339	Meprobamate-D ₇
Norbuprenorphine	$\text{C}_{25}\text{H}_{35}\text{NO}_4$	414.2639	Norbuprenorphine-D ₃
Norfentanyl	$\text{C}_{14}\text{H}_{20}\text{N}_2\text{O}$	233.1648	Norfentanyl-D ₅
Pregabalin	$\text{C}_8\text{H}_{17}\text{NO}_2$	160.1332	Pregabalin-D ₆
Ritalinic Acid	$\text{C}_{13}\text{H}_{17}\text{NO}_2$	220.1332	Ritalinic Acid-D ₁₀
Tapentadol	$\text{C}_{14}\text{H}_{23}\text{NO}$	222.1852	Tapentadol-D ₃



Results and Discussion

This preliminary qualitative screening method was developed to detect the presence of 9 drug analytes in urine samples. The use of the deuterated internal standards (**Table 1**) resulted in a matrix factor very close to 1. The concentration range for all analytes was 20 – 1000 ng/mL, except for Fentanyl and Norfentanyl, where was 2 – 100 ng/mL. Precision and accuracy was assessed at five levels with 3 preps of 6 replicate samples at each concentration: QC Neg (1.2/12.5 ng/mL), QC Pos (2.5/25 ng/mL), QC Mid (300 ng/mL – data not shown), QC High (700 ng/mL), and QC Gluc for analytes Norbuprenorphine and Tapentadol (free 147 and 117 ng/mL, respectively) (**Table 2**).

Calibration curves (**Figure 1**), with a regression model of 1/x showed good linearity with R² values ≥ 0.999 for all analytes (**Table 3**). The specificity assessment (**Figure 2**) included evaluation of interference from matrix/system components and internal standards; all acceptance criteria were met. A small subset of urine samples presented interferences for Gabapentin, Norfentanyl, Pregabalin, Ritalinic Acid, and Norbuprenorphine. These interferences are suspected to be endogenous substances present in urine. This phenomena explains the discrepancies observed between this screening method and the LC-MS/MS method used for confirmation.

Assessments of carryover, dilution linearity, and impact from hematuria were performed as well as stability assessments including freeze/thaw, short-term, long-term, post-preparative, and autosampler stability. In all cases, results met acceptance criteria as per the 2001 Bioanalytical Guidance Document and current industry best practices.

Table 2. Inter-Assay Precision and Accuracy. Acceptance Criteria: ≤ 20 % Deviation from Nominal. (N=18)

Compound	QC Negative			QC Positive			QC High			QC Gluc		
	Nominal Conc. (ng/mL)	Mean % Nominal	%RSD	Nominal Conc. (ng/mL)	Mean % Nominal	%RSD	Nominal Conc. (ng/mL)	Mean % Nominal	%RSD	Nominal Conc. (ng/mL)	Mean % Nominal	%RSD
Carisoprodol	12.0	92.8	4.4	25.0	96.2	4.9	700	97.2	1.6	-	-	-
Fentanyl	1.20	100	3.0	2.50	96.2	5.6	70	98.7	2.3	-	-	-
Gabapentin	12.0	81.5	11.3	25.0	96.0	6.1	700	100	1.2	-	-	-
Meprobamate	12.0	90.0	10.0	25.0	94.0	10.8	700	98.0	2.2	-	-	-
Norbuprenorphine	12.0	93.8	4.0	25.0	99	4.3	700	98.7	1.3	147	98.9	1.4
Norfentanyl	1.20	101	1.8	2.50	93.6	7.7	70	90.3	6.0	-	-	-
Pregabalin	12.0	98.6	6.0	25.0	99.1	6.3	700	97.7	1.8	-	-	-
Ritalinic Acid	12.0	88.4	1.1	25.0	95.6	4.0	700	100	1.1	-	-	-
Tapentadol	12.0	90.8	0.4	25.0	96.7	3.0	700	98.1	1.2	117	106	0.8

Figure 1. Representative Calibration Curve.




Table 3. Summary of Correlation Coefficient Calibration Results.

Compound	Run 1	Run 2	Run 3
Carisoprodol	0.999	1.000	0.999
Fentanyl	0.999	0.999	0.999
Gabapentin	0.999	1.000	1.000
Meprobamate	1.000	1.000	0.999
Norbuprenorphine	1.000	0.999	0.999
Norfentanyl	1.000	0.999	0.999
Pregabalin	1.000	1.000	1.000
Ritalinic Acid	0.999	0.999	0.999
Tapentadol	0.999	1.000	0.999




Figure 2 . Selectivity/Specificity of 100 Samples Confirmed Positive (50) and Negative (50) for the Presence of 9 Analytes Using Tandem Quadrupole Definitive LC-MS/MS Methods and Qualitatively Compared to Q Exactive™ Screen Results.

Carisoprodol			Fentanyl			Gabapentin		
QE Screen	Positive	Negative	QE Screen	Positive	Negative	QE Screen	Positive	Negative
QQQ Positive	50	0	QQQ Positive	50	0	QQQ Positive	50	0
QQQ Negative	0	50	QQQ Negative	0	50	QQQ Negative	1	49

Meprobamate			Norbuprenorphine			Norfentanyl		
QE Screen	Positive	Negative	QE Screen	Positive	Negative	QE Screen	Positive	Negative
QQQ Positive	50	0	QQQ Positive	49	1	QQQ Positive	47	3
QQQ Negative	0	50	QQQ Negative	1	49	QQQ Negative	0	50

Pregabalin			Ritalinic Acid			Tapentadol		
QE Screen	Positive	Negative	QE Screen	Positive	Negative	QE Screen	Positive	Negative
QQQ Positive	48	2	QQQ Positive	50	0	QQQ Positive	50	0
QQQ Negative	0	50	QQQ Negative	2	48	QQQ Negative	0	50

Figure 3 . Representative Chromatograms at LLOQ (2 ng/mL or 20 ng/mL).



Conclusions

This method provides a cost effective, highly selective and specific alternative to EIA screening. The use of a purified recombinant Glucuronidase enzyme and the use of a Kinetex™ 2.6 µm EVO C18 LC column, in conjunction with the high resolution of the Q Exactive™ mass spectrometer, allows for a fast Dilute-and-Shoot sample preparation for 9 drug panel analytes.

Ordering Information

Phases	Kinetex 2.6 µm Minibore Columns (mm)						SecurityGuard™ ULTRA Cartridges*
	20 x 2.1	30 x 2.1	50 x 2.1	75 x 2.1	100 x 2.1	150 x 2.1	3/pk
EVO C18	00M-4725-AN	00A-4725-AN	00B-4725-AN	—	00D-4725-AN	00F-4725-AN	AJ0-9298
PS C18	—	00A-4780-AN	00B-4780-AN	—	00D-4780-AN	00F-4780-AN	AJ0-8951
Polar C18	—	00A-4759-AN	00B-4759-AN	—	00D-4759-AN	00F-4759-AN	AJ0-9532
Biphenyl	00M-4622-AN	00A-4622-AN	00B-4622-AN	—	00D-4622-AN	00F-4622-AN	AJ0-9209
XB-C18	—	00A-4496-AN	00B-4496-AN	00C-4496-AN	00D-4496-AN	00F-4496-AN	AJ0-8782
C18	00M-4462-AN	00A-4462-AN	00B-4462-AN	00C-4462-AN	00D-4462-AN	00F-4462-AN	AJ0-8782
C8	—	00A-4497-AN	00B-4497-AN	00C-4497-AN	00D-4497-AN	00F-4497-AN	AJ0-8784
HILIC	—	00A-4461-AN	00B-4461-AN	00C-4461-AN	00D-4461-AN	00F-4461-AN	AJ0-8786
Phenyl-Hexyl	—	00A-4495-AN	00B-4495-AN	00C-4495-AN	00D-4495-AN	00F-4495-AN	AJ0-8788
F5	—	00A-4723-AN	00B-4723-AN	—	00D-4723-AN	00F-4723-AN	AJ0-9322

*SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#)

for 2.1 mm ID



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5019 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 6559 4364
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

🌐 **All other countries/regions
Corporate Office USA**

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

**BE-HAPPY™
GUARANTEE**

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sale.

Trademarks

Kinetex, SecurityGuard, and BE-HAPPY are trademarks of Phenomenex. Waters and ACQUITY are registered trademarks of Waters Technologies Corporation. Q Exactive and Orbitrap are trademarks of Thermo Fisher Scientific.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is in now way affiliated with Waters Technologies Corporation or Thermo Fisher Scientific.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2023 Phenomenex, Inc. All rights reserved.

