


**NEW** Push  
to Filter

Ready  
to Inject






# Verex Filter Vials Two-Step Vials for Filtration and Analysis

## Table of Contents

The Most Innovative 4-in-1 Technology .....	<b>3</b>
Sample Preparation in Two Steps .....	<b>4</b>
Filter Vial Selection Chart .....	<b>5</b>
Membrane Options for Your Unique Sample .....	<b>6</b>
Applications: Filtration of Samples .....	<b>7-9</b>
Filter Vial Chemical Compatibility .....	<b>10</b>
Ordering Information .....	<b>11</b>

# The Most Innovative 4-in-1 Technology

Verex Filter Vials combines syringe filter and vial technology, eliminating the need for separate syringes, syringe filters, vials, and cap/septa, allowing you to reduce lab waste and simplify your workflow.



## Reduce Sample Loss and Contamination

Eliminate multiple transfers with this all-in-one filtration device

## Internal Plunger

## External Vial

## Particulate-Free Sample

The filter membrane is attached to the internal plunger to ensure removal of particulates

## Save Your Sample

With a low dead-volume of 30  $\mu$ L, more of your sample can be analyzed


### Verex Filter Vial Specifications

- Dimensions: 12 x 32 mm
- Vial material: Polypropylene
- Septa: PTFE/Silicone preSlit
- Filtering capacity: 450  $\mu$ L
- Dead-volume: 30  $\mu$ L
- Cap: 11 mm snap-top cap

# Sample Preparation in Two Steps


Verex Filter Vials are an easy two-step sample preparation device that consists of two parts: an external vial to be filled with sample, and an internal plunger with a filtration membrane and cap with a pre-slit septa.

## Filter Vials



Simply dispense your sample and filter!

## Traditional Syringe Filtration





Compared to the traditional syringe filtration, Verex Filter Vials requires less steps and saves you time!

# Find Your Filter Vial with Two EASY Selections

Verex Filter Vials are offered in a variety of chemically compatible membranes suitable for any application. Proper membrane and pore size selection are the keys to choosing the best product to maintain the integrity of your sample components as well as to protect your system from particulate contamination.


## 1 Select Your HPLC/UHPLC Column Particle Size

$\geq 3\ \mu\text{m}$	$< 3\ \mu\text{m}$
0.45 $\mu\text{m}$ 	0.20 $\mu\text{m}$ 



## 2 Choose Your Sample Type

Aqueous			Solvents	
Biological Samples / Protein Analysis	Tissue, Media, Buffers	Solvent Mixtures	Aqueous Mixtures, Hydrophilic	Non-Aqueous, Hydrophobic / Strong Acids



# Membrane Options for Your Unique Sample



**RC (Regenerated Cellulose)** **Hydrophilic**

**Compound Classes:** Aqueous and organic solutions

**Benefits:** Compatible with a broad range of solvents; fast-flow and ultra-low protein, and non-specific binding characteristics

**Typical Applications:** Broadly recommended as an excellent general purpose/high-performance samples filter for most applications; filtration of aqueous and organic solutions; protein chemistry; and clarification



**PTFE (Polytetrafluoroethylene)** **Hydrophobic**

**Compound Classes:** Organic solvents, acids, alcohols, bases, aromatics

**Benefits:** Compatible with organic solvents, strong acids and bases; chemically and biologically inert; gases or aggressive organic solvents

**Typical Applications:** Filtration of organic-based, highly acidic or basic samples and solvents; drug metabolite studies; clarification of aqueous and organic solvent solutions



**NY (Nylon)** **Hydrophilic**

**Compound Classes:** Bases, solvents, alcohols, aromatic hydrocarbons

**Benefits:** Compatible with organic solvents; not recommended when maximum protein recovery is important

**Typical Applications:** Clarification of aqueous and organic solvent solutions; general filtration or medical assays



**PES (Polyethersulfone)** **Hydrophilic**

**Compound Classes:** Critical biological samples, tissue culture media, additives, and buffers

**Benefits:** Fast-flow and ultra-low protein binding characteristics; removal of particulate matter

**Typical Applications:** Biological studies; dissolution testing; ICP sample preparation



**PVDF (Polyvinylidene Fluoride)** **Hydrophilic**


**Compound Classes:** Alcohols, biomolecules

**Benefits:** Binds less protein than Nylon or PTFE membranes; high flow rates and throughput; low extractables

**Typical Applications:** Biological studies; clarification studies; dissolution testing

# Filtration of Samples

## Monoclonal Antibody (mAbs) Fragments



### HPLC Conditions

**Column:** Synergi™ 4µm Hydro-RP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** [00F-4375-E0](#)  
**Pressure (bar):** 142  
**Mobile Phase:** A: 0.1 % TFA in Water  
 B: 0.1 % TFA in Acetonitrile

Gradient:	Time (min)	% B
	0	25
	10	45
	11	80
	13	80
	13.1	25
	17.1	25


**Flow Rate:** 0.8 mL/min  
**Detector:** UV-Vis @ 280 nm  
**Sample:** Nivolumab, IdeZ Digested

### Sample Preparation

**Reconstitute:** FabRICATOR® Z in 50 µL ddH<sub>2</sub>O to a concentration of 40 units/µL  
**Add:** 1 unit FabRICATOR Z / 1 µg IgG. Final Concentration should be 0.5-10 mg/mL  
**Digest:** Incubate mixture for 16 hrs at 37 °C  
**Load:** Sample into **Verex Filter Vial 0.45 µm, PES** ([ARO-F208-12](#))  
**Inject:** 2 µL final eluate onto HPLC-UV

Minimize sample loss  
and keep consistent  
method performance

## Separation of Chiral Atorvastatin



### LC-UV Conditions

**Column:** Lux® 5µm Cellulose-2  
**Dimension:** 100 x 4.6 mm  
**Part No.:** [00D-4457-E0](#)  
**Pressure (bar):** 108  
**Mobile Phase:** A: 0.1 % TFA in Water  
 B: 0.1 % TFA in Acetonitrile

Gradient:	Time (min)	% B
	25	35

**Flow Rate:** 1000 µL/min  
**Temperature:** Ambient


**Detection:** UV @ 210 nm  
**Injection Volume:** 3 µL  
**Instrument:** Agilent® 1100 HPLC with Quaternary Pump  
**Sample:** 1. Atorvastatin Calcium  
 2. Impurity E

### Sample Preparation

**Make:** 1 mg/mL concentrations of each Atorvastatin Calcium Trihydrate CRS and Atorvastatin Impurity E using methanol  
**Mix:** Standards at a ratio of 1:1 or 5:1, CRS:Impurity E  
**Load:** 400 µL of standard mixes into **Verex Filter Vial 0.45 µm, NY** ([ARO-F207-12](#))  
**Inject:** 3 µL final eluate onto HPLC-UV

# Filtration of Samples

## Multivitamin Tablet Analysis by HPLC-UV



### LC-UV Conditions

**Column:** Synergi™ 4 µm Hydro-RP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** [00F-4375-E0](#)  
**Pressure (bar):** 142  
**Mobile Phase:** A: 20 mM Potassium phosphate + 0.1% Hexane sulfonate, pH 3.0  
 B: Acetonitrile  
**Gradient:**

Time (min)	% B
0	3
3	3
18	50
18.1	3
20	3


**Flow Rate:** 1.5 mL/min  
**Temperature:** 30 °C  
**Detection:** UV @ 210 nm  
**Injection Volume:** 3 µL  
**Instrument:** Agilent® 1100 HPLC with Quaternary Pump  
**Sample:** 1. Thiamine  
 2. Pyridoxine  
 3. Pantothenic Acid  
 4. P-Aminobenzoic Acid  
 5. Riboflavin

### Sample Preparation

**Make:** 1 mg/mL concentrations of each sample with water from the tablet matrix. For Riboflavin, add 1-2 drops of Ammonium Hydroxide to help dissolve  
**Mix:** 50 µL of each sample into a mastermix  
**Load:** 250 µL of mastermix into **Verex Filter Vial 0.45 µm, NY (ARO-F207-12)**  
**Inject:** 3 µL final eluate onto HPLC-UV

No change in peaks with complex tablet and cream formulations!

## Topical Cream Analysis of a Hydrocortisone Extraction by HPLC-UV



### LC-UV Conditions

**Column:** Luna® 5 µm C18  
**Dimension:** 150 x 4.6 µm  
**Part No.:** [00F-4252-E0](#)  
**Pressure (bar):** 284  
**Mobile Phase:** Methanol:Water (50:50)  
**Flow Rate:** 1 mL/min  
**Temperature:** 35 °C  
**Detection:** UV @ 220 nm  
**Injection Volume:** 10 µL  
**Instrument:** Agilent® 1290 with DAD  
**Sample:** 1. Hydrocortisone


### Sample Preparation

**Pre-treatment:** Dissolve 1 g of hydrocortisone cream in 10 mL ethyl acetate  
**Dilute:** 2 mL of pre-treated sample with 8 mL methanol:water (1:1)  
**Load:** 400 µL of diluted cream sample into the **Verex Filter Vial 0.45 µm, RC (ARO-F203-12)**  
**Inject:** 10 µL final eluate onto HPLC-UV



# Filtration of Samples

## Separation of 23 PCBs by GC-MS



### GC-MS Conditions

**Column:** Zebtron™ ZB-Dioxin  
**Dimension:** 60 meter x 0.25 mm x 0.20 µm  
**Part No.:** [7KG-G045-10](#)  
**Injection:** Splitless @ 280 °C, 1 µL  
**Recommended Liner:** Zebtron PLUS Single Taper Z-Liner™  
**Liner Part No.:** [AG2-0A13-05](#) (for Agilent systems)  
**Carrier Gas:** Helium @ 2.0 mL/min (constant flow)  
**Oven Program:** 100 °C for 2 mins, 300 °C @ 6 °C/min for 5 min  
**Detector:** GC-MS  
**Sample:** 23 PCBs.  
 Find the full sample list online at [www.phenomenex.com/AN1005](http://www.phenomenex.com/AN1005)

### Sample Preparation

**Dilute:** 1 mL or 0.5 mL (Verex) of 0.4 µg/mL of PCB in Isooctane  
**Load:** Diluted sample into **Verex Filter Vial** 0.45 µm, PTFE ([ARO-F202-12](#))  
**Inject:** 1 µL filter eluate onto GC-MS

Particulate-free samples ensures no changes in chromatography

## See the Ease of Use Firsthand!



# Filter Vial Chemical Compatibility

Chemical	Filter Media					Housing	
	Regenerated Cellulose	Polytetrafluoroethylene	Nylon	Polyether-sulfone	Polyvinylidene Fluoride	Polypropylene	
	RC	PTFE	NY	PES	PVDF	PP	
<b>ACIDS</b>							
Acetic Acid, 5 %	R	R	R	R	R	R	
Acetic Acid, 10 %	R	R	L	R	R	R	
Acetic Acid, Glacial	R	R	N	R	R	L	
Boric Acid	T	R	L	T	T	R	
Hydrochloric, 6N	N	R	N	R	L	T	
Hydrochloric, Conc.	N	R	N	R	R	T	
Hydrofluoric, 10 %	L	R	N	T	R	R	
Hydrofluoric, 35 %	N	T	N	T	R	T	
Nitric Acid, 6N	N	L	N	N	T	T	
Nitric Acid, Conc.	N	N	N	N	R	T	
Sulfuric Acid, 6N	L	L	N	T	R	T	
Sulfuric Acid, Conc.	N	N	N	N	T	T	
<b>ALCOHOLS</b>							
Amyl Alcohol	R	R	R	N	R	R	
Benzyl Alcohol	R	R	L	N	R	R	
Butyl Alcohol	T	R	R	R	R	R	
Butyl Cellosolve	T	R	R	T	T	T	
Ethyl Alcohol	T	R	R	R	R	T	
Ethylene Glycol	R	R	R	R	R	R	
Glycerin	R	R	R	R	R	R	
Isobutyl Alcohol	T	R	R	T	R	T	
Isopropanol	R	R	R	R	R	T	
Methanol	R	R	T	R	R	T	
Methyl Cellosolve	T	R	R	T	R	T	
Propanol	R	R	R	T	R	R	
<b>BASES</b>							
Ammonium Hydroxide, 6N	L	R	N	R	R	T	
Potassium Hydroxide, 6N	L	R	R	T	R	T	
Sodium Hydroxide, 6N	L	R	N	R	R	T	
<b>SOLVENTS</b>							
Acetone	R	R	R	N	N	R	
Acetonitrile	R	R	T	R	R	R	
Amyl Acetate	R	R	R	L	R	L	
Aniline	R	R	R	R	T	L	
Benzene	R	L	T	R	R	L	
Bromoform	T	R	R	T	T	T	
Butyl Acetate	R	R	R	L	T	L	
Carbon Tetrachloride	R	L	R	R	R	N	
Cellosolve	R	R	R	T	T	T	
Chloroform	R	L	NR	N	R	L	
Cyclohexane	R	R	R	T	T	R	
Cyclohexanone	R	R	T	N	N	R	
Diethyl Acetamide	R	N	R	T	T	T	
Dimethyl Formamide	L	R	R	N	N	R	
Dimethyl Sulfoxide (DMSO)	R	R	R	N	N	T	
Dioxane	R	R	R	L	R	R	
Ethyl Ether	R	R	R	R	R	N	
Ethylene Dichloride	T	R	R	T	T	T	
Formaldehyde	T	R	R	R	R	R	
Freon TF	T	R	R	R	R	T	
Gasoline	R	R	R	T	R	N	
Hexane	R	R	R	T	R	T	
Isopropyl Acetate	R	R	R	T	N	R	
Kerosene	R	R	R	T	R	T	
Methyl Acetate	R	R	R	T	R	R	
Methyl Ethyl Ketone (MEK)	R	R	R	N	NR	T	
Methyl Isobutyl Ketone	R	R	R	T	N	T	
Methylene Chloride	NR	R	T	N	R	N	
Nitrobenzene	NR	R	T	N	R	R	
Pentane	NR	L	R	R	R	T	
Perchloroethylene	R	R	R	N	T	L	
Pyridine	R	R	T	N	N	L	
Tetrahydrofuran	R	L	T	N	N	L	
Toluene	R	L	R	N	R	L	
Trichloroethane	NR	R	T	L	T	T	
Trichlorethylene	R	L	T	R	R	N	
Triethylamine	R	R	R	T	T	T	
Xylene	R	L	T	L	R	R	
<b>MISCELLANEOUS</b>							
Cottonseed Oil	T	R	R	T	T	R	
Hydrogen Peroxide (30 %)	R	R	R	T	R	R	
Kodak KMER FTFR	T	R	R	T	T	T	
Peanut Oil	T	R	R	T	T	T	
Petroleum Oils	R	T	T	L	R	R	
Sesame Oil	T	R	R	T	T	T	
Shipley (AS-111,340,1350)	T	R	R	T	T	T	
Silicone Oils	R	R	R	R	R	R	
Turpentine	T	R	R	T	T	T	
Waycoat 59	T	R	R	T	T	T	











**Key**

- R** Recommended
- L** Limited Resistance (testing before use is recommended)
- N** Not Recommended
- T** Testing Recommended
- NR** Not Resistant

# Verex Filter Vials

## Ordering Information




Description		Pore Size	Part No.	Unit
<b>Verex Filter Vial-RC</b> (Regenerated Cellulose)		0.20 µm	<a href="#">ARO-F103-12</a>	100/pk
		0.45 µm	<a href="#">ARO-F203-12</a>	100/pk
<b>Verex Filter Vial-PTFE</b> (Polytetrafluoroethylene)		0.20 µm	<a href="#">ARO-F102-12</a>	100/pk
		0.45 µm	<a href="#">ARO-F202-12</a>	100/pk
<b>Verex Filter Vial-NY</b> (Nylon)		0.20 µm	<a href="#">ARO-F107-12</a>	100/pk
		0.45 µm	<a href="#">ARO-F207-12</a>	100/pk
<b>Verex Filter Vial-PES</b> (Polyethersulfone)		0.20 µm	<a href="#">ARO-F108-12</a>	100/pk
		0.45 µm	<a href="#">ARO-F208-12</a>	100/pk
<b>Verex Filter Vial-PVDF</b> (Polyvinylidene Fluoride)		0.20 µm	<a href="#">ARO-F106-12</a>	100/pk
		0.45 µm	<a href="#">ARO-F206-12</a>	100/pk

**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)



**NEW** Push to Filter  
Ready to Inject

**Australia**

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 5010 9707  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 800-852-3944  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thainfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
info@phenomenex.com

🌐 **All other countries/regions**  
**Corporate Office USA**

t: +1 (310) 212-0555  
info@phenomenex.com



**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Verex, Synergi, Zebron, Z-Liner, and BE-HAPPY are trademarks and Luna and Lux are registered trademarks of Phenomenex. SCIEIX and QTRAP are registered trademarks of AB SCIEIX Pte. Ltd. Agilent is a registered trademark of Agilent Technologies, Inc. FABRICATOR is a registered trademark of Genovis, Inc.

**Disclaimer**

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Agilent Technologies, Inc. or Genovis, Inc.

**FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.**

© 2021 Phenomenex, Inc. All rights reserved.