

STEP 2: Protocols

Strata-X-Drug B

	1	2	3
	Opiates, 6-MAM, PCP, Amphetamines, Methadone, Healthcare Opiates, and Propoxyphene*	Marijuana Metabolites	Cocaine Metabolites
Condition	Not Required		
Load	Pre-treated urine sample	Pre-treated urine sample	Pre-treated urine sample
Wash 1	600 µL of 100 mM Sodium acetate buffer (pH 5.0)	600 µL of 100 mM Sodium acetate buffer (pH 5.0)	600 µL of 0.1 N Hydrochloric acid
Wash 2	600 µL Methanol	600 µL of Acetonitrile/100 mM Sodium acetate buffer (pH 5.0) (30:70)	600 µL Methanol
Dry	10 minutes under full vacuum	15 minutes under full vacuum	10 minutes under full vacuum
Elute	2x 300 µL of Ethyl acetate/ Isopropanol/ Ammonium hydroxide (70:20:10)	2x 300 µL of Ethyl acetate/Isopropanol (85:15)	2x 300 µL of Ethyl acetate/Isopropanol/ Ammonium hydroxide (70:20:10)

* Opiates, 6-MAM, PCP, Amphetamines, Methadone, Healthcare Opiates, and Propoxyphene can be extracted simultaneously or separately using the same SPE methodology.

Methods are written for 30 mg/well Strata-X-Drug B 96-Well Plate; however they can be scaled to accommodate smaller or larger sample sizes and sorbent masses.

Strata-X-Drug N

	1	2
	Barbiturates	Benzodiazepines
Condition	Not Required	
Load	Pre-treated urine sample	Pre-treated urine sample
Wash 1	600 µL of 0.1 N Hydrochloric acid (HCl)	600 µL of Acetonitrile/Water (20:80)
Wash 2	2x 600 µL of Methanol/ 0.1 N HCl (30:70)	—
Dry	10 minutes under full vacuum	10 minutes under full vacuum
Elute	2x 300 µL of Ethyl acetate/ Isopropanol (85:15)	2x 300 µL of Ethyl acetate/ Isopropanol (85:15)

Methods are written for 30 mg/well Strata-X-Drug N 96-Well Plate; however they can be scaled to accommodate smaller or larger sample sizes and sorbent masses.

Round Well Collection Plates (polypropylene)

Part No.	Well Bottom	Well Volume	Unit
AH0-7279	Round	1 mL	50/pk
AH0-8636	Round	2 mL	50/pk

Square Well Collection Plates (polypropylene)

Part No.	Well Bottom	Well Volume	Unit
AH0-7192	Conical	350 µL	50/pk
AH0-7193	Conical	1 mL	50/pk
AH0-7194	Conical	2 mL	50/pk
AH0-8635	Round-Conical	2 mL	50/pk

Round Well Sealing Mats

Part No.	Description	Material	Unit
AH0-8631	Pierceable, 7 mm diameter	Silicone	50/pk
AH0-8632	Pre-Slit, 7 mm diameter	Silicone	50/pk
AH0-8633	Pierceable, 8 mm diameter	Silicone	50/pk
AH0-8634	Pre-Slit, 8 mm diameter	Silicone	50/pk
AH0-7362	Sealing Tap Pad	—	10/pk

Square Well Sealing Mats


Part No.	Description	Material	Unit
AH0-8597	Pierceable	Silicone	50/pk
AH0-8598	Pre-Slit	Silicone	50/pk
AH0-8199	Pierceable	Santoprene™	100/pk
AH0-7195	Pierceable	Ethylene Vinyl Acetate (EVA)	50/pk
AH0-7362	Sealing Tap Pad	—	10/pk

For additional resources, visit www.phenomenex.com/strataxdrug


STEP 3: Analyze Extracts

Kinetex® UHPLC/HPLC Core Shell LC Columns


Increased resolution, efficiency, and greater sensitivity.



Kinetex EVO C18
pH 1-12 stable C18 that delivers robust methods and improved peak shape for basic analytes.



Kinetex Biphenyl
100% aqueous stable reversed phase chemistry with hydrophobic and aromatic enhanced polar and aromatic selectivity.



Kinetex Polar C18
Combined C18 and polar modified surface that provide polar and non-polar retention alongside 100% aqueous stability.

Go to www.phenomenex.com/Kinetex to learn more.

Zebron™ GC Columns

Zebron ZB-Drug 1

- Optimized separation for drugs of abuse
- Provides fast analysis with better peak shape
- Specifically deactivated to improve quantitation for drug compounds

Zebron ZB-5PLUS™

- Highly inert for improved peak shape of drugs of abuse
- MS certified low bleed levels provide maximum sensitivity
- Traditional bonding chemistry provides the same selectivity as other 5% phenyl columns

To find a Zebron product, go to www.phenomenex.com/Zebron



Trademarks

Strata and Kinetex are registered trademarks and Zebron and 5PLUS are trademarks of Phenomenex. Santoprene is a trademark of Exxon Mobile Corporation. IMCSzyme is a registered trademark of Integrated Micro-Chromatography Systems, LLC.

Strata-X is patented by Phenomenex. U.S. Patent No. 7,119,145

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2019 Phenomenex, Inc. All rights reserved.

strata[®]X-Drug

THANK YOU
for choosing Strata-X-Drug
Solid Phase Extraction (SPE) Products



Phenomenex products are available worldwide.

www.phenomenex.com/mysupport

QUALITY
MANAGEMENT SYSTEM
CERTIFIED BY DNV GL
■ ISO 9001:2015 ■

IL86090319_WW

strata[®]X-Drug

General Methods

phenomenex[®]
...breaking with tradition[™]

If using Strata-X-Drug B Plus, perform hydrolysis directly in-well and turn on vacuum to initiate loading of sample. Verify that pH is between 4.0 and 6.0.

*For analysis of THC-OH and THC, use β -glucuronidase enzyme protocol.

PCP	To each 2 mL urine sample, add 2 mL of 100 mM Sodium acetate buffer (pH 5.0) then vortex. Verify that pH is between 4.0 and 6.0.
Cocaine Metabolites	To each 2 mL urine sample, add 100 μ L of 100 mM Phosphate buffer (pH: 6.0) and 1000 μ L of 0.35 M Sodium periodate. Vortex and incubate at room temperature for 25 minutes. Verify that the pH is between 5.5 and 6.5.
Amphetamines	To each 2 mL urine sample, add 100 μ L of 1.8 N Potassium hydroxide. Vortex then incubate at 60 °C for 20 minutes. Cool and add ~450 μ L glacial acetic acid then vortex. Verify that the pH is between 4.0 and 6.0.
Marijuana	To each 2 mL urine sample, add 100 μ L of 1.8 N Potassium hydroxide. Vortex then incubate at 60 °C for 20 minutes. Cool and add ~450 μ L glacial acetic acid then vortex. Verify that the pH is between 4.0 and 6.0.
Metabolites*	To each 2 mL urine sample, add 100 μ L of 1.8 N Potassium hydroxide. Vortex then incubate at 60 °C for 20 minutes. Cool and add ~450 μ L glacial acetic acid then vortex. Verify that the pH is between 4.0 and 6.0.

General Methods Inside

Enzyme:	1. Aliquot 200 μ L urine 2. Add 100 μ L 0.1 M Ammonium acetate buffer (pH 4, adjusted) 3. Add 100 μ L Campbell Science β -Glucuronidase Enzyme (Campbell Part No. DR2102) 4. Incubate 1 hour at 50 °C
IMCSzyme[®]:	1. Aliquot 140 μ L urine 2. Add 80 μ L IMCS buffer (IMCS Part No. 04-EZ-RHB-20) 3. Add 30 μ L IMCS Enzyme (IMCS Part No. 04-E1F-010) 4. Incubate 30 min at 60 °C
Kura Enzyme:	1. Aliquot 200 μ L urine 2. 100 μ L 0.1 M Ammonium acetate buffer (pH 5, adjusted) 3. 100 μ L Kura Enzyme Solution (Kura Part No. BG100-10ML) 4. Incubate 1 hr at 50 °C

Opiates, 6-MAM, Benzodiazepines, β -Glucuronidase Enzyme Hydrolysis

STEP 1: Sample Pre-treatments