



Ph. Eur. Monograph 2423: Valsartan Related Substances with Ph. Eur. Method Modernization

Krishna Chaitanya Routhu¹, Rajesh Babu Dandamudi, PhD¹, Swetha Kotikalapudi¹, Heiko Behr, PhD², and Bryan Tackett, PhD³
¹India Phenologix Lab, Phenomenex India, Hitech Defence and Aerospace Park Industrial Area, Mahadeva Kodigehalli, Holbi, Jala Taluka, Bengaluru 562149, India
²Phenomenex Ltd. Deutschland, Zeppelinstr. 5, 63741 Aschaffenburg, Germany
³Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

Valsartan is an antihypertensive drug which selectively inhibits angiotensin receptor type II. The European Pharmacopeia (Ph. Eur.) monograph method for the LC-UV determination of related compounds for Valsartan is based on an end-capped octadecylsilyl silica gel column run under isocratic mobile phase conditions. This application note demonstrates the potential method improvements that can be achieved within the allowed adjustments of chromatographic conditions. Different dimensions of a Luna™ Omega C18 column were used for evaluating the related substances method, within the new allowable adjustments published by European Pharmacopeia, Chapter 2.2.46, updated July 2022.

System suitability per Ph. Eur. Monograph 2423 for Valsartan Related Substances is a minimum resolution of 3.0 between the peaks due to impurity C and Valsartan. The results in the three methods show that the system suitability criteria were met. The use of a Luna Omega 3 µm C18, 100 x 3.0 mm column, a Luna Omega 5 µm C18, 100 x 4.6 mm column, and a Luna Omega 1.6 µm C18, 50 x 2.1 mm column is an allowed adjustment to the original column dimension with the flow rates scaled accordingly to accommodate the adjustments to column length (L), internal diameter (ID), and particle size (dp). With the Luna Omega 1.6 µm C18, 50 x 2.1 mm column, we demonstrated a reduction in total

analysis time by 75% (from 20 minutes to 5 minutes with a 0.6 mL/min flow) over the Luna Omega 3 µm C18, 100 x 3.0 mm column, and of 83.3% compared to the original column dimensions (5 µm, 125 x 3.0 mm) specified in the monograph.

All the reference solutions were prepared as indicated in Ph. Eur. monograph 2423 for Valsartan. The following certified reference standards (CRS) were purchased from the European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026 F - 67081 Strasbourg (France):

- Y0001132, Valsartan CRS
- Y0001145, Valsartan for System Suitability CRS

Figure 1. Valsartan Structure




Table 1. Ph. Eur. Adjustments of Chromatographic Conditions and Method Comparison.

| Method Parameter | Allowable Adjustments (Isocratic) | Method 1 | Method 2 | Method 3 |
|--|---|--|--|--|
| Stationary Phase | No change of the identity of the substituent permitted | As Specified | As Specified | As Specified |
| Column Length (L), Column Particle Size (dp), L/dp | The particle size and/or length of the column may be modified provided that the ratio of L/dp remains constant or in the range of -25 % to +50 % of the prescribed ratio. | Length = 100 mm Particle Size = 3 µm L/dp = 33.3 Deviation = +33.3 % (Allowed) | Length = 100 mm Particle Size = 5 µm L/dp = 20.0 Deviation = -20.0 % (Allowed) | Length = 50 mm Particle Size = 1.6 µm L/dp = 31.25 Deviation = +25.0 % (Allowed) |
| Column Internal Diameter (dc) | In the absence of a change in particle size and/or length of the column, the internal diameter of the column may be adjusted. | 3.0 mm (As Specified) | 4.6 mm (Allowed) | 2.1 mm (Allowed) |
| Flow Rate (F) | Flow rate is adjusted for changes in column diameter and particle size using the following equation: $F_2 = F_1 \times \left(\frac{dc_2^2 \times dp_1}{dc_1^2 \times dp_2} \right)$ <small>F₁ = flow rate indicated in the monograph, in mL/min F₂ = adjusted flow rate, in mL/min dc₁ = internal diameter of the column indicated in the monograph, in mm dc₂ = internal diameter of the column used, in mm dp₁ = particle size indicated in the monograph, in µm dp₂ = particle size of the column used, in µm</small> After an adjustment due to a change in column dimensions, an additional change in flow rate of ± 50 per cent is permitted. | 0.6 mL/min Deviation = -10.45 % (Allowed) | 0.9 mL/min Deviation = -4.26 % (Allowed) | 0.6 mL/min Deviation = -1.64 % (Allowed) |
| Column Temperature | ± 10 °C | As Specified | As Specified | As Specified |
| Composition of Mobile Phase | The amount of the minor solvent components may be adjusted by ± 30 % relative; no component is altered by more than 10 % absolute. | As Specified | As Specified | As Specified |
| Detection Wavelength | No adjustment permitted | As Specified | As Specified | As Specified |
| Injection Volume | When column dimensions are changed, it may be adjusted with the equation: $V_{inj2} = V_{inj1} \times \left(\frac{L_2 \times dc_2^2}{L_1 \times dc_1^2} \right)$ <small>V_{inj} = injection volume indicated in the monograph, in µL V_{adj} = adjusted injection volume, in µL dc₁ = internal diameter of the column indicated in the monograph, in mm dc₂ = internal diameter of the column used, in mm L₁ = column length indicated in the monograph, in mm L₂ = new column internal diameter, in mm</small> | 8 µL | 19 µL | 2 µL |



LC-UV Conditions

Columns: Luna™ Omega 3 µm C18, 100 x 3.0 mm (00D-4784-Y0) – Method 1
 Luna Omega 5 µm C18, 100 x 4.6 mm (00D-4785-E0) – Method 2
 Luna Omega 1.6 µm C18, 50 x 2.1 mm (00B-4742-AN) – Method 3

Mobile Phase: Mobile Phase (Table 1)

Flow Rate: 0.6 mL/min (Isocratic) – Method 1 and 3
 0.9 mL/min (Isocratic) – Method 2

Injection: 8 µL – Method 1
 19 µL – Method 2
 2 µL – Method 3

Temperature: 25 °C




Detector: UV @ 225 nm

System: Waters® ACQUITY Arc® HPLC – Method 1 and 2
 Waters ACQUITY® H-Class UHPLC – Method 3

Table 1. Preparation of Test and Reference Solutions

| Solution | Composition |
|-------------------------------|--|
| Mobile Phase | Glacial Acetic Acid / Acetonitrile / Water (1:500:500, v/v/v) |
| Test Solution | Dissolve 50 mg of Valsartan CRS in Mobile Phase and dilute to 100 mL with Mobile Phase. |
| Reference Solution (a) | Dilute 1.0 mL of Test Solution to 100.0 mL with Mobile Phase. Dilute 1 mL of the solution to 10 mL with Mobile Phase. |
| Reference Solution (b) | Dissolve the contents of a vial of Valsartan for System Suitability CRS (containing impurity C) in 1.0 mL of the Mobile Phase. |

Figure 2. System Suitability Test for Related Substances Using Reference Solution (b).



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5019 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 6559 4364
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

🌐 **All other countries/regions
Corporate Office USA**

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

**BE-HAPPY™
GUARANTEE**

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sales.

Trademarks

Luna and BE-HAPPY are trademarks of Phenomenex. Waters, ACQUITY, and ACQUITY Arc are registered trademarks of Waters Technologies Corporation.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Waters Technologies Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2023 Phenomenex, Inc. All rights reserved.

