

APPLICATION

Rapid and Sensitive Analysis of Nicotine and Metabolites Using a Kinetex® 2.6 µm EVO C18 Core-Shell Column

Laura Snow, Xianrong (Jenny) Wei, and Simon Lomas
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

With tobacco usage being linked to several respiratory and circulatory diseases and cancers, the major components and resulting metabolites are heavily researched. One of the main alkaloids in tobacco is nicotine, the stimulant responsible for tobacco's addictive properties. Nicotine along with its metabolites represent important biomarkers for monitoring the uptake of nicotine and related toxicants for tobacco use, long-term environmental tobacco smoke (ETS) exposure, and for optimizing nicotine dependency treatment. Anabasine is a tobacco alkaloid not found in nicotine replacement therapy products such as nicotine gum or patches, and can be used as a biomarker for continued tobacco use during treatment. This group of compounds has poor retention under acidic reversed phase conditions; however, they are neutralized under basic conditions thereby increasing hydrophobic retention. The ethane cross-linking of the Kinetex EVO C18 particles provides the necessary ruggedness and inertness needed for running this method at high pH. Additionally, the high efficiency Kinetex 2.6 µm EVO C18 column provides excellent resolution for HPLC and UHPLC instrumentation, greatly increased sensitivity, and a short chromatographic analysis which can increase lab productivity.


LC Conditions

Column: Kinetex 2.6 µm EVO C18
Dimensions: 100 x 3.0 mm
Part No.: 00D-4725-Y0
Mobile Phase: A: 20 mM Ammonium bicarbonate, pH 8.2
 B: Methanol
Gradient:

Time (min)	% B
0	10
3	90
5	90
5.01	10
6	10

Flow Rate: 0.75 mL/min
Temperature: 25 °C
Detection: MS/MS (SCIEX API 5000™)
Sample:

- 3-OH-Cotinine
- Normicotine
- Cotinine
- Anabasine
- Nicotine



APPLICATION

Ordering Information Kinetex[®] EVO C18 Core-Shell LC Columns

5 µm Minibore Columns (mm)				SecurityGuard [™] ULTRA Cartridges [†]
30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
00A-4633-AN	00B-4633-AN	00D-4633-AN	00F-4633-AN	AJO-9298 for 2.1 mm ID

5 µm MidBore [™] Columns (mm)			SecurityGuard [™] ULTRA Cartridges [†]
50 x 3.0	100 x 3.0	150 x 3.0	3/pk
00B-4633-YO	00D-4633-YO	00F-4633-YO	AJO-9297 for 3.0 mm ID

5 µm Analytical Columns (mm)				SecurityGuard [™] ULTRA Cartridges [†]
50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
00B-4633-E0	00D-4633-E0	00F-4633-E0	00G-4633-E0	AJO-9296 for 4.6 mm ID

2.6 µm Minibore Columns (mm)				SecurityGuard [™] ULTRA Cartridges [†]
30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
00A-4725-AN	00B-4725-AN	00D-4725-AN	00F-4725-AN	AJO-9298 for 2.1 mm ID

2.6 µm MidBore [™] Columns (mm)			SecurityGuard [™] ULTRA Cartridges [†]
50 x 3.0	100 x 3.0	150 x 3.0	3/pk
00B-4725-YO	00D-4725-YO	00F-4725-YO	AJO-9297 for 3.0 mm ID

2.6 µm Analytical Columns (mm)				SecurityGuard [™] ULTRA Cartridges [†]
50 x 4.6	100 x 4.6	150 x 4.6		3/pk
00B-4725-E0	00D-4725-E0	00F-4725-E0		AJO-9296 for 4.6 mm ID

1.7 µm Minibore Columns (mm)			SecurityGuard [™] ULTRA Cartridges [†]
50 x 2.1	100 x 2.1	150 x 2.1	3/pk
00B-4726-AN	00D-4726-AN	00F-4726-AN	AJO-9298 for 2.1 mm ID

[†] SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000

Australia
t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria
t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium
t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada
t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China
t: +86 (0)20 2282-6668
f: +86 (0)20 2809-8130
chinainfo@phenomenex.com

Denmark
t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland
t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France
t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany
t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India
t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland
t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy
t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg
t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico
t: 001-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands
t: +32 (0)9-4780951
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand
t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway
t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com


Puerto Rico
t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain
t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden
t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom
t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA
t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

**All other countries
Corporate Office USA** 
t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide.
For the distributor in your country, contact Phenomenex USA,
International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at <http://www.phenomenex.com/TermsAndConditions>.

Trademarks

Kinetex is a registered trademark and SecurityGuard and MidBore are trademarks of Phenomenex. QTRAP is a registered trademark and API 4000 is a trademark of AB SCIEX Pte. Ltd. AB SCIEX[™] is being used under license.

Kinetex EVO is patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

© 2015 Phenomenex, Inc. All rights reserved.



If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to other products of the same phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.