

# LC-MS/MS Separation of Furosemide and its Metabolite Using the Kinetex® 2.6 µm Biphenyl, Luna® Omega 1.6 µm C18, and Luna Omega 1.6 µm Polar C18 Columns

Roxana Eggleston-Rangel and Dr. Bryan Tackett  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

## Overview

Furosemide belongs to the diuretic drugs and can be administered either orally or by injection. It is mainly used for the treatment of edema. One of the metabolites formed from Furosemide is 4-chloro-5-sulfamoyl anthranilic acid (Saluamine) whose pharmacological activity remains unknown (**Figure 1**). This drives the interest in developing a reversed phase HPLC method able to separate and analyze Furosemide and its metabolites successfully and in a short amount of time. In this application note, we present an LC-MS/MS method to separate Furosemide and its metabolite by utilizing three different columns: Kinetex 2.6 µm Biphenyl, Luna Omega 1.6 µm C18, and Luna Omega 1.6 µm Polar C18.

The Kinetex Biphenyl column provided separation for Furosemide and its metabolite through a mixture of pi-pi and polar interactions; and the higher efficiency provided by the core-shell particle morphology also results in narrow peaks and increased MS sensitivity (**Figure 2a**).

The Luna Omega C18 column is a conventional C18 column with TMS endcapping, making this column a good choice for retention of analytes through hydrophobic interactions. The Luna Omega Polar C18 incorporates a polar surface treatment to offer increased retention for polar analytes, in addition to hydrophobic retention via the C18 bonded phase. The Luna Omega Polar C18 did provide a slight increase in retention for the more polar Saluamine (**Figure 2b, 2c**), while both Luna Omega columns showed an overall reduction in retention compared to the Kinetex Biphenyl column, which is likely due to the multiple interactions between the analytes and the Biphenyl stationary phase (**Figure 2a**).

## LC-MS/MS Conditions

**Column:** Kinetex 2.6 µm Biphenyl ([00B-4622-AN](#))  
Luna Omega 1.6 µm C18 ([00B-4742-AN](#))  
Luna Omega 1.6 µm Polar C18 ([00B-4748-AN](#))

**Dimension:** 50 x 2.1 mm

**Mobile Phase:** A: 0.1 % Formic Acid in Water  
B: See Chromatogram

Gradient: Time (min)	%B
0	5
4	95
5	95
5.1	5
7	5

**Flow Rate:** 400 µL/min

**Injection Volume:** 1 µL

**Temperature:** 40 °C

**LC System:** Agilent® 1200 Series


**Detection:** MS/MS

**Detector:** SCIEX® 4500


## MRM Transitions

Analyte	Q1 (m/z)	Q3 (m/z)
Furosemide-1	329	205
Furosemide-2	329	285
Saluamine-1	248.98	205
Saluamine-2	248.98	126


**Figure 1.** Furosemide and Saluamine




**Figure 2a.** Kinetex® 2.6 µm Biphenyl



**Figure 2b.** Luna® Omega 1.6 µm C18



**Figure 2c.** Luna Omega 1.6 µm Polar C18



**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/LiveChat](http://www.phenomenex.com/LiveChat).

**Australia**

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 5010 9707  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 800-852-3944  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
info@phenomenex.com

🌐 **All other countries/regions  
Corporate Office USA**

t: +1 (310) 212-0555  
info@phenomenex.com

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™  
GUARANTEE**

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Luna and Kinetex are registered trademarks and BE-HAPPY is a trademark of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc. SCIEX is a registered trademarks of AB SCIEX Pte. Ltd.

**Disclaimer**

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Agilent Technologies, Inc.

*FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.*

© 2021 Phenomenex, Inc. All rights reserved.

