



*D I S C O V E R*

*S E L E C T I V I T Y*

Pentafluorophenylpropyl Ligand

explore

**LUNA<sup>®</sup>**  
**PFP**



**phenomenex<sup>®</sup>**  
*...breaking with tradition<sup>SM</sup>*




# Discover New Selectivity

Resolution ( $R_s$ ) is the goal of every chromatographic method. It describes the separation power of the complete chromatographic system relative to the particular components of the mixture.

$$R_s = \frac{\sqrt{N}}{4} \left( \frac{\alpha - 1}{\alpha} \right) \left( \frac{k}{k + 1} \right)$$

Selectivity ( $\alpha$ ) describes the relative retention time ( $k_2 / k_1$ ) of two peaks. It is the most important variable in chromatographic performance.



Notice how even small changes in Selectivity can dramatically affect chromatographic performance.

In HPLC, Selectivity is achieved through **5** mechanisms of interaction

	Bond Type	Relative Strength
1	Ionic Interaction	1000
2	Hydrogen Bonding	100
3	Dipole-Dipole Interactions	10
4	Aromatic Interactions ( $\pi$ - $\pi$ )	5
5	Hydrophobic Interaction	1 Alkyl phases will use hydrophobic selectivity.

# Explore the Advantages of Luna PFP

- Multiple mechanisms for alternative selectivity
- Orthogonal selectivity (to traditional C18 phases) for more hits in your method screening
- Unique polar interactions for trace impurity identification

Luna PFP Selectivity is achieved through **4** of the **5** mechanisms of interaction

1	Hydrogen Bonding
2	Dipole-Dipole Interactions
3	Aromatic $\pi$ - $\pi$ Interactions
4	Hydrophobic

A typical alkyl phase (C18, C8) achieves selectivity through only 1 mechanism of interaction.

## TABLE OF CONTENTS


Discover New Selectivity .....	2-3
Explore the Mechanisms .....	4-5
Applications	
Complex Natural Products .....	6
Highly Polar Compounds .....	7
Conjugated Compounds .....	8
Isomeric Compounds .....	9
Aromatic Compounds .....	10-11
A Wide Variety of Compounds .....	12-13
Explore Luna	
Luna Legacy .....	14-15
Luna Chemistries .....	16-17
Ordering Information .....	18

# Explore the Mechanisms

## 1. Hydrogen bonding interactions

Hydrogen bonding is an extremely effective retention mechanism for polar compounds. The presence of hydrogen atoms bonded to heteroatoms (i.e., oxygen, nitrogen, or sulfur) produces polar groups that exhibit a significant difference in electronegativity between the heteroatom (electron rich) and hydrogen atom (electron poor). These functional groups (i.e. -OH, -NH, or -SH) allow potential hydrogen bonding interactions to occur.



- Hydrogen atoms, when covalently bound in a solute molecule, as they are in protonated amines and hydroxyl groups, will interact electrostatically with the electronegative fluorine (F) atom on the PFP ligand
- Low pH mobile phases, which cause bases to be protonated, increase opportunities for hydrogen bonding interactions with fluorine (F) atoms



## 2. Dipole-Dipole interactions


The Luna PFP ligand can participate in both dipole-dipole interactions and induced dipole-dipole interactions, increasing the potential for solute/ stationary phase interactions. Dipole-dipole interactions act upon polar solutes in much the same way as ionic interaction, only the effect is much weaker due to the partial charge. Induced dipole-dipole interactions may actually accommodate a brief electrostatic interaction between the PFP ligands and neutral solutes.

- Polar interactions are created as the highly electronegative fluorine (F) atoms on the PFP ligands are in a constant partially negative charged state that will attract the partial positive charge on polar solutes
- Induced dipole-dipole can create interactions between the stationary phase and neutral solutes as the highly electronegative fluorine atom can induce an asymmetrical distribution of charge in the neutral solute and create a brief electrostatic interaction





Pentafluorophenylpropyl Ligand



### 3. Aromatic interactions/ $\pi$ - $\pi$ interactions

Many important organic compounds will contain conjugated groups such as aromatic rings or double / triple bonds that may be attracted to the Luna PFP ligand. The PFP ligand contains a benzene ring, an unsaturated aromatic hydrocarbon with alternating double and single bonds, which is highly prone to electrophilic interactions due to the delocalized electrons in p-orbitals above and below the planar ring.

- Solutes containing aromatic rings may participate in a stacking interaction occurring with the benzene ring of the PFP ligand.  $\pi$ - $\pi$  interactions caused by the overlapping p-orbitals of both rings create an attraction that may initiate the solute to arrange itself over the stationary phase group in a discriminating interaction
- Non-aromatic solutes containing double or triple bonds have p-orbital electrons that are ready to interact with the delocalized electron field of the PFP benzene ring. Compounds differing only in the presence or absence of double or triple bonds will receive dissimilar stationary phase interaction

### 4. Hydrophobic interactions

Hydrophobic interactions will occur with any carbon containing solute. These are 'reversed phase' type interactions that are difficult to quantify in terms of strength, but are the mainstay of many HPLC methods in practice. What is commonly referred to as 'hydrophobic bonding' is actually a consequence of water's attempt to maximize polar interactions through its own hydrogen bonding network.

- Hydrophobic interaction is not really bonding, but it is often referred to as such. In an aqueous rich HPLC mobile phase acting upon a hydrophobic stationary phase, water will continuously move, through exclusion, hydrophobic organic groups into that stationary phase. This action is more thermodynamically stable than forming a cavity around the hydrophobic solute
- The propyl linkage and the benzene ring of the PFP ligand both permit an organic rich layer on the silica surface that is a partitioning target for the aqueous excluded carbon solutes

DISCOVER SELECTIVITY

## Complex Natural Products

Paclitaxel (Taxol<sup>®</sup>) and taxane analogs have shown results as anticancer and antitumor drugs. With ever increasing demands for testing, more stable and reproducible columns for new testing methods are required.


New Luna PFP column addresses the need for better and more reliable taxane separations

Luna PFP is an improvement in the analysis of complex taxane analogs

### Taxanes on Luna 3 $\mu$ m PFP

App ID 16300

Column:	Luna 3 $\mu$ m PFP
Dimension:	150 x 4.6 mm
Part No.:	00F-4447-E0
Mobile Phase:	A: Water B: Acetonitrile
Gradient:	A/B (75:25) to (30:70) in 20 min
Flow Rate:	1 mL/min
Temperature:	22 °C
Detection:	UV @ 225 nm
Sample:	1. 10-Deacetylbaccatin III 2. Paccatin III 3. 10-Deacetyl-7-xylosyitaxol B 4. Taxinine M 5. 10-Deacetyl-7-xylosyitaxol 6. 10-Deacetyl-7-xylosyitaxol C 7. 10-Deacetyltaxol 8. 7-Xylosyitaxol 9. Cephalomannine 10. 10-Deacetyl-7-epitaxol 11. Paclitaxel 12. Taxol C 13. 7-Epitaxol



## DISCOVER SELECTIVITY


## Highly Polar Compounds

Compounds containing charged amine groups are highly polar and very difficult to be accurately quantitated on traditional reversed phase columns.

- Selectivity to detect and quantify low level nitrogen containing analytes containing primary, secondary and tertiary amine groups


## Antihistamines on Luna PFP at Low pH App ID 16284

**Column:** Luna 5 µm PFP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4448-EO  
**Mobile Phase:** A: 20mM Potassium phosphate, pH 2.5  
 B: Acetonitrile  
**Gradient:** A/B (90:10) to (20:80) in 15 min  
**Flow Rate:** 1 mL/min  
**Temperature:** 22 °C  
**Detection:** UV @ 220 nm  
**Sample:** 1. Malic acid  
 2. Doxylamine  
 3. Chlorpheniramine  
 4. Bromopheniramine  
 5. Diphenhydramine  
 6. Chlorphenoxamine  
 7. Loratadine



## Antihistamines on Synergi™ Polar-RP™ at Low pH App ID 16283

**Column:** Synergi 4 µm Polar-RP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4336-EO  
**Mobile Phase:** A: 20mM Potassium phosphate, pH 2.5  
 B: Acetonitrile  
**Gradient:** A/B (90:10) to (20:80) in 15 min  
**Flow Rate:** 1 mL/min  
**Temperature:** 22 °C  
**Detection:** UV @ 220 nm  
**Sample:** 1. Malic acid  
 2. Doxylamine  
 3. Chlorpheniramine  
 4. Bromopheniramine  
 5. Diphenhydramine  
 6. Chlorphenoxamine  
 7. Loratadine




Halogens can radically increase the polarity of a compound, thus decreasing typical retention characteristics. Luna PFP retains, discriminates, and separates halogens easily.

- Strong attraction for halogenated compounds


## Positional Isomers of Halogenated Phenols on Luna 3 µm PFP App ID 16296

**Column:** Luna 3 µm PFP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4447-EO  
**Mobile Phase:** A: 0.1 % Formic acid in Water  
 B: 0.1 % Formic acid in Acetonitrile  
**Gradient:** A/B (60:40) to (50:50) in 10 min  
**Flow Rate:** 1 mL/min  
**Temperature:** 22 °C  
**Detection:** UV @ 254 nm  
**Sample:** 1. 2,3-Dimethylphenol 6. 2,5-Dichlorophenol  
 2. 2,5-Dimethylphenol 7. 2,6-Dichlorophenol  
 3. 2,6-Dimethylphenol 8. 3,4-Dichlorophenol  
 4. 3,4-Dimethylphenol 9. 3,5-Dichlorophenol  
 5. 3,5-Dimethylphenol 10. 2,4-Dibromophenol



## Positional Isomers of Halogenated Phenols on Luna 3 µm C18(2) App ID 16297

**Column:** Luna 3 µm C18(2)  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4251-EO  
**Mobile Phase:** A: 0.1% Formic acid in Water  
 B: 0.1% Formic acid in Acetonitrile  
**Gradient:** A/B (60:40) to (50:50) in 10 min  
**Flow Rate:** 1 mL/min  
**Temperature:** 22 °C  
**Detection:** UV @ 254 nm  
**Sample:** 1. 2,3-Dimethylphenol 6. 2,5-Dichlorophenol  
 2. 2,5-Dimethylphenol 7. 2,6-Dichlorophenol  
 3. 2,6-Dimethylphenol 8. 3,4-Dichlorophenol  
 4. 3,4-Dimethylphenol 9. 3,5-Dichlorophenol  
 5. 3,5-Dimethylphenol 10. 2,4-Dibromophenol



# DISCOVER SELECTIVITY


## Conjugated Compounds

Compounds that may differ only in their level of conjugation are easily identified and quantified on Luna PFP by electrophilic interaction.

Enhanced selectivity for compounds with double or triple bonds


### Aflatoxins on Luna 3 µm PFP App ID 16356

**Column:** Luna 3 µm PFP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4447-E0  
**Mobile Phase:** 10 mM Ammonium Acetate, pH 4.5/ Acetonitrile (70:30)  
**Flow Rate:** 1 mL/min  
**Detection:** UV @ 220 nm  
**Sample:**  
 1. Aflatoxin B1  
 2. Aflatoxin B2  
 3. Aflatoxin G1  
 4. Aflatoxin G2




### Aflatoxins on Synergi Max-RP App ID 15523

**Column:** Synergi 4 µm Max-RP  
**Dimension:** 150 x 3.0 mm  
**Part No.:** 00F-4337-Y0  
**Mobile Phase:** 10 mM Ammonium Acetate, pH 4.5 /Acetonitrile (70:30)  
**Flow Rate:** 0.5 mL/min  
**Temperature:** 22 °C  
**Detection:** Mass Spectrometer (MS)  
**Sample:**  
 1. Aflatoxin G2  
 2. Aflatoxin B2  
 3. Aflatoxin G1  
 4. Aflatoxin B1



### Aflatoxins on Gemini® 5 µm C18 App ID 15524

**Column:** Gemini 5 µm C18  
**Dimension:** 150 x 3.0 mm  
**Part No.:** 00F-4435-Y0  
**Mobile Phase:** 10 mM Ammonium Acetate, pH 4.5 /Acetonitrile (70:30)  
**Flow Rate:** 0.5 mL/min  
**Detection:** Mass Spectrometer (MS)  
**Sample:**  
 1. Aflatoxin G2  
 2. Aflatoxin B2  
 3. Aflatoxin G1  
 4. Aflatoxin B1



# DISCOVER SELECTIVITY


## Isomeric Compounds

Positional changes on an analyte of interest may effect the compound's dipole moment. This change can be readily noticed by the highly electronegative **fluorine (F) atom** and other retention mechanisms.

Excellent choice for positional isomers due to the multiple retention characteristics of Luna PFP


### Positional Isomers of Methylacetophenone on Luna 3 $\mu$ m PFP App ID 16298

Column: Luna 3  $\mu$ m PFP  
 Dimension: 150 x 4.6 mm  
 Part No.: 00F-4447-E0  
 Mobile Phase: Water/ Methanol (50:50)  
 Flow Rate: 1 mL/min  
 Temperature: 22  $^{\circ}$ C  
 Detection: UV @ 254 nm  
 Sample: 1. o-Methylacetophenone  
 2. m-Methylacetophenone  
 3. p-Methylacetophenone



### Positional Isomers of Methylacetophenone on Luna 3 $\mu$ m C18(2) App ID 16299

Column: Luna 3  $\mu$ m C18 (2)  
 Dimension: 150 x 4.6 mm  
 Part No.: 00F-4251-E0  
 Mobile Phase: Water/ Methanol (50:50)  
 Flow Rate: 1 mL/min  
 Temperature: 22  $^{\circ}$ C  
 Detection: UV @ 254 nm  
 Column: 1. o-Methylacetophenone  
 2. m-Methylacetophenone  
 3. p-Methylacetophenone




# DISCOVER SELECTIVITY

## Aromatic Compounds

Aromatic compounds show unique retention characteristics on Luna PFP compared to traditional reversed phase columns.


### Beta-Blockers on Luna PFP at Low pH App ID 16286

**Column:** Luna 5 µm PFP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4448-E0  
**Mobile Phase:** A: 20 mM Potassium phosphate, pH 2.5  
 B: Acetonitrile  
**Gradient:** A/B (85:15) to (50:50) in 10 min  
**Flow Rate:** 1 mL/min  
**Temperature:** 22 °C  
**Detection:** UV @ 230 nm  
**Sample:** 1. Atenolol  
 2. Pindolol  
 3. Metoprolol  
 4. Labetalol  
 5. Propranolol  
 6. Alprenolol



### Beta-Blockers on Luna C18(2) at Low pH App ID 16287


**Column:** Luna 5 µm C18(2)  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4252-E0  
**Mobile Phase:** A: 20 mM Potassium phosphate, pH 2.5  
 B: Acetonitrile  
**Gradient:** A/B (85:15) to (50:50) in 10 min  
**Flow Rate:** 1 mL/min  
**Temperature:** 22 °C  
**Detection:** UV @ 230 nm  
**Sample:** 1. Atenolol  
 2. Pindolol  
 3. Metoprolol  
 4. Labetalol  
 5. Propranolol  
 6. Alprenolol



The presence of the aromatic benzene ring in the Luna PFP ligand increases the relative attraction between the stationary phase and aromatic analytes, leading to increased retention for these types of compounds


### Steroid Mix on Luna PFP App ID 16288

**Column:** Luna 5 µm PFP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4448-E0  
**Mobile Phase:** A: Water  
 B: Acetonitrile  
**Gradient:** A/B (60:40) to (30:70) in 8 min, hold for 2 min  
**Flow Rate:** 1 mL/min  
**Temperature:** 22 °C  
**Detection:** UV @ 254 nm  
**Sample:** 1. Prednisone  
 2. Prednisolone  
 3. Hydroxyprogesterone  
 4. Progesterone



### Steroid Mix on Luna C18(2) App ID 16289

**Column:** Luna 5 µm C18(2)  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4252-E0  
**Mobile Phase:** A: Water  
 B: Acetonitrile  
**Gradient:** A/B (60:40) to (30:70) in 8 min, hold for 2 min  
**Flow Rate:** 1 mL/min  
**Temperature:** 22 °C  
**Detection:** UV @ 254 nm  
**Sample:** 1. Prednisone  
 2. Prednisolone  
 3. Hydroxyprogesterone  
 4. Progesterone



## DISCOVER SELECTIVITY


## Aromatic Compounds (con't)

Closely related polyphenolic compounds are readily separated with Luna PFP

Catechins on Luna 3  $\mu$ m PFP


App ID 16294

Column: Luna 3  $\mu$ m PFP  
 Dimension: 150 x 4.6 mm  
 Part No.: 00F-4447-E0  
 Mobile Phase: A: 0.1 % Formic acid in Water  
 B: 0.1 % Formic acid in Acetonitrile  
 Gradient: A/B (80:20) to (55:45) in 10 min  
 Flow Rate: 1 mL/min  
 Temperature: 22 °C  
 Detection: UV @ 280 nm  
 Sample: 1. Gallic acid  
 2. Epigallo catechin  
 3. Catechin  
 4. Epicatechin  
 5. Epigallocatechin gallate  
 6. Gallocatechin gallate  
 7. Epicatechin gallate  
 8. Catechin gallate

Catechins on Luna 3  $\mu$ m C18(2)

App ID 16295

Column: Luna 3  $\mu$ m C18(2)  
 Dimension: 150 x 4.6 mm  
 Part No.: 00F-4251-E0  
 Mobile Phase: A: 0.1 % Formic acid in Water  
 B: 0.1 % Formic acid in Acetonitrile  
 Gradient: A/B (80:20) to (55:45) in 10 min  
 Flow Rate: 1 mL/min  
 Detection: UV @ 280 nm  
 Sample: 1. Gallic acid  
 2. Epigallo catechin  
 3. Catechin  
 4. Epicatechin  
 5. Epigallocatechin gallate  
 6. Gallocatechin gallate  
 7. Epicatechin gallate  
 8. Catechin gallate



# DISCOVER SELECTIVITY

## A Wide Variety of Compounds

Minor functional group changes near pi electrons can be readily exploited with Luna PFP.


- Aromatic acids, aromatic antibiotics, aromatic sulfurs, phenoxy herbicides, structural isomers

### Sulfa Drugs on Luna PFP at low pH

App ID 16285

**Column:** Luna 5 µm PFP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4448-E0  
**Mobile Phase:** A: 0.1 % Formic acid in Water  
 B: 0.1 % Formic acid in Acetonitrile  
**Gradient:** A/B (95:5) to (20:80) in 15 min  
**Flow Rate:** 1 mL/min  
**Temperature:** 22 °C  
**Detection:** UV @ 254 nm  
**Sample:**

1. Sulfaguanidine
2. Sulfanilamide
3. Sulfathiazole
4. Sulfamerazine
5. Sulfamethoxazole
6. Sulfaquinoxaline




### Analgesics on Luna 5 µm PFP

App ID 16303

**Column:** Luna 5 µm PFP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4448-E0  
**Mobile Phase:** 20 mM Potassium phosphate, pH 2.5/ Acetonitrile (50:50)  
**Flow Rate:** 1 mL/min  
**Temperature:** 22 °C  
**Detection:** UV @ 230 nm  
**Sample:**

1. Ethyl paraben
2. Naproxen
3. Indomethacin
4. Ibuprofen




### Herbicides on Luna 3 µm PFP

App ID 16305

**Column:** Luna 3 µm PFP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4447-E0  
**Mobile Phase:** A: 0.1 % Formic acid in Water  
 B: 0.1 % Formic acid in Acetonitrile  
**Gradient:** A/B (45:55) to (25:75) in 6 min  
**Flow Rate:** 1 mL/min  
**Temperature:** 22 °C  
**Detection:** UV @ 280 nm  
**Sample:**

1. 4-Chlorophenoxyacetic acid
2. 2,4-Dichlorophenoxyacetic acid
3. 2,4,5-Trichlorophenoxyacetic acid
4. 2-(2,4,5-Trichlorophenoxy) propionic acid




# DISCOVER SELECTIVITY

## A Wide Variety of Compounds (con't)

- The unique retention mechanisms associated with Luna PFP make it an ideal column to place on all column screens


### Cephalosporins on Luna 3 $\mu$ m PFP App ID 16304

**Column:** Luna 3  $\mu$ m PFP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4447-E0  
**Mobile Phase:** 20 mM Potassium phosphate, pH 2.5/ Acetonitrile (85:15)  
**Flow Rate:** 1 mL/min  
**Temperature:** 22  $^{\circ}$ C  
**Detection:** UV @ 254 nm  
**Sample:** 1. Cefadroxil  
 2. Cefaclor  
 3. Cephalixin  
 4. Cephadrine



### Geometric Isomers of Maleic Acid and Fumaric Acid on Luna 3 $\mu$ m PFP App ID 16302

**Column:** Luna 3  $\mu$ m PFP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4447-E0  
**Mobile Phase:** 20 mM Potassium phosphate, pH 2.5/ Acetonitrile (97:3)  
**Flow Rate:** 1 mL/min  
**Temperature:** 22  $^{\circ}$ C  
**Detection:** UV @ 210 nm  
**Sample:** 1. Maleic acid  
 2. Fumaric acid



# The Luna Legacy


## The Luna PFP Advantage

Luna PFP lives up to the exacting standards of quality and customer satisfaction that has been exemplified in the 10 years of the Luna brand.

- Ultra-pure metal free (99.99 % purity) Luna silica supplies a new level of performance and reproducibility
- Proprietary silica treatment and bonding ensures a robust column
- Chemistry bound to this base sorbent is as stable and reproducible as alkyl phases

Steroids		App ID 16290
Column:	Luna 3 $\mu$ m PFP	
Dimension:	150 x 4.6 mm	
Part No.:	00F-4447-E0	
Mobile Phase:	A: Water B: Acetonitrile	
Gradient:	A/B (70:30) to (35:65) in 25 min	
Flow Rate:	1.0 mL/min	
Injection:	2 $\mu$ L	
Detection:	UV @ 220 nm	
Sample:	1. Estriol 2. Hydrocortisone 3. Prednisone 4. Cortisone 5. Corticosterone 6. Estradiol 7. Cortisone acetate 8. 2,1-Hydroxyprogesterone 9. 17-Hydroxyprogesterone 10. Progesterone 11. Estrone	

## Batch-to-Batch Reproducibility




## The Trouble with other Fluorinated Columns

Traditional PFP phases have failed to maximize user confidence by not generating high accuracy of data, precision, reproducibility, and ruggedness of phase.

- Inconsistencies in the base silica can be amplified in the presence of the sensitive PFP stationary phase!
- Residual metal ions in the silica can increase the activity of residual silanols, affecting the characteristics of the bonded fluorine groups.

## The Luna Legacy (con't)


The ideal Luna surface and quality control ensure robustness and stability even when exposed to harsh TFA conditions of pH 2.0 that, in many other fluorinated phases, may result in considerable loss of stationary phase.



### pH Stability

**Column:** Luna 5 µm PFP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4448-E0  
**Mobile Phase:** A: 0.1 % TFA in Water, pH 2.0  
 B: 0.1 % TFA in Acetonitrile  
**Gradient:** 5 % B to 95 % B in 10 min  
 Re-equilibrate at 5 % B for 5 min  
**Flow Rate:** 1.0 mL/min  
**Temperature:** Room temperature  
**Detection:** UV @ 254 nm  
**Sample:** Reversed Phase 2 Test Mix  
 (Uracil, Acetophenone, Benzene, Toluene, Naphthalene)

**Challenge Tests:** Every cycle of 12 flushes  
**Mobile Phases:** A: Water  
 B: Acetonitrile  
**Gradient:** 40 % to 70 % B in 8 min  
 Hold at 70 % B for 3 min  
**Flow Rate:** 1.0 mL/min  
**Temperature:** Room temperature  
**Detection:** UV @ 254 nm  
**Sample:** Steroid Test Mix (Prednisone, Prednisolone, Hydroprogesterone, Progesterone)



### Low MS Bleed

**Column:** Luna 5 µm PFP  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4448-E0  
**Mobile Phase:** A: 0.1 % Formic acid in Water B: 0.1 % Formic acid in Acetonitrile  
**Gradient:** 5 to 95 % B in 10 min, hold for 2 min, eq. for 5 min  
**Flow Rate:** 1.0 mL/min MS splitting flow: 0.31 mL/min  
**Temperature:** —  
**Detection:** API 3000 MS/MS  
 TurbolonSpray® heater gas flow: 6000 cc/min  
 TurbolonSpray heater temperature: 425 °C, ESI+, Q1 scan




L U

# Luna PFP is a Primary Method Development Tool

For development work in any industry, screening with Luna PFP greatly increases the chance of quick method optimization

- Luna PFP may provide enhanced selectivity for difficult-to-separate compounds and alternate elution orders
- Successful separations will be robust throughout validation
- You may avoid ion pair reagents and complex mobile phase preparations
- Complementary phase can aid in identification, proof of purity and quantitation





# NA

## Luna Chemistries

LUNA PHASES	Particle Size (µm)	Pore Size (Å)	Surface Area (m <sup>2</sup> /g)	Carbon Load (%)	Bonded Phase Coverage (µmole/m <sup>2</sup> )	pH Stability
<b>Silica(2)</b>	3, 5, 10, 10-PREP, 15	100	400	—	—	2.0 - 7.5
<b>C5</b>	5, 10	100	440	12.5	7.85	1.5 - 10
<b>C8</b>	5, 10	100	440	14.75	5.50	1.5 - 10
<b>C8(2)</b>	3, 5, 10, 10-PREP, 15	100	400	13.5	5.50	1.5 - 10
<b>C18</b>	5, 10	100	440	19.0	3.00	1.5 - 10
<b>C18(2)</b>	2.5, 3, 5, 10, 10-PREP, 15	100	400	17.5	3.00	1.5 - 10
<b>CN</b>	3, 5, 10	100	400	7.0	3.80	1.5 - 7.0
<b>NH<sub>2</sub></b>	3, 5, 10	100	400	9.5	5.80	1.5 - 11
<b>Phenyl-Hexyl</b>	3, 5, 10, 10-PREP, 15	100	400	17.5	4.00	1.5 - 10
<b>SCX</b>	5, 10	100	400	Binding Capacity: 0.15 meq/g		2.0 - 7.0
<b>HILIC</b>	3, 5	200	200	5.7	4.30	1.5 - 8.0
<b>PFP</b>	3, 5	100	400	5.7	2.2	1.5 - 8.0

Luna and Gemini are registered trademarks of Phenomenex, Inc. SecurityGuard, Polar-RP and Synergi are trademarks of Phenomenex, Inc. Taxol is a registered trademark of Bristol-Myers Squibb Co. TurbolonSpray is a registered trademark of Applied Biosystems. © 2008 Phenomenex, Inc. All rights reserved. Subject to Phenomenex Standard Terms & Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions)

Description	Application	Description	Application
<b>Silica(2)</b> Unbonded silica	Non-Polar compounds	<b>NH<sub>2</sub></b> Possibly the most rugged and reproducible NH <sub>2</sub> phase for reversed or normal phase separations. Stable from pH 1.5 to 11.0 and under 100 % aqueous conditions.	Excellent for reversed-phase analysis of sugars, sugar alcohols and anionic compounds or for hydrogen bonding compounds under normal phase.
<b>C5</b> Reversed-phase chemistry that offers greater hydrolytic stability and an alternative selectivity over the classic C4 phase.	Highly hydrophobic compounds. Good alternative to C8 or C18 when less retention is desired.	<b>Phenyl-Hexyl</b> Novel (patented) phenyl phase which employs a hexyl linker, as opposed to the traditional propyl chain. Stable under all conditions as well as alternative selectivity to most phenyl phases.	Unique selectivity for very polar and aromatic compounds when used in methanol containing mobile phases.
<b>C8</b> Original C8 phase. Highly retentive; optimized for long lifetime at extremely high pH levels (pH 10).	A “workhorse” for extremely high pH applications (pH 10), or when a more retentive C8 is desired.	<b>SCX</b> Benzene sulfonic acid, strong cation exchange. Guaranteed to provide sharper peak shape and better resolution compared to traditional SCX columns.	Amine and polyamine containing compounds.
<b>C8(2)</b> Optimized for maximum efficiency, peak shape and resolution. Significantly improved performance over traditional C8 phases due to high surface coverage. Slightly lower carbon load than original Luna C8.	Great starting point for methods development of pharmaceuticals, nucleotides and polar compounds. Excellent for highly aqueous and LC/MS applications.	<b>HILIC NEW</b> Unique cross-linked diol offers polar selectivity under HILIC conditions.	Powerful retention and selectivity of polar compounds while offering increased sample throughput and higher MS sensitivity.
<b>C18</b> Original Luna C18 phase. Maximum retention, optimized for long lifetime at extremely high pH levels (pH 10).	A “workhorse” for extremely high pH applications (pH 10), or when maximum retention is desired.	<b>PFP NEW</b> Pentafluorophenyl with a propyl linkage provides multiple retention mechanisms unique to other reversed phase medias.	Unique selectivity for highly polar compounds, complex natural products, isomers and other closely related compounds.
<b>C18(2)</b> Optimized for maximum efficiency, excellent peak shape and resolution; slightly lower carbon load than original Luna C18.	A C18 phase for virtually all HPLC applications. Range of particle sizes offers amazing versatility for capillary-LC/MS, to prep and process scale applications.	<b>CN</b> One of the best performing CN phases available. Can be used in reversed- or normal-phase modes. Based on same 99.99 % pure silica as all Luna phases.	Excellent for improving the retention of polar compounds. Extremely rapid equilibration makes it ideal for rapid screening and gradient applications.

# Ordering Information



If Luna does not provide at least an equivalent separation as compared to a competing column of the same particle size, similar phase and dimensions, send in your comparative data within 45 days and keep the Luna column for FREE.

3 µm Minibore Columns (mm)				
	50 x 2.0	100 x 2.0	150 x 2.0	4 x 2.0*
Phases				/10pk
PFP	00B-4447-B0	00D-4447-B0	00F-4447-B0	AJO-8326
Silica(2)	00B-4162-B0	00D-4162-B0	00F-4162-B0	AJO-4347
C8(2)	00B-4248-B0	00D-4248-B0	00F-4248-B0	AJO-4289
C18(2)	00B-4251-B0	00D-4251-B0	00F-4251-B0	AJO-4286
CN	00B-4254-B0	00D-4254-B0	00F-4254-B0	AJO-4304
Phenyl-Hexyl	00B-4256-B0	00D-4256-B0	00F-4256-B0	AJO-4350
NH <sub>2</sub>	00B-4377-B0	00D-4377-B0	00F-4377-B0	AJO-4301
HILIC	00B-4449-B0	00D-4449-B0	00F-4449-B0	AJO-8328

SecurityGuard™ Cartridges (mm)

for ID: 2.0-3.0 mm

3 µm Analytical Columns (mm)							
	150 x 3.0	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 2.0*	4 x 3.0*
Phases						/10pk	/10pk
PFP	00F-4447-Y0	00B-4447-E0	00D-4447-E0	00F-4447-E0	00G-4447-E0	AJO-8326	AJO-8327
Silica(2)	00F-4162-Y0	00B-4162-E0	00D-4162-E0	00F-4162-E0	00G-4162-E0	AJO-4347	AJO-4348
C8(2)	00F-4248-Y0	00B-4248-E0	00D-4248-E0	00F-4248-E0	—	AJO-4289	AJO-4290
C18(2)	00F-4251-Y0	00B-4251-E0	00D-4251-E0	00F-4251-E0	00G-4251-E0	AJO-4286	AJO-4287
CN	00F-4254-Y0	00B-4254-E0	00D-4254-E0	00F-4254-E0	—	AJO-4304	AJO-4305
Phenyl-Hexyl	00F-4256-Y0	00B-4256-E0	00D-4256-E0	00F-4256-E0	00G-4256-E0	AJO-4350	AJO-4351
NH <sub>2</sub>	00F-4377-Y0	00B-4377-E0	00D-4377-E0	00F-4377-E0	—	AJO-4301	AJO-4302
HILIC	00F-4449-Y0	—	00D-4449-E0	00F-4449-E0	—	AJO-8328	AJO-8329

SecurityGuard™ Cartridges (mm)

for ID: 2.0-3.0 mm 3.2-8.0 mm

5 µm Minibore Columns (mm)			
	50 x 2.0	150 x 2.0	4 x 2.0*
Phases			/10pk
PFP	00B-4448-B0	00F-4448-B0	AJO-8326
Silica(2)	00B-4274-B0	00F-4274-B0	AJO-4347
C5	00B-4043-B0	00F-4043-B0	AJO-4292
C8	00B-4040-B0	00F-4040-B0	AJO-4289
C8 (2)	00B-4249-B0	00F-4249-B0	AJO-4289
C18	00B-4041-B0	00F-4041-B0	AJO-4286
C18 (2)	00B-4252-B0	00F-4252-B0	AJO-4286
CN	00B-4255-B0	00F-4255-B0	AJO-4304
Phenyl-Hexyl	00B-4257-B0	00F-4257-B0	AJO-4350
NH <sub>2</sub>	00B-4378-B0	00F-4378-B0	AJO-4301

SecurityGuard™ Cartridges (mm)

for ID: 2.0-3.0 mm

5 µm Analytical Columns (mm)							
	150 x 3.0	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 2.0*	4 x 3.0*
Phases						/10pk	/10pk
PFP	00F-4448-Y0	00B-4448-E0	00D-4448-E0	00F-4448-E0	00G-4448-E0	AJO-8326	AJO-8327
Silica(2)	—	00B-4274-E0	00D-4274-E0	00F-4274-E0	00G-4274-E0	AJO-4347	AJO-4348
C5	00F-4043-Y0	00B-4043-E0	00D-4043-E0	00F-4043-E0	00G-4043-E0	AJO-4292	AJO-4293
C8	00F-4040-Y0	00B-4040-E0	00D-4040-E0	00F-4040-E0	00G-4040-E0	AJO-4289	AJO-4290
C8(2)	00F-4249-Y0	00B-4249-E0	00D-4249-E0	00F-4249-E0	00G-4249-E0	AJO-4289	AJO-4290
C18	00F-4041-Y0	00B-4041-E0	00D-4041-E0	00F-4041-E0	00G-4041-E0	AJO-4286	AJO-4287
C18(2)	00F-4252-Y0	00B-4252-E0	00D-4252-E0	00F-4252-E0	00G-4252-E0	AJO-4286	AJO-4287
CN	00F-4255-Y0	00B-4255-E0	00D-4255-E0	00F-4255-E0	00G-4255-E0	AJO-4304	AJO-4305
Phenyl-Hexyl	00F-4257-Y0	00B-4257-E0	00D-4257-E0	00F-4257-E0	00G-4257-E0	AJO-4350	AJO-4351
NH <sub>2</sub>	00F-4378-Y0	00B-4378-E0	00D-4378-E0	00F-4378-E0	00G-4378-E0	AJO-4301	AJO-4302
SCX	—	00B-4398-E0	00D-4398-E0	00F-4398-E0	00G-4398-E0	AJO-4307	AJO-4308
HILIC	00F-4450-Y0	—	—	—	—	AJO-8328	AJO-8329

SecurityGuard™ Cartridges (mm)

for ID: 2.0-3.0 mm 3.2-8.0 mm

\* SecurityGuard™ Analytical Cartridges require holder, Part No.: KJ0-4282



## Is your HPLC column protected?



- Easy to use
- Cost-effective
- Universal fit to virtually all manufacturer's columns
- No change in chromatography

*"I consider SecurityGuard a necessity, not an accessory"*

— V. Agarwal / Connecticut, USA



For more information:  
[www.phenomenex.com/SecurityGuard](http://www.phenomenex.com/SecurityGuard)

### additional accessories:

#### Control Temperature


- Improve chromatography
- More reproducible results
- Eliminate varying column temperature



25 - 90 °C

#### Column Selectors

- Quickly screen up to 6 HPLC columns
- Develop, optimize, and validate methods faster



#### Syringe Filters

- Rapid filtration
- High quality
- Available in 6 membrane types





[www.phenomenex.com](http://www.phenomenex.com)

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department by telephone, fax or email: [international@phenomenex.com](mailto:international@phenomenex.com).



<p><b>Australia</b></p> <p>tel.: 02-9428-6444            fax: 02-9428-6445            email: <a href="mailto:auinfo@phenomenex.com">auinfo@phenomenex.com</a></p>	<p><b>Austria</b></p> <p>01-319-1301            01-319-1300            anfrage@phenomenex.com</p>	<p><b>Canada</b></p> <p>(800) 543-3681            (310) 328-7768            info@phenomenex.com</p>	<p><b>Denmark</b></p> <p>4824 8048            4810 6265            dkinfo@phenomenex.com</p>	<p><b>France</b></p> <p>01 30 09 21 10            01 30 09 21 11            franceinfo@phenomenex.com</p>	<p><b>Germany</b></p> <p>06021-58830-0            06021-58830-11            anfrage@phenomenex.com</p>
<p><b>Ireland</b></p> <p>tel.: 01 247 5405            fax: +44 1625-501796            email: <a href="mailto:eireinfo@phenomenex.com">eireinfo@phenomenex.com</a></p>	<p><b>Italy</b></p> <p>051 6327511            051 6327555            italiainfo@phenomenex.com</p>	<p><b>New Zealand</b></p> <p>09-4780951            09-4780952            nzinfo@phenomenex.com</p>	<p><b>Puerto Rico</b></p> <p>(800) 541-HPLC            (310) 328-7768            info@phenomenex.com</p>	<p><b>United Kingdom</b></p> <p>01625-501367            01625-501796            ukinfo@phenomenex.com</p>	<p><b>USA</b></p> <p>(310) 212-0555            (310) 328-7768            info@phenomenex.com</p>