

# The Ultimate Pre-Packed Preparative Column for HPLC and SFC GUARANTEED!

# Axia PREP LC columns offer:

- Increased Performance
- Groundbreaking Lifetimes
- Optimized Loadability
- Increased Reproducibility



# The Axia<sup>™</sup> Advantage

Available in over 40 unique achiral and chiral selectivities, Axia advanced preparative column packing and column hardware design offer several advantages. Unlike traditional column packing methods, the Axia packing method offers increased sorbent bed density for increased performance and eliminates media bed collapse as a source of premature column failure in preparative HPLC/SFC columns.





If Axia packed columns do not provide at least an equivalent separation as compared to a competing preparative column of the same particle size, same phase, and dimensions, return the product with comparative data within 45 days for a FULL REFUND. Only applies to 21.2 mm ID columns.

© 2015 Phenomenex, Inc. All rights reserved.

Axia <sup>™</sup> Technology	pp. 4-7
Award winning column packing technology	
Media Selection Chart	pp. 8-9
Quickly find the media for your purification needs	
Media Selectivity Details	pp. 10-25
Achiral Chemistries	
Core-Shell Media	
Kinetex <sup>®</sup>	pp. 10-13
Aeris <sup>®</sup>	pp. 14-15
Fully Porous Media	
Gemini <sup>®</sup>	pp. 16-17
Synergi <sup>™</sup>	pp. 18-19
Luna <sup>®</sup>	pp. 20-21
Jupite <sup>®</sup>	p. 22
Clarity <sup>®</sup>	p. 23
Chiral Chemistries	
Lux <sup>®</sup>	pp. 24-27
PhenoLogix	pp. 28-29

Screening, Optimization and Scale-up Services

"I find Axia Columns to be very robust and durable. I often use the prep column for much longer than predicted with reproducible peaks. This saves us a significant amount of money."

> David Wisnoski GlaxoSmithKline, USA



# Axia™ Technology vs. Traditional "OBD" Prep Column Packing

# **Traditional Slurry Packing**

Traditional slurry packing processes, like the Waters® OBD™ (Optimum Bed Density) column packing approach, involve the column being removed from the column packing station once it is packed.

Several potential problems with this packing method are:

- Variability in column performance due to increased number of manual operations required for assembly
- · Potential silica media damage during recompression
- Level of process control is based on traditional slurry packing technology

#### Conventional packing process involves: Compression → Decompression → Recompression → Final Column





#### Diagram from Waters Corporation U.S. Patent No. 7,399,410

# **Axia Packing Technology**

Axia packed preparative columns involve a single axial compression step unlike conventional packed preparative columns. The ideal column bed density is custom calculated and automated for each specific media and column size. Computer control of the entire process ensures both proper bed density and column uniformity every time.

During the Axia packing process, the packing piston is locked in place, eliminating any decompression and then recompression of the media sorbent, thus maintaining media and column bed integrity. This solves common lifetime and performance problems associated with conventional packing processes for preparative columns.



# Axia Packing Process Involves:

Compression → Final Column



# Traditional packed preparative columns produce non-uniform media beds with sheared and crushed particles

Decompression and then recompression during packing can damage the media and lead to increased column-to-column variability, flow disturbances, and decreased column lifetimes.



\*The images are believed to be representative, but individual columns may vary.

\*SEM of Waters<sup>®</sup> OBD<sup>™</sup> inlet frit



Crushed media or silica fines at frit surface after packing

# Axia<sup>™</sup> packed columns produce uniform media bed with intact particles

The highly tuned patented process and hardware eliminates potential decompression ensuring bed stability and optimal packing density.

The media found on the inlet frit of the Axia packed column shows no signs of damage unlike the media found on inlet frit of traditionally packed prep columns.



\*SEM of Axia inlet frit

Intact media at frit surface after packing

We are using chromatography media from Phenomenex for GPL/GMP purposes, therefore we audited Phenomenex USA as a manufacturer. From the beginning, we were impressed with Phenomenex and the attitude of their employees. Phenomenex is a unique company in many aspects. Their degree of dedication to customer service, to the organization of the QMS system and last but not least the positive atmosphere in the company is impressive. The outcome of the audit was to our fullest satisfaction.

Major Generic Pharma Company, Europe

View an animated packing process comparison at phenomenex.com/axia

# Axia<sup>™</sup> Technology Outperforms Traditional Packing Processes!

Because of the constant pressure placed on the integrated packing piston, Axia packed columns possess the dynamic capability of maintaining a consistent, homogeneous media bed. This results in superior column performance no matter which media selectivity you choose.

To better understand how much Axia technology improves column performance over traditionally slurry packed preparative columns we scaled-up a 5 µm Lux<sup>®</sup> Cellulose-1 chiral media analytical column and packed the same media into two different 150 x 21.2 mm I.D. columns. One column was packed using Axia technology and the other prep column was packed using the traditional slurry packing process.

The Axia packing technology had a substantial increase in column efficiency resulting in increased resolution over traditionally packed preparative columns. With increased resolution you are able to increase your sample load enabling you to purify more target compound(s) per purification run. This equates to better throughput and economics.

# Warfarin Chiral Purification in Normal Phase Mode



 Resolution calculated with peak width at baseline and center retention time due to the overloaded peaks being off-scale

For more detailed information on this warfarin application, request technical note:

TN-9002: Scaling from Analytical to Preparative Chiral Chromatography While Balancing Purity, Yield, and Throughput under HPLC and SFC Conditions

Tip:

# **Unmatched Column Reproducibility**

The completely automated Axia<sup>™</sup> packing system provides feedback control and infinite tuning of packing density for specific media characteristics such as mechanical strength and porosity. An optimum bed density can be consistently reproduced columnto-column.

This directly translates into consistent efficiency and peak asymmetry measurements and decreases the column variability seen in traditionally packed preparative columns.



Axia columns provide me with first rate quality and engineering. Reliability, reproducibility, and durability are provided with all Axia columns that I use. I can literally purify 2500 samples per column. The time and cost savings are tremendous.

> Derrrick Miyao —Large Biotech Manufacturer, USA

We have used Phenomenex Axia prep-HPLC columns for several years and they consistently provide excellent separation and reproducibility for a variety of different compounds.

> Jeremy R. Wolf ABC Laboratories, USA

# Phase Selection Chart





# First Core-Shell Preparative HPLC/SFC Column Ever!

Kinetex<sup>®</sup> Core-Shell Technology produces increased efficiencies over traditional, fully porous columns, yielding remarkable chromatographic resolution, higher peak capacities, and greater sensitivity, so labs can get even more out of their HPLC analyses!

The benefits of Kinetex Core-Shell Technology include:

- Increased efficiencies over traditional fully porous columns
- Seamless scalability from HPLC/UHPLC to Preparative LC
- Kinetex 5 µm provides better performance than traditional fully porous 5 and 3 µm materials



# **High Column Efficiency**

Combining 5 µm Kinetex core-shell and Axia<sup>™</sup> technology can provide the highest separation efficiency of any pre-packed preparative HPLC column.

#### Waters<sup>®</sup> XBridge<sup>®</sup> 5µm C18 Prep OBD<sup>™</sup>



#### Kinetex 5 µm XB-C18 Axia Packed



							Applicatio	ns		Type of C	ompounds	;	Loading
Key: Best Suited	Very Go	ood				Small Molecules Peptides	Proteins	Chiral	Oligonucle- otides	Polar	Hydro- phobic	Bases	Available Surface Area
Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m²/g)	Carbon Load (%)	pH Range								
Kinetex C18	1.3, 1.7, 2.6, 5	100	200	12	1.5-8.5*								
Kinetex XB-C18	1.7, 2.6, 5	100	200	10	1.5-8.5*								
Kinetex EVO C18	5	100	200	11	1-12								
Kinetex C8	1.7, 2.6, 5	100	200	8	1.5-8.5*								
Kinetex Phenyl-Hexyl	1.7, 2.6, 5	100	200	11	1.5-8.5*								
Kinetex Biphenyl	1.7, 2.6, 5	100	200	11	1.5-8.5*								
Kinetex HILIC	1.7, 2.6, 5	100	200	0	2.0-7.5	$\bullet  \bullet$							

\*Columns are pH stable from 1.5-10 under isocratic conditions. Columns are pH stable 1.5-8.5 under gradient conditions. Comparative separations may not be representative of all applications.



# **Excellent Loadability!**

With narrower peak widths than fully porous columns across every sample load, Axia<sup>™</sup> packed Kinetex 5 µm columns give you the capability of increased sample load and higher throughput for vastly improved purification performance and economics.

#### Waters<sup>®</sup> XBridge<sup>®</sup> 5µm C18 Prep OBD<sup>™</sup>



#### Kinetex 5µm C18 Axia Packed



Kinetex Axia Preparative columns are fantastic! I currently use two Kinetex 5µm C18 150 x 21.2mm columns in parallel for high throughput purifications (<100 mg scale), and Kinetex core shell media delivers significantly improved peak shape and lower back pressure compared to many of the industry. I can also analyze quickly my purified fractions with the same core shell phase on my analytical UPLC<sup>®</sup> system.

Chris DeVore Neurocrine Biosciences, USA

Tip:

If you would like to see a loading study performed with the combination of Axia Packing and Kinetex Core-Shell Particle Technologies, request technical note: TN-1058.

# Seamless Scalability from HPLC/UHPLC to PREP

Kinetex<sup>®</sup> packed with Axia<sup>™</sup> technology makes it the first coreshell sorbent commercially available for small-scale preparative applications. Combine this with the fact that the entire Kinetex core-shell line is fully scalable from  $1.3\,\mu$ m to  $5\,\mu$ m, means that transferring high performance HPLC/UHPLC methods to preparative HPLC and SFC formats is fast and simple.

## Waters<sup>®</sup> XBridge<sup>®</sup> 5µm C18 150 x 4.6mm



# Kinetex 5µm EVO C18 150 x 4.6mm



containens ior an	columna.
Columns:	Kinetex 5 µm EVO C18
	XBridge 5µm C18
Dimensions:	150 x 4.6 mm
	150 x 21.2 mm (Kinetex AXIA Packed)
Mobile Phase:	A: 0.1 % TFA in Water
	B: 0.1 % TFA in Acetonitrile
Gradient:	20% to 70% B over 10 minutes
Flow Rate:	1.5 mL/min
	30 mL/min (Kinetex AXIA)
Temperature:	Ambient
Detection:	UV @ 254 nm
Sample:	Proprietary Pharmaceutical Sample
•	

nditions for all columns

Axia packed column has a great efficiency for the separation of several classes of natural compounds. Due to its low back pressure and therefore high flow work conditions, time for conditioning the columns is sped up greatly!

request technical note: TN-1135.

Tip:

Sylvian Cretton -Europe

6

8

10

12 min

4

For more information on the power of Kinetex core-shell scalability,

2

Comparative separations may not be representative of all applications.

# Kinetex Medio

# **A Broad Spectrum of Column Selectivities**

Kinetex<sup>®</sup> core-shell columns are available in a wide range of stationary phases, allowing you to optimize your separation for maximum resolution and loadability across HPLC, UHPLC, and preparative HPLC and SFC applications.















#### Conditions for all columns: Columns: Kinetex 5 µm C18 Kinetex 5 µm XB-C18 Kinetex 5 µm EVO C18 Kinetex 5 µm C8 Kinetex 5 µm Biphenyl Kinetex 5 µm Phenyl-Hexyl Dimensions: 100 x 4.6 mm Mobile Phase: A: 0.1 % TFA in Water B: 0.1 % TFA in Acetonitrile Gradient: Time (min) % B 0 20 20 22 20 22.5 5 25 5

25 Flow Rate: 1.5 mL/min

# Phenyl-Hexyl



Temperature: 22 °C Detection: UV @ 330 nm Sample: 1. Chlorogenic Acid Others: Antioxidants from green coffee

Sample: 1. Chlorogenic Acid Others: Antioxidants 1



For more information on Chlorogenic Acids from Green Coffee by HPLC, request technical note: TN-1134.



# Increased Performance for Peptide Purifications

Based on core-shell particle technology, Aeris<sup>™</sup> PEPTIDE media is designed with small pores (100Å), an inert XB-C18 surface chemistry, and multiple particle sizes to meet the selectivity, resolution and loading demands of chemists working with synthetic peptides. The benefits of Aeris PEPTIDE columns include:

- Optimized media for peptide purifications
- Multiple particle size options for method development flexibility and peptide impurity analysis
- Seamless scalability from HPLC/UHPLC to preparative HPLC



XB-C18 chemistry best suited for resolving peptides

# Multiple Particle Sizes For Added Versatility



							A	pplication	IS		Type of C	ompounds	5	Loading
						Small				Oligonucle-		Hydro-		
Key: Best Suited	Very Go	od				Molecules	Peptides	Proteins	Chiral	otides	Polar	phobic		Surface Area
Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m²/g)	Carbon Load (%)	pH Range									
Aeris	1.7, 2.6, 3.6, 5	100	200	12	1.5-9									



# Develop, Purify, and Analyze Peptide Fractions with One Media

Aeris PEPTIDE is fully scalable in retention and selectivity with its 4 unique particle sizes (1.7  $\mu$ m, 2.6 µm, 3.6 µm, and 5 µm) for easy transfer from HPLC and UHPLC methods to preparative applications.

# Aeris PEPTIDE 2.6 µm XB-C18



# Aeris PEPTIDE 5µm XB-C18



**Analytical method** 



#### Preparative scale-up and fraction collection



# Aeris PEPTIDE 2.6 µm XB-C18



# Setting the Standard for pH Method Development



Gemini features a pH stability from 1-12, making it optimal for high alkaline washes and high pH purifications of basic drugs.

#### **Optimized parameters include:**

- Innovative surface layer for increased pH stability
- High-surface area for increased loading
- Silica smoothness for stable packing beds
- Bonding density for excellent reproducibility

# Second-Generation TWIN-NX<sup>™</sup> Technology

### Gemini NX-C18

TWIN-NX technology uses an improved patented organo-silica grafting process which incorporates highly stabilizing ethane cross-linking. These organic groups are evenly incorporated into the grafted layers on the silica surface while maintaining a pure silica core. This not only provides resistance to high pH attack, but also maintains the high efficiency and mechanical strength of a silica particle.

\*This bonding technology is also available in Core-Shell media. See Kinetex® EVO on page 10.



App ID 22850

min

Dramatically improve sample resolution, productivity and performance of any preparative column media with Axia<sup>™</sup> column hardware and packing technology. Axia packed prep column offers the opportunity for longer lifetime, higher loading and increased throughput.

#### Gemini 5µm NX-C18 Axia Packed



								Application	ons			Type of C	ompounds	5	Loading
		0				Small				Oligonucle-			Hydro-		Available
Key: Best Suited	U very	Good				Molecules	Peptides	Proteins	Chiral	otides	Acids	Polar	phobic	Bases	Surface Area
Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m²/g)	Carbon Load (%)	pH Range						_				
Gemini C18	3, 5, 10	110	375	14	1.0-12.0										
Gemini C6-Phenyl	3, 5	110	375	12	1.0-12.0										
Gemini NX-C18	3, 5, 10	110	375	14	1.0-12.0										

Comparative separations may not be representative of all applications.

# Gemini Medio

# Flexibility in pH Adjustments Allows for Increased Purification Performance.

Separating basic compounds at higher pH levels produces dramatic changes when compared to low pH conditions. At pH 10.5, the basic compounds become neutralized and are more hydrophobic. The retention for the uncharged basic compounds increases providing an increase in separation along with superior peak shapes.



# Increased Loading with Unique Selectivities



Synergi is available in four unique phases, each offering dramatic differences in:

- Selectivity
- Retention time
- Resolution

The unique selectivity profiles found within the Synergi product line offer complementary selectivity to the standard C18, C8, or silica phases traditionally employed in preparative HPLC.

## Synergi Polar-RP

For Polar and Aromatic Mixtures

#### (100 % Aqueous Stable)

#### Synergi Fusion-RP

Balanced Non-polar and Polar Performance (100 % Aqueous Stable)



# Embedded polar group complements C18 ligand with balanced polar selectivity (Ag) TMS TMS TMS Ultra-pure Silica

#### Synergi Hydro-RP

Strong Non-polar and Polar Retention (100 % Aqueous Stable)



### Synergi Max-RP

**Excellent for Basic Compounds at Neutral pH** 



							A	pplication	IS			Type of C	ompounds	5	Loading
Key: Best Suited	Very Good	d				Small	Pontidos	Protoine	Chiral	Oligonu-		Polar	Hydro-		Available
			Surface			Wolecules	replices	FIOLEINS	Onna	ciedides	Acias	r Ulai	phobic	Dases	
Packing Material	Particle Size (µm)	Pore Size (Å)	Area (m²/g)	Carbon Load (%)	pH Range										
Synergi Fusion-RP	4, 10	80	475	12	1.5-10.0*										
Synergi Max-RP	4, 10	80	475	17	1.5-10.0*										
Synergi Hydro-RP	4, 10	80	475	19	1.5-7.5										
Synergi Polar-RP	4, 10	80	475	11	1.5-7.0										

\*pH range is 1.5-10 under isocratic conditions and 1.5-9.0 under gradient conditions.

# **Selectivity Like No Other**

Offering a balanced combination of hydrophobic and polar selectivity, Synergi<sup>™</sup> Fusion-RP separates compounds exhibiting moderately polar and hydrophobic characteristics.

#### Hydrophobic basic compounds



from other RP phases."

The slightest variations in compound polarity and aromaticity are exploited by Synergi Polar-RP to achieve the greatest separation between polar and/or aromatic compounds.

#### Increased resolution of polar compounds with Synergi Polar-RP compared to traditional C18 phases



CARBOGEN AMCIS, Switzerland

Ordering information on page 31

to often show the desired selectivity, distinguishing this phase



Comparative separations may not be representative of all applications

# Media for One of the World's Leading PREP HPLC Columns

Luna<sup>®</sup> high surface area (400 m<sup>2</sup>/g) silica packing materials provide optimized parameters specifically designed for the purification of small molecules and peptides. This media allows high loading with excellent lifetimes.

#### **Optimized loading parameters include:**

- Silica smoothness for stable packed beds
- Optimum pore size/distribution provide outstanding performance
- High pore volume offers increased surface area
- · Fine tuned bonding density for excellent reproducibility
- Greater loading capacity with an extended pH range of 1.5 to 10.0\*



We routinely use Axia<sup>™</sup> packed columns from Phenomenex for peptide purifications. Among various preparative HPLC columns we have used, the Axia packed Luna columns (5µm) stand out. We have been very satisfied with the increased loading capacity and excellent performance.

> Guangcheng Jiang Ferring Research Institute, Inc., USA

								Application	S			Type of C	ompounds		Loading
Key: Best Suited	Very Good					Small Molecules	Peptides	Proteins	Chiral	Oligonucle- otides	Acids	Polar	Hydro- phobic	Bases	Available Surface Area
Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m²/g)	Carbon Load (%)	pH Range										
Luna C18(2)	3, 5, 10, 10- <i>PREP</i> , 15	100	400	17.5	1.5-10.0*										
Luna C8(2)	3, 5, 10, 10- <i>PREP</i> , 15	100	400	13.5	1.5-10.0*										
Luna C5	5, 10	100	440	12.5	1.5-10.0*										
Luna Phenyl-Hexyl	3, 5, 10, 10- <i>PREP</i> , 15	100	400	17.5	1.5-10.0*										
Luna Silica(2)	3, 5, 10, 10- <i>PREP</i> , 15	100	400	-	-										
Luna CN	3, 5, 10	100	400	7	1.5-7.0										
Luna NH <sub>2</sub>	3, 5, 10	100	400	9.5	1.5-11.0										
Luna SCX	5, 10	100	400		2-7.0										
Luna HILIC	3, 5	200	200	5.7	1.5-8.0										
Luna PFP(2)	5, 10	100	400	11.5	1.5-9.0										
*nH range is 1 5-10 unde	er isocratic conditions and	d 1 5-8 5 under	gradient cond	litions											



# **Simple Scale-Up**

Axia<sup>™</sup> column technology provides the same high efficiency chromatographic performance for preparative scale columns (21.2, 30, and 50 mm ID) as obtained in 4.6 mm ID analytical columns. This improvement in preparative column performance across

all lengths and internal diameters makes it easier to select the appropriate column size to achieve the desired purity and yield without having to compromise on performance.



# **Proven Media for Peptide Purifications**

Optimal compromise between throughput, recovery, yield. Ability to perform high loading (0.74 g on column) and achieve high purity (>98 %) in a single purification run.

# Preparative Purification of Bivalirudin (20 amino acid peptide also know as Angiomax)



#### **Purity Confirmation of Combined Fractions**



#### 11 Combined fractions 27.8 – 29.8 min; Recovery 80.5 % with purity $\geq$ 98.5 %

Peak No.	Time (min)	Area	Area %
1	12.74	73.7	0.35
2	13.83	40.6	0.19
3	14.37	21118.7	98.53
4	15.858	200.5	0.93

Column: Luna 5 µm C18(2) Dimensions: 250 x 4.6 mm

Flow Rate: 1 mL/min Temperature: 25 °C Detection: UV @ 220 nm

Mobile Phase: A: 0.1 % TFA in Water B: 0.1 % TFA in Acetonitrile Injection Volume: 2 µL Gradient: 20% to 50% B in 30 min

Sample: Combined Fractions



Ordering information on page 31

# Media for Biomolecules

The Jupiter HPLC column portfolio, including Jupiter 300 and Jupiter Proteo, offers optimized reversed phase solutions for peptide and protein purification. Identify, purify, and analyze almost any protein with Jupiter columns.



# Jupiter Proteo 90 Å

- For separation of proteins and peptides < 10,000 MW</li>
- C12 bonded onto an ultra-high surface area (475 m<sup>2</sup>/g) silica for increased peak capacity and resolution of peptide separations
- Direct scale-up from analytical to preparative and bulk materials

#### **Resolve Peptides with Similar Hydrophobicity**

Jupiter Proteo is able to fully resolve peptides that differ in hydrophobicity by one methyl group.



- 4. Ac-Arg-Gly-Val-Gly-Gly-Gly-Leu-Gly-Leu-Gly-Leu-Gly-Lys-Amide
- 4. Ac-Arg-Gly-Val-Gly-Gly-Leu-Gly-Leu-Gly-Lys-Amide 5. Ac-Arg-Gly-Val-Val-Gly-Gly-Leu-Gly-Leu-Gly-Lys-Amide

# Jupiter 300Å

- For separation of proteins > 10,000 MW
- Available with C18 and C4 bonded phases
- 1.5 10 pH stability for method ruggedness and easy protein removal
- Direct scale up to preparative and bulk materials

#### **Compare PEGylated vs. Native Forms of Proteins**

Reversed phase separation of PEGylated and native proteins on a Jupiter 300 C4 column. Note the good resolution of multiple PEGylated forms for all proteins tested.



We purchased the Jupiter 300 C18 300 Å column a few months ago and have been quite impressed with its performance. The Jupiter 300 column provides better separation of the proteins. As for reproducibility, the control profiles have not changed since day one of its use.

Major Biotech Company, Europe

							A	pplication	S		1	Type of Co	ompound trais	s	Loading
Key: Best Suite	ed 🌔 Very	y Good				Small Molecules	Peptides		Chiral	Oligonu- cleotides		Polar	Hydro- phobic	Bases	Available Surface Area
Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m²/g)	Carbon Load (%)	pH Range										
Jupiter C18	5, 10, 15	300	170	13.3	1.5-10.0										
Jupiter C4	5, 10, 15	300	170	5	1.5-10.0										
Jupiter Proteo	4, 10	90	475	15	1.5-10.0										

Comparative separations may not be representative of all applications.

# Purification of Synthetic Oligonucleotides



# Clarity Oligo-RP<sup>™</sup>

Unique media specifically designed for reversed phase purification of oligonucleotides with balanced hydrophobicity and polar selectivity. The media is based on composite particle TWIN<sup>™</sup> technology that provides improved selectivity and efficiency for oligonucleotides when compared to competing hybrid, polymer, and silica media.

# **RP-HPLC Preparative Purification**

- Easily separate N-1 failure sequences from target oligo with >90 % purities
- Purify oligos up to 60 nt in length
- Trityl-off purification of DNA, RNA, thioates, and modified/ labeled oligonucleotides
- 3µm, 5µm, 10µm particles for seamless scaling

#### Preparative 20nt DNA Oligo-RP Purification



# Clarity Oligo-WAX™

Clarity Oligo-WAX is a crosslinked weak anion exchanger media designed for successful ion-exchange purification of synthetic DNA/RNA. Oligo-WAX is an advantageous combination of purity, capacity, mechanical strength, cost, and efficiency.

- Excellent efficiency column results in > 90 % purities due to good fractionation of closely eluting compounds
- High loading capacity due to very high density ligand
- Increase productivity by running at higher flow rates and pressures

## Purify Failure Sequences and Contaminants from Target Sequence

Ion-exchange is an excellent separation mode for purifying contaminants and failure sequences from target sequences. Clarity Oligo-WAX, due to its increased efficiency compared to other ion-exchange columns, has the ability to recognize minute charge differences in nucleotide sequences such as failure sequences or base substitutions.

#### **DNA Purification of N-1 Sequence from Target N Sequence**



We have used the Axia<sup>™</sup> prep columns and not had problems with them. I have never had to adjust for retention gaps. This speaks directly to the quality of Phenomenex's phases and the quality of their PREP columns.

-Major Biotech Company, USA

							A	pplication	s		Type of C	ompounds	;	Loading
									-		Neu	itrals		
Key: Best	t Suited	) Very Good	d			Small Molecules	Peptides	Proteins	Chiral	Oligonucle- otides	Polar	Hydro- phobic		Available Surface Area
Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m²/g)	Carbon Load (%)	pH Range									
Oligo-RP	3, 5, 10	110	375	14	1.0-12.0									
Oligo-Wax	10	360	150	-	1.0-11.0									

# Complete Chiral Solutions



Achieving optimal chiral separation is easier than ever with five unique Lux polysaccharide stationary phases to screen. Choose a phase, then transfer the method to lab scale, process, pilot, and commercial scale.

# Lux chiral preparative columns simplify the separation process:

- Unique and traditional phases that increase the success rate of the chiral screen
- Consistent particle size distribution so performance is maintained
- Mechanically strong media for increased stability
- Available in multiple particle sizes for direct scale up (3µm and 5µm packed columns for screening and small scale purifications; 10µm and 20µm bulk media for process scale purifications)



Cellulose tris(3,5-dimethylphenylcarbamate) Guaranteed Alternative to CHIRALCEL® OD®, OD-H®, OD-3, OD-RH, and OD-3R



#### Lux Cellulose-3

Cellulose tris(4-methylbenzoate) Guaranteed Alternative to CHIRALCEL OJ<sup>®</sup> , OJ-H<sup>®</sup>, OJ-3, OJ-RH, and OJ-3R



#### Lux Amylose-1 Amylose tris(3,5-dimethylphenylcarbamate) Guaranteed Alternative to

CHIRALPAK<sup>®</sup> AD<sup>®</sup>, AD-H<sup>®</sup>, AD-3, AD-RH<sup>®</sup>, and AD-3R

# **Resolve Your Enantiomers with Five** Unique Phases

The Lux family of bulk cellulose chiral selectors provides a variety of complementary selectivities.

# Screen for the most effective chiral separation under the following conditions:

- Reversed Phase
- Polar Organic
- Normal Phase
- Supercritical Fluid Chromatography (SFC)



Cellulose tris(3-chloro-4-methylphenylcarbamate) Guaranteed Alternative to CHIRALCEL OZ, OZ-H®, OZ-3, OZ-RH, and OZ-3R



Cellulose tris(4-chloro-3-methylphenylcarbamate) Guaranteed Alternative to CHIRALCEL OX-H, OX-3, OX-RH, and OX-3R



Lux Amylose-2

Amylose tris(5-chloro-2-methylphenylcarbamate) Guaranteed Alternative to CHIRALPAK® AY®, AY-H®, AY-3, AY-RH, and AY-3R

							Chir	al Applicat	tions		Тур	oe of Chir	al Compou	inds	Loading
												Ne	utrals		
Key: Best Suited	Very G	iood				Small Molecules	Pep- tides	Proteins	Chiral	Oligonu- cleotides	Acids	Polar	Hydro- phobic	Bases	Available Surface Area
Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m²/g)	Carbon Load (%)	pH Range										
Lux Cellulose-1	3, 5, 10, 20	1,000	-	-	2-9.0										
Lux Cellulose-2	3, 5, 10, 20	1,000	-	-	2-9.0										
Lux Cellulose-3	3, 5, 10, 20	1,000	-	-	2-9.0										
Lux Cellulose-4	3, 5, 10, 20	1,000	-	-	2-9.0										
Lux Amylose-1	5	1,000	-	-	2-9.0										
Lux Amylose-2	3, 5	1,000	-	-	2-9.0										

CHIRALPAK, CHIRALCEL, AY, AY-H, OD, OD-H, OJ, OJ-H, and OZ-H are registered trademarks of DAICEL Corporation. Phenomenex is in no way affiliated with DAICEL Corporation.

# **Column Screening for Optimal Chiral Resolution**

Being able to utilize differences in selectivity in each of the five Lux<sup>®</sup> columns can help develop methods more efficiently by offering broad and complementary chiral recognition abilities.

In the example below, a simple screen determined which column gave the best separation.

#### Etozolin

Based on a five phase screen under reversed phase conditions, the optimal chiral stationary phase for resolving Etozolin is Lux Cellulose-3.

#### Lux 5 µm Cellulose-3



Conditions for all columns: Column: As noted Dimension: 250 x 4.6 mm Mobile Phase: Acetonitrile / 20 mM Ammonium bicarbonate with 0.1 % Diethylamine (60:40) Flow Rate: 1 mL/min Temperature: Ambient Detection: UV @ 220 nm







#### Lux 3 µm Cellulose-4



# Innovative chiral selector will succeed where others fail.

16

18 min

10

12

14

#### Lux 5µm Cellulose-4



CHIRALPAK, CHIRALCEL, OD-H, and AD are registered trademarks of DAICEL Corporation. Columns used for comparison were manufactured by DAICEL Corporation. Phenomenex is in no way affiliated with DAICEL Corporation. Comparative separations may not be representative of all applications.

# Load More with an Increase in Column Length

Axia<sup>™</sup> column technology allows separation to scale up directly based on column length. With the 100 mm length column a 32 mg/load separation was achieved and an increased sample

load of 80 mg/load was achieved on the longer 250 mm length column. As expected when increasing the load, the peak width and tailing increased but there was no loss of resolution.



Lux Axia preparative column are wonderful! I regularly use Lux chiral stationary phase Cellulose-2 and Cellulose-4 and less frequently, the Lux Amylose-2. In our community of chiral analysis/purification scientists, there are some who use the CC4 column instead of the \*equivalent\* Lux Cellulose-4. On several occasions we've seen separation and good peak shape on the Lux Cellulose-4 that was completely missing from the CC4. Customer support and delivery times are always within a few days.

> Julia G. Christie GlaxoSmithKline, USA

# **Easy SFC Scale-up**

# SFC Purfication of Terfenadine



For SFC column screening, use Lux 150 x 3.0 mm ID columns.



Ordering information on page 31

Tip:

# A New Era of Technical Support Services

# Let Us Do the Work for You

PhenoLogix, our in-house application support lab, saves you time and money by screening multiple scout columns and solvent strategies for new purification methods or revalidating your current methods. We work together to make you successful by minimizing your process purification development time and optimizing your purification method.

#### **Chiral Screening**

- Normal Phase
- Reversed Phase
- Polar Organic
- SFC

## **Method Optimization Services**

- Fast Turnaround
- Easy Method Transfer
- Continued Support

#### **Preparative and Process Scale-Up**

- Media Screening
- Small Scale Purification
- DAC Packing Assistance





# Your Method, Our Scientists

Quality Products, Advanced Performance, Complete Support

For more information or to begin a project today, please contact your local Phenomenex representative or email us at **phenologix@phenomenex.com** 

## You can also visit us online:

#### www.phenomenex.com/phenologix



Our scientists at American Peptide have taken advantage of Phenomenex's column packing services, application development, and project-specific consultation services for some of our most challenging separations.

American Peptide Company, USA



# Ordering Information

# **Achiral Phases**

Finase 5µm PEPTIDE XB-C18 Kinetex® Phase 5µm XB-C18 EVO C18	Length 150 250 Length 50 100 150 150 250	ID 21.2 21.2 ID 21.2 21.2 21.2 21.2 21.2	Part No. 00F-4632-P0-AX 00G-4632-P0-AX Part No. 00B-4627-P0-AX
PEPTIDE XB-C18 Kinetex® Phase 5µm XB-C18 EVO C18	150 250 Length 50 100 150 150 250	21.2 21.2 ID 21.2 21.2 21.2 21.2	00F-4632-P0-AX 00G-4632-P0-AX Part No. 00B-4627-P0-AX
Kinetex® Phase 5µm XB-C18	250 Length 50 100 150 150 250	21.2 21.2 1D 21.2 21.2 21.2 21.2 21.2	00G-4632-P0-AX
Kinetex® Phase 5µm XB-C18 EVO C18	Length 50 100 150 150 250	ID 21.2 21.2 21.2 21.2	Part No.
Kinetex® Phase 5µm XB-C18 EVO C18	Length 50 100 150 150 250	ID 21.2 21.2 21.2	<b>Part No.</b> 00B-4627-P0-AX
Phase 5µm XB-C18 EVO C18	Length 50 100 150 150 250	ID 21.2 21.2 21.2 21.2	Part No. 00B-4627-P0-AX
5µm XB-C18 EVO C18	50 100 150 150 250	21.2 21.2 21.2	00B-4627-P0-AX
XB-C18 EVO C18	50 100 150 150 250	21.2 21.2 21.2	00B-4627-P0-AX
EVO C18	100 150 150 250	21.2 21.2	
EVO C18	150 150 250	21.2	00D-4627-P0-AX
EVO C18	150 250		00F-4627-P0-AX
EVO C18	250	30	00F-4627-U0-AX
EVO C18		21.2	00G-4627-P0-AX
	50	21.2	00B-4633-P0-AX
	100	21.2	00D-4633-P0-AX
	100	30	00D-4633-U0-AX
	150	21.2	00F-4633-P0-AX
	150	30	00F-4633-U0-AX
	250	21.2	00G-4633-P0-AX
	250	30	00G-4633-U0-AX
Biphenyl	50	21.2	00B-4605-P0-AX
	50	30	00B-4605-U0-AX
	100	21.2	00D-4605-P0-AX
	100	30	00D-4605-U0-AX
	150	21.2	00F-4605-P0-AX
	150	30	00F-4605-U0-AX
	250	21.2	00G-4605-P0-AX
	250	30	00G-4605-U0-AX
HILIC	100	21.2	00D-4606-P0-AX
	150	21.2	00F-4606-P0-AX
	250	21.2	00G-4606-P0-AX
C18	50	21.2	00B-4601-P0-AX
	50	30	00B-4601-U0-AX
	100	21.2	00D-4601-P0-AX
	100	30	00D-4601-110-4X
	150	21.2	00E-4601-P0-4X
	150	30	00F-4601-10-4X
	250	21.2	00G-4601-P0-AX
	250	30	00G-4601-U0-AX
C8	50	21.2	00B-4608-P0-AX
	50	30	00B-4608-U0-AX
	100	21.2	00D-4608-P0-4X
	100	30	00D-4608-110-4X
	150	21.2	00F-4608-P0-AY
	150	30	00F-4608-110-AY
	250	21.2	006-4608-20-44
	250	30	00G-4608-110 AV
Phonyl_Hoyyl	50	21.2	008-4000-00-AX
і пенуі-пехуі	50	21.2	00B-4003-FU-AX
	100	21.2	00D-4003-00-AX
	100	21.2	00D-4003-PU-AX
	100	30	00D-4003-00-AX
	150	21.2	00E-4603-P0-4V

бµт NX-C18 С18	50 50 75 100 100 150 150 250 250	21.2 30 30 21.2 30 21.2 30 21.2 30 21.2	00B-4454-P0-AX 00B-4454-U0-AX 00C-4454-U0-AX 00D-4454-P0-AX 00D-4454-U0-AX 00F-4454-P0-AX
NX-C18 C18	50 50 75 100 100 150 150 250 250	21.2 30 30 21.2 30 21.2 30 21.2 30 21.2	00B-4454-P0-AX 00B-4454-U0-AX 00C-4454-U0-AX 00D-4454-P0-AX 00D-4454-U0-AX 00F-4454-P0-AX
C18	50 75 100 100 150 150 250 250	30 30 21.2 30 21.2 30 21.2	00B-4454-U0-AX 00C-4454-U0-AX 00D-4454-P0-AX 00D-4454-U0-AX 00F-4454-P0-AX
C18	75 100 100 150 150 250 250	30 21.2 30 21.2 30 21.2	00C-4454-U0-AX 00D-4454-P0-AX 00D-4454-U0-AX 00F-4454-P0-AX
C18	100 100 150 150 250 250	21.2 30 21.2 30 21.2	00D-4454-P0-AX 00D-4454-U0-AX 00F-4454-P0-AX
C18	100 150 150 250 250	30 21.2 30 21.2	00D-4454-U0-AX 00F-4454-P0-AX
C18	150 150 250 250	21.2 30 21.2	00F-4454-P0-AX
C18	150 250 250	30 21.2	
C18	250 250	21.2	00F-4454-U0-AX
C18	250		00G-4454-P0-AX
C18		30	00G-4454-U0-AX
	50	21.2	00B-4435-P0-AX
	50	30	00B-4435-U0-AX
	75	30	00C-4435-U0-AX
	100	21.2	00D-4435-P0-AX
	100	30	00D-4435-U0-AX
	150	21.2	00F-4435-P0-AX
	150	30	00F-4435-U0-AX
	250	21.2	00G-4435-P0-AX
	250	30	00G-4435-U0-AX
C6-Phenyl	75	30	00C-4444-U0-AX
	100	21.2	00D-4444-P0-AX
	100	30	00D-4444-U0-AX
	150	21.2	00F-4444-P0-AX
	250	21.2	00G-4444-P0-AX
0 µm			
NX-C18	50	21.2	00B-4455-P0-AX
	50	50	00B-4455-V0-AX
	100	21.2	00D-4455-P0-AX
	100	30	00D-4455-U0-AX
	100	50	00D-4455-V0-AX
	150	21.2	00F-4455-P0-AX
	150	30	00F-4455-U0-AX
	150	50	00F-4455-V0-AX
	250	21.2	00G-4455-P0-AX
	250	30	00G-4455-U0-AX
	250	50	00G-4455-V0-AX
C18	50	21.2	00B-4436-P0-AX
	100	21.2	00D-4436-P0-AX
	100	30	00D-4436-U0-AX
	100	50	00D-4436-V0-AX
	150	21.2	00F-4436-P0-AX
	150	30	00F-4436-U0-AX
	150	50	00F-4436-V0-AX
	250	21.2	00G-4436-P0-AX
	250	30	00G-4436-U0-AX
	250	50	00G-4436-V0-AX
Jupiter <sup>®</sup>			

Longth	ID	Bort No
Lengui	עו	Fall NO.
250	21.2	00G-4397-P0-AX
250	30	00G-4397-U0-AX
100	21.2	00D-4055-P0-AX
250	30	00G-4055-U0-AX
100	21.2	00D-4168-P0-AX
150	21.2	00F-4168-P0-AX
250	21.2	00G-4168-P0-AX
	Length 250 250 100 250 100 150 250	Length ID   250 21.2   250 30   100 21.2   250 30   100 21.2   250 30   100 21.2   250 30   100 21.2   250 21.2   250 21.2   250 21.2

# guarantee

21.2

30

00G-4603-P0-AX

00G-4603-U0-AX

250

250

If Axia<sup>™</sup> packed columns do not provide at least an equivalent separation as compared to a competing preparative column of the same particle size, same phase, and dimensions, return the product with comparative data within 45 days for a FULL REFUND. Only applies to 21.2 mm ID columns.

For additional sizes not displayed, please contact your Phenomenex technical consultant or local distributor.

Luna®			
Phase	Length	ID	Part No.
5µm			
C18(2)	50	21.2	00B-4252-P0-AX
	50	30	00B-4252-00-AX
	100	21.2	00D-4252-PU-AX
	150	30 21.2	00D-4232-00-AX
	250	21.2	001-4252-P0-AX
	250	30	000-4252-10-AX
C8(2)	50	21.2	00B-4249-P0-AX
(-)	50	30	00B-4249-U0-AX
	100	21.2	00D-4249-P0-AX
	100	30	00D-4249-U0-AX
	150	21.2	00F-4249-P0-AX
C5	150	21.2	00F-4043-P0-AX
CN	50	21.2	00B-4255-P0-AX
	100	30	00D-4255-U0-AX
	150	21.2	00F-4255-P0-AX
	250	21.2	00G-4255-P0-AX
Phenyl-Hexyl	50	21.2	00B-4257-P0-AX
	100	21.2	00D-4257-P0-AX
	100	30	00D-4257-00-AX
NU	100	21.2	00D 4279 D0 AX
N11 <sub>2</sub>	150	21.2	00D-4378-P0-AX
	250	21.2	00G-4378-P0-AX
HILIC	50	21.2	00B-4450-P0-AX
	100	21.2	00D-4450-P0-AX
	150	21.2	00F-4450-P0-AX
	250	21.2	00G-4450-P0-AX
	250	30	00G-4450-U0-AX
PFP(2)	50	21.2	00B-4448-P0-AX
	100	21.2	00D-4448-P0-AX
	100	30	00D-4448-U0-AX
	150	21.2	00F-4448-P0-AX
	250	21.2	00G-4448-P0-AX
0111	250	30	00G-4448-U0-AX
Silica (2)	50	21.2	00B-42/4-P0-AX
	100	21.2	00E 4274-PU-AX
	250	21.2	00G-4274-PO-AX
	250	30	00G-4274-U0-AX
10 um	200		
C18(2)	50	21.2	00B-4253-P0-AX
	50	30	00B-4253-U0-AX
	100	21.2	00D-4253-P0-AX
	250	21.2	00G-4253-P0-AX
	250	30	00G-4253-U0-AX
00(0)	250	50	00G-4253-V0-AX
C8(2)	50	21.2	00B-4250-P0-AX
	250	21.2	00G-4250-PU-AX
C5	100	21.2	000-4230-00-ΑΧ
00	250	21.2	00G-4092-P0-AX
	250	50	00G-4092-V0-AX
CN	250	21.2	00G-4300-P0-AX
Phenyl-Hexyl	250	21.2	00G-4285-P0-AX
	250	30	00G-4285-U0-AX
NH <sub>2</sub>	250	21.2	00G-4379-P0-AX
Silica (2)	250	21.2	00G-4091-P0-AX
	250	30	00G-4091-U0-AX
	250	50	00G-4091-V0-AX
15 µm			
C18(2)	250	21.2	00G-4273-P0-AX
	250	30	00G-4273-U0-AX
09(2)	250	50	00G-4273-V0-AX
GO(Z)	250	21.2	00G 4272-PU-AX
51116a ( <i>L)</i>	250	30	00G-4271-10-AX

Tip:

## Protect your Axia Prep Columns with SecurityGuard. SecurityGuard

Contact your local Phenomenex representative for ordering information or visit: www.phenomenex.com/guardit

# Synergi<sup>™</sup>

Phase	Length	ID	Part No.
4 µm			
Fusion-RP	50	21.2	00B-4424-P0-AX
	100	21.2	00D-4424-P0-AX
	100	30	00D-4424-U0-AX
	150	21.2	00F-4424-P0-AX
	250	21.2	00G-4424-P0-AX
Hydro-RP	50	21.2	00B-4375-P0-AX
	75	30	00C-4375-U0-AX
	100	21.2	00D-4375-P0-AX
	100	30	00D-4375-U0-AX
	150	21.2	00F-4375-P0-AX
	250	21.2	00G-4375-P0-AX
	250	30	00G-4375-U0-AX
Max-RP	50	21.2	00B-4337-P0-AX
	100	21.2	00D-4337-P0-AX
	100	30	00D-4337-U0-AX
	150	21.2	00F-4337-P0-AX
	250	21.2	00G-4337-P0-AX
	250	30	00G-4337-U0-AX
Polar-RP	50	21.2	00B-4336-P0-AX
	50	30	00B-4336-U0-AX
	75	30	00C-4336-U0-AX
	100	21.2	00D-4336-P0-AX
	100	30	00D-4336-U0-AX
	150	21.2	00F-4336-P0-AX
	250	21.2	00G-4336-P0-AX
	250	30	00G-4336-U0-AX
10 µm			
Fusion-RP	250	21.2	00G-4425-P0-AX
Hydro-RP	250	21.2	00G-4376-P0-AX
-	250	30	00G-4376-U0-AX
Max-RP	100	21.2	00D-4350-P0-AX
	250	21.2	00G-4350-P0-AX
	250	30	00G-4350-U0-AX
Polar-RP	250	21.2	00G-4351-P0-AX

### **Clarity**<sup>®</sup>

Phase	Length	ID	Part No.
5µm			
Oligo-RP <sup>™</sup>	100	21.2	00D-4442-P0-AX
	250	21.2	00G-4442-P0-AX
10 µm			
Oligo-RP	150	21.2	00F-4445-P0-AX
	150	30	00F-4445-U0-AX
	250	21.2	00G-4445-P0-AX
<b>Oligo-WAX</b> ™	150	21.2	00F-4451-P0-AX
	150	30	00F-4451-U0-AX
	250	21.2	00G-4451-P0-AX

# **Chiral Phases**

Lux®			
Phase	Length	ID	Part No.
5µm			
Amylose-1	150	21.2	00F-4732-P0-AX
NEW	150	30	00F-4732-U0-AX
	250	21.2	00G-4732-P0-AX
	250	30	00G-4732-U0-AX
	250	50	00G-4732-V0-AX
Amylose-2	150	21.2	00F-4472-P0-AX
	250	21.2	00G-4472-P0-AX
	250	30	00G-4472-U0-AX
Cellulose-1	150	21.2	00F-4459-P0-AX
	250	21.2	00G-4459-P0-AX
	250	30	00G-4459-U0-AX
	250	50	00G-4459-V0-AX
Cellulose-2	150	21.2	00F-4457-P0-AX
	250	21.2	00G-4457-P0-AX
	250	30	00G-4457-U0-AX
	250	50	00G-4457-V0-AX
Cellulose-3	150	21.2	00F-4493-P0-AX
	250	21.2	00G-4493-P0-AX
	250	30	00G-4493-U0-AX
	250	50	00G-4493-V0-AX
Cellulose-4	150	21.2	00F-4491-P0-AX
	250	21.2	00G-4491-P0-AX
	250	30	00G-4491-U0-AX
	250	50	00G-4491-V0-AX

# The Ultimate Pre-packed Preparative HPLC/SFC Column

Australia t: +61 (0)2-9428-6444 f: +61 (0)2-9428-6445 auinfo@phenomenex.com

Austria t: +43 (0)1-319-1301 f: +43 (0)1-319-1300 anfrage@phenomenex.com

**Belgium** t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) f: +31 (0)30-2383749

beinfo@phenomenex.com

**Canada** t: +1 (800) 543-3681 f: +1 (310) 328-7768

info@phenomenex.com

China t: +86 (0)20 2282-6668 f: +86 (0)20 2809-8130

chinainfo@phenomenex.com

**Denmark** t: +45 4824 8048 f: +45 4810 6265 nordicinfo@phenomenex.com

**Finland** t: +358 (0)9 4789 0063 f: +45 4810 6265

nordicinfo@phenomenex.com

France t: +33 (0)1 30 09 21 10 f: +33 (0)1 30 09 21 11 franceinfo@phenomenex.com

**Germany** t: +49 (0)6021-58830-0 f: +49 (0)6021-58830-11 anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400 f: +91 (0)40-3012 2411 indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405 f: +44 1625-501796 eireinfo@phenomenex.com

Italy

t: +39 051 6327511 f: +39 051 6327555 italiainfo@phenomenex.com

# phenomenex ..breaking with tradition

#### www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at www.Phenomenex.com/TermsAndConditions

Trademarks

Aeris, Clarity, Gemini, Kinetex, Luna, and Lux are registered trademarks, and Oligo-RP, Oligo-WAX, SecurityGuard, Synergi, TWIN-NX, TWIN Technology, and Two-In-One Technology are trademarks of Phenomenex. UPLC, XBridge, Symmetry, XTerra, and Waters are registered trademarks and OBD and SymmetryShield are trademarks of Vater Corporation. ZORBAX is a registered trademark of Agilent Technologies, Inc. Vydac is a registered trademark of Alltech Associates, Inc. CHIRALPAK, CHIRALCEL, AD, AY, AY-H, OD, OD-H, OJ, OJ-H, and OZ-H are registered trademarks of DAICEL Corporation.

#### Disclaimer

7890815\_W

BR17

Comparative separations may not be representative of all applications. Phenomenex is not affiliated with Waters Corporation. Agilent Technologies, Inc., Alltech Associates, Inc. or DAICEL Corporation, Columns used for comparison studies were manufactured by and purchased from Waters Corp., Agilent Technologies, Inc., Alitech Associates, Inc. and DAICEL Corporation.

The opinions stated herein are solely those of the speaker and not necessarily those of any company or organization. Axia column and packing technology is patented by Phenomenex. U.S. Patent No. 7, 674, 383

Gemini and Kinetex EVO are patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign

counterparts © 2015 Phenomenex, Inc. All rights reserved.

- Luxembourg t: +31 (0)30-2418700 f: +31 (0)30-2383749 nlinfo@phenomenex.com

# Mexico t: 001-800-844-5226 f: 001-310-328-7768

tecnicomx@phenomenex.com

**The Netherlands** t: +31 (0)30-2418700 f: +31 (0)30-2383749 nlinfo@phenomenex.com

**New Zealand** t: +64 (0)9-4780951 f: +64 (0)9-4780952 nzinfo@phenomenex.com

# **Norway** t: +47 810 02 005

- f: +45 4810 6265 nordicinfo@phenomenex.com

- **Puerto Rico** t: +1 (800) 541-HPLC f: +1 (310) 328-7768 info@phenomenex.com
- Spain
- t: +34 91-413-8613 f: +34 91-413-2290
- espinfo@phenomenex.com
- **Sweden** t: +46 (0)8 611 6950
- f: +45 4810 6265 nordicinfo@phenomenex.com

#### **United Kingdom** t: +44 (0)1625-501367

f: +44 (0)1625-501796 ukinfo@phenomenex.com

#### USA

- t: +1 (310) 212-0555 f: +1 (310) 328-7768
- info@phenomenex.com

- f: +1 (310) 328-7768
- All other countries Corporate Office USA t: +1 (310) 212-0555 info@phenomenex.com