

APPLICATIONS

A Fast and Effective Quantitation Method for Uracil, 5,6-Dihydrouracil, and 5-Fluorouracil from Human Serum by LC-MS/MS

Shahana Huq
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

Uracil, 5-Fluorouracil (5FU), and 5,6-Dihydrouracil (UH₂) offer unique markers for enzymatic activity. In order to determine the 5FU catabolic rate, uracil and its homologues are tested in correlation to each other to ultimately determine if proper metabolism is occurring. The primary goal of this study was to develop a sample preparation and LC-MS/MS method for quantitative analysis of uracil and its two homologues from human serum. This application focuses on developing a sample preparation method and LC-MS/MS analysis using a Strata[®]-X PRO solid phase extraction (SPE) and a Kinetex[®] PS C18 column respectively. The Kinetex PS C18 Column is a positively charged, surface modified C18 phase that caters its unique selectivity to separating the very polar uracil and its homologues in this analysis, while the novel Strata-X PRO greatly reduces the phospholipids in the sample and provides cleaner extracts.

SPE Protocol

Pretreatment: Add 100 µL Human serum (doubly stripped) to a tube. Dispense 900 µL of chilled (~0°C) Acetonitrile to the serum while vortexing. Centrifuge samples at 6000 rpm for 5 minutes.

96-Well Plate: Strata-X PRO, 30 mg/well

Part No.: [8B-S536-TGA](#)

Condition: 1 mL Acetonitrile

Load: Pass the supernatant (from pre-treatment) and apply vacuum to collect eluted extract

Dry Down: Under Nitrogen and heat around 40-45°C

Reconstitute: 100 µL initial mobile phase

Note: For phospholipid analysis human plasma EDTA was utilized and direct injection (bypass dry down and reconstitution steps) of the eluted sample was made.

LC-MS/MS Conditions

Quantitative Analysis for Uracil Homologues

Column:	Kinetex 2.6 µm PS C18	
Dimensions:	150 x 3.0 mm	
Part No.:	00F-4780-Y0	
Mobile Phase:	A: 0.1% Formic acid in Water B: Methanol	
Isocratic:	Time (min)	% B
	0	7
	12	7
Flow Rate:	0.2 mL/min	
Injection Volume:	10 µL	
Column Temperature:	25 °C	
Instrument:	Agilent [®] 1260	
Detection:	MS/MS SCIEX API Triple Quad [™] 4500, Dual Polarity (Positive for U and UH ₂ , Negative for 5 FU)	

Qualitative Analysis for Phospholipids

Column:	Kinetex 2.6 µm C18	
Dimensions:	50 x 2.1 mm	
Part No.:	00B-4462-AN	
Mobile Phase:	A: 0.1% Formic acid in Water B: 0.1% Formic Acid in Methanol	
Gradient:	Time (min)	% B
	0	40
	0.5	95
	11.5	95
	11.51	40
	13.5	40
Flow Rate:	0.4 mL/min	
Injection Volume:	2 µL	
Column Temperature:	40 °C	
Instrument:	Agilent 1260	
Detection:	MS/MS SCIEX API Triple Quad 4500, ESI Source (Positive)	

Figure 1. Structure of Uracil and Homologues





Table 1. Retention Time (RT), MRM Transition and % Recovery for Analytes

Analyte	RT	Q1	Q3	Spike Concentration (ng/mL)	% Recovery	% CV
UH ₂	4.31	114.9	55.04	100	90	5.4
Uracil	4.65	112.9	69.8	100	84	3.2
5FU	7.07	128.8	41.9 86.1 58.9	100	89	3.2
Uracil 1,3-15N ₂ (+Ve IS)	4.65	114.8	96.9	200	N/A	N/A
5 Cl Uracil (-Ve IS)	10.82	145.1	42.1	200	N/A	N/A

Figure 2. Representative Chromatogram of Extracted Human Serum Analyzed by a Kinetex 2.6 μ m PS C18 LC Column Under ESI Positive Polarity



App ID 25777

Figure 3. Representative Chromatogram of Extracted Human Serum Analyzed by a Kinetex[®] 2.6 μ m PS C18 LC column Under ESI Negative Polarity





Figure 4. Representative Chromatogram for Qualitative Evaluation of Phospholipid in Extracted Human Plasma Samples, Analyzed by a Kinetex 2.6 μ m C18 LC Column (A) Protein Precipitation (B) Strata[®]-X PRO Extraction



References

1. Barbara Buchel, Peter Rhyn, Stefan Schurch, Biomed. Chromatogr. 2013; **27**: 7-16
2. Ruta Svobaitė, Isabella Solassol, Frederic Pinguet, Clinical Chemistry 54:9, 1463-1472 (2008)
3. Jenny P. Dai, Amra Tabakovic, Welley Loc, Current Trends in Mass Spectrometry October 2013

APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/ChatNow.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 (12) 881 0121
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

☎ All other countries/regions Corporate Office USA

t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Strata and Kinetex are registered trademarks of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc. SCIEX is a registered trademark and Triple Quad is a trademarks of AB SCIEX Pte. Ltd. AB SCIEX™ is being used under license.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is not affiliated with Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.