

APPLICATIONS

LC-MS/MS Quantitative Analysis of NDMA in Ranitidine Active Pharmaceutical Ingredient (API) and Drug Product using the SCIEX[®] 4500 QTRAP[™]

Sandeep Choudhary¹, Aman Sharma¹, Chandrasekar¹, Manoj Pillai¹, and Swetha Kotikalapudi²

¹SCIEX, Gurgaon, India

²Phenomenex India

Introduction

The U.S. Food and Drug Administration has learned that some ranitidine medicines, including some products commonly known as the brand-name drug Zantac[®], contain low levels of the genotoxic nitrosamine impurity N-nitrosodimethylamine (NDMA).

Nitrosamines such as NDMA are classified as probable human carcinogens based on results from laboratory tests. Nitrosamines have been found to be present in water and foods, including cured meats, dairy products, and vegetables. Nitrosamines have the potential to be intermediates in organic synthesis and due to their potent genotoxicity, it has become a requirement to accurately quantitate nitrosamines in pharmaceuticals during drug development and manufacturing.

NDMA is a hydrophilic impurity and such compounds present retention challenges for the C18 columns typically used as they offer little interaction and retention for polar molecules. The accurate quantitative analysis of a polar compound like NDMA requires good retention and separation from potential matrix interferences, in order to minimize potential matrix generated signal suppression. In order to achieve enhanced retention and separation, a high efficiency Kinetex[®] core-shell column containing a Biphenyl stationary phase was used for the determination of NDMA in ranitidine samples. This column is different from traditional C18 columns as it offers enhanced polar selectivity, which should be useful for improved retention of the polar NDMA compound.

The LC-MS/MS method described here is similar to the LC-MS/MS based analytical method recommended by USFDA¹ as other methods such as GC-MS, which use elevated temperatures have been reported to generate NDMA via reaction of ranitidine with itself, and this will result in an overestimation of NDMA present in ranitidine API and drug products².

Materials

NDMA standard was obtained as a kind gift from a customer. LC-MS grade Methanol was obtained from Biosolve, and LC-MS grade Water was obtained from Merck. Sample diluent was 10% Methanol in Water, by volume.

Standards Preparation

A 1 mg/mL standard stock solution of NDMA was prepared in diluent and stored at -20 °C.

The working standards were prepared from the stock solution by serial dilution in diluent to obtain the concentrations of 500, 100, 50, 25, 12.5, 6.25, 4.5, 1.5, 0.75, and 0.5 ng/mL.

Sample Preparation for API

The API sample was accurately weighed (100 mg) and the volume was made up to 2 mL with diluent and vortexed until completely dissolved. This solution was used for analysis.

Sample Preparation for Drug Product

The tablets equivalent to 600 mg ranitidine were accurately weighed and crushed. The powdered tablets were taken, weighed and dissolved in appropriate volume of diluent to achieve a target concentration of 50 mg/mL of ranitidine. The suspension was then vortexed for 15 minutes followed by sonication for another 10 minutes. Finally, this solution was centrifuged at 4000 rpm for 5 minutes and the resulting supernatant filtered through a PVDF syringe filter. The first 1 mL of filtrate was discarded, and the remainder of the filtered solution was used for analysis.

The HPLC and MS/MS conditions are noted below. Data acquisition and processing was performed using Analyst[®] 1.6.3 and MultiQuant[™] 3.0.3 was used for data analysis. A 1/x2 weighted linear regression was used to generate the calibration curve.

LC-MS/MS Conditions

Analytical Column: Kinetex 2.6 µm Biphenyl			
Dimensions: 100 x 4.6 mm			
Part No.: 00D-4622-E0			
Mobile Phase: A: 0.05% Formic acid in Water B: 0.05% Formic acid in Methanol			
Gradient:	Time (min)	% B	Flow Rate (mL/min)
	0	5	0.4
	6	5	0.4
	6.1	80	0.4
	6.5	80	0.4
	6.6	80	1.0
	7.6	80	1.2
	13.5	80	1.0
	13.6	5	0.4
	17	5	0.4
Valco Valve:	Time (min)	Position	
	6.9	B (out)	
	13.5	A (in)	
Temperature: 35 °C			
Injection Volume: 20 µL			
HPLC System: Shimadzu® NEXERA® (Shimadzu Corporation)			
Mass Spectrometer: SCIEX® 4500 QTRAP™ (AB SCIEX Pte. Ltd.)			
MS Mode: APCI Positive			
Scan Type: MRM			
Nebulizer Current: 3 µA			
Curtain Gas: 30 psi			
GS1: 40 psi			
CAD: 40 °C			
Source Temperature: 350 °C			
Dwell Time: 100 msec			
UV Detection: 254 nm			

Compound Retention Times and Dependent Mass Spectrometer Values

Name	Q1	Q3	DP (V)	EP (V)	CE (V)	CXP (V)	RT (min)
NDMA-1	75.05	58.1	45	10	25	10	5.6
NDMA-2	75.05	43	45	10	15	10	5.6
Ranitidine	315	229.1	45	10	25	10	7.3

Results and Discussion

Figure 1 shows the MRM chromatogram of the impurity NDMA (75.05/58.1 & 75.05/43) and ranitidine (315.1/229.9). The corresponding UV chromatogram is shown at the bottom of **Figure 1**. The separation between the NDMA and ranitidine peaks was optimized using the gradient mobile phase conditions detailed in the LC-MS/MS conditions. The MS/MS parameters and analyte retention times are also noted above.

Linearity Range, LOD & LOQ

The calibration curve (**Figure 2**) consisted of 10 calibration points covering the concentration range from 0.5 ng/mL to 500 ng/mL. The calibration curves were constructed by fitting the analyte concentrations versus the peak area of the analyte with a regression line with $1/x^2$ weighting. The regression coefficient, r^2 , determined from the calibration curves was 0.99. The linearity was obtained over the concentration range of 0.5 ng/mL to 500 ng/mL, corresponding to 0.01 – 10 ppm in 50 mg API. The LOQ and LOD for NDMA were 1.5 ng/mL and 0.5 ng/mL, corresponding to 0.03 and 0.01 ppm, respectively, with respect to 50 mg/mL API (**Figure 3**).

Accuracy, precision and recovery

The accuracy and precision were demonstrated by spiking a known concentration of NDMA at three different concentrations (LOD, LOQ and specification limit). The accuracy and precision of the method was evaluated based on recovery studies using standard addition to API. Recovery of NDMA was determined in fortified samples. The overall recoveries observed at all three levels were between 80 to 120% as summarized in **Table 1**. The results demonstrate that this method is suitable for the accurate analysis of NDMA in API and drug product.

Figure 1.

Chromatogram of Ranitidine sample containing NDMA contamination, upper chromatogram showing 3 transition of MS/MS for NDMA (75/43 - red, 75/58.1 - blue) and Ranitidine (315.1/229.9 - green), and the lower chromatogram showing UV detection at 254 nm with ranitidine eluting at 7.01 – 8.7 in Valco waste line.





Figure 2.

Calibration Curve for NDMA over the calibration range 0.5 to 500 ng/mL for m/z 75.0 ▶ 58.1 (top) and 75.0 ▶ 43.0 (bottom).

**Figure 3.**

Extracted Ion Chromatograms of NDMA for diluent, specification, LOQ and LOD level.




Table 1.
Accuracy, precision and recovery in spiked sample

Concentration Sample (ppm)	NDMA (in 50 mg/mL API)			
	Concentration relative to 50 mg/mL ranitidine	Accuracy (n=6)	Precision (% RSD, n=6) (Spiked Standard)	% Recovery (n=6)
Control sample (50 mg/mL ranitidine)	5 ng/mL	101.18%	3.19	(80-120)
LOD (0.01 ppm)	0.5 ng/mL	101.04%	1.49	(80-120)
LLOQ (0.03 ppm)	1.5 ng/mL	102.77%	2.5	(80-120)
Spec Level (0.09 ppm)	4.5 ng/mL	104.07%	1.43	(80-120)

Conclusion

The LC-MS/MS method for analysis for NDMA in ranitidine was successfully developed and applied to the analysis of NDMA impurities in pharmaceutical API and drug product. The method accurately quantifies the levels of NDMA below the levels required by USFDA and EMA, ensuring the safety of the API and pharmaceutical drug product, while avoiding the use of high temperatures (as in the GC method) which can result in overestimation of NDMA in ranitidine samples. The choice of the Kinetex Biphenyl column was demonstrated to give retention for the polar NDMA impurity, while providing sufficient resolution between NDMA and Ranitidine to minimize any potential for matrix interference.

References

1. Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS) Method for the Determination of NDMA in Ranitidine Drug Substance and Solid Dosage Drug Product; <https://www.fda.gov/media/131868/download>.
2. <https://www.valisure.com/wp-content/uploads/Valisure-Ranitidine-FDA-Citizen-Petition-v4.12.pdf>

Kinetex® LC Column Ordering Information

2.6 µm Analytical Columns (mm)			SecurityGuard ULTRA Cartridges [‡]
Phases	100 x 4.6	150 x 4.6	3/pk
Biphenyl	00D-4622-E0	00F-4622-E0	AJ0-9207 for 4.6 mm ID

[‡]SecurityGuard ULTRA Cartridges require Holder, Part No.: AJ0-9000

BE-HAPPY™
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.
www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex is a registered trademark and SecurityGuard and BE-HAPPY are trademarks of Phenomenex. Shimadzu and Nexera are registered trademarks of Shimadzu Corporation. SCIEX and Analyst are registered trademarks, QTRAP and MultiQuant are trademarks of AB SCIEX Pte. Ltd. AB SCIEX™ is being used under license. Zantac is a registered trademark of Boehringer Ingelheim Pharmaceuticals, Inc., used under license.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.

Australia
t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria
t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium
t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada
t: +1 (800) 543-3681
info@phenomenex.com

China
t: +86 400-606-8099
cninfo@phenomenex.com

Denmark
t: +45 4824 8048
nordicinfo@phenomenex.com

Finland
t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France
t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany
t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India
t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland
t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy
t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg
t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico
t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands
t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand
t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway
t: +47 810 02 005
nordicinfo@phenomenex.com

Poland
t: +48 (12) 881 0121
pl-info@phenomenex.com

Portugal
t: +351 221 450 488
ptinfo@phenomenex.com

Singapore
t: +65 800-852-3944
sginfo@phenomenex.com

Spain
t: +34 91-413-8613
esinfo@phenomenex.com

Sweden
t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland
t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan
t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom
t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA
t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com