

Multi-Residue Pesticide Screening Method using GC/MS

Sky Countryman and Kory Kelly
Phenomenex, Inc., Torrance, CA, USA

Pesticides are widely used by farmers to control pests, weeds and molds that would otherwise decrease crop production. While this has significantly increased worldwide food productions, these same pesticides pose significant health and environmental risks.


The restrictions for specific pesticides differ from one country to the next. As world trade increases, the potential threat to other countries' populations increases. This is especially true in the European Union, where produce can be transported from one country to another quite easily.

For this reason, pesticides are the subjects of increasing regulation. Since many different types of pesticides can be used on the same food product, Multi-Residue screening approaches are used to look for more than 300 compounds at a time. Gas Chromatography (GC) is still the most commonly used method for the majority of the pesticide classes. While ECD or NPD may

be used for screening, Mass Spectrometer (MS) detection must be employed to provide positive confirmation.

Zebron MultiResidue™ columns represent a solution for all classes of pesticides analysis. The columns were developed using two new stationary phases that are unlike any commercially available today. Each phase was optimized to resolve a different set of analytes. However, both are good for a wide variety of pesticides.

Zebron MultiResidue columns are well suited for use on all types of detectors. Both columns are MS certified, so they can also be used with GC/MS for multi-residue pesticide methods. The unique selectivity offered by each phase allows for increased resolution of critical compounds vs. standard 5ms type phases. Retention time data is available for over 300 different pesticides on GC/MS, further simplifying new method development.



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department by telephone, fax or email: international@phenomenex.com.

Australia

tel.: 02-9428-6444
 fax: 02-9428-6445
 email: info@phenomenex.com.au

Ireland

tel.: 01 247 5405
 fax: +44 1625-501796
 email: eireinfo@phenomenex.com

Austria

01-319-1301
 01-319-1300
 anfrage@phenomenex.com

Italy

051 736176
 051 735302
 italiainfo@phenomenex.com

Canada

(800) 543-3681
 (310) 328-7768
 info@phenomenex.com

New Zealand

09-4780951
 09-4780952
 info@phenomenex.co.nz

Denmark

4824 8048
 4810 6265
 dkinfo@phenomenex.com

Puerto Rico

(800) 541-HPLC
 (310) 328-7768
 info@phenomenex.com

France

01 30 09 21 10
 01 30 09 21 11
 franceinfo@phenomenex.com

United Kingdom

01625-501367
 01625-501796
 ukinfo@phenomenex.com


Germany

06021-58830-0
 06021-58830-11
 anfrage@phenomenex.com

USA

(310) 212-0555
 (310) 328-7768
 info@phenomenex.com

phenomenex
 ...breaking with tradition™



Peak No.	Sample Analyte
1	Dichlorvos
2	EPTC
3	Butylate
4	3,5-Dichlorobenzoic acid (methyl ester)
5	Vernolate
6	Pebulate
7	Mevinphos
8	4-Nitrophenol (methyl ester)
9	Mevinphos Isomer
10	Trichlorfon
11	Dicamba (methyl ester)
12	MCPP (methyl ester)
13	Molinate
14	Tebuthiuron
15	MCPA (methyl ester)
16	DEET
17	Tetraethyl pyrophosphate (methyl ester)
18	Demeton
19	Dichloroprop (methyl ester)
20	Trifluralin
21	Thionazin
22	Cycloate
23	Benefin
24	Propachlor
25	Ethoprop
26	Chlorpropham
27	2,4-D (methyl ester)
28	Sulfotep
29	Naled

Peak No.	Sample Analyte
30	Phorate
31	Dicrotophos
32	Pentachlorophenol (methyl ester)
33	Profluralin
34	Demeton Isomer
35	Prometon
36	Atraton
37	Monocrotophos
38	Atraton Isomer
39	Silvex (methyl ester)
40	Terbufos
41	Propazine
42	Diazinon
43	Pronamide
44	Atrazine
45	Simazine
46	Terbutylazine
47	Dioxathion
48	Fonofos
49	Dimethoate
50	2,4,5-T Methyl ester
51	Disulfoton
52	Chloramben (methyl ester)
53	Phosphamidon Isomer
54	Secbumeton
55	Dichlofenthion
56	2,4-DB (methyl ester)
57	Terbacil
58	Dinoseb (methyl ester)
59	Alachlor

Peak No.	Sample Analyte
60	Chlorpyrifos methyl
61	Phosphamidon
62	Ronnel
63	Prometryn
64	Ametryn
65	Bentazon (methyl ester)
66	Aspon
67	Simetryn
68	Metribuzin
69	Methyl parathion
70	Terbutryn
71	Metolachlor
72	Malathion
73	DCPA
74	Fenitrothion
75	Chlorpyrifos
76	Trichloronate
77	Triadimeton
78	Pichloram (methyl)
79	Isopropalin
80	Fenthion
81	MGK-264 Isomer
82	Parathion
83	Merphos
84	Bromacil
85	Clofenvinfos Isomer
86	MGK-624
87	Pendimethalin
88	Diphenamid
89	Clofenvinfos
90	Butachlor

Peak No.	Sample Analyte
91	Crotoxyphos
92	Stirofos
93	Tokuthion
94	Oxadiazon
95	Merphos Oxide
96	Napropamide
97	Fenamiphos
98	Oxyfluorfen
99	Acifluorfen
100	Carboxin
101	Ethion
102	Tricyclazole
103	Fensulfthion
104	Carbofenotion
105	Famfur
106	Norflurazon
107	Hexazinone
108	EPN
109	Phosmet
110	Leptophos
111	Azinphos-Methyl
112	Fenarimol
113	Azinphos-ethyl
114	Coumaphos

For a digital copy of this Technical Note, please visit www.Phenomenex.com/TechNotes/2021

Zebron and MultiResidue are trademarks of Phenomenex Inc. Phenomenex is a registered trademark of Phenomenex Inc. ©2007 Phenomenex, Inc. All rights reserved. Copying or re-use of this information is not allowed without written permissions from Phenomenex.

Retention times for additional pesticides using Zebron MultiResidue-1

Column: Zebron MultiResidue-1
Dimensions: 30 meter x 0.25 mm x 0.25 µm
Part No.: 7HG-G016-11
Injection: Splitless at 280 °C, 1 µL
Carrier Gas: Constant Flow Helium 1.2 mL/min
Oven Program: 60 °C (Hold 1 min) to 130 °C at 10 °C/min to 230 °C at 4 °C/min to 300 °C (Hold @ 300 °C for 10 min)
Detector: MSD; 300 °C

Retention Time	Compound
6.34	Metaldehyde
7.14	Phenmedipham (frag-1)
7.81	Monolinuron IP (frag)
8.28	Propoxur - C ₂ H ₃ NO
8.61	Ethiolate
9.07	Isoproturon (frag)
10.21	Bendiocarb - C ₂ H ₃ NO
10.70	Carbofuran - C ₂ H ₃ NO
10.85	Linuron IP (frag)
11.43	EPTC/Eptam
11.85	Dichlobenil
12.52	Biphenyl
12.78	Butylate
13.36	Vernolate
13.49	Dioxacarb - C ₂ H ₃ NO
13.73	Pebulate
13.88	Etridiazole
14.00	3,4-Dichloroaniline
14.36	Propham
14.68	Metolcarb
14.96	d10-Acena-phthalene
15.30	Folpet IP (frag)
15.44	Chloroneb
15.71	Pentachlorobenzene
15.91	o-Phenylphenol
15.92	Captan IP frag (THPI)
15.95	Carbaryl - C ₂ H ₃ NO
16.14	Molinate
16.24	Tebuthiuron
16.38	Isoprocab
17.80	Tecnazene
18.08	Fenobucarb
18.14	Propachlor
18.31	Cycloate/Ro-Neet
18.40	Propoxur - C ₂ H ₃ NO
18.50	Ethilfluralin
18.53	2,3,5,6-Tetrachloroaniline

Retention Time	Compound
18.82	Trifluralin
18.86	Diphenylamine
18.99	Benfluralin
19.27	Chlordimeform
19.50	CIPC
19.57	2,3,5-Trimethacarb
19.63	Tebutam
19.67	Diallate-1
19.83	Penmediapham (frag-2)
20.11	Bendiocarb - C ₂ H ₃ NO
20.16	MBTZ
20.19	Diallate-2
20.31	Promecarb
20.36	Alpha-BHC
20.55	2,6-DIPN
20.66	Hexachlorobenzene
21.30	Ethoxyquin
21.56	Prometon
21.58	Profluralin
21.65	3,4,5-Trimethacarb
21.68	Bufencarb-1
21.77	Desmedipham
21.79	Dichloran
21.86	Command/Clomazone
21.97	Carbofuran
21.98	Simazine
22.01	Monolinuron
22.02	Atrazine
22.02	Terbumeton
22.23	Quintozene
22.38	gamma-BHC (Lindane)
22.45	Cycluron
22.46	Terbutylazine
22.47	Aminocarb
22.63	Fluchloralin
22.74	Pronamide
22.76	Tefluthrin
22.84	Isocarbamid
22.98	Triallate
23.11	Dinitramine
23.19	Pyrimethanil
23.19	Pentachlorobenzonitrile
23.25	Bufencarb-2
23.58	Secbumeton
23.59	Mexacarbate
23.93	Pirimacarb
23.99	beta-BHC
23.99	Fenfuram
24.19	Octhilinone
24.22	beta-Spiroxamine
24.48	Pentachloroaniline

Retention Time	Compound
24.83	Ethiofencarb
24.84	Dimethachlor
24.98	delta-BHC
25.08	Desmetryn
25.23	Heptachlor
25.29	Vinclozolin
25.51	Dioxacarb - C ₂ H ₃ NO
25.55	Chlorothalonil
25.59	a-Spiroxamine
25.66	Propanil
25.73	Cymiazole
25.91	Metalaxyl
25.94	Prometryn
25.98	Metribuzen
25.98	Ametryn
25.98	Simetryn
26.20	Kresoxim-Methyl
26.35	Carbaryl - C ₂ H ₃ NO
26.37	3-OH Carbofuran
26.55	Terbutryn
26.57	Fuberidazole
26.70	Pentachlorothioanisole
26.84	Aldrin
26.86	Ethofumesate
26.95	Fenpropimorph
26.95	Dichlofluanid
27.04	Metolachlor
27.05	Methiocarb
27.14	Dursban
27.15	Norea
27.25	Dacthal/DCPA
27.49	Naphthalene Acetamide
27.66	Butralin
27.74	Diethofencarb
27.97	Triadimefon
28.09	Isopropalin
28.10	Nitrothal isopropyl
28.28	MGK-264 (peak-1)
28.39	Methfuroxan
28.44	Tetraconazole
28.50	p,p'-Dicofol (frag)
28.51	Cyanazine
28.81	Pendimethalin

Continued on next page

Retention times for additional pesticides using Zebron MultiResidue - 1 (cont.)

Retention Time	Compound	Retention Time	Compound	Retention Time	Compound
28.86	Diphenamid	34.51	Endosulfan-II	39.72	Pyridaben
28.89	Desmethyl diphenamid	34.62	Pyrethrin-2	39.87	Prochloraz
28.91	Oxychlorane	34.81	Oxadixyl	39.96	Fluquinconazole
29.01	MGK-264 (peak-2)	35.05	Benalaxyl	40.42	Cypermethrins (4 peaks)
29.20	Heptachlor Epoxide	35.11	Mepronil	40.61	Fenbuconazole
29.24	Cyprodinil	35.14	Trifloxystrobin	40.73/40.99	Flucythrinate (2 peaks)
29.25	Pyracarbolid	35.53	Propiconazole/Tilt-1	40.88	Etofenprox
29.34 / 29.40	Allethrin	35.55	Quinoxifen	41.00	Boscalid
29.48	Tolyfluanid	35.68	p,p'-DDT	41.78/41.91	Fluvalinate (2 peaks)
29.60	Pyrifenoxy-1	35.73	o,p'-Methoxychlor	41.80	Fenvalerate-1
29.64	Penaconazole	35.86	Propiconazole/Tilt-2	42.12	Fenvalerate-2
29.85	Fipronil	35.92 / 35.97	Propargite	42.24/44.98	Dimethomorph (2 peaks)
30.00	Furalaxyl	35.95	Lenacil	42.41	Pyraclostrobin
30.06	Triflumizole	36.00	Endosulfan Sulfate	42.75/43.12	Decamethrin (2 peaks)
30.18	Procymidone	36.04	Piperonyl Butoxide	43.06/43.18	Difenoconazole (2 peaks)
30.33/30.76	Triadimenol (2 Peaks)	36.09	Resmethrin	43.87	Azoxystrobin
30.51	trans-Chlordane	36.17	Nuarimol	44.63	Famoxadone
30.55	Aniten/Flurecol Butyl E.	36.27	Hexazinone		
30.60	Captan	36.30	Carbosulfan		
30.81	Pyrifenoxy-2	36.34	Tebuconazole		
30.82	Folpet	36.44	Nitralin		
30.86	cis-Chlordane	36.45	Sethoxydim		
30.88	Thiabendazole	36.65	Bifenthrin		
30.97	Endosulfan-I	36.70	Epoxiconazole		
30.98	trans-Nonachlor	36.91	Tetramethrin-1		
31.62	Napropamide	37.05	Bromopropylate		
31.94	Oxydiazon	37.13	Fenpropathrin		
32.02	Hexaconazole	37.16	Tetramethrin-2		
32.04	Isoprothiolane	37.36	p,p'-Methoxychlor		
32.10	Flutolanil	37.39	Iprodione		
32.18	Fludioxinil	37.40	Tebufenpyrad		
32.24	p,p'-DDE	37.45	Fenoxycarb		
32.25	Imazalil	37.47	Bifenazate		
32.36	Pyrethrin-1	37.48	Fenhexamid		
32.39	Dieldrin	37.55	Fenazaquin		
32.42	Buprofezin	37.80	d-Phenothrin		
32.49	Bupirimate	37.93	Acetamiprid		
33.13	Methoprotryne	38.08	Tetradifon		
33.20	Flusilazole	38.22	Pyriproxifen		
33.30	Carboxin	38.24	lambda-Cyhalothrin		
33.36	Myclobutanil	38.29	Acrinathrin		
33.36	Tricyclazole	38.34	Triticonazole		
33.39	Endrin	38.37	Mirex		
33.85/33.95	Cyproconazole (2 peaks)	38.42	Amitraz		
33.88	Chlorobenzilate	38.66	Tralkoxydim		
34.01	Nitrofen	38.93	Fenarimol		
34.06	o,p'-DDT	38.98	Spirodiclofen		
34.10	cis-Nonachlor	38.99	Naproanalid		
34.20	Diniconazole	39.29	cis-Permethrin		
34.32	Etaconazole-1	39.49	trans-Permethrin		
34.49	Etaconazole-2	39.56	Biteranol-1		
		39.72	Biteranol-2		

We would like to thank **Greg Mercer** at the **US Food and Drug Administration** for supplying the retention time data for the pesticides.

ORDERING INFORMATION

Part No.
7HG-G016-11-TN
Description
MultiResidue-1, 30 m x 0.25 mm x 0.25 µm
Part No.
7HG-G017-10-TN
Description
MultiResidue-2, 30 m x 0.25 mm x 0.20 µm