

# APPLICATIONS

## Reproducible FAME Analysis on 100 meter Zebron™ ZB-FAME GC Column

Ramkumar Dhandapani, Zara Jalali, Rola Elabaji and Richard Jack  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

### Overview

Fatty Acid Methyl Esters (FAMES) are derivatized esters of long chain fatty acids and are monitored in food for finger printing authenticity of expensive olive oils, determining trans fat levels, and ascertaining omega-3 and omega-6 content. The separation is extremely challenging because of the presence of cis and trans isomers of the FAMES that are structurally similar to each other and exhibit similar boiling points. In this study, a highly crosslinked GC stationary phase ZB-FAME was explored for the separation of cis/trans FAME compounds. ZB-FAME GC columns are made of unique, high cyano stationary phase selectivity that can recognize the difference between FAME compounds and isomers that are very similar in physical properties. In addition to providing the required selectivity, ZB-FAME offers a highly crosslinked stationary phase that provides low bleed and high temperature limits of 280 °C . ZB-FAME is available in a traditional 100-meter dimension as well as shorter column dimensions that can provide up to 6 times faster analysis.

### GC-FID Conditions

**Column (Figure 1 & 2):** Zebron ZB-FAME

**Dimension:** 100 meter x 0.25 mm x 0.20 µm

**Part No.:** [7MG-G033-10](#)

**Column (Figure 3):** Agilent® J&W HP-88

**Dimension:** 100 meter x 0.25 mm x 0.20 µm

**Recommended Z-Guard™:** [7AG-G000-00-GZK](#)

**Injection:** Split (50:1) @ 250 °C, 1 µL

**Recommended Liner:** Zebron PLUS Liner Compatible with Agilent GC Instrument, Single Taper Z-Liner with Wool

**Liner Part No.:** [AG2-0A13-05](#)


**Carrier Gas:** Hydrogen @ 1.2 mL/min (constant flow)

**Oven Program:** 100 °C for 5 min, 140 °C @ 3 °C/min for 5 min, 190 °C @ 1.5 °C/min, 260 °C @ 6 °C/min for 5 min

**Detector:** GC-FID


**Temperature:** 280 °C

**Figure 1: Separation of 37 Component FAME Mix on a 100-meter Zebron ZB-FAME GC Column**




Analyte details in Table 1

**Figure 2: Precision of 6 Injections of 37 Component FAME Mix on a 100-meter ZB-FAME GC Column**



Same GC Method Parameters as in page 1

**Figure 3: Separation of 37 Component FAME Mix on a 100-Meter Agilent® J&W HP-88 GC Column**



App ID 26281

Same GC Method Parameters as in page 1

Comparative separations may not be representative of all applications.

**Table 1.**  
**Peak Area and Retention Time of 37 Component FAME Mix on a 100-meter ZB-FAME GC Column**


Peak	Analyte Name	Alternate Nomenclature	Resolution	RT (min)	RT %RSD (n=6)	Area %RSD (n=6)
1	Butanoic Acid Methyl Ester	C4:0	11.0	6.87	0.02	1.32
2	Hexanoic Acid Methyl Ester	C6:0	43.5	8.30	0.01	1.51
3	Octanoic Acid Methyl Ester	C8:0	79.9	11.20	0.01	1.96
4	Decanoic Acid Methyl Ester	C10:0	113.1	15.78	0.01	2.38
5	Undecanoic Acid Methyl Ester	C11:0	63.1	18.45	0.01	2.59
6	Dodecanoic Acid Methyl Ester	C12:0	64.4	21.49	0.01	2.72
7	Tridecanoic Acid Methyl Ester	C13:0	67.3	25.16	0.02	2.58
8	Myristic Acid Methyl Ester	C14:0	69.0	29.27	0.02	2.74
9	Myristoleic Acid Methyl Ester	C14:1 cis 9	40.0	31.77	0.03	2.62
10	Pentadecanoic Acid Methyl Ester	C15:0	29.3	33.58	0.03	2.74
11	cis-10-Pentadecenoic Acid Methyl Ester	C15:1 cis 10	40.2	36.14	0.03	2.62
12	Hexadecanoic Acid Methyl Ester	C16:0	27.7	38.09	0.03	2.86
13	Palmitoleic Acid Methyl Ester	C16:1 cis 9	27.2	39.99	0.03	2.80
14	Heptadecanoic Acid Methyl Ester	C17:0	38.6	42.53	0.02	2.86
15	cis-10-Heptadecenoic Acid Methyl Ester	C17:1 cis 10	28.3	44.41	0.02	2.95
16	Stearic Acid Methyl Ester	C18:0	36.5	47.02	0.02	2.95
17	Elaidic Acid Methyl Ester	C18:1 trans 9	12.6	47.95	0.03	2.95
18	Oleic Acid Methyl Ester	C18:1 cis 9	7.2	48.45	0.03	3.05
19	Linolelaidic Acid Methyl Ester	C18:2 trans 9,12	24.1	50.10	0.03	3.30
20	Linoleic Acid Methyl Ester	C18:2 cis 9,12	17.7	51.31	0.03	3.26
21	γ-Linolenic Acid Methyl Ester	C18:3 cis 6,9,12	29.7	53.34	0.03	3.29
22	α-Linolenic Acid Methyl Ester	C18:3 cis 9,12,15	22.8	54.86	0.03	3.24
23	Arachidic Acid Methyl Ester	C20:0	11.0	55.66	0.03	3.11
24	cis-11-Eicosenoic Acid Methyl Ester	C20:1 cis 11	15.9	56.84	0.03	3.08
25	cis-11,14-Eicosadienoic Acid Methyl Ester	C20:2 cis 11,14	37.1	59.11	0.02	3.23
26	Heneicosanoic Acid Methyl Ester	C21:0	3.3	59.29	0.03	3.29
27	cis-8,11,14-Eicosatrienoic Acid Methyl Ester	C20:3 cis 8,11,14	20.4	60.39	0.02	3.15
28	Arachidonic Acid Methyl Ester	C20:4 cis 5,8,11,14	18.3	61.27	0.02	3.19
29	cis-11,14,17-Eicosatrienoic Acid Methyl Ester	C20:3 cis 11,14,17	1.6	61.34	0.02	3.64
30	Behenic Acid Methyl Ester	C22:0	12.0	61.94	0.02	3.39
31	Erucic Acid Methyl Ester	C22:1 cis 13	11.8	62.52	0.02	3.42
32	cis-5,8,11,14,17-Eicosapentaenoic Acid Methyl Ester	C20:5 cis 5,8,11,14,17	17.2	63.27	0.02	3.41
33	cis-13,16-Docosadienoic Acid Methyl Ester	C22:2 cis 13,16	13.3	63.81	0.02	3.46
34	Tricosanoic Acid Methyl Ester	C23:0	4.0	63.98	0.02	3.48
35	Lignoceric Acid Methyl Ester	C24:0	37.1	65.74	0.02	3.79
36	Nervonic Acid Methyl Ester	C24:1 cis 15	8.6	66.13	0.02	3.98
37	cis-4,7,10,13,16,19-Docosahexaenoic Acid Methyl Ester	C22:6 cis 4,7,10,13,16,19	29.7	67.25	0.02	4.28

### Summary of Results

This application note demonstrates reproducible retention of FAME components on traditional 100 meter selectivity. **Figure 1** represents complete separation of 37 component FAME mixture on a traditional selectivity after an hour bakeout at 280 °C. **Table 1** presents the precision data for both peak area and retention time on a ZB-FAME 100 meter column. The %RSD for retention time for six replicate injection is less than 0.5% which demonstrates intact stationary phase in ZB-FAME. With conventional high cyano coated FAME analysis columns, this type of reproducibility is not feasible. **Figure 2** represents the overlapped chromatograms of 6 injections with negligible difference in peak area and retention time. On contrary, **Figure 3** shows separation of FAME on an Agilent® J&W HP-88 GC column that lost resolution of peak 33 and 34 after an hour bakeout at 280 °C.

For fast FAME analysis, ZB-FAME is also available in 60, 30, and 20-meter formats. Due to similar phase volume ratio among dimensions, any FAME method can be seamlessly transferred to achieve same resolution, yet a shorter run time. **Figure 4** presents separation of 37 component FAME on a 30 meter column with run time as short as 30 min. Further, with a 20 meter dimension (**Figure 5**), a run time of 11 min was achievable with baseline resolution of all the peaks.

**Figure 4: Separation of 37 Component FAME Mix on a 30-meter ZB-FAME GC Column**




App ID 23821

### GC-FID Conditions

<b>Column</b>	Zebtron ZB-FAME
<b>Dimension:</b>	30 meter x 0.25 mm x 0.20 µm
<b>Part No.:</b>	<a href="#">7HG-G033-10</a>
<b>Recommended Z-Guard™:</b>	<a href="#">7AG-G000-00-GZK</a>
<b>Injection:</b>	Split (50:1) @ 240 °C, 1 µL
<b>Recommended Liner:</b>	Zebtron PLUS Liner Compatible with Agilent GC Instrument, Single Taper with Wool
<b>Liner Part No.:</b>	<a href="#">AG2-OA11-05</a>
<b>Carrier Gas:</b>	Hydrogen @ 1.2 mL/min (constant flow)
<b>Oven Program:</b>	100 °C for 2 min, 140 °C @ 10 °C/min, 190 °C @ 3 °C/min, 260 °C @ 30°C/min for 2 min
<b>Detector:</b>	FID
<b>Temperature:</b>	260 °C
<b>Analytes:</b>	see details in Table 1

**Figure 5: Separation of 37 Component FAME Mix on a 20-meter ZB-FAME GC Column**



App ID 23839

<b>Column</b>	Zebtron ZB-FAME
<b>Dimension:</b>	20 meter x 0.18 mm x 0.15 µm
<b>Part No.:</b>	<a href="#">7FD-G033-05</a>
<b>Recommended Z-Guard™:</b>	<a href="#">7AG-G000-00-GZK</a>
<b>Injection:</b>	Split (100:1) @ 250 °C, 1 µL
<b>Recommended Liner:</b>	Zebtron PLUS Liner Compatible with Agilent GC Instrument, Single Taper with Wool
<b>Liner Part No.:</b>	<a href="#">AG2-0A11-05</a>
<b>Carrier Gas:</b>	Helium @ 1.0 mL/min (constant flow)
<b>Oven Program:</b>	80 °C for 1.5 min, 160 °C @ 40 °C/min, 185 °C @ 5 °C/min, 260 °C @ 30°C/min
<b>Detector:</b>	FID
<b>Temperature:</b>	260 °C
<b>Analytes:</b>	see details in Table 1

## Conclusion

Zebtron™ ZB-FAME GC column provides optimal separation of cis/trans FAME components. Due to extensive crosslinking through ESC™ process, this column can withstand higher temperature limit of 280 °C and provide reproducible retention of FAME compounds. ZB-FAME is available in traditional 100 meter dimension as well as in fast analysis dimensions for reproducible FAME profiling.

# APPLICATIONS

## Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/LiveChat](http://www.phenomenex.com/LiveChat).

<p><b>Australia</b> t: +61 (0)2-9428-6444 auiinfo@phenomenex.com</p>	<p><b>India</b> t: +91 (0)40-3012 2400 indiainfo@phenomenex.com</p>	<p><b>Portugal</b> t: +351 221 450 488 ptinfo@phenomenex.com</p>
<p><b>Austria</b> t: +43 (0)1-319-1301 anfrage@phenomenex.com</p>	<p><b>Ireland</b> t: +353 (0)1 247 5405 eireinfo@phenomenex.com</p>	<p><b>Singapore</b> t: +65 800-852-3944 sginfo@phenomenex.com</p>
<p><b>Belgium</b> t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com</p>	<p><b>Italy</b> t: +39 051 6327511 italiainfo@phenomenex.com</p>	<p><b>Spain</b> t: +34 91-413-8613 espinfo@phenomenex.com</p>
<p><b>Canada</b> t: +1 (800) 543-3681 info@phenomenex.com</p>	<p><b>Luxembourg</b> t: +31 (0)30-2418700 nlinfo@phenomenex.com</p>	<p><b>Sweden</b> t: +46 (0)8 611 6950 nordicinfo@phenomenex.com</p>
<p><b>China</b> t: +86 400-606-8099 cninfo@phenomenex.com</p>	<p><b>Mexico</b> t: 01-800-844-5226 tecnicomx@phenomenex.com</p>	<p><b>Switzerland</b> t: +41 (0)61 692 20 20 swissinfo@phenomenex.com</p>
<p><b>Denmark</b> t: +45 4824 8048 nordicinfo@phenomenex.com</p>	<p><b>The Netherlands</b> t: +31 (0)30-2418700 nlinfo@phenomenex.com</p>	<p><b>Taiwan</b> t: +886 (0) 0801-49-1246 twinfo@phenomenex.com</p>
<p><b>Finland</b> t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com</p>	<p><b>New Zealand</b> t: +64 (0)9-4780951 nzinfo@phenomenex.com</p>	<p><b>United Kingdom</b> t: +44 (0)1625-501367 ukinfo@phenomenex.com</p>
<p><b>France</b> t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com</p>	<p><b>Norway</b> t: +47 810 02 005 nordicinfo@phenomenex.com</p>	<p><b>USA</b> t: +1 (310) 212-0555 info@phenomenex.com</p>
<p><b>Germany</b> t: +49 (0)6021-58830-0 anfrage@phenomenex.com</p>	<p><b>Poland</b> t: +48 22 104 21 72 pl-info@phenomenex.com</p>	<p>🌐 <b>All other countries/regions</b> <b>Corporate Office USA</b> t: +1 (310) 212-0555 info@phenomenex.com</p>

## www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™**  
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.  
[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

### Trademarks

Zebron, Z-Guard, ESC (Engineered Self Cross-linking), and BE-HAPPY are trademarks of Phenomenex.. Agilent is a registered trademark of Agilent Technologies, Inc.

### Disclaimer

Comparative separations may not be representative of all applications.

Phenomenex is in no way affiliated with Agilent.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.