

## LC-MS/MS Separation of Diclofenac and its Metabolite 4'-Hydroxydiclofenac using Luna® Omega Polar C18 and Kinetex® Biphenyl

Roxana Eggleston-Rangel, Dr. Bryan Tackett, and Anania Assefa Amenu  
Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501-1430, USA


### Overview

Diclofenac is a phenylacetic acid derivative and non-steroidal anti-inflammatory drug. Diclofenac is completely absorbed from the GI tract but undergoes significant metabolism. The primary and prevalent metabolite is 4'-Hydroxydiclofenac, which has significantly lower pharmacologic activity than Diclofenac. This drives the interest in developing a reversed phase HPLC method able to separate and analyze Diclofenac and its metabolites successfully and in a short amount of time.

This application illustrates the rapid separation of Diclofenac and 4'-Hydroxydiclofenac using Luna Omega Polar C18 and Kinetex Biphenyl.

In this application note we utilized these different stationary phases to highlight the selectivity differences and the benefits to be realized in selecting an appropriate stationary phase chemistry for the separation challenge.

Under the same mobile phase conditions, we can see that the Kinetex Biphenyl provided increased retention for the Diclofenac and its hydroxy metabolite through a mixture of pi-pi and polar interactions (**Figure 1b**); and the higher efficiency provided by the core-shell particle morphology also results in narrower peaks and increased MS sensitivity.



### LC-UV Conditions

**Columns:** Luna Omega 1.6 µm Polar C18 ([OOB-4748-AN](#))  
Kinetex 2.6 µm Biphenyl ([OOB-4622-AN](#))

**Dimensions:** 50 x 2.1 mm

**Mobile Phase:** A = 0.1 % Formic Acid in Water  
B = 0.1 % Formic Acid in Acetonitrile

**Injection:** 1 µL

**Gradient:**

Time (min)	%B
0	5
4	95
5	95
5.1	5
7	5

**Pressure (bar):** 400 (Luna Omega 1.6 µm Polar C18)  
300 (Kinetex 2.6 µm Biphenyl)

**Flow Rate:** 0.4 mL/min

**Column Temp.:** 40 °C

**Detection:** MS/MS

**Detector:** SCIEX® 4500

**System:** Agilent® 1260 Infinity UHPLC



### MRM Transitions

Peak #	Analyte	Q1	Q3	CXP
1	Diclofenac-1	294	250	-10
1	Diclofenac-2	294	214	-8
2	4-Hydroxydiclofenac-1	310	230	-13
2	4-Hydroxydiclofenac-2	310	194	-13




Figure 1. Results and Observations

(a) Luna® Omega 1.6 µm Polar C18



The Luna Omega Polar C18 incorporates a polar surface treatment to offer increased retention for polar analytes, in addition to hydrophobic retention via the C18 bonded phase.

(b) Kinetex® 2.6 µm Biphenyl



The Kinetex Biphenyl stationary phase offers multiple potential interactions to increase retention for polar aromatic compounds. These include pi-pi interactions and increased hydrogen bond accepting capacity.



**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/LiveChat](http://www.phenomenex.com/LiveChat).

**Australia**

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 5010 9707  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 800-852-3944  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
info@phenomenex.com

**🌐 All other countries/regions  
Corporate Office USA**

t: +1 (310) 212-0555  
info@phenomenex.com

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™  
GUARANTEE**

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Luna and Kinetex are registered trademarks and BE-HAPPY is a trademark of Phenomenex.

Agilent is a registered trademark of Agilent Technologies, Inc. SCIEX is a registered trademark of AB SCIEX Pte. Ltd.

**Disclaimer**

Comparative separations may not be representative of all applications.

Phenomenex is in no way affiliated with Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.

