

APPLICATIONS

LC-HRMS Analysis of Haloacetic Acids (HAAs) and Dalapon in Water using a Luna[®] Omega Polar C18 3 μm 100 x 4.6 mm Column

Cristina Postigo¹ and Dirk Hansen²

¹Institute of Environmental Assessment and Water Research (IDAEA-CSIC), C/Jordi Girona 18-26, 08034 Barcelona, Spain

²Phenomenex Ltd, Zeppelinstraße 5, 63741 Aschaffenburg, Germany

Overview

Haloacetic Acids (HAAs) are small polar contaminants in disinfected water. One major source for these molecules is the chlorination of drinking water. The respective synthetic pathways are supported by low pH conditions and high organic content in the water. Because of the negative health effects of HAAs when consumed over an extended period at elevated concentrations, their concentration in drinking water is regulated in various countries, and consequently must be monitored strictly.

The use of LC-MS/MS to determine HAAs in water has gained attention in recent years. This is linked to advances in LC column technology, that allowed the development of silica-based stationary phases with increased affinity for polar compounds like the Luna Omega Polar C18 used for this application.

In this application note we are presenting a LC-HRMS method for the effective monitoring of 13 HAAs and dalapon in drinking water.

Experimental Conditions

The analyses were performed on a SCIEX[®] ExionLC[™] system equipped with a Shimadzu[®] FCV -11AL reservoir selection valve and 0.5 mL injection loop connected to a SCIEX X500R QTOF mass spectrometer.

HPLC Parameters

Column: Luna Omega Polar C18 3 μm

Dimensions: 100 x 4.6 mm

Part No.: [00D-4760-E0](#)

Guard Column: SecurityGuard[™] Polar C18 4 x 3 mm

Part No.: [AJ0-7601](#)

Mobile Phase: A: 0.1 % Formic acid in water

B: 0.1 % Formic acid in ACN

Gradient: Time [min] %B

0	5
3	5
5	20
7	100
10	100
11.5	5
16	5

Flow Rate: 1 mL/min

Injection: 200 μL acidified water sample (0.1% FA)

Detection: MSD

Mass Spectrometer Parameters

Mass acquisition was performed using negative electrospray ionization and a high resolution multiple reaction monitoring (MRM^{HR}) workflow.

Table 1. Ion Transitions and Compound Dependent Parameters

Analyte	TOF MS	Parent Ion	Fragment Ion	CE
Chloroacetic Acid (CAA)	92.9743	92.97	34.9693	-15
Bromoacetic Acid (BAA)	136.9238	136.92	78.9189	-15
Iodoacetic Acid (IAA)	184.9099	184.91	126.9050	-15
Dichloroacetic Acid (DCAA)	126.9353	126.93	82.9461	-10
Dibromoacetic Acid (DBAA)	216.8330	216.83	172.8430	-15
Diiodoacetic Acid (DIAA)	310.8070	310.81	266.8172	-10
Bromo-chloroacetic Acid (BCAA)	170.8848	170.88	78.9189	-25
Bromo-iodoacetic Acid (BIAA)	262.8204	262.82	126.9050	-35
Chloro-iodoacetic Acid (CIAA)	218.8709	218.87	126.9050	-20
Trichloroacetic Acid (TCAA)	116.9067	116.91	34.9693	-10
Tribromoacetic Acid (TBAA)	250.7536	250.75	78.9189	-20
Dichloro-bromoacetic Acid (DCBAA)	162.8540	162.85	78.9189	-10
Dibromo-chloroacetic Acid (DBCBA)	206.8034	206.80	78.9189	-15
2,2 Dichloropropionic Acid (DPN)	140.9511	150.95	34.9693	-25
2,3-Dibromopropanoic Acid (ISTD)		78.92	78.9192	-20

Table 2. MS Source Parameters

Source Parameter	Setting
Source Temp [°C]	650
Curtain Gas	Nitrogen
Curtain Gas Pressure [psi]	45
Ion Spray Voltage [V]	-3,500

Chromatograms and Retention Times

Table 3. Retention Times

Analyte	t _R	Analyte	t _R
CAA	2.4	CIAA	3.1
BAA	2.9	TCAA	4.2
IAA	4.6	TBAA	6.5
DCAA	2.0	DCBAA	4.8
DBAA	2.7	DBCBA	5.7
DIAA	6.5	DPN	4.3
BCAA	2.3	ISTD	7.7
BIAA	3.9		

Figure 1a. Tap water fortified with only Cl containing HAAs




Figure 1b. Tap water fortified with only Br containing HAAs




Figure 1c. Tap water fortified with I containing HAAs





Figure 1d. Tap water fortified with dalapon and Br and Cl containing HAAs



Results and Discussion

The presented reversed phase method allows a fast and reliable quantification of HAAs in drinking water samples under gradient conditions. The use of the Luna Omega Polar C18 phase provides the retention required to minimize ion suppression effects. Therefore, it is possible to inject the water samples directly with only a minimal sample pretreatment (internal standard addition and acidification). No sample extraction process is needed to monitor these compounds at concentrations below the regulated levels. The column used for this study has been used for the analysis of 480 injections of various waters including fortified tap water, groundwater and surface water without any loss of separation performance.

Acknowledgements

Special thanks to Cristina Postigo at the Institute of Environmental Assessment and Water Research in Barcelona for contributing this method.

APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

<p>Australia t: +61 (0)2-9428-6444 auiinfo@phenomenex.com</p>	<p>India t: +91 (0)40-3012 2400 indiainfo@phenomenex.com</p>	<p>Portugal t: +351 221 450 488 ptinfo@phenomenex.com</p>
<p>Austria t: +43 (0)1-319-1301 anfrage@phenomenex.com</p>	<p>Ireland t: +353 (0)1 247 5405 eirinfo@phenomenex.com</p>	<p>Singapore t: +65 800-852-3944 sginfo@phenomenex.com</p>
<p>Belgium t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com</p>	<p>Italy t: +39 051 6327511 italiainfo@phenomenex.com</p>	<p>Spain t: +34 91-413-8613 espinfo@phenomenex.com</p>
<p>Canada t: +1 (800) 543-3681 info@phenomenex.com</p>	<p>Luxembourg t: +31 (0)30-2418700 nlinfo@phenomenex.com</p>	<p>Sweden t: +46 (0)8 611 6950 nordicinfo@phenomenex.com</p>
<p>China t: +86 400-606-8099 cninfo@phenomenex.com</p>	<p>Mexico t: 01-800-844-5226 tecnicomx@phenomenex.com</p>	<p>Switzerland t: +41 (0)61 692 20 20 swissinfo@phenomenex.com</p>
<p>Denmark t: +45 4824 8048 nordicinfo@phenomenex.com</p>	<p>The Netherlands t: +31 (0)30-2418700 nlinfo@phenomenex.com</p>	<p>Taiwan t: +886 (0) 0801-49-1246 twinfo@phenomenex.com</p>
<p>Finland t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com</p>	<p>New Zealand t: +64 (0)9-4780951 nzinfo@phenomenex.com</p>	<p>United Kingdom t: +44 (0)1625-501367 ukinfo@phenomenex.com</p>
<p>France t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com</p>	<p>Norway t: +47 810 02 005 nordicinfo@phenomenex.com</p>	<p>USA t: +1 (310) 212-0555 info@phenomenex.com</p>
<p>Germany t: +49 (0)6021-58830-0 anfrage@phenomenex.com</p>	<p>Poland t: +48 22 104 21 72 pl-info@phenomenex.com</p>	<p>🌐 All other countries/regions Corporate Office USA t: +1 (310) 212-0555 info@phenomenex.com</p>

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Luna is a registered trademark and SecurityGuard is a trademark of Phenomenex. ExionLC is a trademark and SCIEX is a registered trademark of AB SCIEX Pte. Ltd. AB SCIEX is used under license. Shimadzu is a registered trademark of Shimadzu Corp.

Disclaimer

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362.

CAUTION: this patent only applies to the analytical-sized guard cartridge holder,

and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridge

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.