

# APPLICATIONS

## Analysis of Antibody Drug Conjugates using High Efficiency Gel Filtration Columns

Michael Klein, Brian Rivera, and Michael McGinley  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Gel filtration chromatography is a common method used for the analysis of macromolecules, including antibody drug conjugates (ADCs). This presents unique challenges since ADCs are more hydrophobic than their parent monoclonal antibodies (mAbs). In this study, we investigate the feasibility of using a non-toxic ADC model with a Yarra<sup>TM</sup> 3  $\mu$ m SEC-3000 column for aggregate analysis.

### Introduction

Antibody drug conjugates (ADCs) are immunoconjugates designed for delivery of a cytotoxic agent, or payload, directly to tumor cells. ADCs are composed of a payload covalently bonded to the monoclonal antibody via lysine, cysteine, and unnatural amino acid residues. Site-specific linker chemistries generally produce homogenous conjugate species. However, depending on the hydrophobicity of the cytotoxic drug and the linker chemistry, heterogeneous species can still form and may have more of a propensity for aggregation<sup>1</sup>.

The most common analytical method for aggregate analysis of ADCs and other monoclonal antibody therapeutics is Gel Filtration Chromatography (GFC). Adjustments to existing GFC methods, such as the minor addition of IPA to the mobile phase, are typically made for acceptable peak shape and chromatography<sup>2</sup>.

As reported by Wagner-Rousset and colleagues<sup>3</sup>, an antibody fluorophore conjugate can be used as an alternative to ADCs for demonstrating analytical method development, including GFC. In this study, we used a non-toxic ADC model and investigated the feasibility of GFC method development using Yarra 3  $\mu$ m SEC-3000 HPLC columns.

### Materials and Methods

FITC-IgG was purchased from Fitzgerald Industries International (Acton, MA, USA). The instrument used for all GFC separations was an Agilent<sup>®</sup> 1290 LC system (Agilent Technologies, Palo Alto, CA, USA) with an upper pressure limit of 1200 bar, equipped with a binary pump, autosampler, and UV-Vis Detector. All separations used a Yarra 3  $\mu$ m SEC-3000 300 x 7.8 mm HPLC column (p/n 00H-4513-K0) obtained from Phenomenex (Torrance, CA, USA). Mobile phases consisted of low salt buffer solutions with or without an IPA modifier (noted in each chromatogram).

### Results and Discussion


Using a low salt buffer (100 mM Sodium Phosphate, pH 6.8), Yarra 3  $\mu$ m SEC-3000 gave acceptable peak shape of monomer ADC and separation of dimer and HMW aggregate is observed (Figure 1).



**Michael Klein**


Michael Klein is a brand manager for HPLC and UHPLC products. He enjoys spending time on the beach pretending to be a decent beach volleyball player.

**Figure 1.** hIgG-FITC in Low Salt Buffer



The addition of 10% IPA did not significantly change the profile, and with this non-toxic ADC model it was not necessary for acceptable chromatography (Figure 2). This also indicates that the surface of the Yarra particles is highly inert and the addition of IPA is not necessary to reduce adsorption or obtain sharper peak shape of the larger aggregates.

**Figure 2.** hIgG-FITC in Low Salt Buffer + 10% IPA



# APPLICATIONS

## Conclusion

GFC for aggregate analysis is a common application for monoclonal antibodies and is particularly necessary for ADC analysis since conjugation can cause a propensity for aggregation. Typically, modifications must be made to existing GFC methods because of the hydrophobic nature of ADCs.

In this study, we use a non-toxic ADC model to demonstrate that aggregate analysis could be performed on Yarra™ 3 µm SEC-3000 HPLC column. No modification to the existing GFC method running conditions was necessary to obtain acceptable results. This might prove useful for applications requiring screening of various ADC isoforms.

## References

1. Junutula, et al. "Site-specific Conjugation of a Cytotoxic Drug to an Antibody Improves the Therapeutic Index." *Nature Biotechnology* 26.8 (2008): 925-32
2. Wakankar, Aditya, Yan Chen, Yatin Gokarn, and Fredric S. Jacobson. "Analytical Methods for Physicochemical Characterization of Antibody Drug Conjugates." *MAbs* 3.2 (2011): 161-72
3. Wagner-Rousset, et al. "Antibody-drug Conjugate Model Fast Characterization by LC-MS following IdeS Proteolytic Digestion." *MAbs* 6.1 (2014)

## Yarra Column Ordering Information

Yarra 3 µm SEC Columns (mm)	Narrow Bore	Analytical	Analytical	SecurityGuard™ Cartridges (mm)
Phases	300 x 4.6	150 x 7.8	300 x 7.8	4 x 3.0*
Yarra 3 µm SEC-2000	00H-4512-E0	00F-4512-K0	00H-4512-K0	AJ0-4487
Yarra 3 µm SEC-3000	00H-4513-E0	00F-4513-K0	00H-4513-K0	AJ0-4488
Yarra 3 µm SEC-4000	00H-4514-E0	00F-4514-K0	00H-4514-K0	AJ0-4489

\*SecurityGuard™ Analytical Cartridges require holder, Part No.: KJ0-4282

for ID: 4.6 - 7.8 mm

### Australia

t: +61 (0)2-9428-6444  
f: +61 (0)2-9428-6445  
auinfo@phenomenex.com

### Austria

t: +43 (0)1-319-1301  
f: +43 (0)1-319-1300  
anfrage@phenomenex.com

### Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
f: +31 (0)30-2383749  
beinfo@phenomenex.com

### Canada

t: +1 (800) 543-3681  
f: +1 (310) 328-7768  
info@phenomenex.com

### Denmark

t: +45 4824 8048  
f: +45 4810 6265  
nordicinfo@phenomenex.com

### Finland

t: +358 (0)9 4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com

### France

t: +33 (0)1 30 09 21 10  
f: +33 (0)1 30 09 21 11  
franceinfo@phenomenex.com

### Germany

t: +49 (0)6021-58830-0  
f: +49 (0)6021-58830-11  
anfrage@phenomenex.com

### India

t: +91 (0)40-3012 2400  
f: +91 (0)40-3012 2411  
indiainfo@phenomenex.com

### Ireland

t: +353 (0)1 247 5405  
f: +44 1625-501796  
eireinfo@phenomenex.com

### Italy

t: +39 051 6327511  
f: +39 051 6327555  
italiainfo@phenomenex.com

## www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

### Luxembourg

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

### Mexico

t: 001-800-844-5226  
f: 001-310-328-7768  
tecnicomx@phenomenex.com

### The Netherlands

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

### New Zealand

t: +64 (0)9-4780951  
f: +64 (0)9-4780952  
nzinfo@phenomenex.com

### Norway

t: +47 810 02 005  
f: +45 4810 6265  
nordicinfo@phenomenex.com

### Puerto Rico

t: +1 (800) 541-HPLC  
f: +1 (310) 328-7768  
info@phenomenex.com

### Sweden

t: +46 (0)8 611 6950  
f: +45 4810 6265  
nordicinfo@phenomenex.com

### United Kingdom

t: +44 (0)1625-501367  
f: +44 (0)1625-501796  
ukinfo@phenomenex.com

### USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com

### All other countries Corporate Office USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com



If Yarra analytical columns do not provide at least an equivalent or better separation as compared to competing column with similar dimension, phase, and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

#### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

#### Trademarks

Agilent is a registered trademark of Agilent Technologies, Inc. Yarra and SecurityGuard are trademarks of Phenomenex.

#### Disclaimer

Comparative separations may not be representative of all applications.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362

**CAUTION:** this patent only applies to the analytical-sized guard cartridge holder, and does apply to SemiPrep, Prep, or ULTRA holders, or to any cartridges.

© 2015 Phenomenex, Inc. All rights reserved.