


# APPLICATIONS

## Deep Proteome Coverage on HeLa Lysate using a bioZen™ Fractionation Column for High pH (Basic Reversed Phase) Fractionation in Combination with a bioZen Nano Peptide XB-C18 Column for Nano Flow LC-MS

Roxana Eggleston-Rangel, Dr. Jason Anspach, Brian Rivera, and Dr. Bryan Tackett  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



**Roxana Eggleston-Rangel**  
Application Scientist

Roxana likes to spend time with her dog and family. She has a German Shepperd mix named Cobi who was named after the official mascot of the 1992 Summer Olympics and not the basketball player. Besides the lab, you might find Roxana in old episodes of the Ghost Whisperer, MAD TV, The L word and others as she used to be a TV extra during her school years.

### Introduction


Basic reversed phase (pH > 9) or high pH fractionation is a strategy commonly used to increase proteome coverage and its main purpose is to maximize protein identification in complex sample matrices. This methodology generally includes a high pH reversed phase fractionation step together with a concatenation strategy followed by nano low pH reversed phase LC-MS for maximum sensitivity (**Figure 1**). The purpose of the concatenated strategy is to produce evenly distributed ion chromatograms that reduce analysis time while increasing proteome identification<sup>1</sup>. High pH fractionation can be problematic due to silica instability at high pH (>8) that results in degradation by dissolution and by the nonspecific adsorption of peptides into the column's hardware.

bioZen Fractionation columns offer a unique combination of titanium column hardware to minimize nonspecific adsorption and TWIN-NX™ Technology, a patented organo-silica grafting process which incorporates highly stabilizing ethylene crosslinking. These organic groups are evenly incorporated into the grafted layers on the silica surface while maintaining a pure silica core. This not only provides resistance to high pH attack but also maintains the high efficiency and mechanical strength of a silica particle<sup>2</sup> (**Figures 2 and 3**). In this application note we show how to achieve deep proteome coverage on HeLa cell lysate using a high pH reversed phase fractionation strategy in combination with nano flow LC-MS.


### Figure 1.

Example of sensitivity increase on a 20 peptide mixture (SCIEX® PepCalMix) when using a bioZen nano LC column (250 nL/min) (**Left**) as opposed to micro flow (10 µL/min) on a Kinetex® micro LC column (**Right**).

#### bioZen Nano Peptide XB-C18 Column



#### Kinetex XB-C18 Micro Column



**Column:** Kinetex 2.6 µm XB-C18  
bioZen 2.6 µm Peptide XB-C18  
**Dimensions:** 150 x 0.3 mm  
150 x 0.075 mm  
**Part No.:** [OOF-4496-AC](#)  
[OOF-4768-AW-21](#)  
**Pressure (bar):** 200 bar  
**Mobile Phase:** A: 0.1 % Formic Acid in Water  
B: 0.1 % Formic Acid in Acetonitrile  
**Gradient:**


Time (min)	%B
0	3
30	40
35	70
40	80
41	3
45	3

**Flow Rate:** Kinetex: 4 µL/min  
bioZen: 250 nL/min  
**Temperature:** 25 °C  
**LC System:** NanoLC™ 425 (SCIEX)  
**Detection:** MS/MS  
**Detector:** 6500 QTRAP® (SCIEX)  
**Injection Volume:** 2.5 µL  
**µg on Column:** 1 µL


**Figure 2.**

The TWIN-NX™ Technology in the bioZen™ Fractionation column contains crosslinking that creates a high pH stable surface.

### Multi-Point Ligand Attachment Resists Low pH Ligand Cleavage



### Ethylene Crosslinking Resists High pH Attack





## Results and Discussion

High pH reversed phase fractionation UV response shows a mostly even separation of peptides through the chromatographic gradient (**Figure 4, Table 1**). Similarly, concatenated pooled fractions resulted in even distribution of peptides in the ion chromatogram, with normalized intensities higher than 1e9 throughout the nano-LC gradient, which results in higher identification of peptides and proteins (**Figure 5**). Comparing the number of identified peptides in the unfractionated HeLa lysate digest to the fractionated samples shows a significant increase in the number of peptides identified by fractionation/concatenation (**Figure 6**). This significant increase is attributed to the fact that fractionation reduces the complexity of a sample while increasing peak capacity. Additionally, the bioZen fractionation column with its titanium hardware is expected to minimize nonspecific adsorption resulting in higher peptide identifications than regular stainless steel hardware (**Figure 3**). Normally, a higher number of peptide identifications translates into a higher number of protein identifications. As shown in **Figure 7**, a significantly higher number of proteins were identified in the fractionated vs. the unfractionated sample. Besides the significant increase of protein and peptide identifications, using the bioZen 2.6 μm Peptide XB-C18 Nano LC column can be beneficial when injecting complex samples for this type of workflow. The core-shell technology in the bioZen 2.6 μm Peptide XB-C18 Nano LC column results in sharper peaks in comparison to columns packed with fully porous particles, while having a lower backpressure. Sharper peaks mean higher intensity and consequently higher peptide/protein identifications<sup>3</sup>. The lower backpressure helps to increase the column lifetime and allows for method development flexibility in terms of flow rate and temperature, aspects that are nearly impossible when using 1.6-1.9 μm particle sizes or smaller.


**Figure 3.**

**Left.** Peak area vs. injection number comparing titanium vs. stainless steel hardware for a SEC separation of γ-globulin and ovalbumin. In this case, SEC is similar to basic reversed phase chromatography in the fact that the basic environment prevents the displacement of metal ions in the hardware, which can result in peak tailing and analyte adsorption. **Right.** BioTi™ technology minimizes the need for priming with the titanium infused biocompatible hardware and frit that does not interfere with protein or peptide integrity.



**Figure 4.**

UV trace of a 30 minute basic pH reversed phase separation using 100 µg of HeLa tryptic digest on a bioZen™ High pH Fractionation column.



**Fractionation LC Conditions**

**Column:** bioZen 3 µm High pH Fractionation Column  
**Dimensions:** 150 x 2.1 mm  
**Part No.:** 00F-4793-AN  
**Pressure (bar):** 150 bar  
**Mobile Phase:** A: 10 mM Ammonium formate in Water  
 B: 10 mM Ammonium formate in 90 % Acetonitrile and 10% Water

Gradient:	Time (min)	%B
	0	1
	1	1
	25	25
	27	60
	28	70
	33	70
	34	1

**Flow Rate:** 300 µL/min  
**Temperature:** 50 °C  
**LC System:** Vanquish™ Flex UHPLC  
**Detection:** UV @ 280 nm  
**Injection Volume:** 100 µL  
**µg on Column:** 100 µg

**Fractionation Concatenation Strategy**

A total of 32 fractions resulted from the high pH reversed phase fractionation, concatenation of fractions followed as shown in the table below, resulting in a total of 8 fractions to be analyzed by nano flow LC-MS.


Concatenated Fraction Number
1 (pooled fractions 1, 9, 17, 25)
2 (pooled fractions 2, 10, 18, 26)
3 (pooled fractions 3, 11, 19, 27)
4 (pooled fractions 4, 12, 20, 28)
5 (pooled fractions 5, 13, 21, 29)
6 (pooled fractions 6, 14, 22, 30)
7 (pooled fractions 7, 15, 23, 31)
8 (pooled fractions 8, 16, 24, 32)

**Table 1.**

Amount of HeLa tryptic digest in each of the pooled fractions. Amount was determined using BCA assay. Each fraction was dissolved to a concentration of 400 ng/µL in LC-MS solvent A and 2.5 µL was injected onto the bioZen 2.6 µm Peptide XB-C18 Nano LC column for nano flow LC-MS analysis.

Pooled Fraction Number	1	2	3	4	5	6	7	8
Total µg	15	19	32	31	14	15	16	15

**Figure 5.**  
Representative total ion chromatograms for fractions 1, 3 and 7



### LC Conditions

**Column:** bioZen™ 2.6 µm Peptide XB-C18 Nano LC  
**Dimensions:** 250 x 0.075 mm  
**Part No.:** [00G-4768-AW-21](#)  
**Pressure (bar):** 350 bar  
**Mobile Phase:** A: 0.2% Formic Acid in 98 % Water, 2% Acetonitrile  
 B: 0.2% Formic Acid in 80 % Acetonitrile, 20 % Water

Gradient:	Time (min)	%B
	0	2
	7.5	6
	90	25
	120	40
	121	98
	133	98

**Flow Rate:** 500 nL/min  
**Temperature:** 50 °C  
**LC System:** EASY-nLC™ 1200  
**Detection:** nanoESI  
**Detector:** Orbitrap™ HF  
**Injection Volume:** 2.5 µL  
**µg on Column:** 1 µg


### MS Conditions

**Ion Source:** Nanospray Flex™ Ion Source (NG)  
**Ion Source Parameters:**  
**Voltage:** 2.5kV  
**Ion Mode:** Positive  
**MS:** Q Exactive™ HF


**MS Parameters:**  
**Acquisition Mode:** Data Dependent  
**MS1 Resolution:** 60,000  
**MS1 AGC Target:** 3e6  
**MS1 Injection Time:** 15ms  
**Scan Range:** [375-1500](#) m/z  
**MS2 Resolution:** 30,000  
**MS2 AGC Target:** 1e5  
**MS2 INJECTION TIME:** 45ms  
**Collision Dissociation Type:** HCD  
**Collision Energy:** 28(N)

**Figure 6.**

Number of peptide identifications for each fraction, unique peptides found in all fractions (combined), and in unfractionated HeLa tryptic digest.

**Figure 7.**

Number of protein identifications for each fraction, unique proteins found in all fractions (combined), and in unfractionated HeLa tryptic digest.



## Conclusions

The current application shows that deep proteome coverage can be achieved by combining two strategies that are commonly used to improve peptide identification. By combining offline high pH fractionation, using the bioZen™ 3 µm High pH Fractionation column, with the highly sensitive bioZen 2.6 µm Peptide XB-C18 Nano LC column for nano LC-MS, a 70% increase in peptide identifications and a close to 50% increase in protein identifications was accomplished. For users in need of obtaining high proteome coverage, this strategy is strongly encouraged.

## References

- 1) Bath, T. S. (2017). An Optimized Shotgun Strategy for the Rapid Generation of Comprehensive Human Proteomes. *Cell Syst.*, 587-599.
- 2) Anspach J, Rivera B, Srinivasa Rao. (2018, June 1st). LCGC. Retrieved from [www.chromatographyonline.com](http://www.chromatographyonline.com): <https://www.chromatographyonline.com/view/bioinert-versus-biocompatible-benefits-different-column-materials-liquid-chromatography-separations>
- 3) Joshua J. Coon1, A. S. (2016). Now, More Than Ever, Proteomics Needs Better Chromatography. *Cell Systems*, 321-324.

## Ordering Information

### bioZen Nano LC Columns with Integrated SecurityLINK™ Fitting

Phases	150 x .075	250 x .075	500 x .075
bioZen 3 µm Peptide PS-C18	<a href="#">00F-4771-AW-21</a>	<a href="#">00G-4771-AW-21</a>	–
bioZen 2.6 µm Peptide XB-C18	<a href="#">00F-4768-AW-21</a>	<a href="#">00G-4768-AW-21</a>	–
bioZen 3 µm Polar C18	<a href="#">00F-4782-AW-21</a>	<a href="#">00G-4782-AW-21</a>	–
bioZen 5 µm Peptide XB-C18	–	–	<a href="#">00J-4792-AW-11</a>

### bioZen Nano LC Columns with Open Fused-Silica Inlet/Outlet Fitting

Phases	150 x .075	250 x .075	500 x .075
bioZen 3 µm Peptide PS-C18	<a href="#">00F-4771-AW-11</a>	<a href="#">00G-4771-AW-11</a>	–
bioZen 2.6 µm Peptide XB-C18	<a href="#">00F-4768-AW-11</a>	<a href="#">00G-4768-AW-11</a>	–
bioZen 3 µm Polar C18	<a href="#">00F-4782-AW-11</a>	<a href="#">00G-4782-AW-11</a>	–
bioZen 5 µm Peptide XB-C18	–	–	<a href="#">00J-4792-AW-21</a>

### bioZen High pH Fractionation Column

Fractionation Column		
Part No.	Description	Dimension
<a href="#">00F-4793-AN</a>	bioZen 3 µm High pH Fractionation Column	150 x 2.1 mm

### Nano Trap Columns

Trap Columns		
	10 x 0.075 mm	Unit
RP-1	<a href="#">05N-4252-AW</a>	3/pk
RP-2	<a href="#">05N-4754-AW</a>	3/pk

### Trap Fittings

Trap Fittings		
Part No.	Description	Unit
<a href="#">AQO-7602</a>	PEEKlok™ fittings with 6-40 thread for 1/32" OD tubing (2 x fittings, 6 x ferrules and 1 x tightening tool)	ea
<a href="#">AQO-7603</a>	PEEKlok fittings with 6-32 thread for 1/32" OD tubing (2 x fittings, 6 x ferrules and 1 x tightening tool)	ea
<a href="#">AQO-7600</a>	PEEKlok fittings with 10-32 thread for 1/32" OD tubing with low profile hex head (2 x fittings, 6 x ferrules and 1 x wrench)	ea

**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

### Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/ChatNow](http://www.phenomenex.com/ChatNow).

#### Australia

t: +61 (0)2-9428-6444  
auserinfo@phenomenex.com

#### Austria

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

#### Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

#### Canada

t: +1 (800) 543-3681  
info@phenomenex.com

#### China

t: +86 400-606-8099  
cninfo@phenomenex.com

#### Czech Republic

t: +420 272 017 077  
cz-info@phenomenex.com

#### Denmark

t: +45 4824 8048  
nordicinfo@phenomenex.com

#### Finland

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

#### France

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

#### Germany

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

#### India

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

#### Ireland

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

#### Italy

t: +39 051 6327511  
italiainfo@phenomenex.com

#### Luxembourg

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

#### Mexico

t: 01-800-844-5226  
tecnicomx@phenomenex.com

#### The Netherlands

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

#### New Zealand

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

#### Norway

t: +47 810 02 005  
nordicinfo@phenomenex.com

#### Poland

t: +48 22 104 21 72  
pl-info@phenomenex.com

#### Portugal

t: +351 221 450 488  
ptinfo@phenomenex.com

#### Singapore

t: +65 800-852-3944  
sginfo@phenomenex.com

#### Slovakia

t: +420 272 017 077  
sk-info@phenomenex.com

#### Spain

t: +34 91-413-8613  
espinfo@phenomenex.com

#### Sweden

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

#### Switzerland

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

#### Taiwan

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

#### United Kingdom

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

#### USA

t: +1 (310) 212-0555  
info@phenomenex.com

#### 🌐 All other countries/regions

**Corporate Office USA**  
t: +1 (310) 212-0555  
info@phenomenex.com



### www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

#### Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

#### Trademarks

Kinetex is a registered trademark and bioZen, TWIN-NX, SecurityLINK, and BioTi are trademarks of Phenomenex. PeekLok is a trademark of Trajan Scientific Australia Pty Ltd. Vanquish, EASY-nLC, Orbitrap, Nanospray Flex, and Q Exactive are trademarks of Thermo Fisher Scientific. SCIEX and QTRAP are registered trademarks and NanoLC is a trademark of AB SCIEX Pte. Ltd. AB SCIEX™ is being used under license.

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.