

APPLICATIONS

Analysis of IdeS Digested Monoclonal Antibody Fragments

Ivan Lebedev and Brian Rivera

Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501, USA

Overview

Monoclonal antibodies (mAbs) are well-established therapeutics, with a variety of different analytical methods associated with purity analysis. An analytical technique for mAb purity is reversed phase LC (RPLC). This is a primary technique utilized for mAbs because of its relatively short analysis time, with high resolution, and ability to separate hydrophobic variants such as oxidation, glycoforms, and lysine variants. Another utility of reversed phase is to analyze fragments, with traditional methods using a reducing agent like dithiothreitol (DTT) to reduce interchain disulfides, yielding heavy chain and light chain.

Although the analysis of heavy and light chains is useful, the use of site-specific cysteine proteases is also another approach that allows for better characterization of the antibody. Further, this can be combined with a reducing agent to further identify and quantitate mAb heterogeneity. One particularly useful partial digestion of mAbs is using IdeS, a cysteine protease that cleaves below the hinge region of the IgG. This enzyme generates a (Fab')₂ fragment and two Fc fragments (Fc/2). This allows for better characterization of the Fc fragment.

Figure 1 shows an overlay of intact and IdeS digested infliximab. Again, this digest yields the (Fab')₂ fragment and two single chain Fc fragments. At the intact level, there is not discrimination of lysine variants, which are more hydrophilic thus earlier eluting by reversed phase. After digestion, the lysine variant on the Fc of infliximab is clearly identified, as shown in the **Figure 1** inset. Interestingly, the (Fab')₂ fragment elutes later than the full-length antibody, which is somewhat counterintuitive as the fragment has a lower molar mass. Additionally, the chromatographic profile looks similar to the intact, with a partial post peak. This post peak could be deamidated variants in the heavy chain.¹

In instances where IdeS digested fragments are not sufficient, one can perform a subsequent reduction after digestion, yielding smaller fragments. This process reduces (Fab')₂ into an Fd' and light chain. Again, these peaks might correspond to deamidated variants, though further analytical characterization by orthogonal methods, such as peptide mapping, would need to be performed to confirm peak identity.

In summary, the analysis of antibody fragments can provide key insights into sample heterogeneity. Separation of variants, such as lysine and deamidated variants, is possible by digestion using a site-specific cysteine protease, such as IdeS, and by using a bioZen™ 2.6 µm WidePore C4 LC column. By combining the enzymatic digestion with a subsequent reduction, further characterization of both Fc/2 and Fd' fragments allows for better characterization of both variable and conserved regions of the antibody.


LC Conditions

Column:	bioZen 2.6 µm WidePore C4
Dimensions:	100 x 2.1 mm
Part No.:	OOD-4786-AN
Mobile Phase A:	0.1 % TFA in Water
Mobile Phase B:	0.1 % TFA in Acetonitrile
Gradient Program:	30-40% B in 5 minutes (Fig 1-2)
Flow Rate:	0.8 mL/min
Temperature:	80°C
Detection:	UV @ 214 nm
Injection:	Infliximab, as indicated (0.5 mg/mL), Figures 1



Figure 1. Infliximab, Intact and IdeS Digested


Overlay of intact infliximab (black trace) and fragments generated by IdeS digestion. Lysine variant (*) is shown as pre-peak in the Fc/2 fragment. Although smaller in molecular weight, the (Fab')₂ fragment elutes later than full-length mAb.



App ID 26136

Figure 2. Infliximab, IdeS Digested and Reduced

After subsequent reduction of infliximab, fragments generated include Fc/2 as well as light chain (LC) and Fd', or the hinge containing variable region of the heavy chain. This provides insight on heterogeneity of the Fab region of the antibody that the protease digestion itself could not determine.



App ID 26137

APPLICATIONS

Reference

¹Pisupati, Karthik et al. "Biosimilarity under stress: A forced degradation study of Remicade® and Remsima™." mAbs vol. 9,7 (2017): 1197-1209. doi:10.1080/19420862.2017.1347741

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

<p>Australia t: +61 (0)2-9428-6444 auinfo@phenomenex.com</p>	<p>India t: +91 (0)40-3012 2400 indiainfo@phenomenex.com</p>	<p>Portugal t: +351 221 450 488 ptinfo@phenomenex.com</p>
<p>Austria t: +43 (0)1-319-1301 anfrage@phenomenex.com</p>	<p>Ireland t: +353 (0)1 247 5405 eirinfo@phenomenex.com</p>	<p>Singapore t: +65 800-852-3944 sginfo@phenomenex.com</p>
<p>Belgium t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com</p>	<p>Italy t: +39 051 6327511 italiainfo@phenomenex.com</p>	<p>Spain t: +34 91-413-8613 espinfo@phenomenex.com</p>
<p>Canada t: +1 (800) 543-3681 info@phenomenex.com</p>	<p>Luxembourg t: +31 (0)30-2418700 nlinfo@phenomenex.com</p>	<p>Sweden t: +46 (0)8 611 6950 nordicinfo@phenomenex.com</p>
<p>China t: +86 400-606-8099 cninfo@phenomenex.com</p>	<p>Mexico t: 01-800-844-5226 tecnicomx@phenomenex.com</p>	<p>Switzerland t: +41 (0)61 692 20 20 swissinfo@phenomenex.com</p>
<p>Denmark t: +45 4824 8048 nordicinfo@phenomenex.com</p>	<p>The Netherlands t: +31 (0)30-2418700 nlinfo@phenomenex.com</p>	<p>Taiwan t: +886 (0) 0801-49-1246 twinfo@phenomenex.com</p>
<p>Finland t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com</p>	<p>New Zealand t: +64 (0)9-4780951 nzinfo@phenomenex.com</p>	<p>United Kingdom t: +44 (0)1625-501367 ukinfo@phenomenex.com</p>
<p>France t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com</p>	<p>Norway t: +47 810 02 005 nordicinfo@phenomenex.com</p>	<p>USA t: +1 (310) 212-0555 info@phenomenex.com</p>
<p>Germany t: +49 (0)6021-58830-0 anfrage@phenomenex.com</p>	<p>Poland t: +48 22 104 21 72 pl-info@phenomenex.com</p>	<p>🌐 All other countries/regions Corporate Office USA t: +1 (310) 212-0555 info@phenomenex.com</p>

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

bioZen and Be-HAPPY are trademarks of Phenomenex..

Comparative separations may not be representative of all applications.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.