

TN-1306

Demonstrating Utility of the Phenomenex Online UHPLC/HPLC Method Transfer Tool Using the Luna® Omega C18 Columns

Zeshan Aqeel, Phil Koerner, PhD, and Bryan Tackett, PhD
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Introduction

The analysis of an ABN7 test mix was used to assess the scalability for Luna Omega C18 columns across a range of particle sizes (1.6 µm, 3 µm, and 5 µm) in three different dimensions. In this technical note, we show how the Phenomenex online UHPLC/HPLC Method Transfer Tool (<https://www.phenomenex.com/Tools/KinetexMethodTranslation>) was utilized to determine method parameters when scaling from one particle and column size to another.

LC-UV Conditions

- Column:** Luna Omega 1.6 µm C18
Luna Omega 3 µm C18
Luna Omega 5 µm C18
- Dimension:** 50 x 2.1 mm
150 x 4.6 mm
250 x 4.6 mm
- Part No.:** [00B-4742-AN](#)
[00F-4784-E0](#)
[00G-4785-E0](#)
- Mobile Phase:** A: Water
B: Acetonitrile
- Gradient:** See Figures
- Flow Rate:** See Figures
- Injection Volume:** 1 µL, 50 x 2.1 mm column
5 µL, 150 x 4.6 mm column
5 µL, 250 x 4.6 mm column
- Temperature:** 30 °C
- Detection:** UV @ 254 nm
- LC System:** Agilent® 1260 Infinity
- Analytes:**
 1. Uracil (1.0 mg/mL in Methanol)
 2. Pindolol (5.0 mg/mL in Methanol)
 3. Chlorpheniramine (5.0 mg/mL in Methanol)
 4. Nortriptyline (5.0 mg/mL in Methanol)
 5. 3-Methyl, 4-Nitrobenzoic acid (5.0 mg/mL in Methanol)
 6. 2-Hydroxy, 5-Methylbenzaldehyde (5.0 mg/mL in Methanol)
 7. Hexanophenone (5.0 mg/mL in Methanol)

Results and Discussion

Figure 1 shows the results of ABN7 run on a Luna Omega 1.6 µm C18, 50 x 2.1 mm column. The benefits of the smaller particle size and shorter column are apparent; baseline separation (resolution >2) of all analytes is achieved in less than 6 minutes. This highlights one of the advantages of UHPLC - the separation is achieved in a shorter time.

The Method Transfer Tool can be used to properly scale the gradient times and flow rate when scaling this separation from a Luna Omega 1.6 µm C18, 50 x 2.1 mm column to a Luna Omega 3 µm C18, 150 x 4.6 mm column (**Table 1**). While the analysis time has increased, the chromatographic resolution between the analytes is maintained due to the gradient and flow rate being scaled with the Method Transfer Tool. These new parameters generated the chromatogram displayed in **Figure 2**.

The Method Transfer Tool can now be used to further scale the gradient times and flow rate for a Luna Omega 5 µm C18, 250 x 4.6 mm column (**Table 2**). The resulting chromatogram is shown in **Figure 3**. Again, while the overall analysis time has increased, the chromatographic resolution between the analytes has been maintained due to the scaling of the gradient and flow rate.

The Method Transfer Tool can also be used for making parameter predictions when using the same column on two different instrument types. **Table 3** is an example of using a Luna Omega 5 µm C18, 250 x 4.6 mm column on a HPLC instrument with a larger delay volume but wanting to move it to a UHPLC instrument with a smaller delay volume. The red circle indicates the differences in the gradient table that results.

This tool can aid in the method transfer not only between different particle sizes and column dimensions, but also between new and older generation instruments having different system volumes. **Table 4** is an example of the same scaling that occurred for **Figure 2**; however, the instrument type was changed from a UHPLC to a HPLC. The red circle indicates the differences in the gradient table that results.




Zeshan Aqeel
Senior Application Scientist

Aside from the lab being his favorite place to be, Zeshan enjoys playing vintage videogames with his twin boys and loves every minute of reliving parenthood with his baby girl.



Figure 1. Luna® Omega 1.6 µm C18, 50 x 2.1 mm



LC Column 1				
Length (mm)	Column ID (mm)	Flow Rate (mL/min)	Particle Type / Size (µm)	Porosity
50	2.1	0.4	1.6 FP	0.660

Original Gradient Table				
Time (min)	%A	%B	#CV	%B/CV
0	95	5	3.94	
0.5	95	5	1.75	
5.5	5	95	17.50	5.14
6.6	5	95	3.85	
7	95	5	1.40	
10	95	5		

App ID 26802

Peak No.	Analyte	Retention Time (min)	k'	Resolution
1	Uracil	0.39	-	-
2	Pindolol	2.56	5.55	62.1
3	Chlorpheniramine	3.02	6.73	14.5
4	Nortriptyline	3.74	8.57	21.4
5	3-Methyl, 4-Nitrobenzoic acid	3.84	8.81	2.7
6	2-Hydroxy, 5-Methylbenzaldehyde	4.14	9.59	8.6
7	Hexanophenone	5.42	12.86	33.6

Table 1. Scaling Calculator Predictions – 1.6 µm to 3 µm


LC Instrument 1		Predicted Total Pressure (bar)	LC Instrument 2		Predicted Total Pressure (bar)
Pressure Limit (bar) =	600	392	Pressure Limit (bar) =	600	519
Delay Volume (mL) =	0.75		Delay Volume (mL) =	0.75	
Instrument Type =	UHPLC		Instrument Type =	UHPLC	
Maximum Viscosity =	1.01E-03	ACN / H2O 25C	Separation Goal =	Geometric Scaling	

LC Column 1					LC Column 2				
Length (mm)	Column ID (mm)	Flow Rate (mL/min)	Particle Type / Size (µm)	Porosity	Length (mm)	Column ID (mm)	Flow Rate (mL/min)	Particle Type / Size (µm)	Porosity
50	2.1	0.4	1.6 FP	0.660	150	4.6	1.92	3 FF	0.660

Original Gradient Table					Scaled Gradient Table				
Time (min)	%A	%B	#CV	%B/CV	Time	%A	%B	#CV	%B/CV
0	95	5	6.56		0.00	95	5	0.46	
0.5	95	5	1.75		6.73	95	5	7.86	
5.5	5	95	17.50	5.14	21.73	5	95	17.50	5.14
6.6	5	95	3.85		25.03	5	95	3.85	
7	95	5	1.40		26.23	95	5	1.40	
10	95	5			30.05	95	5	4.46	



Figure 2. Luna® Omega 3 µm C18, 150 x 4.6 mm



LC Column 1				
Length (mm)	Column ID (mm)	Flow Rate (mL/min)	Particle Type / Size (µm)	Porosity
150	4.6	1.92	3 FP	0.660
Original Gradient Table				
Time (min)	%A	%B	#CV	%B/CV
0	95	5	0.46	
6.73	95	5	7.85	
21.73	5	95		5.14
25.03	5	95	3.85	
26.23	95	5	1.40	
30.05	95	5	4.46	

App ID 26803

Peak No.	Analyte	Retention Time (min)	k'	Resolution
1	Uracil	1.15	-	-
2	Pindolol	10.04	7.72	200
3	Chlorpheniramine	11.36	8.87	19.89
4	Nortriptyline	13.66	10.87	34.68
5	3-Methyl, 4-Nitrobenzoic acid	13.90	11.07	3.37
6	2-Hydroxy, 5-Methylbenzaldehyde	14.70	11.77	10.45
7	Hexanophenone	18.49	15.06	49.25

Table 2. Scaling Calculator Predictions – 3 µm to 5 µm


LC Instrument 1		Predicted Total Pressure (bar)	LC Instrument 2		Predicted Total Pressure (bar)				
Pressure Limit (bar) =	600	520	Pressure Limit (bar) =	1000	402				
Delay Volume (mL) =	0.75		Delay Volume (mL) =	0.75					
Instrument Type =	UHPLC		Instrument Type =	UHPLC					
Maximum Viscosity =	1.01E-03	ACN / H2O 25C	Separation Goal =	Geometric Scaling					
LC Column 1					LC Column 2				
Length (mm)	Column ID (mm)	Flow Rate (mL/min)	Particle Type / Size (µm)	Porosity	Length (mm)	Column ID (mm)	Flow Rate (mL/min)	Particle Type / Size (µm)	Porosity
150	4.6	1.92	3 FP	0.660	250	4.6	1.92	5 FF	0.660
Original Gradient Table					Scaled Gradient Table				
Time (min)	%A	%B	#CV	%B/CV	Time	%A	%B	#CV	%B/CV
0	95	5	0.46		0.00	95	5	0.27	
6.73	95	5	7.85		11.48	95	5	8.04	
21.73	5	95		5.14	36.48	5	95	17.50	5.14
25.03	5	95	3.85		41.98	5	95	3.85	
26.23	95	5	1.40		43.98	95	5	1.40	
30.05	95	5	4.46		50.08	95	5	4.27	



LC Column 1				
Length (mm)	Column ID (mm)	Flow Rate (mL/min)	Particle Type / Size (µm)	Porosity
250	4.6	1.92	5 FP	0.660

Original Gradient Table				
Time (min)	%A	%B	#CV	%B/CV
0	95	5	0.27	
11.48	95	5	8.04	
36.48	5	95	17.50	5.14
41.98	5	95	3.85	
43.98	95	5	1.40	
50.08	95	5	4.27	

Figure 3. Luna® Omega 5 µm C18, 250 x 4.6 mm



App ID 26804

Peak No.	Analyte	Retention Time (min)	k'	Resolution
1	Uracil	1.87	-	-
2	Pindolol	16.96	8.08	200
3	Chlorpheniramine	19.11	9.24	21.01
4	Nortriptyline	23.02	11.33	38.12
5	3-Methyl, 4-Nitrobenzoic acid	23.35	11.51	3.02
6	2-Hydroxy, 5-Methylbenzaldehyde	24.71	12.24	10.83
7	Hexanophenone	30.94	15.57	47.35

Table 3. Scaling Calculator Predictions – HPLC to UHPLC, Same Column

LC Instrument 1		Predicted Total Pressure (bar)	LC Instrument 2		Predicted Total Pressure (bar)
Pressure Limit (bar) =	600	401	Pressure Limit (bar) =	600	520
Delay Volume (mL) =	1.45		Delay Volume (mL) =	0.75	
Instrument Type =	HPLC		Instrument Type =	UHPLC	
Maximum Viscosity =	1.01E-03	ACN / H2O 25C	Separation Goal =	Geometric Scaling	

LC Column 1					LC Column 2				
Length (mm)	Column ID (mm)	Flow Rate (mL/min)	Particle Type / Size (µm)	Porosity	Length (mm)	Column ID (mm)	Flow Rate (mL/min)	Particle Type / Size (µm)	Porosity
150	4.6	1.92	3 FP	0.660	150	4.6	1.92	3 FF	0.660

Original Gradient Table					Scaled Gradient Table				
Time (min)	%A	%B	#CV	%B/CV	Time	%A	%B	#CV	%B/CV
0	95	5	0.88		0.00	95	5	0.46	
6.73	95	5	7.85		7.09	95	5	8.28	
21.73	5	95		5.14	22.09	5	95	17.50	5.14
25.03	5	95	3.85		25.39	5	95	3.85	
26.23	95	5	1.40		26.59	95	5	1.40	
30.05	95	5	4.46		30.41	95	5	4.46	



Table 4. Scaling Calculator Predictions – Scaling UHPLC to HPLC, 1.6 μm to 3 μm

LC Instrument 1		Predicted Total Pressure (bar)	LC Instrument 2		Predicted Total Pressure (bar)
Pressure Limit (bar) =	600	392	Pressure Limit (bar) =	600	401
Delay Volume (mL) =	0.75		Delay Volume (mL) =	1.45	
Instrument Type =	UHPLC		Instrument Type =	HPLC	
Maximum Viscosity =	1.01E-03	ACN / H2O 25C	Separation Goal =	Geometric Scaling	

LC Column 1					LC Column 2				
Length (mm)	Column ID (mm)	Flow Rate (mL/min)	Particle Type / Size (μm)	Porosity	Length (mm)	Column ID (mm)	Flow Rate (mL/min)	Particle Type / Size (μm)	Porosity
50	2.1	0.4	1.6 FP	0.660	150	4.6	1.92	3 FF	0.660

Original Gradient Table					Scaled Gradient Table				
Time (min)	%A	%B	#CV	%B/CV	Time	%A	%B	#CV	%B/CV
0	95	5	6.56		0.00	95	5	0.88	
0.5	95	5	1.75		6.37	95	5	7.43	
5.5	5	95	17.50	5.14	21.37	5	95	17.50	5.14
6.6	5	95	3.85		24.67	5	95	3.85	
7	95	5	1.40		25.87	95	5	1.40	
10	95	5			30.05	95	5	4.88	

Table 5. Summary of Data From all Three Columns

Peak No.	Analyte	Luna® Omega 1.6 μm C18, 50 x 2.1 mm			Luna Omega 3 μm C18, 150 x 4.6 mm			Luna Omega 5 μm C18, 250 x 4.6 mm		
		Retention Time (min)	k'	Resolution	Retention Time (min)	k'	Resolution	Retention Time (min)	k'	Resolution
1	Uracil	0.39	-	-	1.15	-	-	1.87	-	-
2	Pindolol	2.56	5.55	62.1	10.04	7.72	200	16.96	8.08	200
3	Chlorpheniramine	3.02	6.73	14.5	11.36	8.87	19.89	19.11	9.24	21.01
4	Nortriptyline	3.74	8.57	21.4	13.66	10.87	34.68	23.02	11.33	38.12
5	3-Methyl, 4-Nitrobenzoic acid	3.84	8.81	2.7	13.90	11.07	3.37	23.35	11.51	3.02
6	2-Hydroxy, 5-Methylbenzaldehyde	4.14	9.59	8.6	14.70	11.77	10.45	24.71	12.24	10.83
7	Hexanophenone	5.42	12.86	33.6	18.49	15.06	49.25	30.94	15.57	47.35



Conclusion

The Method Transfer Tool illustrated in this technical note allows for scaling methods for transfer from HPLC to UHPLC, or vice versa, while preserving chromatographic resolution. The data summarized in **Table 5** demonstrate that scaling the gradient times and flow rate will result in an equivalent separation when there is a need or desire to change the particle size and dimensions of the column used. This tool can also be used when it is necessary to move from a UHPLC to a HPLC instrument. For example, when a method is developed on UHPLC but transferred to

routine use on HPLC due to instrument or column constraints. Alternatively, when an existing HPLC method is desired to transfer to UHPLC to take advantage of reduced analysis times while maintaining overall chromatographic resolution. The data here represents the accuracy of geometric scaling between Luna® Omega C18 columns with different particle sizes and column dimensions using the UHPLC/HPLC Method Transfer Tool.

Luna Omega Ordering Information

1.6 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges† (mm)
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
Polar C18	00A-4748-AN	00B-4748-AN	00D-4748-AN	00F-4748-AN	AJ0-9505
PS C18	00A-4752-AN	00B-4752-AN	00D-4752-AN	00F-4752-AN	AJ0-9508
C18	00A-4742-AN	00B-4742-AN	00D-4742-AN	00F-4742-AN	AJ0-9502

for 2.1 mm ID

3 µm Analytical Columns (mm)					SecurityGuard Cartridges (mm)
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0* /10 pk
Polar C18	00B-4760-E0	00D-4760-E0	00F-4760-E0	00G-4760-E0	AJ0-7601
PS C18	00B-4758-E0	00D-4758-E0	00F-4758-E0	00G-4758-E0	AJ0-7606
C18	00B-4784-E0	00D-4784-E0	00F-4784-E0	00G-4784-E0	AJ0-7612
SUGAR	—	00D-4775-E0	00F-4775-E0	00G-4775-E0	AJ0-4495

for ID: 3.2-8.0 mm

5 µm Analytical Columns (mm)					SecurityGuard Cartridges (mm)
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0* /10 pk
Polar C18	00B-4754-E0	00D-4754-E0	00F-4754-E0	00G-4754-E0	AJ0-7601
PS C18	00B-4753-E0	00D-4753-E0	00F-4753-E0	00G-4753-E0	AJ0-7606
C18	00B-4785-E0	00D-4785-E0	00F-4785-E0	00G-4785-E0	AJ0-7612

for ID: 3.2-8.0 mm

† SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#)
 *SecurityGuard™ Analytical Cartridges require holder, Part No.: [KJ0-4282](#)



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: 0120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 All other countries/regions

Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Luna is a registered trademark and SecurityGuard and BE-HAPPY are trademarks of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Agilent Technologies, Inc.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362

CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP, or ULTRA holders, or to any cartridges.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.

