

LC-MS/MS Separation of MDMA and its Metabolites Using the Kinetex® 2.6 µm Biphenyl Column


Roxana Eggleston-Rangel and Dr. Bryan Tackett
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

MDMA (3,4-Methylenedioxyamphetamine) is a ring-substituted amphetamine analog with reinforcing psychoactive effects, in addition to neurotoxic potential toward central 5-HT neurons. It is metabolized into MDA (3,4-methylenedioxyamphetamine), HMMA (4-hydroxy-3-methoxymethamphetamine), and HMA (4-hydroxy-3-methoxyamphetamine) (**Figure 1**). Being able to successfully separate these drug substances is a useful tool for forensic labs. In this application note, we present an LC-MS/MS method to separate MDMA and its metabolites by utilizing the Kinetex 2.6 µm Biphenyl column.

The Kinetex Biphenyl column provided separation for MDMA and its metabolites through a mixture of pi-pi and polar interactions; and the higher efficiency provided by the core-shell particle morphology also results in narrow peaks and increased MS sensitivity (**Figure 2**).

Figure 1. MDMA Metabolism



LC-MS/MS Conditions

Column: Kinetex 2.6 µm Biphenyl
Dimension: 50 x 2.1 mm
Part No.: [00B-4622-AN](#)
Mobile Phase: A: 0.1 % Formic Acid in Water
 B: 0.1 % Formic Acid in Methanol

Gradient: Time (min)	%B
0	5
4	95
5	95
5.1	5
7	5


Flow Rate: 400 µL/min
Injection Volume: 1 µL
Temperature: 40 °C
LC System: Agilent® 1260 Infinity
Detection: MS/MS
Detector: SCIEX® 4500

MRM Transitions

Peak No.	Analyte	Q1 (m/z)	Q3 (m/z)
1	HMA-1	182	165
1	HMA-2	182	105
2	HMMA-1	196	165
2	HMMA-2	196	105
3	MDA-1	180	163.1
3	MDA-2	180	135
4	MDMA-1	194.1	163.2
4	MDMA-2	194.1	135.2



Figure 2. Kinetex® 2.6 µm Biphenyl



App ID 26780



The Kinetex Biphenyl stationary phase offers multiple potential interactions to increase retention for polar aromatic compounds. These include pi-pi interactions and increased hydrogen bond accepting capacity.

Peak No.	Analyte	Retention Time (min)	Concentration (mg/mL)
1	HMA	1.78	100
2	HMMA	2.09	100
3	MDA	2.77	100
4	MDMA	3.01	100



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions
Corporate Office USA**

t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

**BE-HAPPY™
GUARANTEE**

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex is a registered trademark and BE-HAPPY is a trademark of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc. SCIEX is a registered trademarks of AB SCIEX Pte. Ltd.

Disclaimer

Phenomenex is in no way affiliated with Agilent Technologies, Inc.
FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.
© 2021 Phenomenex, Inc. All rights reserved.

