

APPLICATIONS

Effect of Temperature on Reversed Phase Separations of Monoclonal Antibodies

Cody Warhurst, Dawn Chen, Brian Rivera
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501

The analysis of intact monoclonal antibodies (mAbs) by reversed phase HPLC is a common technique for assessing protein variation, such as clipping and heterogeneity due to post-translational modifications. Whether for LC-UV impurity analysis or for confirming primary sequence using LC-MS, the analysis of intact mAbs can provide useful insight on a protein therapeutic.

As with small molecule analysis, proteins run by reversed phase will elute in order of hydrophobicity; i.e. more hydrophobic proteins will retain longer. However, unlike small molecules, the secondary and tertiary structure of the protein impact the separation and the so-called “native” structure of proteins is disrupted by typical reversed phase conditions. Acidic modifier and organic solvent partially dissociate the protein, leading to the exposure of hydrophobic side chains, and potentially contributing to peak broadening. To overcome peak broadening, a [bioZen™ 2.6 μm WidePore C4](#) LC Column was utilized in combination with method parameter enhancements.

Although these denaturing conditions are inherent to the detect, one particularly useful method parameter to adjust is column temperature. Modulation in temperature modulates the protein secondary and tertiary structure, which can be utilized to ultimately effect peak shape and even selectivity; presumably, slight changes in protein conformation affect the interaction of hydrophobic amino acid residues in the primary sequence with the moderately hydrophobic stationary phase.

Figure 1 shows improvements in peak shape with for intact infliximab, when comparing 70 °C to 80 °C column temperature. Interestingly, we observe a decrease in performance with further increases to 90 °C. Upon further investigation, this increase in peak width is due to partial separation of another variant, observed as a later eluting shoulder in **Figure 2**.

The effect of temperature is more clearly observed in subunit analysis of mAbs (**Figure 3**). Increase in temperature improves peak width at half height when comparing 70 °C and 80 °C, with improvements in peak widths for light chain. However, peak widths are slightly worse when comparing 80 °C and 90 °C.

More interesting are the impurity profiles for the heavy chain, which vary between each method column temperature (**Figure 4**). Depending on the importance in characterizing each impurity, one profile may be more desirable than another.

In summary, increases in temperature can lead to improvements in peak shape in large molecule reversed phase LC methods, as seen with the [bioZen WidePore C4](#). Further, selectivity needs to be assessed separately as the increases in temperature can be detrimental to the separation. At the intact level, one may observe partial separation of another impurity, which could lead to overall peak broadening but a more detailed impurity profile. At the subunit level, differences in impurity profile can be observed and as such, temperature should be implemented as part of the design of experiment.

LC Conditions

Column:	bioZen 2.6 μm WidePore C4
Catalog No.:	00D-4786-AN
Dimensions:	100 x 2.1 mm
Mobile Phase A :	0.1 % TFA in Water
Mobile Phase B:	0.1 % TFA in Acetonitrile
Gradient:	20-55 % B in 6 minutes (<i>Fig 1,2</i>) 25-45 % B in 6 minutes (<i>Fig 3,4</i>)
Flow Rate:	0.5 mL/min
Detection:	280 nm (<i>Fig 1,2</i>) 214 nm (<i>Fig 3,4</i>)
Temperature:	<i>As Indicated</i>
Sample:	Infliximab (1 mg/mL) NIST mAb (1 mg/mL)

Figure 1:

Stacked chromatogram showing the effect of temperature on peak widths at half height for intact infliximab. Optimal peak widths are observed at 80 °C.




Figure 2:

Detailed overlay of intact infliximab, showing a nominal improvement in impurity profile at 90 °C.




Figure 3:

Chromatogram stack of reduced NIST mAb. Light chain show improved peak width at half height run at the 80 °C when compared to 70 °C, though no improvement as temperature is increased to 90 °C.





Figure 4:

Chromatogram overlay of heavy chain of NIST mAb, showing differences in impurity profile. Both 70 °C and 80 °C show improved separation of later eluting shoulder impurities.



APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eirinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY[™]
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

bioZen is a trademark of Phenomenex.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.