

Determination of Volatile Organic Compounds in Environmental Samples as per HJ639 using Zebron™ ZB-624PLUS™ GC Column

Anne Jurek¹, Bryan Tackett², Ramkumar Dhandapani², and Richard Jack²

¹EST Analytical Inc., 503 Commercial Drive, Fairfield, OH 45014 USA

²Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

A common technique for the sampling of Volatile Organic Compounds (VOCs) in water is purge and trap technique. This technique is technically referred to as dynamic headspace extraction as it is an exhaustive extraction and preconcentration technique that allows low level detection of volatile compounds. Purge and Trap in conjunction with a GC-MS is a very powerful analytical technique for volatile organic compounds (VOCs) and BTEX compounds in environmental samples. Due to the expansive list of analytes, the method requires a GC-MS for accurate identification and quantification of analytes that may have chromatographic coelutions.

A 624 selectivity is commonly used for the GC separation. However, traditional GC columns with 624 selectivity are not mass spectrometer compatible and can cause excessive bleed that can contaminate the mass spec source. Presented in this application note are the separation of VOCs using Zebron ZB-624PLUS, which is a GC-MS compatible low bleed stationary phase (**Figure 1**). This methodology is adopted from regulation HJ639-2012 and retention data is presented in **Table 1**. The method demonstrated in this application note is precise with less than 4% RSD for peak area, linearity, and accurate with % recovery ranging from 93-110% at concentration levels 5 µg/L and 50 µg/L (**Table 2**). As represented in **Figure 2**, ZB-624PLUS provides optimal separation of over 60 VOCs within 27 minutes analysis time. Due to low bleed and MS compatibility of ZB-624PLUS stationary phase through Engineered Self Cross-linking™ (ESC™), the column provides steady baseline that facilitates detection of VOC as low as 1 ppb (**Figure 3**).

Figure 1. Benefits of Zebron™ ZB-624PLUS for VOC Analysis





Figure 1. Volatile Organic Compounds by Purge and Trap GC-MS on Zebron™ ZB-624PLUS™ GC Column at 100 ppb



Purge and Trap Conditions

Purge and Trap Concentrator	EST Analytical Evolution 2
Trap Type	K
Valve Oven Temp.	140 °C
Transfer Line Temp.	140 °C
Trap Temp.	35 °C
Moisture Reduction Trap (MoRT) Temp.	39 °C
Purge Time	11 min
Purge Flow	40 mL/min
Dry Purge Temp.	Off
Dry Purge Flow	40 mL/min
Dry Purge Time	1.0 min
Desorb Pressure Control	On
Desorb Pressure	5 psi
Desorb Time	2.0 min
Desorb Preheat Delay	10 sec
Desorb Temp.	250 °C
Moisture Reduction Trap (MoRT) Bake Temp.	250 °C
Bake Temp	210 °C
Sparge Vessel Bake Temp.	110 °C
Bake Time	6 min
Bake Flow	40 mL/min
Purge and Trap Auto-Sampler	EST Analytical Centurion WS
Sample Type	Water
Internal Standard Vol.	5 µL
Sample Volume	5 mL

GC-MS Conditions


Column: Zebron ZB-624PLUS™
Dimension: 30 meter x 0.25 mm x 1.4 µm
Part No.: [7HG-G040-27](#)
Injection: Split 30:1 @ 220 °C
Carrier Gas: Helium @ 1.0 mL/min (constant flow)
Oven Program: 40 °C for 2 min to 120 °C @ 5 °C/min for 2 min, to 220 °C @ 10 °C/min for 2 min
Detector: GC-MS
Transfer Line Temperature: 180 °C
Mode: Scan (35-300 m/z)
Source Temperature: 230 °C
Quad Temperature: 150 °C
Solvent delay: 0.7 min

Analyte details and retention time of analytes on page 4

Data Courtesy : EST Analytical Inc. Phenomenex is not associated with EST Analytical Inc.



Figure 3. Volatile Organic Compounds by Purge and Trap GC-MS on Zebron™ ZB-624PLUS™ GC Column at 1 ppb



App ID 26319

Purge and Trap Conditions

Purge and Trap Concentrator	EST Analytical Evolution 2
Trap Type	K
Valve Oven Temp.	140 °C
Transfer Line Temp.	140 °C
Trap Temp.	35 °C
Moisture Reduction Trap (MoRT) Temp.	39 °C
Purge Time	11 min
Purge Flow	40 mL/min
Dry Purge Temp.	Off
Dry Purge Flow	40 mL/min
Dry Purge Time	1.0 min
Desorb Pressure Control	On
Desorb Pressure	5 psi
Desorb Time	2.0 min
Desorb Preheat Delay	10 sec
Desorb Temp.	250 °C
Moisture Reduction Trap (MoRT) Bake Temp.	250 °C
Bake Temp	210 °C
Sparge Vessel Bake Temp.	110 °C
Bake Time	6 min
Bake Flow	40 mL/min
Purge and Trap Auto-Sampler	EST Analytical Centurion WS
Sample Type	Water
Internal Standard Vol.	5 µL
Sample Volume	5 mL

GC-MS Conditions

Column: Zebron ZB-624PLUS™
Dimension: 30 meter x 0.25 mm x 1.4 µm
Part No.: [7HG-G040-27](#)
Injection: Split 30:1 @ 220 °C
Carrier Gas: Helium @ 1.0 mL/min (constant flow)
Oven Program: 40 °C for 2 min to 120 °C @ 5 °C/min for 2 min, to 220 °C @ 10 °C/min for 2 min
Detector: GC-MS
Transfer Line Temperature: 180 °C
Mode: Scan (35-300 m/z)
Source Temperature: 230 °C
Quad Temperature: 150 °C
Solvent delay: 0.7 min

Analyte details and retention time of analytes on page 4

Data Courtesy : EST Analytical Inc. Phenomenex is not associated with EST Analytical Inc.



Table 1. Retention time peak identification of Volatile Organic Compounds

Compound	Average Retention Time (n=7)	Compound	Average Retention Time (n=7)
Vinyl Chloride	2.02	Chlorobenzene-d5 (IS)	13.95
1,1-Dichloroethene	3.36	Chlorobenzene	14.02
Methylene Chloride	3.96	1,1,1,2-Tetrachloroethane	14.28
trans-1,2-Dichloroethene	4.32	Ethylbenzene	14.32
1,1-Dichloroethane	4.9	Xylene (m+p)	14.67
2-chloro-1,3-butadiene	5.01	Xylene (o)	15.61
cis-1,2-Dichloroethene	5.72	Styrene	15.66
2,2-Dichloropropane	5.73	Bromoform	16.04
Bromochloromethane	6.07	Isopropylbenzene	16.6
Chloroform	6.24	BFB SUR	17
Dibromofluoromethane SUR	6.48	Bromobenzene	17.29
1,1,1-Trichloroethane	6.51	1,1,2,2-Tetrachloroethane	17.44
1,1-Dichloropropene	6.77	1,2,3-Trichloropropane	17.52
Carbon Tetrachloride	6.78	n-Propylbenzene	17.71
Benzene	7.08	2-Chlorotoluene	17.87
1,2-Dichloroethane	7.15	4-Chlorotoluene	18.2
Fluorobenzene (IS)	7.59	1,3,5-Trimethylbenzene	18.25
Trichloroethene	8.26	tert-Butylbenzene	19.13
1,2-Dichloropropane	8.69	1,2,4-Trimethylbenzene	19.33
Dibromomethane	8.85	sec-Butylbenzene	19.87
Bromodichloromethane	9.22	1,3-Dichlorobenzene	20.12
1-chloro-2,3-epoxypropane	9.86	1,4-Dichlorobenzene-d4 (IS)	20.4
cis-1,3-Dichloropropene	10.11	Isopropyltoluene	20.45
Toluene-d8 SUR	10.66	1,4-Dichlorobenzene	20.47
Toluene	10.8	1,2,-Dichlorobenzene	21.49
trans-1,3-Dichloropropene	11.33	n-Butylbenzene	21.66
1,1,2-Trichloroethane	11.74	1,2-Dibromo-3-chloropropane	23.52
Tetrachloroethene	11.96	1,2,4-Trichlorobenzene	25.35
1,3-Dichloropropane	12.09	Hexachlorobutadiene	25.73
Dibromochloromethane	12.58	Naphthalene	25.81
1,2-Dibromoethane	12.8	1,2,3-Trichlorobenzene	26.27

Data Courtesy : EST Analytical Inc. Phenomenex is not associated with EST Analytical Inc.



Table 2. Precision, Accuracy and Linearity of Volatile Organic Compounds on Zebron™ ZB-624PLUS™

Compound	Calibration Curve RF	Calibration Curve Linearity (%RSD)	Calibration Curve Linearity (R ²)	Precision %RSD (5µg/L)	% Recovery (5µg/L)	Precision %RSD (50µg/L)	% Recovery (50µg/L)
Vinyl Chloride	0.469	6.10	0.999	4.30	93.34	3.49	97.97
1,1-Dichloroethene	0.348	7.69	0.999	3.20	99.11	4.46	105.19
Methylene Chloride	0.433	4.86	1.000	2.30	103.23	3.31	106.19
1,1-dichloroethane	0.772	7.32	1.000	3.60	100.49	3.52	107.51
2-chloro-1,3-butadiene	0.715	7.85	0.999	3.74	98.17	4.24	107.62
trans-1,2-Dichloroethene	0.402	6.13	0.999	3.14	101.40	3.64	107.02
cis-1,2-Dichloroethene	0.468	7.80	1.000	2.50	98.71	3.54	108.30
2,2-Dichloropropane	0.689	5.55	0.999	6.60	98.00	5.44	108.68
Bromochloromethane	0.273	6.19	1.000	2.81	100.57	2.91	104.87
Chloroform	0.779	6.94	1.000	2.27	102.49	3.49	108.26
Dibromofluoromethane SUR	0.409	8.16	1.000	2.44	100.17	3.38	108.55
1,1,1-Trichloroethane	0.683	7.01	1.000	4.26	98.69	4.21	107.11
Carbon Tetrachloride	0.572	7.88	0.999	5.49	94.34	4.55	106.46
1,1-Dichloropropene	0.601	6.66	0.999	4.62	98.40	4.59	106.80
Benzene	1.717	6.48	1.000	2.96	100.31	3.48	107.84
1,2-Dichloroethane	0.657	7.24	1.000	3.65	102.83	3.41	109.36
Trichloroethene	0.444	8.04	1.000	2.61	98.06	3.62	107.34
1-chloro-2,3-epoxypropane	0.143	6.10	1.000	5.98	89.00	3.06	97.02
1,2-Dichloropropane	0.459	7.47	1.000	2.59	100.14	3.54	109.03
Dibromomethane	0.302	8.12	1.000	2.56	98.46	3.00	106.53
Bromodichloromethane	0.625	8.49	1.000	2.36	97.74	3.40	109.88
cis-1,3-Dichloropropene	0.778	8.65	1.000	2.38	97.34	3.65	109.78
Toluene-d8 SUR	1.494	7.67	1.000	2.75	99.20	3.63	108.95
Toluene	1.095	7.56	1.000	2.50	99.23	3.63	108.98
trans-1,3-Dichloropropene	0.732	9.10	1.000	2.31	97.69	3.69	109.99
1,1,2-Trichloroethane	0.418	8.50	1.000	2.52	100.94	3.19	108.20
Tetrachloroethene	0.373	9.98	0.999	3.60	95.77	3.28	105.99
1,3-Dichloropropane	0.758	7.70	1.000	2.47	102.00	3.29	108.67
Dibromochloromethane	0.481	10.72	1.000	2.66	93.31	3.30	109.50
1,2-Dibromoethane	0.466	8.89	1.000	1.86	99.57	3.02	108.86
Chlorobenzene	1.451	7.26	1.000	2.43	100.91	3.21	107.36
1,1,1,2-Tetrachloroethane	0.542	7.74	1.000	2.54	97.69	3.06	107.12

Data Courtesy : EST Analytical Inc. Phenomenex is not associated with EST Analytical Inc.



Table 2. Precision, Accuracy, and Linearity of Volatile Organic Compounds on Zebron™ ZB-624PLUS™

Compound	Calibration Curve RF	Calibration Curve Linearity (%RSD)	Calibration Curve Linearity (R ²)	Precision %RSD (5µg/L)	% Recovery (5µg/L)	Precision %RSD (50µg/L)	% Recovery (50µg/L)
Ethylbenzene	2.545	6.56	0.999	3.15	100.34	3.44	107.37
Xylene (m+p)	3.919	6.53	0.999	3.00	101.00	3.35	107.10
Xylene (o)	2.044	7.07	1.000	2.48	100.34	3.12	107.75
Styrene	1.651	7.85	1.000	2.72	99.06	3.26	108.43
Bromoform	0.458	11.04	1.000	3.81	89.63	2.63	107.35
Isopropylbenzene	2.471	7.04	0.999	2.96	98.63	3.65	107.25
BFB SUR	0.726	7.27	1.000	2.66	100.91	3.25	107.43
Bromobenzene	1.026	6.87	1.000	2.13	102.34	3.29	107.78
1,2,3-trichloropropane	1.202	7.16	1.000	3.87	99.14	2.34	106.11
1,1,2,2-Tetrachloroethane	0.894	6.95	1.000	3.46	101.54	2.53	106.02
n-Propylbenzene	2.918	6.42	0.999	3.05	99.43	3.84	107.30
2-Chlorotoluene	0.569	6.94	1.000	2.72	100.57	3.16	107.27
4-Chlorotoluene	0.593	7.90	1.000	2.65	99.31	3.24	107.87
1,3,5-Trimethylbenzene	2.037	6.52	0.999	2.82	99.74	3.34	106.62
tert-Butylbenzene	2.018	6.85	0.999	3.43	98.20	4.01	107.02
sec-Butylbenzene	0.517	8.13	0.999	3.17	97.20	3.93	107.49
1,2,4-Trimethylbenzene	2.074	7.24	0.999	2.42	99.54	3.18	107.00
1,3-Dichlorobenzene	1.138	7.22	1.000	2.50	100.17	3.01	106.68
Isopropyltoluene	2.189	7.47	0.999	3.23	97.97	3.91	107.14
1,4,-Dichlorobenzene	2.202	7.55	1.000	2.85	99.57	3.24	106.05
1,2,-Dichlorobenzene	2.123	7.45	1.000	3.31	98.97	3.04	105.53
n-Butylbenzene	3.912	7.70	0.999	3.32	96.14	4.31	107.41
1,2-Dibromo-3-chloropropane	0.552	8.49	1.000	5.69	93.69	2.67	103.26
1,2,4-Trichlorobenzene	1.463	9.11	1.000	3.22	95.57	3.14	106.07
Naphthalene	5.006	7.20	1.000	3.15	97.97	2.80	104.39
Hexachlorobutadiene	0.604	8.45	0.999	2.99	96.14	4.58	104.33
1,2,3-Trichlorobenzene	1.389	8.64	1.000	3.46	95.89	2.94	104.75
Average (page 6 & 7)	1.166	7.55	1.000	3.16	98.58	3.46	106.94

Data Courtesy : EST Analytical Inc. Phenomenex is not associated with EST Analytical Inc.

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
 auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
 anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
 t: +32 (0)2 511 8666 (Dutch)
 beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
 info@phenomenex.com

China

t: +86 400-606-8099
 cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
 cz-info@phenomenex.com

Denmark

t: +45 4824 8048
 nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
 nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
 franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
 anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
 hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
 indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
 eireinfo@phenomenex.com

Italy

t: +39 051 6327511
 italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
 nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
 tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
 nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
 nzinfo@phenomenex.com

Norway

t: +47 810 02 005
 nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
 pl-info@phenomenex.com

Portugal

t: +351 221 450 488
 ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
 sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
 sk-info@phenomenex.com

Spain

t: +34 91-413-8613
 espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
 nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
 swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
 twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
 thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
 ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
 info@phenomenex.com

All other countries/regions

Corporate Office USA
 t: +1 (310) 212-0555
 info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPYTM

GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Zebtron, Engineered Self Cross-linking, ESC, 624PLUS, and BE-HAPPY are trademarks of Phenomenex.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.

AN22150321_W

