

# Heat and High pH Forced Degradation of Antibodies on a Biozen™ dSEC-2 Column

Tanya Napolitano, PhD, Ivan Lebedev, PhD, and Bryan Tackett, PhD  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

## Overview

The shelf life and storage considerations of various monoclonal antibodies (mAbs) heavily depend on their inherent biophysical properties. Size exclusion chromatography (SEC) provides an opportunity to elucidate the behaviors of proteinaceous material with respect to time. Temperature and alkalinity can affect the longevity of an antibody in solution. In this application note, we demonstrate the aggregation propensity of multiple antibodies in response to their environmental conditions using a Biozen dSEC-2 column.

To test high temperature conditions for mAb degradation, we heated two antibodies, Herceptin and Cetuximab, at 60°C while taking aliquots every 24 hours across 168 hours. High temperature stress is one of the most relevant conditions for providing information about potential long-term degradation at the intended storage condition. High temperature was found to accelerate the formation of aggregates (**Figures 1 and 2**) which were successfully separated from the monomer using the Biozen dSEC-2 column.

To evaluate pH, the storage solution was adjusted to pH 9 from the initial storage pH of 6-7. Infliximab and Cetuximab were each incubated under alkaline conditions at 37°C for 72 hours with aliquots sampled at 24-hour intervals. It was observed that the combination of physiological temperature as well as high pH leads to the degradation of the antibodies through cleavage of their disulfide bonds at the hinge and domain interfaces (**Figures 3 and 4**).


This application note demonstrates the applicability of the Biozen dSEC-2 column in understanding the biophysical integrity monoclonal antibodies and related proteins. Stability studies have been shown to generate impurities that feature a degree of denaturation thus exposing some hydrophobic regions of the folded protein's interior. The hydrophilic nature of the Biozen dSEC-2 stationary phase shows favorable aggregate profiles for forced degradation work under high salt conditions which is a known challenge when working with more hydrophobic hybrid particles for these studies.

## LC Conditions

- Column:** Biozen 1.8 µm dSEC-2  
**Dimension:** 300 x 4.6 mm  
**Part No.:** [00H-4787-E0](#)  
**Mobile Phase:** 200 mM Potassium Phosphate + 250 mM Potassium Chloride, pH 6.2  
**Flow Rate:** 0.35 mL/min  
**Injection** 2 µL (Herceptin)  
**Volume:** 10 µL (Cetuximab, Infliximab)  
**Temperature:** 25 °C  
**Detector:** UV @ 280 nm  
**System:** Waters® ACQUITY® UPLC H-Class  
**Samples:** 1. Herceptin (10 mg/mL)  
2. Cetuximab (2 mg/mL)  
3. Infliximab (2 mg/mL)




**Figure 1.** Herceptin Day 0, 3, and 7 Forced Degradation by Heat Analysis on a Biozen™ dSEC-2 Column.



App ID Z7178


**Figure 2.** Cetuximab Day 0, 3, and 7 Forced Degradation by Heat Analysis on a Biozen dSEC-2 Column.



App ID Z7177




**Figure 3.** Cetuximab Day 0 and 3 Forced Degradation by High pH Analysis on a Biozen™ dSEC-2 Column.



App ID 27179

**Figure 4.** Infliximab Day 0 and 3 Forced Degradation by High pH Analysis on a Biozen dSEC-2 Column.



App ID 27180



**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/Chat](http://www.phenomenex.com/Chat).

**Australia**

t: +61 (0)2-9428-6444  
auiinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 5019 9707  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 6559 4364  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

🌐 **All other countries/regions  
Corporate Office USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Biozen and BE-HAPPY are trademarks of Phenomenex. Waters and ACQUITY are registered trademarks of Waters Technologies Corporation.

**Disclaimer**

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Waters Technologies Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.

