



Fast and Optimal Extraction of PFAS Compounds Using Strata™-XL-AW

Agustin Pierri, PhD¹, Ramkumar Dhandapani, PhD², Richard Jack, PhD², and Bryan Tackett, PhD²

¹Weck Laboratories, Inc., 14859 Clark Ave., Industry, CA 91745, USA

²Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

EPA 533 recommends SPE cartridges containing weak anion exchange, mixed-mode polymeric sorbent (polymeric backbone and a diamino ligand), with a particle size of approximately 33 µm. The Solid Phase Extraction (SPE) sorbent must have a pKa above 8 for optimal extraction. In this application note we compared the recovery of PFAS compounds using a Strata-XL-AW SPE cartridge, a similar functionality but with 100 µm particle size, pore size of 300 Å, and a 500 mg/6 mL format. The EPA 533 method flexibility allows alternative SPE phases to be used. Six sets of extractions were performed on drinking water using a Strata X-AW and a Strata XL-AW SPE cartridge. The extractions performed using the Strata-XL-AW met EPA recovery requirements, while providing the advantage of faster sample loading. This reduces the sample loading time by half for clean samples (like drinking water) compared to traditional extraction methods using smaller pore and particle sizes, while maintaining % recovery, as presented in **Table 1**.

LC Conditions

Column: Gemini™ 3 µm C18

Dimensions: 50 x 2.0 mm

Part No.: [00B-4439-B0](#)

Mobile Phase: A: 20 mM Ammonium Acetate
B: Methanol

Gradient:	Time (min)	%B
	0	5
	0.5	5
	3	40
	16	80
	18	80
	20	95
	22	95
	25	5
	35	5

Flow Rate: 0.25 mL/min

Injection Volume: 2 µL

Detection: ESI-MS/MS

Sample Preparation

Step	Description
Sample Pre-treatment:	250 mL sample is fortified with isotopically labeled analogues of the method analytes.
Load:	250 mL pre-treated sample through a Strata-XL-AW 500 mg/6 mL (8B-S051-HCH) cartridge or through a Strata-X-AW 500 mg/6 mL (8B-S038-HCH) cartridge.
Wash:	Cartridges with aqueous Ammonium Acetate, followed by Methanol.
Elute:	Sample with Ammonium Hydroxide in Methanol.
Dry Down:	Under a gentle stream of Nitrogen in a heated water bath.
Reconstitute:	Sample to a final volume to 1 mL with 20 % Water in Methanol (v/v) before analyzing by LC-MS/MS.

“The Strata-XL-AW tubes result in a significant time savings during the sample loading step. The Strata-XL-AW tubes take about half the time to load clean samples, and about 20 % of the time required for dirty samples. This reduces the impact of the SPE bottleneck, allowing us to process more samples per shift.”

- Agustin Pierri, PhD,
Technical Director, Weck Laboratories



Table 1. Comparison of % Recovery of PFAS Compounds Using Strata™-XL-AW versus Strata-X-AW SPE Cartridges.

Target Analyte	Strata XL-AW			Target Analyte	Strata-X-AW		
	Avg % Recovery (n=6)	STD	%RSD		Avg % Recovery (n=6)	STD	%RSD
¹³ C ₄ -PFBA	89.23	2.47	3%	¹³ C ₄ -PFBA	91.28	3.37	4%
PFBA	97.45	3.12	3%	PFBA	97.2	4.13	4%
PFMPA	101.54	3.07	3%	PFMPA	100.62	3.40	3%
¹³ C ₅ -PFPeA	90.55	2.88	3%	¹³ C ₅ -PFPeA	92.54	3.15	3%
PFPeA	104.59	3.97	4%	PFPeA	101.93	4.51	4%
¹³ C ₃ -PFBS	88.56	2.48	3%	¹³ C ₃ -PFBS	87.13	2.03	2%
PFBS	103.54	2.87	3%	PFBS	104.89	3.93	4%
PFMBA	104.9	4.68	4%	PFMBA	104.09	4.66	4%
PFEESA	104.71	4.06	4%	PFEESA	104.92	3.91	4%
NFDHA	108.74	7.44	7%	NFDHA	112.38	7.87	7%
¹³ C ₂ -4:2 FTS	83.73	3.79	5%	¹³ C ₂ -4:2 FTS	82.47	2.88	3%
4:2 FTS	92.23	3.41	4%	4:2 FTS	91.46	3.04	3%
¹³ C ₅ -PFHxA	90.99	3.43	4%	¹³ C ₅ -PFHxA	91.35	3.67	4%
PFHxA	96.42	4.63	5%	PFHxA	95.57	4.76	5%
PFPeS	101.69	4.70	5%	PFPeS	101.71	2.16	2%
¹³ C ₃ -HFPO-DA	89.75	4.07	5%	¹³ C ₃ -HFPO-DA	93.37	2.23	2%
HFPO-DA	106.86	4.11	4%	HFPO-DA	103.16	3.27	3%
¹³ C ₄ -PFHpA	90.66	4.29	5%	¹³ C ₄ -PFHpA	90.34	3.42	4%
PFHpA	99.51	3.76	4%	PFHpA	99.45	4.33	4%
¹³ C ₃ -PFHxS	91.27	3.82	4%	¹³ C ₃ -PFHxS	91.22	4.46	5%
PFHxS	103.11	4.96	5%	PFHxS	102.05	2.73	3%
ADONA	102	3.31	3%	ADONA	102.8	3.80	4%
¹³ C ₂ -6:2 FTS	86.01	5.28	6%	¹³ C ₂ -6:2 FTS	84.55	2.58	3%
6:2 FTS	93.71	4.85	5%	6:2 FTS	94.44	3.89	4%
¹³ C ₈ -PFOA	91.4	4.12	5%	¹³ C ₈ -PFOA	91.08	3.54	4%
PFOA	102.97	3.67	4%	PFOA	102.9	4.60	4%
PFHpS	105.16	4.49	4%	PFHpS	103.84	5.91	6%
¹³ C ₉ -PFNA	89.91	4.04	4%	¹³ C ₉ -PFNA	89.96	3.05	3%
PFNA	102.76	1.98	2%	PFNA	101.27	5.45	5%
¹³ C ₈ -PFOS	88.85	4.72	5%	¹³ C ₈ -PFOS	88.02	3.13	4%
PFOS	105.93	4.61	4%	PFOS	102.97	5.48	5%
9Cl-PF3OUdS	99.65	2.93	3%	9Cl-PF3OUdS	95.64	3.81	4%
¹³ C ₂ -8:2 FTS	85.24	4.22	5%	¹³ C ₂ -8:2 FTS	83.6	2.37	3%
8:2 FTS	100.36	3.75	4%	8:2 FTS	99.39	4.86	5%
¹³ C ₆ -PFDA	89.83	5.43	6%	¹³ C ₆ -PFDA	89.62	2.55	3%
PFDA	100.45	4.10	4%	PFDA	97.12	4.22	4%
¹³ C ₇ -PFUnA	87.57	4.32	5%	¹³ C ₇ -PFUnA	84.69	3.25	4%
PFUnA	100.51	2.95	3%	PFUnA	99.46	4.25	4%
11Cl-PF3OUdS	103.01	2.32	2%	11Cl-PF3OUdS	100.1	6.08	6%
¹³ C ₂ -PFDoA	83.61	4.49	5%	¹³ C ₂ -PFDoA	82.21	2.36	3%
PFDoA	103.67	4.07	4%	PFDoA	101.45	4.30	4%



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5019 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 6559 4364
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

**🌐 All other countries/regions
Corporate Office USA**

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Strata, Gemini, and BE-HAPPY are trademarks of Phenomenex.

Disclaimer

Strata-X is patented by Phenomenex. U.S. Patent No. 7,119,145.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.

