

PAIN MANAGEMENT TESTING UPGRADED

Combine Phenomenex Sample Preparation
and Core-Shell HPLC/UHPLC Solutions for:

- Fast Results
- Reproducible High Performance
- Increased Sensitivity
- Excellent Productivity



Guaranteed Solutions

Effective and highly reproducible methodology helps you ensure results and reduce costs. With guaranteed high quality products from Phenomenex, you can easily push your research to the next level.

guarantee

We guarantee that if Phenomenex products in this guide do not provide at least an equivalent separation as compared to other products of the same phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.

Solid Phase Extraction (SPE)

Strata™-X-Drug B

- Designed and quality controlled for basic drugs of abuse
- Time and solvent saving polymeric strong cation-exchange SPE sorbent
- No equilibration or condition steps needed



Protein Precipitation (PPT)

Impact™

- Complete in plate removal of proteins (β -glucuronidase)
- Increase throughput 3x
- Leak-free oleophobic membrane
- Easily automated




New long drip plate

HPLC/UHPLC Analysis

Kinetex® Biphenyl Core-Shell Technology

- Rugged and reliable core-shell columns
- Incredible performance and efficiencies
- Enhanced polar and basic retention
- Quality controlled for basic drugs of abuse



For Research Use Only. Not for use in diagnostic procedures.

© 2015 Phenomenex, Inc. All rights reserved.

Phenomenex | WEB: www.phenomenex.com

Rapid, Simple Cleanup

for Pain Management Drugs and Drugs of Abuse

Skip the method development, significantly reduce your sample preparation time, and save solvent using Strata™-X-Drug B SPE tubes and 96-well plates. One method is all you will need to extract over 41 pain management drugs. Or easily upgrade the throughput of your existing protein precipitation methods with the new Impact™ PPT long drip plates.

Sample Pretreatment		
Urine:		Serum / Plasma / Saliva:
Acid Hydrolysis	Enzymatic Hydrolysis	
1. To 2 mL urine sample, add 500 µL concentrated hydrochloric acid (HCl) 2. Heat at 90 °C for 2 hours 3. Add 2 mL of 200 mM Sodium acetate buffer (pH 4.0) 4. Add 1 mL of 6 N Potassium hydroxide (KOH). Vortex. 5. Centrifuge for 5-6 minutes at 5,000-5,500 rpm (20-22 °C) 6. Verify pH of sample is between 4.0 – 6.0 7. Load pretreated sample directly onto Strata-X-Drug B sorbent (see procedure below)	1. Dilute 500 µL urine with 100 µL buffer* and 20 µL of 1,000,000 units/mL beta-glucuronidase in a 96-well collection plate 2. Vortex for 5-6 seconds 3. Incubate in a water bath at 63 °C for 30 minutes 4. If performing SPE, transfer sample to a 96-well collection plate or autosampler vial, seal and centrifuge for 10 minutes at 2,000 rpm. Load supernatant onto Strata-X-Drug B sorbent (see procedure below). If performing rapid removal of β-glucuronidase, load sample directly onto Impact Protein Precipitation plate (see procedure below).	1. Dilute plasma/serum 3:1 with 1 % Acetic acid 2. Vortex 3. Centrifuge for 5-6 minutes at 3,000-3,500 rpm (20-22 °C) 4. Load pretreated sample directly onto Strata-X-Drug B sorbent (see procedure below)

*Buffer prepared by adding 800 mL deionized water and 111 mL glacial acetic acid to a 1 L volumetric flask, adjust final volume to the line with a 50 % KOH solution and mix by inversion.

For Ultra Clean and Concentrated Samples

Strata-X-Drug B SPE Steps

Procedure designed for 60 mg bed mass Strata-X-Drug B and can be adjusted for smaller or larger bed masses.

Condition:	Not required
Equilibrate:	Not required
Load:	Pretreated sample (see above for pretreatment recommendations)
Wash:	2 mL of 100 mM Sodium acetate buffer (pH 5.0)
Wash 2:	2 mL Methanol
Elute:	2 mL Ethyl acetate/Isopropanol/Ammonium hydroxide (pH 5.0) (7:2:1)
Evaporate:	to dryness under a stream of Nitrogen at 50 °C
Reconstitute:	1 mL 15 % Methanol
Inject:	5 µL



For Rapid Removal of Proteins (Ex. β-Glucuronidase)


Impact Protein Precipitation Steps

1. Load 300 µL Acetonitrile to the wells of an Impact Plate
2. Load 100 µL hydrolyzed sample directly onto the Impact Plate
3. Vortex for 2 minutes at maximum possible speed
4. Apply vacuum at 2-7 inches Hg for 2-3 minutes and collect filtrate in a collection plate



Rapid, Reliable, and High Performance Pain Panel Testing

The new Kinetex® Biphenyl ensures that your comprehensive screens display the needed performance with a ruggedness and reliability that will keep you at ease.



Conditions for all separations:

Column: Kinetex 2.6 μ m Biphenyl (50 x 2.1 mm)
Kinetex 2.6 μ m Biphenyl (50 x 3.0 mm)
Kinetex 5 μ m Biphenyl (50 x 3.0 mm)

Mobile Phase: A: Water with 0.1 % Formic Acid
B: Methanol with 0.1 % Formic Acid



Gradient	Time (min)	%B
	0	10
	0.5	10
	2	25
	4.5	80
	4.51	85
	5.5	85
	5.51	10
	7	10

Flow Rate: 0.6 mL/min

Temperature: 40°C

Detection: MS/MS (AB SCIEX API 5000™)


- Samples:**
1. Morphine
 2. Oxycodone
 3. Hydromorphone
 4. Amphetamine
 5. Methamphetamine
 6. MDA
 7. Naloxone
 8. Codeine
 9. 6-MAM
 10. Oxycodone
 11. MDMA
 12. Naltrexone
 13. Hydrocodone
 14. MDEA
 15. Norfentanyl
 16. Tramadol
 17. Benzoylcegonine
 18. Meprobamate
 19. Normeperidine
 20. Meperidine
 21. Norbuprenorphine
 22. Fentanyl
 23. Buprenorphine
 24. Flurazepam
 25. PCP
 26. Carisoprodol
 27. Propoxyphene
 28. Midazolam
 29. Sufentanil
 30. Norpropoxyphene
 31. EDDP
 32. Lorazepam
 33. Methadone
 34. Clonazepam
 35. Oxazepam
 36. Hydroxylalprazolam
 37. Nordiazepam
 38. Flunitrazepam
 39. Temazepam
 40. Alprazolam
 41. Diazepam



Go Even Faster!

Using the high performance and versatility of the Kinetex® 2.6µm Biphenyl with a ballistic gradient allows you to cut down even more time from your cycles and significantly improve productivity. Imagine what you can do with this solution on your multiplex or dual-stream systems.

Kinetex 2.6µm Biphenyl 50 x 3.0 mm



Conditions:

Column: Kinetex 2.6µm Biphenyl
Dimension: 50 x 3.0 mm
Part No.: 00B-4622-Y0
Mobile Phase: A: Water with 0.1 % Formic Acid
 B: Methanol with 0.1 % Formic Acid
Gradient:

Time (min)	%B
0	10
2.5	100
3.5	100
3.51	10
5	10

Flow Rate: 0.7 mL/min
Temperature: 40°C
Detection: MS/MS (AB SCIEX API 5000™)


Analyte	Q1	Q3	Retention Time (min)
Morphine	286.2	152	1.38
Oxymorphone	302.1	227	1.44
Hydromorphone	286.2	185	1.53
Amphetamine	136.1	91	1.68
Naloxone	328.1	165.2	1.78
Methamphetamine	150.1	91.2	1.84
Codeine	300.1	152.1	1.85
MDA	180.1	133	1.86
Oxycodone	316.2	241	1.92
Naltrexone	342.2	267.1	1.94
Hydrocodone	300.1	199.1	1.96
MDMA	194.1	163.2	1.99
MDEA	208.2	163.2	2.11
Norfentanyl	233.2	84.1	2.2

Analyte	Q1	Q3	Retention Time (min)
Tramadol	264.2	58.1	2.24
Benzoylcegonine	290.2	168.1	2.3
Meperidine	248.2	220	2.33
Meprobamate	219.2	158.2	2.35
Norbuprenorphine	414.3	83.2	2.41
Fentanyl	337.2	105.1	2.63
Buprenorphine	468.3	55.2	2.63
Flurazepam	388.2	315.1	2.7
Carisoprodol	261.2	176.2	2.73
PCP	244.3	91	2.74
Propoxyphene	340.2	266.3	2.75
Sufentanil	387.2	238.2	2.77
6-MAM	328.1	165.2	2.8
Midazolam	326.1	291.1	2.83

Analyte	Q1	Q3	Retention Time (min)
Normeperidine	234.1	160.2	2.85
EDDP	278.3	234.2	2.95
Methadone	310.2	265.2	2.98
Lorazepam	321	275.1	2.98
Clonazepam	316	270.1	3.03
Norpropoxyphene	308.2	100.1	3.05
Oxazepam	287.1	241	3.05
Hydroxylalprazolam	325.1	297	3.12
Nordiazepam	271.1	140	3.12
Flunitrazepam	314.1	268.2	3.21
Temazepam	301.1	255.1	3.24
Alprazolam	309.1	281	3.26
Diazepam	285.1	193.1	3.37

Extended Lifetime with Biological Samples

Kinetex Biphenyl Stability in Urine (First 1000 injections)




- 7-Aminoclonazepam
- 7-Aminoflunitrazepam
- Alprazolam
- Amphetamine
- Benzoylcegonine
- Buprenorphine
- Clonazepam
- Codeine
- Estazolam
- Flunitrazepam
- Hydrocodone
- Hydromorphone
- Lorazepam
- MDMA
- Methadone
- Midazolam
- Morphine
- Norbuprenorphine
- Nordiazepam
- Norhydrocodone
- Oxazepam
- Temazepam
- Triazolam


It's Time to Upgrade to Kinetex® Biphenyl

Easily adjust your current fully porous methodology to get performance benefits:

Restek® 3 µm Ultra Biphenyl




The Kinetex Upgrade



**Faster run time,
higher sensitivity**

50 x 2.1 mm



**Excellent retention,
sensitivity, and
low pressure**

50 x 3.0 mm

Conditions for all separations:

Column: Restek 3 µm Ultra Biphenyl (50 x 2.1 mm)
Kinetex 2.6 µm Biphenyl (50 x 2.1 mm)
Kinetex 2.6 µm Biphenyl (50 x 3.0 mm)

Mobile Phase: A: Water with 0.1 % Formic Acid
B: Methanol with 0.1 % Formic Acid

Gradient:	Time (min)	%B
	0	10
	0.5	10
	2	25
	4.5	80
	4.51	85
	5.5	85
	5.51	10
	7	10

Flow Rate: 0.6 mL/min

Temperature: 40°C

Detection: MS/MS (AB SCIEX API 5000™)

Samples:

- | | | |
|--------------------|----------------------|------------------------|
| 1. Morphine | 17. Benzoylcegonine | 33. Methadone |
| 2. Oxycodone | 18. Meprobamate | 34. Clonazepam |
| 3. Hydromorphone | 19. Normeperidine | 35. Oxazepam |
| 4. Amphetamine | 20. Meperidine | 36. Hydroxylalprazolam |
| 5. Methamphetamine | 21. Norbuprenorphine | 37. Nordiazepam |
| 6. MDA | 22. Fentanyl | 38. Flunitrazepam |
| 7. Naloxone | 23. Buprenorphine | 39. Temazepam |
| 8. Codeine | 24. Flurazepam | 40. Alprazolam |
| 9. 6-MAM | 25. PCP | 41. Diazepam |
| 10. Oxycodone | 26. Carisoprodol | |
| 11. MDMA | 27. Propoxyphene | |
| 12. Naltrexone | 28. Midazolam | |
| 13. Hydrocodone | 29. Sufentanil | |
| 14. MDEA | 30. Norpropoxyphene | |
| 15. Norfentanyl | 31. EDDP | |
| 16. Tramadol | 32. Lorazepam | |

Restek is a registered trademark of Restek Corporation. Phenomenex is not affiliated with Restek Corporation.
Comparative separations may not be representative of all applications.

Ordering Information



STRATA-X-DRUG B

Format	Sorbent Mass	Part Number	Unit	Price
Tube				
	10 mg	8B-S128-AAK	1 mL (100/box)	
	30 mg	8B-S128-TAK	1 mL (100/box)	
	30 mg	8L-S128-TAK†	1 mL (100/box)	
	30 mg	8B-S128-TBJ	3 mL (50/box)	
	60 mg	8B-S128-UBJ	3 mL (50/box)	
	60 mg	8B-S128-UCH	6 mL (30/box)	
	60 mg	8B-S128-UCL	6 mL (200/box)	

96-Well Plate



† Tab-less tube

10 mg	8E-S128-AGB	2 Plates/box
30 mg	8E-S128-TGB	2 Plates/box
60 mg	8E-S128-UGB	2 Plates/box



IMPACT PROTEIN PRECIPITATION


96-Well Plate

Part No.	Description	Unit	Price
CE0-7565	Impact Protein Precipitation 2 mL Square Well Filter Plate	2/Box	
CE0-7566	Impact Protein Precipitation Long Drip 2 mL Square Well Filter Plate	2/Box	



Kinetex core-shell columns come in a variety of phases, dimensions, and particle sizes.

Visit www.phenomenex.com/kinetex for more info.



guarantee

We guarantee that if Phenomenex products in this guide do not provide at least an equivalent separation as compared to other products of the same phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.

KINETEX 2.6µm COLUMNS (mm)

SecurityGuard™ ULTRA Cartridges†									
Phases	50 x 2.1	100 x 2.1	50 x 3.0	50 x 4.6	100 x 4.6	150 x 4.6	3/pk	3/pk	3/pk
Biphenyl	00B-4622-AN	00D-4622-AN	00B-4622-Y0	00B-4622-E0	00D-4622-E0	00F-4622-E0	AJ0-9209	AJ0-9208	AJ0-9207
							for 2.1 mm ID	for 3.0 mm ID	for 4.6 mm ID

KINETEX 5µm COLUMNS (mm)

SecurityGuard™ ULTRA Cartridges†									
Phases	50 x 2.1	50 x 3.0	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk	3/pk	3/pk
Biphenyl	00B-4627-AN	00B-4627-Y0	00B-4627-E0	00D-4627-E0	00F-4627-E0	00G-4627-E0	AJ0-9209	AJ0-9208	AJ0-9207
							for 2.1 mm ID	for 3.0 mm ID	for 4.6 mm ID

†SecurityGuard ULTRA cartridges require holder, Part No. AJ0-9000.

Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex is a registered trademark, Strata-X, Impact, and SecurityGuard are trademarks of Phenomenex. Restek is a registered trademark of Restek Corporation. API 5000 is a trademark of AB SCIEX Pte. Ltd. AB SCIEX is being used under license.


Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is not affiliated with Restek Corporation.

Strata-X is patented by Phenomenex. U.S. Patent No. 7,119,145

© 2015 Phenomenex, Inc. All rights reserved.

Download guides, access on demand webinars,
and search our application database at:



www.phenomenex.com/clinicalapps

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 (0)20 2282-6668
f: +86 (0)20 2809-8130
chinainfo@phenomenex.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 001-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain

t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country,
contact Phenomenex USA, International Department at international@phenomenex.com