

# APPLICATIONS

## GC-MS/MS analysis of Tetra through Octa Dioxins and Furans on Zebron™ ZB-Dioxin GC column

Ramkumar Dhandapani<sup>1</sup>, Rola El abaji<sup>1</sup>, Eric Cull<sup>2</sup> and Agustin Pieri<sup>2</sup>  
<sup>1</sup>Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA  
<sup>2</sup>Weck Laboratories, Inc., 14859 Clark Avenue, Industry, CA 91745 USA



**Ramkumar Dhandapani, Ph.D.**  
 Product Manager - Gas Chromatography

*He has a PhD in Analytical Chemistry and over 15 years experience in chromatographic method development and troubleshooting. Ramkumar loves to write poems, read Shakespeare, and attend Shakespeare plays.*

### Introduction

Dioxins are toxic compounds and exhibit stable structure due to the presence of electron withdrawing group Chlorine attached to aromatic ring. Due to the stability of the bond, these analytes remain stable in the environment and falls under Persistent Organic Pollutants (POPs) category. They are continuously monitored in environmental samples. In addition to environmental samples, dioxins also enter into food chain. So, there is a necessity for highly efficient method to monitor these toxic compounds and separate them from matrix. Common method of analysis utilizes a GC-HRMS or a GC-MS/MS to accurately quantitate and separate toxic Dioxins and furans. EPA-1613 and EPA-8290 recommends analysis by GC-HRMS on 2 different GC columns. Presented here is a GC-MS/MS method that utilizes single column to exceed EPA-1613 and EPA-8290 requirements. In addition, Food samples have a lot of interfering components that needs monitoring of parent to daughter ion transition in order to authentically quantitate toxic dioxins by GC-MS/MS.


The method presented utilizes specific isotopically labeled dioxin and furan compounds as internal standard to accurately quantify dioxins and furans. Detailed evaluation of analytical figure of merits for the GC-MS/MS method was performed. The selectivity of ZB-Dioxin GC column provides increased resolution for 2,3,7,8-TCDD and 2,3,7,8-TCDF from its isomer components, thereby making GC-MS/MS method optimization easier and as an alternative technique to traditional HRMS analysis.

### GC-MS/MS Method Parameters for all chromatograms

**Column:** Zebron ZB-Dioxin  
**Dimensions:** 60 meter x 0.25 mm x 0.20 µm with 5 meter Guardian™  
**Part No.:** 7KG-G045-10-GGA  
**Injection:** Pulse Splitless (1.0 min, 30 psi) @ 110 °C for 0.31 min to 330 °C at 600 °C/min  
**Recommended Liner:** Zebron PLUS 2 mm ID Straight Wool in Middle  
**Liner Part No:** AG2-0E00-05  
**Carrier Gas:** Helium @ 1.056 mL/min (constant flow)  
**Oven Program:** 100 °C for 0 min to 220 °C @ 30 °C/min for 15.5 min to 240 °C @ 2 °C/min for 5 min to 270 °C @ 5 °C/min for 5.0 min to 330 °C @ 15 °C/min for 7.0 min  
**Detector:** GC-MS/MS  
**Transfer Line Temperature:** 340 °C


**Figure 1**

GC-MS/MS Separation of 2,3,7,8-TCDD and its isomer on a 60 meter Zebron ZB-Dioxin GC Column with 5 meter Guardian.



**Table1**  
MS/MS Parameters for Determination of Tetra Through Octa Dioxins and Furans.


Analyte	Start Time (min)	Precursorion	Production	Dwell Time (ms)	CE (V)	
13C-TCDD	22.0	333.9	269.9	50	26	
		331.9	267.9			
37Cl-TCDD		328	263.5	10	26	
		326	261.5			
TCDD		321.9	258.9	100	26	
		319.9	256.9			
13C-TCDF		317.9	253.9	50	40	
		315.9	251.9			
TCDF		305.9	242.9	100	40	
		303.9	240.9			
13C-PeCDF		31.2	351.9	287.9	25	40
PeCDF			349.9	285.9		
			13C-TCDD	339.9	276.9	75
337.9				274.9		
13C-TCDD	333.9		269.9	25	26	
	331.9		267.9			
13C-TCDD	321.9		258.9	75	26	
	319.9		256.9			
13C-TCDF	317.9		253.9	25	40	
	315.9		251.9			
TCDF	305.9		242.9	75	40	
	303.9		240.9			
13C-PeCDD	33.0		367.9	303.9	50	26
PeCDD			365.9	301.9		
		13C-PeCDF	355.9	292.9	100	26
353.9			290.9			
PeCD		351.9	287.9	50	40	
		349.9	285.9			
13C-HxCDDs		339.9	276.9	100	40	
		337.9	274.9			
HxCDDs		403.9	339.9	50	25	
		401.9	337.9			
13C-HxCDFs	391.8	328.8	100	25		
	389.8	326.8				
HxCDFs	387.9	323.9	50	40		
	385.9	321.9				
13C-HpCDD	375.8	312.8	100	40		
	373.8	310.8				
HpCDD	437.8	373.8	50	24		
	435.8	371.8				
13C-HpCDF	425.8	362.8	100	24		
	423.8	360.8				
HpCDF	421.8	357.8	50	40		
	419.8	355.8				
13C-OCDD	409.8	346.8	100	40		
	407.8	344.8				
OCDD	471.8	407.8	50	24		
	469.8	405.8				
13C-OCDF	459.7	396.7	100	24		
	457.7	394.7				
OCDF	455.8	391.8	50	40		
	453.8	389.8				
OCDF	443.7	380.7	100	40		
	441.7	378.7				

**Figure 2**GC-MS/MS Separation of 2,3,7,8-TCDF and its isomer on a 60 meter Zebron<sup>TM</sup> ZB-Dioxin GC Column with 5 meter Guardian<sup>TM</sup>.

Analyte Name	RT (min)	
1	1,3,6,8-TCDF	23.7
2	2,3,7,8-TCDF	29.4
3	1,2,8,9-TCDF	32.4

**Figure 3**


GC-MS/MS Separation of HxCDD and its isomer on a 60 meter Zebron ZB-Dioxin GC Column with 5 meter Guardian.



Analyte Name	RT (min)	
1	1,2,4,6,7,9-HxCDD	42.5
2	1,2,3,4,7,8-HxCDD	45.8
3	1,2,6,7,8-HxCDD	46.2
4	1,2,3,7,8,9-HxCDD	45.6

**Figure 4**


GC-MS/MS Separation of HxCDF and its isomer on a 60 meter Zebron ZB-Dioxin GC Column with 5 meter Guardian.



Analyte Name	RT (min)	
1	1,2,3,4,7,8-HxCDF	41.1
2	1,2,3,4,6,8-HxCDF	43.9
3	1,2,3,6,7,8-HxCDF	44.3
4	2,3,4,6,7,8-HxCDF	45.6
5	1,2,3,7,8,9-HxCDF	47.2

**Figure 5**


GC-MS/MS Separation of PeCDD and its isomer on a 60 meter Zebron™ ZB-Dioxin GC Column with 5 meter Guardian™.



Analyte Name	RT (min)
1 1,2,4,7,9-PeCDD	34.8
2 1,2,3,7,8-PeCDD	38.9
3 1,2,3,8,9-PeCDD	40.1

**Figure 6**


GC-MS/MS Separation of PeCDF and its isomer on a 60 meter Zebron ZB-Dioxin GC Column with 5 meter Guardian.



Analyte Name	RT (min)
1 1,3,4,6,8-PeCDF	31.6
2 1,2,3,7,8-PeCDF	37.0
3 2,3,4,7,8-PeCDF	39.0
4 1,2,3,8,9-PeCDF	40.6

**Figure 7**


GC-MS/MS Separation of OCDD on a 60 meter Zebron ZB-Dioxin GC Column with 5 meter Guardian.



Analyte Name	RT (min)
1 OCDD	52.8

**Figure 8**


GC-MS/MS Separation of OCDF on a 60 meter Zebron™ ZB-Dioxin GC Column with 5 meter Guardian™.

**Figure 9**

Zebron ZB-Dioxin GC column with integrated Guardian

The ZB-Dioxin column with integrated Guardian, provides the following benefits:

- Protects against non-volatile contamination by acting as a guard
- Prevents excessive trimming of analytical column thereby by extending the resolution of critical pairs
- Provides solvent focusing by acting as a retention gap for large volume samples and for low boiling solvent when used as the extracting solvent
- Seamless union of analytical and guard column leak-free connection that reduces Mass spec repair and downtime



**Table 2**

Linearity and Accuracy of Tetra through Octa Dioxins and Furans on ZB-Dioxin by GC-MS/MS.

Analyte	Retention Time (min)	Internal Standard	Range (ppb)	R <sup>2</sup>
13C-1,2,3,4-TCDD	29.41	13C-2,3,7,8-TCDD	100	-
2,3,7,8-TCDF	29.37	13C-2,3,7,8-TCD	0.05-219	0.98996
37Cl-2,3,7,8-TCDD	30.364	13C-2,3,7,8-TCDD	0.2-222	0.98562
2,3,7,8-TCDD	30.38	13C-TCD	0.05-224	0.98996
1,2,3,7,8-PeCDF	36.95	13C-1,2,3,7,8-PeCDF	0.3-1045	0.99778
2,3,4,7,8-PeCDF	38.97	13C-2,3,4,7,8-PeCDF	0.3-1042	0.99789
1,2,3,7,8-PeCDD	39.39	13C-1,2,3,7,8-PeCDD	0.2-1055	0.99638
1,2,3,4,7,8-HxCDF	43.93	13C-1,2,3,4,7,8-HxCDF	0.2-1109	0.98727
1,2,3,6,7,8-HxCDF	44.20	13C-1,2,3,6,7,8-HxCDF	0.25-1073	0.99391
2,3,4,6,7,8-HxCDF	45.53	13C-2,3,4,6,7,8-HxCDF	0.22-1121	0.98465
2,3,4,6,7,8-HxCDD	45.81	13C-1,2,3,4,7,8-HxCDD	0.24-1114	0.98631
1,2,3,7,8,9-HxCDD	46.53	13C-1,2,3,7,8,9-HxCDD	0.22-1120	0.98552
1,2,3,4,6,7,8-HpCDF	48.96	13C-1,2,3,4,6,7,8-HpCDF	0.26-1003	0.99995
1,2,3,4,6,7,8-HpCDD	50.11	13C-1,2,3,4,6,7,8-HpCDD	0.25-1021	0.99998
1,2,3,4,7,8,9-HpCDF	50.66	13C-1,2,3,4,7,8,9-HpCDF	0.25-1021	0.99942
OCDD	52.83	13C-OCDD	0.6-1950	0.99917
OCDF	53.08	13C-OCDF	0.5 to 2043	0.99941

**Table 3**

Detection Limit and Quantitation Limit Calculation for 2,3,7,8-TCDF from Linearity.

Concentration in ppb	Relative Response
0.05	0.0005
0.45	0.005
1.90	0.0212
10.26	0.1145
42.45	0.4745
219.21	2.4468
Slope	0.01
Residual Standard Deviation (RSD)	0.000316
Detection Limit for 2,3,7,8-TCDF	<b>3.3*RSD/Slope</b>
	0.09 ppb
Quantitation Limit for 2,3,7,8-TCDD	<b>10*RSD/Slope</b>
	0.28 ppb

## Results and Discussion

Dioxin analysis is complex due to numerous dioxin and furan compounds that are close in structure. Under EI mass spec source, the fragmentation in the source are similar for closely resembling compounds. There is a necessity for GC stationary phase to resolve critical pairs chromatographically, irrespective of high end detectors like GC-HRMS or GC-MS/MS. While a GC-MS/MS method can certainly resolve coelutions spectrally in many cases, however with enhanced chromatographic resolution, optimization of a GC-MS/MS method becomes easier. Presented in **Figure 1** is the separation of toxic dioxin isomer 2,3,7,8-TCDD from its isomer. Regulated methods like EPA-1613 mandate that the resolution is less than 25% valley for these critical pairs. With Zebron™ ZB-Dioxin GC column, a baseline resolution of 2,3,7,8-TCDD from its neighboring isomer 1,2,3,9-TCDD is achieved. **Figure 1-8** shows separation of Tetra through Octa Dioxins and Furans on a single method. Traditionally, a confirmation column like a 225 phase would be used to resolve coelution of TCDF isomers. With Zebron ZB-Dioxin GC column, the selectivity of the stationary phase allows for the complete resolution of 2,3,7,8-TCDF from its isomers, thereby cancelling out the necessity for a secondary confirmation column.

Presented in **Table 2** are the analytical figures of merits for the developed GC-MS/MS method on Zebron ZB-Dioxin GC column. The method provided a wide linear range from low to high parts per billion for various analytes. The exact range for individual dioxin is listed in the same table. In order to eliminate discrepancy in accuracy due to difference in physico-chemical properties of the analytes, it is essential to use surrogates and internal standards that are close in structure to the analyte of interest. For the same reason, multiple isotopically labeled dioxins and furans with <sup>13</sup>C label and <sup>37</sup>Cl were chosen. A six point calibration curve was plotted for individual dioxins and furans to establish the linearity and range of the method. The R<sup>2</sup> values are displayed in the table demonstrating linearity of the developed method. Since analysis of samples by GC-MS/MS eliminates the noise level by being selective in the first and third quadrupole, it is prudent to consider an alternate technique to determine the detection and quantitation limit of analytes instead of adopting traditional add hyphens signal-to-noise method. Presented in **Table 3** is the detection limit and quantitation limit calculation for 2,3,7,8-TCDF calculated using linearity data. This method is instrument independent and method specific and considers the error in the regression line to calculate accurate detection and quantification limits.

This method can be considered for both environmental and food based samples. Specifically, when the samples are expected to have non-volatile contaminants, it is very essential to use a guard column to trap non-volatile matrix contaminants. The chosen column dimension for this method utilizes a 5 meter integrated Guardian™, which is a seamless union of guard column to the analytical column (**Figure 9**). This offers leak-free connection which is extremely useful for MS based detection. With traditional guard columns, there is a high probability of leaks from the connections that can lead to expensive Mass spec repair and maintenance.

The Zebron ZB-Dioxin column with integrated Guardian, provides the following benefits:

- Protects against non-volatile contamination by acting as a guard
- Prevents excessive trimming of analytical column thereby extending the resolution of critical pairs
- Provides solvent focusing by acting as a retention gap for large volume samples and for low boiling solvent when used as the extracting solvent
- Seamless union of analytical and guard column provides leak-free connection that reduces Mass spec repair and downtime

These proven benefits of Guardian ensures less frequent column change, reduces instrument downtime and improves lab productivity.

In addition to providing enhanced resolution of TCDD and TCDF isomers, Zebron ZB-Dioxin stationary phase undergoes extensive cross-linking through Engineered Self Cross-Linking™ (ESC™) process to provide mass spec compatible low bleed column. This reduces potential column bleed to reduce instrument downtime and maintenance. Thus, the proprietary selectivity, enhanced resolution of TCDD and TCDF, ESC, integrated Guardian protection and high efficiency makes ZB-Dioxin, the ideal choice for the separation of dioxin and furan compounds by GC-MS/MS.

## Conclusion

The Zebron ZB-Dioxin GC column provides the highest resolution for 2,3,7,8-TCDD and 2,3,7,8-TCDF from its isomer compounds and serves as a single GC column solution for dioxin separation by GC-MS/MS. The integrated Guardian option prevents excessive trimming of analytical column thereby increasing lab throughput.

## Acknowledgement

*This study was performed at Weck Laboratories, Inc. and we thank them for their collaboration on this project. Phenomenex is not affiliated with Weck Laboratories, Inc.*

## Ordering Information

### Zebron™ ZB-Dioxin GC Column

ID(mm)	df(μm)	Temp. Limits °C	Part No.
<b>60-Meter</b>			
0.25	0.2	40 to 320/340	<a href="#">7KG-G045-10</a>
<b>60-Meter with 5-Meter Guardian™</b>			
0.25	0.2	40 to 320/340	<a href="#">7KG-G045-10-GGA</a>
<b>40-Meter</b>			
0.18	0.14	40 to 320/340	<a href="#">7PD-G045-47</a>

### Zebron PLUS Liners

Description	Inlet Style	Dimensions ID x L (mm)	Deactivation	Part No.	Unit
<b>For Agilent® and Thermo®</b>					
Single taper	S/SL	4 x 78.5	PLUS Inert	<a href="#">AG2-0A10-05</a>	5/pk
<b>For Agilent and Thermo</b>					
Single taper	S/SL	2 x 78.5	PLUS Inert	<a href="#">AG2-0E00-05</a>	5/pk

### Standard Z-Guard™ Columns and Kits

ID (mm)	Description	Part No.	
		5-Meter	10-Meter
0.25	Guard Column	<a href="#">7AG-G000-00-GZ0</a>	<a href="#">7CG-G000-00-GZ0</a>
	Guard Column Kit	<a href="#">7AG-G000-00-GZK</a>	<a href="#">7CG-G000-00-GZK</a>

### Easy Seals™ Inlet Base Seals

				10/pk
Description	Injection Type	Groove Style	Inlet Hole Diameter (mm)	Part No.
Easy Seals Gold Inlet Seal	Splitless	Single	0.8	<a href="#">AG0-8620</a>



**BE-HAPPY™**  
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)



**Australia**

t: +61 (0)2-9428-6444  
 auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
 anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
 t: +32 (0)2 511 8666 (Dutch)  
 beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
 info@phenomenex.com

**China**

t: +86 400-606-8099  
 cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
 cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
 nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
 nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
 franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
 anfrage@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
 indiainfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
 eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
 italiainfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
 nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
 tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
 nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
 nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
 nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
 pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
 ptinfo@phenomenex.com

**Singapore**

t: +65 800-852-3944  
 sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
 sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
 espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
 nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
 swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
 twinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
 ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
 info@phenomenex.com

🌐 **All other countries/regions  
 Corporate Office USA**

t: +1 (310) 212-0555  
 info@phenomenex.com

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

**Terms and Conditions**

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Zebron, Guardian, Engineered Self Cross-linking, ESC, Z-Guard, Easy Seals, and BE-HAPPY are trademarks of Phenomenex.

Agilent is a registered trademark of Agilent Technologies, Inc. Thermo Scientific is a registered trademark of Thermo Fisher Scientific Inc.

**Disclaimer**

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Agilent or Thermo.

This study was performed at Weck Laboratories Inc. and we thank them for their collaboration on this project. Phenomenex is not affiliated with Weck Laboratories Inc.

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.