

# APPLICATIONS

## Stationary Phases for the Process Scale Purification of Peptides and Insulin Analogs

Matt Trass, J Preston, Dani Xing, Genevieve Hodson  
Phenomenex Inc., 411 Madrid Avenue, Torrance, CA 90501, USA

### Introduction

The crude peptide analog used in this purification study is crude Insulin Glargine first brought to market by Sanofi-Aventis, Germany in 2000 mainly for the purpose of managing type 1 diabetes. Since coming off patent in 2015, there has been a large interest in the purification process for generic alternatives.

Well established methods have been used in industry for the purification of insulin while improvements are being made continuously. One of the limitations to high-throughput purification is the cost and time prohibitive nature of the media used. Here we examine two silicas, both C8 phases, for the purification of an insulin analog as well as evaluating their advantage for the removal of aggregate (build-up) post purification. While some insulin-like peptides such as Liraglutide have shown not to have issue with aggregate buildup, others do require significantly longer rewashes to remove build-up.

The crude material was obtained from a major insulin and insulin analogs manufacturer. Insulin Glargine is a 53-amino acid peptide and its chemical structure is represented in Figure 1. The analog has the base structure of A chain and B chain of human insulin but the asparagine residue at A21 on the A chain has been replaced with glycine. Additionally, two arginine were added to the C-terminus of the B-chain at B31 and B32.

The development of a multi-step purification process on two types of silicas is presented here for commercially significant crude insulin analog, Insulin Glargine. The focus of this tech note is to demonstrate the efficacy of both prep C8 phases on a well-established USP method that has been further optimized on the analytical format and then scaled to the prep format to produce material of a suitable purity. The investigative parameters include base silica and gradient conditions.

### Materials and Methods

The crude material was provided from a major insulin and insulin analogs manufacturer. Sodium phosphate, sodium chloride,

1-propanol were obtained from Fisher Scientific (Waltham, MA, USA); acetonitrile was obtained from Honeywell (Morris Plains, NJ, USA); and acetic acid, ammonium acetate and ammonium chloride were obtained from Sigma Chemical (St. Louis, MO, USA).


### Results and Discussion

The purity of the crude insulin glargine material was obtained with a slight modification of the USP insulin methodology. The arginine modification of insulin Glargine shifts the PI of human insulin from 5.4 to 6.7. As such the pH of the mobile phase was modified from pH 2.3 to pH 4. This modification was made in the analytical format. Prep conditions were then used to show efficacy of both C8 media for purification. While the 20 minute gradient was sufficient to demonstrate purification, a longer run at 40 minutes, with a prolonged gradient, proved to be more effective for larger scale purification.

The main intent was to show the Gemini C8(3) hybrid silica has comparable selectivity to a typical C8 media found in industry. As proof of concept, we only took a representative single step purification approach. Fractions were collected every 30 secs and all fractions were assayed. The analytical results were used to calculate purities and yields for theoretical pools. On both Luna and Gemini media, the chromatographic profile were similar and the purity and yield results were comparable. The purity obtained from these one step purifications were at least 98%.

An additional study was conducted to verify the applicability of the Gemini C8(3) material for insulin purification. To remove aggregate buildup, a caustic wash is often employed. Figure 2 shows the sustainability of the silica material chromatographically even after 20 washes using 50:50 1 N NaOH/1-propanol (3.5 CV), then 55:45 0.02 M Acetic acid/1-propanol (12 CV).

### Chemical Structure of Insulin Glargine



Peptide chemical formula  
C<sub>267</sub>H<sub>404</sub>N<sub>72</sub>O<sub>78</sub>S<sub>6</sub>  
Peptide average weight 6063.0 Da

#### Sequences

>A chain  
GIVEQCCTSICSLYQLENYCG

>B chain  
FVNQHLCGSHLVEALYLVCGERGFFYTPKTRR


<https://www.drugbank.ca/drugs/DB00047>



# APPLICATIONS

**Figure 1.**

**Comparison of Luna C8(3) and Gemini C8(3) for Different Gradient Times and Formats**



Conditions same for both columns:

**Columns:** As described  
**Dimensions:** 250 x 4.6 mm  
**Mobile Phase:** A: 20 mM Ammonium acetate, 200 mM Ammonium chloride, 10 % Acetonitrile, pH 4  
 B: 20 mM Ammonium acetate, 50 % Acetonitrile, pH 4  
**Gradient:**

| Time (min) | % B |
|------------|-----|
| 0          | 35  |
| 20         | 55  |

  
**Flow Rate:** 1.2 mL/min  
**Injection Volume:** 20 µL  
**Temperature:** 35 °C  
**Detection:** UV @ 214 nm (ambient)

Conditions same for both columns:

**Columns:** As described  
**Dimensions:** 250 x 4.6 mm  
**Mobile Phase:** A: 20 mM Ammonium acetate, 200 mM Ammonium chloride, 10 % Acetonitrile, pH 4  
 B: 20 mM Ammonium acetate, 50 % Acetonitrile, pH 4  
**Gradient:**

| Time (min) | % B |
|------------|-----|
| 0          | 35  |
| 40         | 50  |
| 45         | 50  |

  
**Flow Rate:** 1.2 mL/min  
**Injection Volume:** 200 µL (10 mg)  
**Temperature:** 35 °C  
**Detection:** UV @ 214 nm (ambient)


Conditions same for both columns:

**Columns:** As described  
**Dimensions:** 250 x 4.6 mm  
**Mobile Phase:** A: 20 mM Ammonium acetate, 200 mM Ammonium chloride, 10 % Acetonitrile, pH 4  
 B: 20 mM Ammonium acetate, 50 % Acetonitrile, pH 4  
**Gradient:**

| Time (min) | % B |
|------------|-----|
| 0          | 35  |
| 40         | 50  |
| 45         | 50  |

  
**Flow Rate:** 1.2 mL/min  
**Injection Volume:** 500 µL (25 mg)  
**Temperature:** 35 °C  
**Detection:** UV @ 214 nm (ambient)

# APPLICATIONS

**Figure 2.**
**Collected Fractions from Gemini C8(3) vs Luna C8(3)**


Column: Kinetex<sup>®</sup> 5 $\mu$ m C18  
 Dimensions: 250 x 4.6 mm  
 Mobile Phase: A: 58 mM Sodium phosphate, 417 mM Sodium chloride, 25% Acetonitrile, pH2.5  
 B: 125 mM Sodium phosphate, 893 mM Sodium chloride, 65% Acetonitrile, pH2.5  
 Gradient:
 

| Time (min) | % B |
|------------|-----|
| 0          | 4   |
| 20         | 17  |
| 30         | 37  |

  
 Flow Rate: 1.2 mL/min  
 Injection Volume: 5  $\mu$ L  
 Temperature: 35  $^{\circ}$ C  
 Detection: UV @ 214 nm (ambient)

**Conclusion**

Chromatographic media is often a significant portion of the cost for large scale peptide purification. The cost is not solely due to the media, but also the time spent regenerating the column back to initial conditions. Therefore, it is often more cost and time effective to use a durable silica media that can be used under caustic conditions.

This work used the analog Insulin Glargine to demonstrate the comparability of the high pH resistant silica to that of a typical silica used in industry chromatographically. Multiple fractions were assayed from both Luna C8(3) and high pH stable Gemini C8(3). These two media were found to be comparable and yield insulin glargine at greater than 98% purity.

Gemini C8(3) and Luna C8(3) were used as the stationary phase for the purification of Insulin Glargine. Luna C8(3) was introduced in 2014 and Gemini C8(3) in 2017.

For both columns, 25 mg of crude (1% load) was loaded and fractions were collected every 30 seconds from before the main peak eluted. These fractions were assayed by the analytical USP methodology using a Kinetex C18 column.

From the fraction analysis data, purity and yield results for different possible pools were calculated.

**Luna 10  $\mu$ m-PREP C8(3)\***

| Pool | % purity | % yield |
|------|----------|---------|
| 1    | 98.41    | 90.9    |
| 2    | 98.91    | 87.4    |
| 3    | 99.30    | 82.2    |
| 4    | 99.60    | 75.4    |

**Gemini 10  $\mu$ m C8(3)\***

| Pool | % purity | % yield |
|------|----------|---------|
| 1    | 98.24    | 92.8    |
| 2    | 98.91    | 84.3    |
| 3    | 99.44    | 71.9    |
| 4    | 99.71    | 61.5    |

**\* Table note:**

% yield is the peak area in the pool divided by the total peak area available.  
It is not pool divided by crude.

# APPLICATIONS

**Australia**  
t: +61 (0)2-9428-6444  
auiinfo@phenomenex.com

**Austria**  
t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**  
t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**  
t: +1 (800) 543-3681  
info@phenomenex.com

**China**  
t: +86 400-606-8099  
cninfo@phenomenex.com

**Denmark**  
t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**  
t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**  
t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**  
t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**India**  
t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Ireland**  
t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**  
t: +39 051 6327511  
italiainfo@phenomenex.com

**Luxembourg**  
t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**  
t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**  
t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**  
t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**  
t: +47 810 02 005  
nordicinfo@phenomenex.com

**Portugal**  
t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**  
t: +65 800-852-3944  
sginfo@phenomenex.com

**Spain**  
t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**  
t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**  
t: +41 61 692 20 20  
swissinfo@phenomenex.com

**United Kingdom**  
t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**  
t: +1 (310) 212-0555  
info@phenomenex.com

**All other countries  
Corporate Office USA**   
t: +1 (310) 212-0555  
info@phenomenex.com

## Luna<sup>®</sup> Ordering Information

| Luna 10 µm- <i>PREP</i> | Columns     | Bulk Media   |          |          |          |          |
|-------------------------|-------------|--------------|----------|----------|----------|----------|
|                         |             | 250 x 4.6 mm | 1 kg     | 5 kg     | 10 kg    | 50 kg    |
| C18(3)                  | 00G-4616-E0 | 04K-4616     | 04L-4616 | 04M-4616 | 04N-4616 | 04P-4616 |
| C8(3)                   | 00G-4623-E0 | 04K-4623     | 04L-4623 | 04M-4623 | 04N-4623 | 04P-4623 |
| Silica(3)               | 00G-4617-E0 | 04K-4617     | 04L-4617 | 04M-4617 | 04N-4617 | 04P-4617 |

## Gemini<sup>®</sup> C8(3) Ordering Information

| 10 µm Analytical, Semi-Prep and Axia <sup>™</sup> Packed Preparative Columns (mm) |             |             |                |                |                |
|---|-------------|-------------|----------------|----------------|----------------|
| Phase   | 250 x 4.6   | 250 x 10    | 250 x 21.2     | 250 x 30       | 250 x 50       |
| C8(3)   | 00G-4763-E0 | 00G-4763-N0 | 00G-4763-P0-AX | 00G-4763-U0-AX | 00G-4763-V0-AX |

| Bulk Media |          |          |          |          |
|------------|----------|----------|----------|----------|
| Phase      | 100 g    | 1 kg     | 5 kg     | 10 kg    |
| C8(3)      | 04G-4763 | 04K-4763 | 04L-4763 | 04M-4763 |

## Kinetex<sup>®</sup> Ordering Information

| 5 µm Analytical Columns (mm) |             |             |             |             | SecurityGuard <sup>™</sup><br>ULTRA Cartridges <sup>‡</sup> |
|------------------------------|-------------|-------------|-------------|-------------|---|
| Phase                        | 50 x 4,6    | 100 x 4,6   | 150 x 4,6   | 250 x 4,6   | 3/pz  |
| EVO C18                      | 00B-4633-E0 | 00D-4633-E0 | 00F-4633-E0 | 00G-4633-E0 | AJ0-9296  |
| F5                           | 00B-4724-E0 | 00D-4724-E0 | 00F-4724-E0 | 00G-4724-E0 | AJ0-9320  |
| Biphenyl                     | 00B-4627-E0 | 00D-4627-E0 | 00F-4627-E0 | 00G-4627-E0 | AJ0-9207  |
| XB-C18                       | 00B-4605-E0 | 00D-4605-E0 | 00F-4605-E0 | 00G-4605-E0 | AJ0-8768  |
| C18                          | 00B-4601-E0 | 00D-4601-E0 | 00F-4601-E0 | 00G-4601-E0 | AJ0-8768  |
| C8                           | 00B-4608-E0 | 00D-4608-E0 | 00F-4608-E0 | 00G-4608-E0 | AJ0-8770  |
| Phenyl-Hexyl                 | 00B-4603-E0 | 00D-4603-E0 | 00F-4603-E0 | 00G-4603-E0 | AJ0-8774  |

per ID 4,6 mm



If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to other products of the same phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.

### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

### Trademarks

Luna, Kinetex, and Gemini are registered trademarks of Phenomenex.

Gemini is patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

© 2018 Phenomenex, Inc. All rights reserved

## www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)