



Simplified Chiral HPLC/SFC Column Screening Strategies




HPLC Screen

Tip We suggest screening all five Lux phases to identify the optimal chiral separation.


Normal Phase (NP)



Polar Organic (PO)



Reversed Phase (RP)




Notes: This screening strategy can be started at any step depending on the properties of the racemates. A common dimension used in chiral screening is 250 x 4.6 mm. For faster screening, use shorter columns.
* Use 0.1 % DEA with basic and neutral compounds and 0.1 % HCOOH with acidic and neutral compounds.
† Changing % IPA in methanol can be occasionally beneficial.

Key: IPA: Isopropanol; DEA: Diethylamine; MeOH: Methanol; CH₃CN: Acetonitrile; EtOH: Ethanol; CH₃COONH₄: Ammonium acetate; HCOOH: Formic acid; NH₄HCO₃: Ammonium bicarbonate; CO₂: Carbon Dioxide

SFC Screen

Tip We suggest screening all five Lux phases to identify the optimal chiral separation.




*Can also use 0.1 % Formic Acid or 0.1% TFA

Solvent Considerations

Solvent Switching

Lux columns are shipped in 90 % Hexane : 10 % IPA



COMPATIBLE	
Polar Organic	<ul style="list-style-type: none"> Methanol Acetonitrile IPA Mixtures of above
Normal Phase	<ul style="list-style-type: none"> Alkane/alcohol mixtures
Reversed Phase	<ul style="list-style-type: none"> Aqueous methanol/acetonitrile Buffer and methanol/acetonitrile mixtures
SFC	<ul style="list-style-type: none"> Supercritical CO₂
AVOID	
<ul style="list-style-type: none"> Tetrahydrofuran Acetone Chlorinated hydrocarbons Ethylacetate 	<ul style="list-style-type: none"> Dimethylsulfoxide Dimethylformamide N-methylformamide Pyridine

**Once column is in reversed phase mode, it is not recommended to solvent switch.
See column care use notes at www.phenomenex.com/lux for more information.

Why Choose Lux Chiral Columns?

- Stable in normal phase, polar organic, SFC, and reversed phase conditions
- 3 µm and 5 µm packed columns, as well as, 10 µm and 20 µm bulk media for scale up
- Pressure stable up to 300 bar
- High efficiency and loading capacity



guarantee

If Lux analytical columns (≤ 4.6 mm ID) do not provide at least an equivalent or better separation as compared to a competing column of the same particle size, similar phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

Trademarks
Lux is a registered trademark of Phenomenex.
CHIRALCEL®, OJ-H®, OZ-H®, OD-H®, OX-H®, and AY-H® are registered trademarks of DAICEL Chemical Industries, Ltd.

© 2012 Phenomenex, Inc. All rights reserved.



FREE Chiral Screening!

For more information or to begin a project today, please contact your local Phenomenex representative.

You can also visit us online:
www.phenomenex.com/phenologix

or email us at:
phenologix@phenomenex.com

phenoLogix
Your Method. Our Scientists.

phenomenex
...breaking with tradition™

WEB: www.phenomenex.com/lux