

Arginine as a Mobile Phase Co-solvent to Improve High Molecular Weight Aggregate Recovery for Size Exclusion Chromatography

James Song, Brian Rivera, and Dr. Bryan Tackett
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

Size Exclusion Chromatography (SEC) is a technique for the separation of biomolecules based on their size in solution and is particularly useful for quantitating size variants and high molecular weight aggregates. Silica-based SEC columns are modified with a hydrophilic stationary phase, typically to minimize electrostatic interactions of positive moieties on proteins and other analytes. However, depending on the physicochemical properties of the analyte, it may be prudent to investigate different mobile phase co-solvents to improve separation and sample recovery. One co-solvent that has been studied extensively is arginine, a common reagent used for protein refolding.¹ In this application note, we explore the impact of sample (i.e. aggregate) recovery on NIST mAb when modulating arginine concentration in mobile phase.

The mobile phase consisted of 200 mM Potassium Phosphate, pH 6.2, with arginine concentration assessed at 100 mM, 200 mM, and 250 mM. For NIST mAb, aggregate recovery was significantly less for 100 mM (80851 mAU) when compared to 200 and 250 mM (96818 and 97146, respectively). Further, percent high molecular weight (HMW) aggregate by peak area was 2.16% for 100 mM, which was again significantly lower than higher concentrations. Visually, post-peak fragment separation was slightly worse with 100 mM. HMW aggregate percentages and resolution values of monomer and dimer are summarized in **Table 1**.

In summary, arginine may be a useful co-solvent for SEC applications for large molecules. Particularly, arginine may improve high molecular weight aggregate recovery for monoclonal antibodies.

LC Conditions

Column: Biozen™ 1.8 µm dSEC-2, 200 Å
Part No.: [00H-4787-E0](#)
Dimensions: 300 x 4.6 mm
Mobile Phase: 200 mM Potassium Phosphate, Arginine as indicated, pH 6.2
Flow Rate: 0.35 mL/min
Detection: UV @ 280 nm
Temperature: 25 °C
Sample: NIST mAb (30 µg)

1. Arakawa, Tsutomu et al. "Suppression of protein interactions by arginine: a proposed mechanism of the arginine effects." *Biophysical chemistry* vol. 127,1-2 (2007): 1-8. doi:10.1016/j.bpc.2006.12.007



Figure 1. SEC chromatographic overlays for ADC mimic, demonstrating the effect of organic solvent to peak shape and recovery.




Table 1: Assessment of chromatographic parameters

Arginine Concentration	HMW Area Count (mAU)	Monomer Peak Height	% HMW	Rs 1,2
100 mM	80851	453324	2.16	3.02
200 mM	96818	493990	2.56	3.26
250 mM	97146	473197	2.58	3.25



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Biozen and BE-HAPPY are trademarks of Phenomenex.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.

