


APPLICATIONS

Robust Separation of Hydrochlorothiazide and Chlorothiazide in Hydrochlorothiazide Tablets Using the Kinetex® 5 µm C18 column

Zeshan Aqeel, J Preston, Ryan Splitstone, and Phil Koerner
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Monograph: Hydrochlorothiazide
pK_{a1}: 7.9, pK_{a2}: 9.2
LogP: -0.07

Overview

This application highlights the robust and reproducible separation of Hydrochlorothiazide and its related impurity Chlorothiazide under the USP monograph conditions for assay of Hydrochlorothiazide Tablets. The application demonstrates the potential method improvements that can be achieved per the allowable adjustments outlined in USP General Chapter <621> relative to the original column and conditions referenced in the monograph.

USP Monograph: Hydrochlorothiazide Tablet Assay

Standard Solution 0.15 mg/mL of USP Hydrochlorothiazide RS in *Mobile Phase*

System Suitability Solution 0.015 mg/mL of Chlorothiazide RS and 0.015 mg/mL Hydrochlorothiazide RS in *Mobile Phase**

* Note: A volume of Acetonitrile not exceeding 10 % of the total volume of solution may be used to dissolve the USP ReferenceStandard (RS)

Column

Size Method 1: 250 x 4.6 mm, Method 2: 250 x 4.6 mm

Stationary Phase Method 1: Symmetry® 5 µm C18, Method 2: Kinetex 5 µm C18

Temperature 30 °C

Mobile Phase Acetonitrile and 0.1 M Monobasic Sodium Phosphate (1:9). Adjust with Phosphoric Acid to a pH of 3.0 ± 0.1

Isocratic Isocratic: (1:9, A:B)
Total Run Time: 30 min

Flow Rate 2.0 mL/min

Detector UV @ 254 nm

Injection Volume 20 µL of System Suitability solution and Standard solution

System Suitability – System Suitability solution and Standard solution

Sample: Standard solution and System Suitability solution:


- Resolution (Rs): NLT 2.0 between Chlorothiazide and Hydrochlorothiazide for System Suitability solution
- Relative Standard Deviation: NMT 1.5 % for Standard solution (5 replicate injections)

Method Overlay

Kinetex 5 µm C18 and Symmetry 5 µm C18 Overlay Standard solution

Standard Solution	Symmetry C18	Kinetex C18	Percent Difference
Retention Time	6.31	4.27	-32%
Peak Area	702.8	705.8	0%
Peak Height	76.6	127.1	66%
Peak Width	0.1413	0.0855	-39%

Column: Symmetry 5 µm C18
Kinetex 5 µm C18
Dimensions: 250 x 4.6 mm
Flow Rate: 2.0 mL/min
Sample: 1. Hydrochlorothiazide



APPLICATIONS


Method 1

System Suitability solution: Hydrochlorothiazide Symmetry[®] 5 μm C18

Column: Symmetry 5 μm C18
Dimensions: 250 x 4.6 mm
Flow Rate: 2.0 mL/min
Sample: 1. Chlorothiazide
2. Hydrochlorothiazide

Chlorothiazide Average Height: 55.94
Peak Height RSD: 0.84 %

Hydrochlorothiazide Average Height: 78.86
Peak Height RSD: 1.45 %



App ID: 25877


Method 2

System Suitability solution: Hydrochlorothiazide Kinetex[®] 5 μm C18

Column: Kinetex 5 μm C18 Core Shell
Dimensions: 250 x 4.6 mm
Part No.: [00G-4601-E0](#)
Flow Rate: 2.0 mL/min
Sample: 1. Chlorothiazide
2. Hydrochlorothiazide

Chlorothiazide Average Height: 87.54
Peak Height RSD: 0.06 %

Hydrochlorothiazide Average Height: 128.84
Peak Height RSD: 0.69 %



App ID: 25879

APPLICATIONS

Adjustments for Meeting System Suitability

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	(As specified)	(As specified)
Concentration of Salts in Buffer	± 10 %	(As specified)	(As specified)
Composition of the Mobile Phase	± 30 % Relative; cannot exceed ± 10 % Absolute adjustment; cannot be reduced to zero	(As specified)	(As specified)
Wavelength of Detector	No deviations permitted	254 nm (as specified)	(As specified)
Injection Volume	Can be adjusted as much as needed; must be consistent with linearity, precision, and detection requirements	20 µL (as specified)	10 µL (Allowed)
Column Temperature	± 10 °C	30 °C (Allowed)	30 °C (Allowed)
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	L1 (as specified)	(As specified)
Column Length	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25 % and +50 %*	250 mm (as specified)	250 mm (As specified)
Column Internal Diameter	Can be adjusted so long as linear velocity is maintained	4.6 mm (as specified)	4.6 mm (As specified)
Particle Size	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25 % and +50 %*	5 µm (as specified)	5 µm (As specified)
Flow Rate	± 50 % (at given ID)	2.0 mL/min (as specified)	2.0 mL/min (As specified)

*Alternatively (as for the application of particle size adjustment to superficially porous particles), other L/dp combinations can be used provided that the number of theoretical plates (N) is within -25% to +50%.

Allowable Column Adjustments: L/dp Ratio -25 % to 50 %

Column	Length (mm)	ID (mm)	dp (µm)	L/dp	Allowable Range 37,500-75,000
Original	250	4.6	5	50,000	37,500 – 75,000
Alternative	250	4.6	5	50,000	ALLOWED

Method Summary and Comparison

	Method 1	Method 2
Column	Symmetry® 5 µm C18	Kinetex 5 µm C18
System Suitability Solution Chlorothiazide Average Rt	5.1 min	3.64 min
System Suitability Solution Hydrochlorothiazide Average Rt	6.27 min	4.36 min
Relative Retention Times for Chlorothiazide and Hydrochlorothiazide are about 0.8 and 1.0, Respectively		
System Suitability Solution Hydrochlorothiazide Average Rs	5.4	5.51
System Suitability (Standard Solution) Hydrochlorothiazide Average Peak Height	74	135
System Suitability (Standard Solution) Hydrochlorothiazide Peak Height RSD (n=5)	0.9 %	1.48 %
Backpressure (Bar)	176	190

APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/ChatNow.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 (12) 881 0121
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

● All other countries/regions
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

BE-HAPPY™
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex is a registered trademark and BE-HAPPY is a trademark of Phenomenex. Symmetry is a registered trademark of Water Technologies Corporation.

Disclaimer

Phenomenex is not affiliated with Waters Corporation. Comparative separations may not be representative of all applications.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.