

# APPLICATIONS

## Comparison of Trypsin Digestion pH for Peptide Mapping Workflows

M. Christina Malinao and Brian Rivera  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501, USA

### Overview

Peptide mapping is a common method for protein characterization. The general workflow includes the isolation of a protein, followed by in-solution digest using a serine protease to yield peptides, which are subsequently separated by LC and analyzed by UV and/or MS techniques. Because of its specificity and the general size of peptides generated, trypsin is most commonly used. Trypsin has an optimal activity at around pH 8. However, under even moderately basic conditions, deamidation can occur. In this application note, we investigate the differences in digestion between a standard ammonium bicarbonate at pH 8.2 and a 1X PBS, pH 7.4.

As we can see in **Figure 1**, there is a clear difference visually when looking at TICs for the standard digestion in ammonium bicarbonate when compared to PBS. Appreciable peaks before 15 minutes are not observed when using PBS, indicating the absence of smaller tryptic peptides. This result may indicate missed cleavages due to the suboptimal pH by which the digestion is occurring.

However, in observing peptides that may be prone to deamidation, namely the "VVSV" peptide, we can see that there is significantly more deamidated peptides in the ammonium bicarbonate buffers. The two later eluting peaks indicate the deamidated VVSV peptides, with the bottom trace indicating significantly higher relative peak heights (**Figure 2**).

One interesting note is that sequence coverage between both samples prepared was relatively similar (91.7% for bicarbonate, 86.6% for PBS). Missed cleavages do not necessarily mean indicate a decrease in sequence coverage, as longer peptides may still be identified by tandem MS. However, missed cleavages may introduce variability in the digestion which should be assessed depending on the scope of the method.

In summary, there may be some variation when digestion is performed at suboptimal pH. Therefore, it is prudent to investigate the trypsin digestion pH to assess the impact on a particular protein of interest.

### Digestion Procedure:

| Step                   | Details   |
|------------------------|---|
| <b>Denaturation</b>    | To sample, add 1:1 (v:v) 5 M Guanidine  |
| <b>Reduction</b>       | 1:10 (v:v) 200 mM DTT:Protein<br>Incubate at 57 °C for 30 min, shaking at 1000 rpm                            |
| <b>Alkylation</b>      | 1:2 (v:v) 400 mM iodoacetamide (IAM): DTT<br>Incubate in the dark 45 min<br>Quench, 1:2 (v:v) 200 mM DTT: IAM |
| <b>Buffer Exchange</b> | 100 mM Ammonium Bicarbonate, overnight or<br>1X Phosphate Buffered Saline, pH 7.4, overnight                  |
| <b>Digestion</b>       | 1:20 (w/w) Trypsin:Sample<br>Incubate 37 °C for 6 h, shaking at 1000 rpm                                      |
| <b>Reaction Quench</b> | Formic acid<br>SpeedVac to dryness, resuspend in mobile phase prior to analysis                               |

### LC Conditions

**Column:** bioZen™ 2.6 µm Peptide XB-C18

**Dimension:** 150 x 2.1 mm

**Part No.:** [00F-4768-AN](#)

**Recommended Guard:** SecurityGuard™ ULTRA

**Guard Cartridge Part No.:** [AJ0-9806](#)

**Guard Holder Part No.:** [AJ0-9000](#)

**Mobile Phase:** A: 0.1 % Formic Acid in Water

B: 0.1 % Formic Acid in Acetonitrile

**Flow Rate:** 0.3 mL/min


**Gradient:** 1-50% B in 50 minutes

**Temperature:** 40 °C

**Detector:** Q-TOF (SCIEX® X500B)

**Sample:** Tryptic digest, NIST mAb

**Figure 1. Comparison of TICs, pH 8.2 vs pH 7.4**




**Comparison of Sequence Coverage**

**Ammonium Bicarbonate Digested Trastuzumab**  
Heavy Chain Sequence Coverage 91.7%

```

EVQLVESGGGLVQPGGSLRLSCAASGFNIKDTYIHWVRQAPGKGLWVARIYPTNGYTRYADSVKGRFTISADTSKNTAYLQMNSLRAEDTAVYYCSRWGGDGFYAMDYWGQGTLLVTVSSASTKGPSVFPLAPSSKSTSGTAALGCLVKDYFPEPVTVSWNSGALTSGVHTFPAVLQSSGLYSLSSVTVPSSSLGTQTYICNVNHKPSNTKVDKKVEPPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMI SRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSRDELTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNV FSCSV MHEALHNHYTQKSLSLSPGK
    
```

**Figure 2. XIC Comparison, VVSV Peptide**



**Phosphate Buffer (PBS) Digested Trastuzumab**  
Heavy Chain Sequence Coverage 86.6%


```

EVQLVESGGGLVQPGGSLRLSCAASGFNIKDTYIHWVRQAPGKGLWVAR IYPTNGYTRYADSVKGRFTISADTSKNTAYLQMNSLRAEDTAVYYCSRWGGDGFYAMDYWGQGTLLVTVSSASTKGPSVFPLAPSSKSTSGTAALGCLVKDYFPEPVTVSWNSGALTSGVHTFPAVLQSSGLYSLSSVTVPSSSLGTQTYICNVNHKPSNTKVDKKVEPPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMI SRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSRDELTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNV FSCSV MHEALHNHYTQKSLSLSPGK
    
```

# APPLICATIONS

## Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/LiveChat](http://www.phenomenex.com/LiveChat).

|  |  |   |
|--|--|---|
| <p><b>Australia</b><br/>t: +61 (0)2-9428-6444<br/>auiinfo@phenomenex.com</p>                                       | <p><b>India</b><br/>t: +91 (0)40-3012 2400<br/>indiainfo@phenomenex.com</p>      | <p><b>Singapore</b><br/>t: +65 800-852-3944<br/>sginfo@phenomenex.com</p>   |
| <p><b>Austria</b><br/>t: +43 (0)1-319-1301<br/>anfrage@phenomenex.com</p>  | <p><b>Ireland</b><br/>t: +353 (0)1 247 5405<br/>eirinfo@phenomenex.com</p>       | <p><b>Spain</b><br/>t: +34 91-413-8613<br/>espinfo@phenomenex.com</p>   |
| <p><b>Belgium</b><br/>t: +32 (0)2 503 4015 (French)<br/>t: +32 (0)2 511 8666 (Dutch)<br/>beinfo@phenomenex.com</p> | <p><b>Italy</b><br/>t: +39 051 6327511<br/>italiainfo@phenomenex.com</p>         | <p><b>Sweden</b><br/>t: +46 (0)8 611 6950<br/>nordicinfo@phenomenex.com</p>   |
| <p><b>Canada</b><br/>t: +1 (800) 543-3681<br/>info@phenomenex.com</p>  | <p><b>Luxembourg</b><br/>t: +31 (0)30-2418700<br/>nlinfo@phenomenex.com</p>      | <p><b>Switzerland</b><br/>t: +41 (0)61 692 20 20<br/>swissinfo@phenomenex.com</p>   |
| <p><b>China</b><br/>t: +86 400-606-8099<br/>cninfo@phenomenex.com</p>  | <p><b>Mexico</b><br/>t: 01-800-844-5226<br/>tecnicomx@phenomenex.com</p>         | <p><b>Taiwan</b><br/>t: +886 (0) 0801-49-1246<br/>twinfo@phenomenex.com</p>   |
| <p><b>Denmark</b><br/>t: +45 4824 8048<br/>nordicinfo@phenomenex.com</p>   | <p><b>The Netherlands</b><br/>t: +31 (0)30-2418700<br/>nlinfo@phenomenex.com</p> | <p><b>United Kingdom</b><br/>t: +44 (0)1625-501367<br/>ukinfo@phenomenex.com</p>  |
| <p><b>Finland</b><br/>t: +358 (0)9 4789 0063<br/>nordicinfo@phenomenex.com</p>                                     | <p><b>New Zealand</b><br/>t: +64 (0)9-4780951<br/>nzinfo@phenomenex.com</p>      | <p><b>USA</b><br/>t: +1 (310) 212-0555<br/>info@phenomenex.com</p>  |
| <p><b>France</b><br/>t: +33 (0)1 30 09 21 10<br/>franceinfo@phenomenex.com</p>                                     | <p><b>Norway</b><br/>t: +47 810 02 005<br/>nordicinfo@phenomenex.com</p>         | <p><b>All other countries<br/>Corporate Office USA</b> <br/>t: +1 (310) 212-0555<br/>info@phenomenex.com</p> |
| <p><b>Germany</b><br/>t: +49 (0)6021-58830-0<br/>anfrage@phenomenex.com</p>  | <p><b>Portugal</b><br/>t: +351 221 450 488<br/>ptinfo@phenomenex.com</p>         |   |

## www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™**

guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.  
[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

### Trademarks

bioZen, SecurityGuard, and Be-Happy are trademarks of Phenomenex. Sciex is a registered trademark of AB Sciex Pte. Ltd. AB SCIEX™ is being used under license.  
*FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.*

© 2020 Phenomenex, Inc. All rights reserved.