

# APPLICATION



## PFAS Analysis in Water Samples using LC/MS/MS Large-Volume Direct Injection

Eric Redman<sup>1</sup>, Simon Roberts<sup>2</sup>, KC Hyland<sup>2</sup>, Christopher Borton<sup>2</sup>, Craig Butt<sup>3</sup>, and Scott Krepich<sup>4</sup>

<sup>1</sup>Test America, Laboratories Sacramento, CA

<sup>2</sup>SCIEX, Redwood City, CA

<sup>3</sup>SCIEX, Framingham, MA

<sup>4</sup>Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

### Overview<sup>1</sup>

Human exposure to PFAS residues has been implicated in the incidence of cancer, obesity, endocrine system disruption, and other adverse health effects. In recognition of these potential risks, sources of human exposure to these chemicals (e.g., via drinking water) are receiving public and scientific attention.

This application note presents a method for the quantitation of per- and polyfluorinated alkyl substances (PFASs) in water samples. Presented here utilizes dilution of a water sample in methanol and direct injection of 950  $\mu$ L of the diluted sample using a 17.5 minute HPLC gradient. The method achieved accurate quantitation at levels of approximately 1-10 ng/L for more than 17 PFASs. Water samples were obtained anonymously from various sources in the United States. Samples were stored in the dark at 4 °C in 250 mL high density polyethylene bottles until analysis. A Phenomenex Luna<sup>®</sup> 5  $\mu$ m C18(2) 30 x 2 mm column (00A-4252-B0) was installed between the pump mixing chamber and the column, outside of a column oven. This column served as a delay or hold-up column to isolate PFAS contamination originating from the pumps and eluents.

### Acknowledgement

Phenomenex acknowledges Test America (Sacramento, CA) for collaborating with SCIEX and Phenomenex to contribute this application.

### Reference


1. For complete application including Mass Spec parameters and sample prep steps, please refer to: "Quantitation of PFASs in Water Samples using LC/MS/MS: Large-Volume Direct Injection and Solid Phase Extraction" at [www.Sciex.com](http://www.Sciex.com)

### LC/MS/MS Conditions

<b>Column:</b>	Gemini <sup>®</sup> 3 $\mu$ m C18	
<b>Dimensions:</b>	100 x 3.0 mm	
<b>Part No.:</b>	00D-4439-Y0	
<b>Mobile Phase:</b>	A: 20 mM Ammonium Acetate in Water B: Methanol	
<b>Gradient:</b>	<b>Time (min)</b>	<b>% B</b>
	0	10
	1.5	65
	8	95
	8.1	99
	12	99
	12.5	10
<b>Injection:</b>	0.950 $\mu$ L	
<b>Flow Rate:</b>	0.6 mL/min	
<b>Temperature:</b>	40 °C	
<b>Detection:</b>	SCIEX Triple Quad <sup>™</sup> 5500 with a Turbo V <sup>™</sup> source	



10 ng/L spike into groundwater matrix diluted with methanol



App ID 23901



# APPLICATION

## Ordering Information

3µm Microbore, Minibore and MidBore™ Columns (mm)										SecurityGuard™ Cartridges (mm)
Phases	50 x 1.0	20 x 2.0	30 x 2.0	50 x 2.0	100 x 2.0	150 x 2.0	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0*
Gemini® C18	00B-4439-A0	00M-4439-B0	00A-4439-B0	00B-4439-B0	00D-4439-B0	00F-4439-B0	00B-4439-Y0	00D-4439-Y0	00F-4439-Y0	AJO-7596

for ID: 2.0-3.0mm

5µm Microbore and Minibore Columns (mm)								SecurityGuard™ Cartridges (mm)
Phases	50 x 1.0	150 x 1.0	250 x 1.0	30 x 2.0	50 x 2.0	150 x 2.0	250 x 2.0	4 x 2.0*
Luna® C18(2)	00B-4252-A0	00F-4252-A0	00G-4252-A0	00A-4252-B0	00B-4252-B0	00F-4252-B0	00G-4252-B0	AJO-4286

for ID: 2.0-3.0mm

### Australia

t: +61 (0)2-9428-6444  
t: +61 (0)2-9428-6445  
a: auinfo@phenomenex.com

### Austria

t: +43 (0)1-319-1301  
t: +43 (0)1-319-1300  
a: anfrage@phenomenex.com

### Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
f: +31 (0)30-2383749  
a: beinfo@phenomenex.com

### Canada

t: +1 (800) 543-3681  
f: +1 (310) 328-7768  
a: info@phenomenex.com

### China

t: +86 (0)20 2282-6668  
f: +86 (0)20 2809-8130  
a: chinainfo@phenomenex.com

### Denmark

t: +45 4824 8048  
f: +45 4810 6265  
a: nordicinfo@phenomenex.com

### Finland

t: +358 (0)9 4789 0063  
f: +45 4810 6265  
a: nordicinfo@phenomenex.com

### France

t: +33 (0)1 30 09 21 10  
f: +33 (0)1 30 09 21 11  
a: franceinfo@phenomenex.com

### Germany

t: +49 (0)6021-58830-0  
f: +49 (0)6021-58830-11  
a: anfrage@phenomenex.com

### India

t: +91 (0)40-3012 2400  
f: +91 (0)40-3012 2411  
a: indiainfo@phenomenex.com

### Ireland

t: +353 (0)1 247 5405  
f: +44 1625-501796  
a: eireinfo@phenomenex.com

### Italy

t: +39 051 6327511  
f: +39 051 6327555  
a: italiainfo@phenomenex.com

## www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

### Luxembourg

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
a: nlinfo@phenomenex.com

### Mexico

t: 01-800-844-5226  
f: 001-310-328-7768  
a: tecnicomx@phenomenex.com

### The Netherlands

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
a: nlinfo@phenomenex.com

### New Zealand

t: +64 (0)9-4780951  
f: +64 (0)9-4780952  
a: nzinfo@phenomenex.com

### Norway

t: +47 810 02 005  
f: +45 4810 6265  
a: nordicinfo@phenomenex.com

### Puerto Rico

t: +1 (800) 541-HPLC  
f: +1 (310) 328-7768  
a: info@phenomenex.com

### Spain

t: +34 91-413-8613  
f: +34 91-413-2290  
a: espinfo@phenomenex.com

### Sweden

t: +46 (0)8 611 6950  
f: +45 4810 6265  
a: nordicinfo@phenomenex.com

### United Kingdom

t: +44 (0)1625-501367  
f: +44 (0)1625-501796  
a: ukinfo@phenomenex.com

### USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
a: info@phenomenex.com

### All other countries Corporate Office USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
a: info@phenomenex.com



If Luna analytical columns do not provide at least an equivalent separation as compared to a competing column of the same particle size, similar phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at <http://www.phenomenex.com/TermsAndConditions>.

### Trademarks

Luna and Gemini are registered trademarks, MidBore and SecurityGuard are trademarks of Phenomenex. Triple Quad is a trademark of AB SCIEX Pte. Ltd. AB SCIEX™ is being used under license.

© 2016 Phenomenex, Inc. All rights reserved.