

## TN-1309

# Nano-flow Lipidomics for Characterization of Lipid Mediators at Trace Levels Using a Biozen™ Peptide XB-C18 Column


Roxana Eggleston-Rangel, Jason Anspach, PhD, Namrata Saxena, and Bryan Tackett, PhD  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

### Introduction

The initial signs of inflammation are triggered by specialized eicosanoids, such as prostaglandins and leukotrienes, which in turn stimulate recruitment of neutrophils. Arachidonic acid (AA) derived lipoxins (LX), and E & D-series resolvins (RvD) that regulate the migration of immune cells, are responsible for the release of cytokines and are also involved in antibody generation. As a result of years of research involving metabolomics and lipidomics, it is now recognized that a novel class of lipid mediators derived mainly from dietary polyunsaturated fatty acids (PUFAs) are involved in signal transduction pathways that are crucial for the regulation and termination of inflammation. These mediators are conjointly called as Specialized Pro-resolving lipid Mediators (SPMs), which include classes of compounds such as “resolvins,” “lipoxins,” “maresins,” and “protectins.”

LC-MS is commonly used to detect and quantify such lipid mediators. However, their detection is challenging due to their physical properties and their bioactive endogenous nanogram to picogram concentrations. In this technical note, we provide a possible solution to these known challenges by using a unique reversed phase core-shell based column chemistry in a miniaturized column format (nano).

**Figure 1.** Structures of Positional Isomers Lipoxin A4 (Top) and Lipoxin B4 (Bottom).



### LC Conditions

<b>Column:</b>	Biozen 2.6 μm Peptide XB-C18 Biozen 3 μm Peptide Polar C18	
<b>Dimension:</b>	250 x 0.075 mm	
<b>Part No.:</b>	<a href="#">00G-4768-AW-21</a> (XB-C18) <a href="#">00G-4782-AW-21</a> (Polar C18)	
<b>Mobile Phase:</b>	A: 0.1 % Formic Acid in Water B: 0.1 % Formic Acid in Acetonitrile	
<b>Gradient:</b>	<b>Time (min)</b>	<b>%B</b>
	0.00	45
	2.00	45
	16.5	80
	16.6	98
	18.5	98
	20.5	10
	30.0	10
<b>Flow Rate:</b>	350 nL/min	
<b>Injection:</b>	1 μL	
<b>Temperature:</b>	40 °C	
<b>Detector:</b>	Q Exactive™ Plus Orbitrap™	
<b>System:</b>	NanoLC™ 425 (SCIEX®)	
<b>Detection:</b>	MS	

### MS Conditions

<b>Scan Type:</b>	Full MS SIM
<b>Polarity:</b>	Negative
<b>Resolution:</b>	140,000
<b>Scan Range:</b>	250 to 500 m/z
<b>AGC Target:</b>	3e6
<b>Maximum IT:</b>	200 ms



## Results and Discussion

Typically, lipids are quite challenging to separate via liquid chromatography due to variations in their geometries, fatty acyl chains, linkages, modifications, and existence of isomeric and isobaric species. A panel of 35 lipid mediators from 5 different lipid mediator classes (**Table 1**), were assayed on a Biozen™ 2.6 µm Peptide XB-C18 column to demonstrate the chromatographic separation performance of the core-shell particle.

The Biozen Peptide XB-C18 stationary phase employs core-shell technology that reduces the analyte's diffusion path, thus, delivers narrower peaks and in turn significantly increases and improves the ion intensities for the analytes. This stationary phase also has protective isobutyl side chains which reduce secondary interactions between analyte molecules and the residual silanol groups present on the stationary phase, thus giving improved peak shapes.


**Figure 2** shows the Extracted Ion Chromatogram (XIC) for the 5 classes of lipid mediators on the Biozen Peptide XB-C18 column. To better understand the effect and benefit of core-shell particle, peak widths and peak capacities were calculated. As expected, the Biozen Peptide XB-C18 core-shell column delivered narrower peak widths, increased peak capacities, and peak ion intensities that assists in achieving the trace detection levels of nanogram and picogram concentrations, as compared to the fully porous Luna™ Omega 3 µm Polar C18 column (**Table 2 & Figure 3**).

Interestingly, the Biozen Peptide XB-C18 column chemistry provided resolution of isobaric species of 3 Eicosanoids ((15)-S-HETE, (12)-S-HETE and 5(S)-HETE) as well as Lipoxins LXA<sub>4</sub> and LXB<sub>4</sub> (**Figures 4 and 5**). These two groups of lipid mediators belonging to the Lipoxin and COX & LOX lipid family, respectively, and share the same molecular formula and m/z making it difficult to separate by mass spectrometry alone. By using core-shell XB-C18 chemistry, these isobaric compounds were fully resolved.

**Table 1.** Lipid Mediator Analytes Used and their Respective M-H Species.

Lipid Type	Analyte - m/z
<b>Polyunsaturated Fatty Acids</b>	Arachidonic Acid - 303.2330[M-H]
	Dihomo-Ω-Linolenic Acid - 305.2486[M-H]
	Docosahexaenoic Acid - 327.2330[M-H]
	Docosapentaenoic Acid - 329.2486[M-H]
	Eicosapentaenoic Acid - 301.2173[M-H]
	Linoleic Acid - 279.2330[M-H]
	α-Linolenic Acid - 279.2330[M-H]
	Ω-Linolenic Acid - 277.2173[M-H]
Stearidonic Acid - 275.2017[M-H]	
<b>COX &amp; LOX</b>	6-keto Prostaglandin F1a - 369.2283[M-H]
	Thromboxane B2 - 369.2283[M-H]
	Prostaglandin F2a - 353.2333[M-H]
	Prostaglandin E2 - 351.2177[M-H]
	Prostaglandin D2 - 351.2177[M-H]
	12(S)-HHTrE - 279.1966[M-H]
	15(S)-HETE - 319.2279[M-H]
12(S)-HETE - 319.2279[M-H]	
5(S)-HETE - 319.2279[M-H]	
<b>SPM-E</b>	Resolvin E1 - 349.2020[M-H]
	(±)18-HEPE - 317.2122[M-H]
	Eicosapentaenoic Acid - 301.2173[M-H]
<b>Lipoxin</b>	LXB4 - 351.2177[M-H]
	LXA4 - 351.2177[M-H]
	15(R)-Lipoxin A4 - 351.2177[M-H]
	Arachidonic Acid - 303.2330[M-H]
<b>SPM-D</b>	Resolvin D3 - 375.2177[M-H]
	17(R)-Resolvin D1 - 375.2177[M-H]
	Resolvin D1 - 375.2177[M-H]
	Resolvin D5 - 359.2228[M-H]
	Docosahexaenoic Acid - 327.2330[M-H]
<b>Labeled Standards</b>	(±)5(6)-EET - 319.2279[M-H]
	(±)11,12-EpETrE - 319.2279[M-H]
	(±)14,15-EET - 319.2279[M-H]
	(±)11,12-EET-d11 - 330.2969[M11D-H] (heavy)
	PGE2-d4 - 355.2428[M4D-H] (heavy)
	LTB4-d4 - 339.2479[M4D-H] (heavy)
	RvE1-d4 - 353.2272[M4D-H] (heavy)
	Lipoxin A4-d5 - 356.2491[M5D-H] (heavy)


**Figure 2.** Extracted Ion Chromatogram (XIC) of 35 Lipids Lipid Mediators and their Metabolites on a Biozen 2.6 µm Peptide XB-C18 column.




**Table 2.** Average Peak Widths and Peak Capacities for 35 Lipid Molecules.

Column	Avg Peak Width	Peak Capacity
Biozen 2.6 μm Peptide XB-C18	0.43	39
Biozen 3 μm Peptide Polar C18	0.59	29


**Figure 3.** Ion Intensity of 10 Polyunsaturated Fatty Acids when Injecting 1 μL of Sample.



**Figure 4.** XIC of Fully Resolved Isobaric Species 3 Eicosanoids: (15)-S-HETE (Left) , (12)-S-HETE (Middle) and (5(S))-HETE (Right) on a Biozen 2.6 μm Peptide XB-C18 Column. Retention Times and Mass Errors are Shown.



**Figure 5.** XIC of Two Lipoxin Isobaric Species, LXB<sub>4</sub> (Left) and LXA<sub>4</sub> (Right) on a Biozen 2.6 μm Peptide XB-C18 Column. Respective Retention Times and Mass Errors are shown.



## Conclusions

The selected core-shell column chemistry of the Biozen Peptide XB-C18, with its ability to reduce analyte diffusion path and silanol activity, resulted in improved peak shapes and significantly increased intensities at low nanogram and picogram concentration levels. Additionally, its ability to fully resolve lipid isobaric species makes Biozen Peptide XB-C18 an appropriate selection for analyte characterization for lipidomics studies.



**Biozen™ Nano LC Columns with Integrated SecurityLINK™ Fitting**

1.6 µm Minibore Columns (mm)			
Phases	150 x 0.075 mm	250 x 0.075 mm	500 x 0.075 mm
Biozen 3 µm Peptide PS C18	<a href="#">00F-4771-AW-21</a>	<a href="#">00G-4771-AW-21</a>	–
Biozen 2.6 µm Peptide XB-C18	<a href="#">00F-4768-AW-21</a>	<a href="#">00G-4768-AW-21</a>	–
Biozen 3 µm Polar C18	<a href="#">00F-4782-AW-21</a>	<a href="#">00G-4782-AW-21</a>	–
Biozen 5 µm Peptide XB-C18	–	–	<a href="#">00J-4792-AW-21</a>



**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/Chat](http://www.phenomenex.com/Chat).

**Australia**

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 5019 9707  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 6559 4364  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

🌐 **All other countries/regions  
Corporate Office USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Kinetex, Luna, SecurityLINK, and BE-HAPPY are trademarks of Phenomenex. Orbitrap and Q. Exactive are trademarks of Thermo Fisher Scientific. SCIEX is a registered trademark and NanoLC is a trademark of AB SCIEX Pte. Ltd.

**Disclaimer**

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Thermo Fisher Scientific.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.

