

APPLICATION

Choosing the Right UHPLC Column for Intact Peptide Analysis

Zeshan Aqeel and Jeff Layne
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Overview

For peptide analysis, including both intact peptides and peptide fragments arising from tryptic digests, reversed phase methods using C18 columns, followed by either UV or MS detection, are an integral component of any thorough characterization workflow. Under classical LC-UV method conditions using water/acetonitrile/TFA gradients, secondary silanol interactions are not problematic because TFA (trifluoroacetic acid) prevents strong secondary interactions between basic amino acids and residual silanols on the silica.

However, analysts are increasingly relying upon LC-MS to perform their peptide analyses, and TFA is rarely used in conjunction with LC-MS due to its strong signal-suppression properties. Thus, when performing peptide analysis by LC-MS using MS-compatible modifiers such as 0.1% formic acid, secondary interaction between basic amino acids in the peptide and the underlying silica can lead to broad peaks, peak tailing, and interfere with resolution between close-eluting peaks. To obtain optimal performance for peptide analysis using LC-MS, the ideal media should be highly efficient (narrow peaks) and highly inert (minimal peak tailing when using MS-compatible modifiers). In addition, a highly retentive media is often favorable, particularly for peptide mapping, to maximize retention of small, polar peptide fragments that may not have a high affinity for a typical C18 phase, particularly in the absence of TFA.

Luna® Omega HPLC and UHPLC products, ranging in size from 1.6 µm for UHPLC work up to 5 µm for conventional HPLC as well as preparative HPLC, are the ideal choice for the analysis of small peptide and tryptic mapping. The high surface area and dense C18 bonding provide Luna Omega with excellent retention for small, more polar peptides and fragments. The thermally modified silica (using a proprietary technology) gives Luna Omega an efficiency advantage over similarly-sized fully porous silica products, which can translate into narrower peaks, improved sensitivity, and improved resolution in many cases. This thermal treatment process also renders the silica highly inert for minimal secondary interactions, even when using LC-MS mobile phase modifiers such as 0.1 % formic acid. Lastly, the unique Luna Omega PS C18 chemistry contains a proprietary, positively-charged functional group on the surface that makes the surface even less likely to exhibit peak tailing when analyzing basic peptides, and can also provide a selectivity quite distinct from a standard C18 phase, which may be useful for separating a target peptide from a matrix interference that co-elutes on a traditional C18.

The three figures inside show the analysis of a basic peptide (first large peak; MW = 1046 g/mol) as well as two isobaric interferences (the two smaller peaks). Under these conditions, you can see that both traditional C18 phases (**Figure 1** = Waters® ACQUITY® BEH C18 1.7 µm, **Figure 2** = Luna Omega

1.6 µm C18) provide essentially identical performance, with similar peak widths, peak shape, and resolution of the isobaric interferences from the main peptide. Clearly, either product is well-suited to this analysis, and other factors such as column lifetime, pricing, and customer service and support would need to be factored into a final column decision.

In contrast, when using the Luna Omega 1.6 µm PS C18 column (**Figure 3**), we get a significant improvement in the separation of the two isobaric interferences from the main peptide peak. As the underlying silica is the same as the Omega C18 phase, we must attribute this shift in retention and selectivity to the reduced bonding density of the PS C18 phase and the presence of the proprietary positive surface functional group. Clearly, under these conditions, the Luna Omega PS C18 provides the best performance of the three columns.

When developing methods for the analysis of peptides using either LC-UV or LC-MS conditions, Luna Omega C18 and Luna Omega PS C18 should be evaluated as part of any thorough method development process.





Figure 1.
Waters® ACQUITY® BEH 1.7µm C18



Conditions for all columns

Column: As specified
Dimension: 50 x 2.1 mm
Mobile Phase: A: 0.1 % Formic acid in Water
 B: 0.1 % Formic acid in Acetonitrile
Gradient:

Time (min)	% B
0	5
4	95

Flow Rate: 0.5 mL/min
Temperature: 50 °C
Detection: MS/MS (SCIEX API 4000™)
Sample: 1. Angiotensin II
 2. Isobaric Interference 1
 3. Isobaric Interference 2

Figure 2.
Luna® Omega 1.6µm C18





Figure 3.
Luna Omega 1.6µm PS C18



Phenomenex is not affiliated with Waters Technologies Corporation. Comparative separations may not be representative of all applications.

Luna[®] Omega Ordering Information

1.6 µm Microbore Columns (mm)			
Phases	50 x 1.0	100 x 1.0	150 x 1.0
Polar C18	00B-4748-A0	00D-4748-A0	00F-4748-A0
C18	00B-4742-A0	00D-4742-A0	00F-4742-A0

1.6 µm Minibore Columns (mm)					SecurityGuard [™] ULTRA Cartridges [‡]
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
Polar C18	00A-4748-AN	00B-4748-AN	00D-4748-AN	00F-4748-AN	AJO-9505
PS C18	00A-4752-AN	00B-4752-AN	00D-4752-AN	00F-4752-AN	AJO-9508
C18	00A-4742-AN	00B-4742-AN	00D-4742-AN	00F-4742-AN	AJO-9502

for 2.1 mm ID

3 µm Minibore Columns (mm)					SecurityGuard [™] Cartridges (mm)
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	4 x 2.0*
Polar C18	00A-4760-AN	00B-4760-AN	00D-4760-AN	00F-4760-AN	AJO-7600
PS C18	00A-4758-AN	00B-4758-AN	00D-4758-AN	00F-4758-AN	AJO-7605

for ID: 2.0 - 3.0 mm

3 µm MidBore [™] Columns (mm)				SecurityGuard [™] Cartridges (mm)
Phases	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0*
Polar C18	00B-4760-Y0	00D-4760-Y0	00F-4760-Y0	AJO-7600
PS C18	00B-4758-Y0	00D-4758-Y0	00F-4758-Y0	AJO-7605

for ID: 2.0 - 3.0 mm

3 µm Analytical Columns (mm)					SecurityGuard [™] Cartridges (mm)
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*
Polar C18	00B-4760-E0	00D-4760-E0	00F-4760-E0	00G-4760-E0	AJO-7601
PS C18	00B-4758-E0	00D-4758-E0	00F-4758-E0	00G-4758-E0	AJO-7606

for ID: 3.2-8.0 mm

5 µm Minibore Columns (mm)					SecurityGuard [™] Cartridges (mm)
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	4 x 2.0*
Polar C18	00A-4754-AN	00B-4754-AN	00D-4754-AN	00F-4754-AN	AJO-7600
PS C18	00A-4753-AN	00B-4753-AN	00D-4753-AN	00F-4753-AN	AJO-7605

for ID: 2.0 - 3.0 mm

5 µm MidBore [™] Columns (mm)				SecurityGuard [™] Cartridges (mm)
Phases	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0*
Polar C18	00B-4754-Y0	00D-4754-Y0	00F-4754-Y0	AJO-7600
PS C18	00B-4753-Y0	00D-4753-Y0	00F-4753-Y0	AJO-7605

for ID: 2.0 - 3.0 mm

5 µm Analytical Columns (mm)					SecurityGuard [™] Cartridges (mm)
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*
Polar C18	00B-4754-E0	00D-4754-E0	00F-4754-E0	00G-4754-E0	AJO-7601
PS C18	00B-4753-E0	00D-4753-E0	00F-4753-E0	00G-4753-E0	AJO-7606

for ID: 3.2-8.0 mm

5 µm Axia [™] Packed Preparative Columns (mm)			SecurityGuard [™] Cartridges (mm)
Phases	150 x 21.2	250 x 21.2	15 x 21.2**
Polar C18	00F-4754-P0-AX	00G-4754-P0-AX	AJO-7603
PS C18	00F-4753-P0-AX	00G-4753-P0-AX	AJO-7608

for ID: 18-29 mm

5 µm Axia [™] Packed Preparative Columns (mm)				SecurityGuard [™] Cartridges (mm)
Phases	150 x 30	250 x 30	250 x 50	15 x 30.0*
Polar C18	00F-4754-U0-AX	00G-4754-U0-AX	00G-4754-V0-AX	AJO-7604
PS C18	00F-4753-U0-AX	00G-4753-U0-AX	00G-4753-V0-AX	AJO-7609

for ID: 30-49 mm

* SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000
 * SecurityGuard Analytical Cartridges require holder, Part No.: KJO-4282
 ** PREP SecurityGuard Cartridges require holder, Part No.: AJO-8223
 † PREP SecurityGuard Cartridges require holder, Part No.: AJO-8277



If Luna analytical columns do not provide at least an equivalent separation as compared to a competing column of the same particle size, similar phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND.



APPLICATION

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 400-606-8099
f: +86 (0)22 2532-1033
phen@agela.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain

t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at <http://www.phenomenex.com/TermsAndConditions>.

Trademarks

Luna is a registered trademark and SecurityGuard, Axia, and MidBore are trademarks of Phenomenex. Waters and ACQUITY are registered trademarks of Waters Technologies Corporation. API 4000 is a trademark of AB SCIEX Pte. Ltd. AB SCIEX is being used under license.

Disclaimer

Comparative separations may not be representative of all applications Phenomenex is not affiliated with Waters Technologies Corporation. Axia column and packing technology is patented by Phenomenex. U.S. Patent No. 7, 674, 383 SecurityGuard is patented by Phenomenex U.S. Patent No. 6,162,362

CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridges.

© 2017 Phenomenex, Inc. All rights reserved.