

# APPLICATIONS

## Validation of an Automated Solution for Post Hydrolysis Enzyme Removal and LC/MS/MS Analysis of Illicit and Pain Related Compounds from Urine

Shahana Wahab Huq<sup>1</sup>, Richard Thomas<sup>2</sup>, Agnes Cua<sup>2</sup>, Seyed Sadjadi<sup>1</sup>, Sean Orlowicz<sup>1</sup>

<sup>1</sup>Phenomenex, 411 Madrid Avenue, Torrance, CA 90501

<sup>2</sup>Precision Toxicology, 3030 Bunker Hill Street, San Diego, CA 92109

The escalating abuse of pain medicines has mandated that health care providers and practitioners routinely perform urine drug testing. Most laboratories opt for a dilute-and-shoot approach, along with enzymatic hydrolysis that requires no sample cleanup. The addition of extra protein as a result of hydrolysis using  $\beta$ -glucuronidase to the urine can result in fouling of the LC column and loss of productivity. In this work, we present an automated method that uses an Impact<sup>TM</sup> Protein Precipitation 96-well plate to remove a majority of the proteins from hydrolyzed urine.

### Introduction

The gradual dependence and addiction of pain medication has made pain management and monitoring one of the fastest growing clinical market segments today. The escalating abuse of pain medicines has mandated that health care providers and practitioners routinely perform urine drug testing using a reliable approach that is fast and cost effective. To meet the demand for a large accrual of samples, most laboratories adopt a dilute-and-shoot approach that requires virtually no sample cleanup. However, the samples must undergo some form of hydrolysis procedure to de-conjugate the metabolized compounds back to their native form. An enzymatic procedure using  $\beta$ -glucuronidase is the most readily accepted form of hydrolysis. However, the addition of extra protein to the urine can result in plugging or otherwise fouling (Figure 1) of the LC column<sup>1</sup>. The continuous increase in system pressure resulting from denatured enzyme within the column will reduce the column and assay performance. In this work, we present a validated procedure that uses an Impact Protein Precipitation 96-well plate to remove the enzyme (and other proteins) from urine (Figure 2). To increase productivity and accuracy, we utilized a Tecan Freedom EVO<sup>®</sup> 100 Liquid Handling System (or similar) in conjunction with a Kinetex<sup>®</sup> 2.6  $\mu$ m Phenyl-Hexyl 50 x 4.6 mm core-shell HPLC column. The chromatographic conditions were adequately efficient to accommodate fifty one (51) pain panel compounds in less than six (6) minutes.

### Experimental Conditions

#### HPLC Conditions

<b>Column:</b>	Kinetex 2.6 $\mu$ m Phenyl-Hexyl
<b>Dimensions:</b>	50 x 4.6 mm
<b>Part No.:</b>	00B-4495-E0
<b>Mobile Phase:</b>	A: 0.1 % Formic acid in Water B: 0.1 % Formic acid in Methanol
<b>Gradient:</b>	<b>Time (min)</b> <b>% B</b>
	0.0            5
	5.5            95
	6.0            95
	6.01          5
	6.5            5
<b>Flow Rate:</b>	1.0 mL/min
<b>Injection Volume:</b>	5 $\mu$ L
<b>Instrument:</b>	Shimadzu <sup>®</sup> XR
<b>Detection:</b>	MS/MS, SCIEX QTRAP <sup>®</sup> 6500 (Positive Ionization)

Sample Preparation Automation: Tecan Freedom EVO<sup>®</sup> 100 with MultiChannel Arm<sup>™</sup> 96 (MCA96)

Note: Additional equilibration time is included in the acquisition method



**Sean Orlowicz**  
 Manager, PhenoLogix  
 When not in the lab, Sean enjoys just about anything involving the outdoors; hiking, climbing, surfing, etc. He is especially at home in the mountains, being an avid skier and motorcyclist.

### Sample Preparation

#### Sample Hydrolysis Procedure

A 500  $\mu$ L sample of urine was diluted with 100  $\mu$ L acetate buffer (pH 4.5-4.8) and 20  $\mu$ L  $\beta$ -glucuronidase, 106 units per mL (DR2100, www.campbellscience.com) in a 96-well collection plate. The samples were vortexed for 10-15 seconds and then incubated in a water bath at 63  $^{\circ}$ C for 30 minutes.


#### Dilute-and-Shoot Protocol

The hydrolyzed samples were sealed and centrifuged for 10 minutes at 2000 rpm (or the maximum possible speed by the centrifuge). The supernatant was then transferred to a LC/MS/MS for analysis.


#### Protein Precipitation

A 100  $\mu$ L volume of the hydrolyzed sample was loaded directly to an Impact (2 mL Square Well Filter Plate, Part No. CE0-7565) Protein Precipitation 96-well plate that had been pre-loaded with 400  $\mu$ L methanol. The plate was sealed and then vortexed for 2 minutes at the maximum possible speed. A vacuum of 2-7" of Hg was applied for 2-3 minutes until filtrate was collected. The resulting extract was then evaporated to dryness and reconstituted in starting mobile phase before being transferred for LC/MS/MS analysis.


**Figure 1.**  
 Number of injections vs. backpressure for dilute-and-shoot samples




**Figure 2.**  
Number of injections vs. backpressure for protein precipitated samples




**Figure 3.**  
LC/MS/MS analysis of 51 compounds in less than 6 minutes




**Figure 4.**  
Calibration curve for benzodiazepines




**Figure 5.**  
Calibration curve for opiates



**Figure 6.**  
Calibration curve for amphetamines



**Figure 7.**  
Calibration curve for analgesics



**Figure 8.**  
Matrix effect summary

Compound	% Matrix Effect	Compound	% Matrix Effect	Compound	% Matrix Effect
6-MAM	97.39	MDMA	105.72	Oxymorphone	106.48
$\alpha$ -OH-Alprazolam	105.88	Meperidine	100.88	PCP	96.62
Alprazolam	103.71	Meprobamate	109.15	Propoxyphene	94.88
Amitriptyline	101.26	Methadone	104.77	Temazepam	102.66
Amphetamine	103.66	Methamphetamine	96.82	Tramadol	102.9
Benzoyllecgonine	107.37	Morphine	98.97	Butalbital	95.53
Buprenorphine	118.52	Naloxone	106.81	Phenobarbital	97.9
Carisoprodol	113.79	Norbuprenorphine	102.84	Methylphenidate	101.19
Clonazepam	88.85	Nordiazepam	101.99	Desipramine	100.38
Codeine	110.88	Norfentanyl	104.61	Imipramine	99.75
EDDP	98.41	Norhydrocodone	105.19	7-Aminoclonazepam	120.89
Fentanyl	102.21	Nortriptyline	102.04	Zolpidem	109.61
Hydrocodone	98.5	Oxazepam	103.37	Gabapentin	98.58
Hydromorphone	103.56	Oxycodone	100.4	THCA	90.26

**Figure 9. Continued**  
Lower and upper limits of quantitation (ng/mL)

Drug Name	LLOQ	ULOQ
Norfentanyl	0.5	100
Norhydrocodone	10	2000
Noroxycodone	5	1000
Nortriptyline	10	2000
Oxazepam	5	1000
Oxycodone	5	1000
Oxymorphone	5	1000
PCP	2.5	500
Propoxyphene	10	2000
Temazepam	5	1000
Tramadol	5	1000
Butalbital	12.5	2500
Phenobarbital	12.5	2500
Tapentadol	5	1000
Methylphenidate	5	1000
Desipramine	5	1000
Imipramine	5	1000
7-Aminoclonazepam	5	1000
Duloxetine	25	1000
Zolpidem	5	1000
MDPV	5	1000
Mephedrone	5	1000
Methylone	5	1000
Gabapentin	80	16000
Sovaldi	10	2000
THCA	12.5	500

**Figure 9.**  
Lower and upper limits of quantitation (ng/mL)

Drug Name	LLOQ	ULOQ
6-MAM	2.5	500
$\alpha$ -OH-Alprazolam	5	1000
Alprazolam	5	1000
Amitriptyline	10	2000
Amphetamine	10	2000
Benzoyllecgonine	5	1000
Buprenorphine	1.5	300
Carisoprodol	10	2000
Clonazepam	5	1000
Codeine	10	2000
Cyclobenzaprine	10	2000
EDDP	10	2000
Fentanyl	0.5	100
Hydrocodone	5	1000
Hydromorphone	5	1000
Lorazepam	5	1000
MDMA	10	2000
Meperidine	5	1000
Meprobamate	10	2000
Methadone	10	2000
Methamphetamine	10	2000
Morphine	5	1000
Naloxone	10	2000
Norbuprenorphine	2.5	500
Nordiazepam	5	1000

**Figure 10.**  
Method precision for high conc. (replicates = 12)

Analyte	Conc. (ng/mL)	Within-Run Precision (%)	Within-Laboratory Precision (%)	Interday Precision (%)	Intraday Precision (%)
6-MAM	250	1.67	2.15	1.15	0.73
$\alpha$ -OH-Alprazolam	500	2.68	3.53	2.33	0.4
Alprazolam	500	2.62	2.41	1.2	1.57
Amitriptyline	1000	7.29	8.23	2.64	2.75
Amphetamine	1000	3.2	2.79	1.49	2.16
Benzoylcegonine	500	2.17	3.47	2.95	1.17
Buprenorphine	150	12.23	12.23	6.12	6.12
Carisoprodol	1000	6.19	15.09	12.3	6.17
Clonazepam	500	3.37	3.95	2.55	1.49
Codeine	1000	3.19	4.84	3.78	1.05
EDDP	1000	5.79	9.73	5.34	5.72
Fentanyl	50	3.16	4.46	2.88	1.26
Hydrocodone	500	2.36	2.31	0.82	0.96
Hydromorphone	500	1.99	3.52	3.18	1.3
Lorazepam	500	3.4	3.59	1.22	0.46
MDMA	1000	1.91	2.99	2.27	0.39
Meperidine	500	3.04	3.52	2.52	1.79
Meprobamate	1000	6.25	6.23	3.38	3.34
Methadone	1000	5.29	9.87	4.8	9.61
Methamphetamine	1000	2.84	6.29	4.94	2.66
Morphine	500	1.63	2.78	2.28	0.38
Naloxone	1000	2.92	2.61	0.6	1.45
Norbuprenorphine	250	5.89	5.48	3.14	3.81
Nordiazepam	500	2.63	2.41	0.99	0.41
Norfentanyl	50	2.49	2.67	1.93	1.67
Norhydrocodone	1000	1.75	1.58	0.68	1.02
Nortriptyline	1000	3.54	3.38	2.09	2.35
Oxazepam	500	2.92	2.76	0.38	0.88
Oxycodone	500	3.29	3.47	2.4	2.11
Oxymorphone	500	1.73	2.4	1.92	0.96
PCP	250	5.04	8.56	1.27	6.8
Propoxyphene	1000	8.76	14.67	10.16	5.95
Temazepam	500	2.91	3.87	1.8	1.81
Tramadol	500	1.74	2.77	2.39	1.03
Butalbital	1250	2.12	4.31	1.36	3.99
Phenobarbital	1250	4.57	6.79	4.88	1.18
Desipramine	500	2.95	2.98	1.71	1.65
Imipramine	500	2.64	4.68	3.79	0.78
7-Aminoclonazepam	500	2.92	3.58	2.38	1.16
Zolpidem	500	3.06	3.73	2.61	1.5
MDPV	500	3.01	3.48	2.43	1.69
Gabapentin	8000	2.84	3.07	2.23	1.9
THCA	250	1.37	5.56	5.19	1.46

# APPLICATIONS

**Figure 11.**  
Method precision for low conc. (replicates = 12)

Analyte	Conc. (ng/mL)	Within-Run Precision (%)	Within-Laboratory Precision (%)	Interday Precision (%)	Intraday Precision (%)
6-MAM	25	4.31	6.15	4.21	1.24
$\alpha$ -OH-Alprazolam	50	2.32	5.58	4.95	1.12
Alprazolam	50	2.92	5.62	4.96	1.23
Amitriptyline	100	5.8	8.08	3.62	4.3
Amphetamine	100	3.27	13.93	13.38	2.08
Benzoylcegonine	50	3.26	6.27	5.37	0.45
Buprenorphine	15	16.06	15.38	4.33	6.33
Carisoprodol	100	10.31	17.58	13.72	3.8
Clonazepam	50	4.92	5.15	3.6	3.26
Codeine	100	3.14	4.28	3.29	1.53
EDDP	100	5.64	5.94	4.01	3.56
Fentanyl	5	4.47	4.47	0	0
Hydrocodone	50	2.57	11.19	10.85	1
Hydromorphone	50	2.93	5.04	4.38	1.55
Lorazepam	50	3	4.27	2.66	1.48
MDMA	100	1.8	2.95	2.6	1.14
Meperidine	50	2.61	5.02	4.49	1.34
Meprobamate	100	8.75	15.31	3.13	12.16
Methadone	100	5.76	6.93	3.38	5.13
Methamphetamine	100	2.66	3.99	3.27	1.34
Morphine	50	3.29	5.01	4.21	1.84
Naloxone	100	3.57	3.51	1.85	1.97
Norbuprenorphine	25	14.14	14.38	7.62	7.16
Nordiazepam	50	3.49	5.03	3.63	0
Norfentanyl	5	2.16	3.28	2.68	1.02
Norhydrocodone	100	3.75	3.56	2.22	2.5
Nortriptyline	100	4.92	5.7	4.46	3.41
Oxazepam	50	2.79	9.64	8.95	2.28
Oxycodone	50	3.1	4.54	3.55	1.26
Oxymorphone	50	3.55	4.95	4.14	2.28
PCP	25	9.8	8.29	3.35	4
Propoxyphene	100	12.18	13.32	4.36	6.94
Temazepam	50	2.32	4.63	3.53	1.9
Tramadol	50	2.97	4.75	4.14	1.84
Butalbital	125	5.08	6.28	3.68	5.22
Phenobarbital	125	0	6.79	4.88	1.18
Methylphenidate	50	4.77	5.01	3.42	3.07
Desipramine	50	3.66	4.96	3.77	1.73
Imipramine	50	3	7.05	6.46	1
7-Aminoclonazepam	50	2.72	11.62	11.31	0.45
Zolpidem	50	3.16	5.2	4.5	1.79
MDPV	50	2.97	4.33	3.54	1.61
Gabapentin	800	1.98	5.19	4.97	1.29
THCA	25	2.68	7.67	6.8	2.37

## Results and Discussion

1. The presence of additional protein ( $\beta$ -glucuronidase) in a urine sample can denature in the column during a normal chromatographic run. The resulting fouling will quickly plug the HPLC column and render it useless (**Figure 1**). However, removal of the protein from the samples by a simple precipitation procedure can enhance the longevity of the column (**Figure 2**).
2. The protein precipitation method evaluated here produced good chromatographic separation of the compounds (**Figure 3**) with acceptable calibration parameters (**Figures 4-7**).
3. The method produced minimal matrix effect despite a simple sample preparation procedure from an inherently dirty sample matrix (**Figure 8**).
4. In addition, great method precision was achieved at both low and high quantitation levels (**Figures 10-11**).

## Conclusion

- Performing protein precipitation after enzymatic hydrolysis provides cleaner extracts and extends the HPLC column lifetime.
- Fifty one (51) compounds and fourteen (14) deuterated analogs were efficiently and quickly resolved in 6 minutes using a Kinetex<sup>®</sup> Phenyl-Hexyl core-shell HPLC column.
- The sample preparation proved to be robust requiring no method development.
- This method produced excellent calibration range.
- The calculated matrix effect results for all the evaluated analytes fell within acceptable limits.
- This bioanalytical method is proven to be suitable for a wide range of acidic, basic, and neutral compounds.

## References

1. S Huq, S Sadjadi, and S Countryman, "Removal of Beta-Glucuronidase Enzyme from Urine Post-Hydrolysis to Improve Assay Performance and Column Lifetime." Mass Spec Application for Clinical Laboratory Conference, 2013

# APPLICATIONS

## Ordering Information

### Impact™ Precipitation Products

Part No.	Description	Unit
<b>Impact Precipitation Products</b>		
<b>CE0-7565</b>	Impact Protein Precipitation, Square Well, Filter Plate, 2 mL	2/pk
<b>CE0-7566</b>	Impact Protein Precipitation, Square Well, Long Drip, Filter Plate, 2 mL	2/pk
<b>Impact Starter Kit for Protein Precipitation</b>		
<b>CE0-8201</b>	Impact Protein Precipitation Plate (CE0-7565) (2 ea) Collection Plate 2 mL (2 ea) Sealing Mat, Santoprene™ (AH0-8199) (2 ea)	ea

### Kinetex® HPLC Columns

Phases	2.6 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges†
	30 x 2.1	50 x 2.1	75 x 2.1	100 x 2.1	150 x 2.1	3/pk
<b>EVO C18</b>	00A-4725-AN	00B-4725-AN	—	00D-4725-AN	00F-4725-AN	AJO-9298
<b>F5</b>	00A-4723-AN	00B-4723-AN	—	00D-4723-AN	00F-4723-AN	AJO-9322
<b>Biphenyl</b>	00A-4622-AN	00B-4622-AN	—	00D-4622-AN	00F-4622-AN	AJO-9209
<b>XB-C18</b>	00A-4496-AN	00B-4496-AN	00C-4496-AN	00D-4496-AN	00F-4496-AN	AJO-8782
<b>C18</b>	00A-4462-AN	00B-4462-AN	00C-4462-AN	00D-4462-AN	00F-4462-AN	AJO-8782
<b>C8</b>	00A-4497-AN	00B-4497-AN	00C-4497-AN	00D-4497-AN	00F-4497-AN	AJO-8784
<b>HILIC</b>	00A-4461-AN	00B-4461-AN	00C-4461-AN	00D-4461-AN	00F-4461-AN	AJO-8786
<b>Phenyl-Hexyl</b>	00A-4495-AN	00B-4495-AN	00C-4495-AN	00D-4495-AN	00F-4495-AN	AJO-8788

for 2.1 mm ID

Phases	2.6 µm MidBore™ Columns (mm)					SecurityGuard™ ULTRA Cartridges†
	30 x 3.0	50 x 3.0	75 x 3.0	100 x 3.0	150 x 3.0	3/pk
<b>EVO C18</b>	—	00B-4725-YO	—	00D-4725-YO	00F-4725-YO	AJO-9298
<b>F5</b>	—	00B-4723-YO	—	00D-4723-YO	00F-4723-YO	AJO-9321
<b>Biphenyl</b>	—	00B-4622-YO	—	00D-4622-YO	00F-4622-YO	AJO-9208
<b>XB-C18</b>	00A-4496-YO	00B-4496-YO	00C-4496-YO	00D-4496-YO	00F-4496-YO	AJO-8775
<b>C18</b>	00A-4462-YO	00B-4462-YO	00C-4462-YO	00D-4462-YO	00F-4462-YO	AJO-8775
<b>C8</b>	00A-4497-YO	00B-4497-YO	00C-4497-YO	00D-4497-YO	00F-4497-YO	AJO-8777
<b>HILIC</b>	00A-4461-YO	—	—	—	00F-4461-YO	AJO-8779
<b>Phenyl-Hexyl</b>	—	00B-4495-YO	—	00D-4495-YO	00F-4495-YO	AJO-8781

for 3.0 mm ID

Phases	2.6 µm Analytical Columns (mm)					SecurityGuard™ ULTRA Cartridges†
	30 x 4.6	50 x 4.6	75 x 4.6	100 x 4.6	150 x 4.6	3/pk
<b>EVO C18</b>	—	00B-4725-E0	—	00D-4725-E0	00F-4725-E0	AJO-9296
<b>F5</b>	—	00B-4723-E0	—	00D-4723-E0	00F-4723-E0	AJO-9320
<b>Biphenyl</b>	—	00B-4622-E0	—	00D-4622-E0	00F-4622-E0	AJO-9207
<b>XB-C18</b>	—	00B-4496-E0	00C-4496-E0	00D-4496-E0	00F-4496-E0	AJO-8768
<b>C18</b>	00A-4462-E0	00B-4462-E0	00C-4462-E0	00D-4462-E0	00F-4462-E0	AJO-8768
<b>C8</b>	—	00B-4497-E0	00C-4497-E0	00D-4497-E0	00F-4497-E0	AJO-8770
<b>HILIC</b>	—	00B-4461-E0	00C-4461-E0	00D-4461-E0	00F-4461-E0	AJO-8772
<b>Phenyl-Hexyl</b>	—	00B-4495-E0	00C-4495-E0	00D-4495-E0	00F-4495-E0	AJO-8774

for 4.6 mm ID

† SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000



If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to other products of the same phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.



# APPLICATIONS

**Australia**

t: +61 (0)2-9428-6444  
f: +61 (0)2-9428-6445  
auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
f: +43 (0)1-319-1300  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
f: +31 (0)30-2383749  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
f: +1 (310) 328-7768  
info@phenomenex.com

**China**

t: +86 (0)20 2282-6668  
f: +86 (0)20 2809-8130  
chinainfo@phenomenex.com

**Denmark**

t: +45 4824 8048  
f: +45 4810 6265  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
f: +33 (0)1 30 09 21 11  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
f: +49 (0)6021-58830-11  
anfrage@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
f: +91 (0)40-3012 2411  
indiainfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
f: +44 1625-501796  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
f: +39 051 6327555  
italiainfo@phenomenex.com

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

**Mexico**

t: 001-800-844-5226  
f: 001-310-328-7768  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
f: +64 (0)9-4780952  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
f: +45 4810 6265  
nordicinfo@phenomenex.com

**Puerto Rico**

t: +1 (800) 541-HPLC  
f: +1 (310) 328-7768  
info@phenomenex.com

**Spain**

t: +34 91-413-8613  
f: +34 91-413-2290  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
f: +45 4810 6265  
nordicinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
f: +44 (0)1625-501796  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com

**All other countries  
Corporate Office USA **

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Kinetex is a registered trademark and Impact, SecurityGuard, and MidBore are trademarks of Phenomenex. Shimadzu is a registered trademark of Shimadzu Corp. Freedom EVO is a registered trademark and MultiChannel Arm is a trademark of Tecan Group Ltd. QTRAP is a registered trademark of AB SCIEX Pte. Ltd. AB SCIEX is being used under license.

**Disclaimer**

The products mentioned are not intended for clinical use. Because they are not intended for clinical use, no claim or representation is made or intended for their clinical use (including, but not limited to diagnostic, prognostic, therapeutic or blood banking). It is the user's responsibility to validate the performance of Phenomenex products for any particular use, since the performance characteristics are not established. Phenomenex products may be used in clinical diagnostic laboratory systems after the laboratory has validated their complete system as required by the Clinical Laboratory Improvements Amendments of 1988 (CLIA '88) regulation in the U.S. or equivalent in other countries.

© 2015 Phenomenex, Inc. All rights reserved.