

The Chromatographer's Guide
**To Improving UHPLC
Column Selectivity**

 **phenomenex**[®]
...breaking with traditionSM

 www.phenomenex.com/UHPLC

Guide to Choosing the Most Effective Selectivity

Step 1

Understanding Solid Support and Efficiency

- p. 3 Selecting the Right Solid Support
- p. 4 How Resolution is Affected by Selectivity and Efficiency
- p. 5 Efficiency Gains with Luna® Omega and Kinetex®

Step 2

How to Use Selectivity

- p. 6 Profiling the Mechanisms of Selectivity
- p. 7 Relating Selectivity to UHPLC Stationary Phases

Step 3

Our UHPLC Selectivities

- p. 9 Alkyl Phases
- p. 10 Phenyl Phases
- p. 11 Polar Retentive Phases

Step 4

Applying Column Selectivities to Your UHPLC Analysis

- pp. 13-15..... Thinking Outside the Traditional C18 Phases
- pp. 16-17..... Improve Polar Separations with 100 % Aqueous Stable Phases
- p. 18 Utilizing Unique Selectivities
- p. 19 Aromatic Based Selectivity
- p. 20 Protect your Column's Selectivity

Step 5

Order Your New Selectivity Tools

- pp. 21-23

Select the Right Solid Support



Phenomenex offers a full range of solid supports including core-shell and thermally modified fully porous. The morphology of the solid support has a significant impact on the resulting material characteristics and column performance.

Core-Shell and Organo-Silica Core-Shell

Unique solid silica core and porous shell that results in faster chromatography and higher efficiencies than conventional fully porous particles.

Well suited for:

- Performance gains on ANY LC system
- Easy system-to-system and lab-to-lab method transfer
- Methods where increased sensitivity is required
- Significantly improving the productivity of older, established methods




Fully Porous – Thermally Modified Silica

Unique high efficiency and extremely robust fully porous silica that offers astounding performance and inertness alongside versatile selectivities.


Well suited for:

- Astounding UHPLC, HPLC, and Preparative HPLC performance and efficiencies
- Greater separation muscle
- Better peak shape through an inert foundation
- Extreme ruggedness and dependability



Thermal Modified Pore Structure

Most importantly, through our proprietary process, we eliminate micropores, further improving column efficiency, inertness, and reproducibility.





How Resolution is Affected by Selectivity and Efficiency

Selectivity (α) and efficiency (N) have the greatest impact on observed resolution (R), when compared to other chromatographic parameters. Often the simplest and most effective way to improve your chromatographic results is to change your column's phase or solid support. Phenomenex offers a wide breadth of phase chemistries across multiple solid supports for simplified method development and optimization.

Selectivity is the most important parameter for increasing resolution.


The Impact of Selectivity on Resolution



Efficiency Gains with Luna Omega and Kinetex


The undeniably high efficiency levels found in each Luna® Omega and Kinetex® column provide you with the potential of huge gains in method performance. While traditional silica and hybrid fully porous particles may claim high performance, when compared to Luna Omega or Kinetex, they fall short and prevent HPLC/UHPLC scientists from reaching their goals.

Efficiency Levels (plates/m)




UHPLC Performance – Cyclosporine Impurity Profile


Kinetex 2.6µm Polar C18



Luna Omega 1.6µm Polar C18



Conventional 1.7µm C18



Conditions for all columns same except where noted:


- Columns:** Kinetex 2.6 µm Polar C18
Luna Omega 1.6 µm Polar C18
Conventional Fully Porous 1.7 µm C18
- Dimensions:** 50 x 2.1 mm
- Mobile Phase:** Acetonitrile/Tert-butyl methyl ether/Water/
Phosphoric acid (430:50:520:1)
- Flow Rate:** 0.30 mL/min
- Temperature:** 80 °C
- Detection:** UV @ 210 nm
- Sample:** Cyclosporine

Comparative separations may not be representative of all applications.




Profiling the Mechanisms of Selectivity

Observed selectivity is dictated by several primary molecular interactions. Below you'll find selectivity parameters used to help characterize reversed phase selectivity mechanisms.




Hydrophobicity

The ability of a phase to hydrophobically interact with carbon groups




Steric Interaction

The ability of a phase to separate compounds based on structural differences




Hydrogen Bond Donating Capacity

The ability of a phase to hydrogen-bond with proton accepting groups



Hydrogen Bond Accepting Capacity

The ability of a phase to hydrogen-bond with proton donating groups



Cation Selectivity at pH 2.8

The ability of a phase to interact with cation groups at acidic pH


Cation Selectivity at pH 7.0

The ability of a phase to interact with cation groups at neutral pH

Relating Selectivity to UHPLC Stationary Phases

Follow the colors to connect functional groups to selectivity profiles! This color coordination demonstrates the relationship between atomic fragments of compounds and how to relate them to column selectivity profiles.


Define the Characteristics of Your Target Compounds



Match Your Target Compounds to the 5 Primary Selectivity Parameters

Three primary mechanisms of selectivity relate to the analyte of interest above. By referencing the column selectivity bar chart below, we can see the correlation between selectivity mechanisms and column selectivity profiles.

Example: Luna® Omega 1.6 μ m Polar C18 Selectivity Bar Chart



Selectivity Overview

p. 9 Alkyl Phases

UHPLC selectivities that are best suited for the analysis of hydrocarbons

p. 10 Phenyl Phases

UHPLC selectivities that are best suited for aromatic compounds

p. 11 Polar Retentive Phases

UHPLC selectivities that are best suited for more polar compounds





Alkyl Phases

Application Spotlight



- Cannabinoids
- Analgesics
- Pharmaceuticals (USP: L1)

Find the right amount of hydrophobic selectivity for your separation. Below are UHPLC columns that are recommended for the separation of hydrocarbon-based compounds.





Luna® Omega C18

Rugged and highly efficient C18 with strong focus on hydrophobic retention of non-polar and polar compounds.





Kinetex® XB-C18

Di-isobutyl side chains differentiate this C18 column. Low ligand density and an inactive surface make this column a great hydrogen acceptor. This phase will demonstrate improved peak shape for basic compounds and increased retention of acids.





Kinetex C18

Very well balanced column providing some selectivity through steric, hydrogen, and cationic pathways. This is a great starting point for ultra-high efficiency separations.




Kinetex C8

Brings the benefits of core-shell technology to USP L7 methods. The phase will provide moderate hydrophobicity and good steric and hydrogen donating selectivity.



Kinetex EVO C18

Novel pH 1-12 stable C18 that delivers robust methods and improved peak shape for bases.



Denotes stationary phases that are 100% aqueous stable

Material Characteristics

Phase	Particle Sizes (µm)	Pore Size (Å)	Effective Surface Area (m ² /g)	Effective Carbon Load %	pH Stability	Pressure Stability
Luna Omega C18	1.6	100	260	11	1.5 - 8.5*	1,000/600† bar
Kinetex XB-C18	1.7, 2.6, 3.5, 5	100	200	10	1.5-8.5*	
Kinetex C18	1.3, 1.7, 2.6, 5	100	200	12	1.5-8.5*	
Kinetex C8	1.7, 2.6, 5	100	200	8	1.5-8.5*	
Kinetex EVO C18	1.7, 2.6, 5	100	200	11	1.0-12.0	

* pH stability under gradient conditions. pH stability is 1.5 - 10 under isocratic conditions.

† 2.1 mm ID Kinetex columns are pressure stable up to 1000 bar.

When using Kinetex 1.3µm or 1.7µm, increased performance can be achieved, however high pressure-capable instrumentation is required.

Phenyl Phases



Application Spotlight

- Taxanes
- Mycotoxins
- Opiates

Aromatic interactions greatly promote steric interactions. Below are UHPLC columns that have the highest potential for pi-pi bond interaction.



Aromaticity

Column chemistries that contain ring structures interact with aromatic or phenyl containing compounds via pi-pi interactions (π stacking).





Kinetex® Biphenyl

100% aqueous stable reversed phase chemistry with hydrophobic, aromatic, and enhanced polar selectivity.




Kinetex Phenyl-Hexyl

Aromatic and moderate hydrophobic selectivity result in the great retention and separation of aromatic hydrocarbons.



Kinetex F5

This pentafluorophenyl propyl column provides a very high degree of steric selectivity to separate structural isomers. The electronegative fluorine groups offer high selectivity for cationic compounds.



Denotes stationary phases that are 100% aqueous stable

Material Characteristics

Phase	Particle Sizes (μm)	Pore Size (\AA)	Effective Surface Area (m^2/g)	Effective Carbon Load %	pH Stability	Pressure Stability
Kinetex Biphenyl	1.7, 2.6, 5	100	200	11	1.5-8.5*	1,000/600 [†] bar
Kinetex Phenyl-Hexyl	1.7, 2.6, 5	100	200	11	1.5-8.5*	
Kinetex F5	1.7, 2.6, 5	100	200	9	1.5-8.5	

* pH stability under gradient conditions. pH stability is 1.5 - 10 under isocratic conditions.

[†] 2.1 mm ID Kinetex columns are pressure stable up to 1000 bar.

When using Kinetex 1.3 μm or 1.7 μm , increased performance can be achieved, however high pressure-capable instrumentation is required.

Method Development Tip!



Try using methanol for the organic portion of the mobile phase. It can help promote pi-pi bond interaction!

Polar Retentive Phases

Application Spotlight



- Peptide Mapping
- Pesticides
- Nucleosides

Column phases with hydrogen bond accepting/donating capacity have increased polar selectivity and retention of polar compounds. Below are UHPLC columns that offer increased polar selectivity/retention.





Luna® Omega Polar C18

100 % aqueous stability and enhanced selectivity/retention for polar analytes without diminishing useful non-polar retention. The C18 ligand provides general hydrophobic interactions while a polar modified particle surface provides enhanced polar compound retention.




Luna Omega PS C18

Unique, 100 % aqueous stable mixed-mode phase that provides both polar and non-polar retention. The surface contains a positive charged ligand which aids in the retention of acidic compounds through ionic interactions, while the C18 ligand promotes general reversed phase hydrophobic retention. The positively charged surface also improves basic compound peak shape through ionic repulsion.



Kinetex® Polar C18

Combined C18 and polar modified surface that provide polar and non-polar retention alongside 100 % aqueous stability.



Denotes stationary phases that are 100 % aqueous stable

Material Characteristics

Phase	Particle Sizes (µm)	Pore Size (Å)	Surface Area (m ² /g)	Carbon Load (%)	pH Stability	Pressure Stability
Kinetex Polar C18	2.6	100	200	9	1.5 - 8.5*	1,000/600 [†] bar
Luna Omega Polar C18	1.6, 3, 5	100	260	9	1.5 - 8.5*	
Luna Omega PS C18	1.6, 3, 5	100	260	9	1.5 - 8.5*	

* pH stability under gradient conditions. pH stability is 1.5 - 10.0 under isocratic conditions.

[†] 2.1 mm ID Kinetex columns are pressure stable up to 1000 bar.

Applying Column Selectivities to Your UHPLC Analysis

- pp. 13-15 **Thinking Outside the Traditional C18 Phases**
Advances over traditional C18 selectivity
- pp. 16-17 **Improve Polar Separations with 100 % Aqueous Stable Phases**
The selectivity power of using aqueous stable UHPLC columns
- p. 18 **Utilizing Unique Selectivities**
Comparison of unique stationary phases
- p. 19 **Aromatic Based Selectivity**
Aromatic (π - π interaction)
- p. 20 **Protect the Column's Selectivity**
Extend your column's lifetime





Browse thousands of applications at:
www.phenomenex.com/applications

Thinking Outside the Traditional C18 Phase

The extent to which resolution depends on selectivity becomes evident when a chromatographer takes advantage of technologic advances in polar stationary phases in order to improve their analysis. Below is an example of resolution improvements obtained by switching from a traditional C18 to a column with additional selectivity mechanisms.

Polar Selectivity of Catecholamines



Conditions for both columns:

Columns: Luna Omega 1.6 µm Polar C18
Kinetex 1.7 µm C18

Dimensions: 50 x 2.1 mm

Mobile Phase: A: 10 mM Ammonium Formate with 0.1% Formic Acid
B: Acetonitrile with 0.1% Formic Acid

Gradient Time (min)	% B
0	0
3	90

Flow Rate: 0.4 mL/min

Injection Volume: 1 µL


Temperature: 22°C

Detection: MS/MS (SCIEX API 4000™)

Sample: 1. Norepinephrine
2. Epinephrine
3. Normetanephrine
4. Dopamine
5. Metanephrine
6. Serotonin

Improved Separation with Polar C18

The increased hydrogen bond accepting capacity of the Luna Omega Polar C18 column improves polar selectivity for analytes such as dopamine.





Comparative separations may not be representative of all applications.



Thinking Outside the Traditional C18 Phase

Traditional C18 phases may not always be the best option. A UHPLC column that has both polar and hydrophobic versatility allows for great method development flexibility.

Acids/Bases/Neutrals



Conditions for both columns:

Columns: Luna Omega 1.6µm PS C18
ACQUITY BEH 1.7µm C18

Dimensions: 50 x 2.1 mm

Mobile Phase: A: Water with 0.1% Formic Acid
B: Acetonitrile with 0.1% Formic Acid

Gradient:	Time (min)	% B
	0	5
	5	95
	5.01	5
	8	5

Flow Rate: 0.4 mL/min


Temperature: 22 °C

Detection: UV @ 254 nm

- Sample:**
1. Uracil
 2. Pindolol
 3. Chlorpheniramine
 4. Nortriptyline
 5. 3-Methyl-4-nitrobenzoic acid
 6. 5-Methyl salicylaldehyde
 7. Hexanophenone

Improved Separation with PS C18

The increase in steric interaction of the Luna Omega PS C18 column improves resolution for aromatic analytes such as nortriptyline.



Comparative separations may not be representative of all applications.


Thinking Outside the Traditional C18 Phase

Traditional UHPLC C18 phases can be prone to peak tailing for highly basic compounds and this tailing can become exacerbated when higher loadability is required. The Luna® Omega PS C18's combination of a high surface area and novel surface chemistry allows for better peak shape as loading increases.


Why is Loadability Important?

- Analytical conformation parallel to purification
- Dealing with high API concentrations when conducting stability tests
- High sample loading to visualize low-level analytes of interest

Amitriptyline Loading Study



Better Peak Shape Over All Loads




Conditions for both columns:
Columns: Luna Omega 1.6 µm PS C18
 ACQUITY BEH 1.7 µm C18
Dimensions: 50 x 2.1 mm
Mobile Phase: A: Water with 0.1 % Formic Acid
 B: Acetonitrile with 0.1 % Formic Acid
Gradient:

Time (min)	% B
0	5
5	80

Flow Rate: 0.4 mL/min
Temperature: 22 °C
Detection: UV @ 254 nm
Sample: Amitriptyline

Luna Omega PS C18 50 x 2.1 mm vs. Other Columns



Better Peak Shape Over All Sample Loadings

Comparative separations may not be representative of all applications.


Phenomenex | WEB: www.phenomenex.com




Improve Polar Separations with 100% Aqueous Stable Phases

A powerful tool in the chromatographer's toolbox is the ability to use 100% aqueous conditions to promote polar selectivity and increased retention. Traditional C18 phases are known to collapse under 100% aqueous conditions, causing a loss of retention and method development headaches.

Catecholamines



Greater Retention and Resolution Under 100% Aqueous Conditions




Conditions for both separations:

Column: Luna Omega 1.6µm Polar C18
Dimension: 50 x 2.1 mm
Part No.: 00B-4748-AN
Mobile Phase: A: Water with 0.1% Formic Acid
 B: Acetonitrile with 0.1% Formic Acid
Gradient: Time (min) % B
 0 3 (except where noted)
 3 100
Flow Rate: 0.4 mL/min


Temperature: 40 °C
Detection: MS/MS (SCIEX API 4000™) (ambient)
Sample: 1. Norepinephrine
 2. Epinephrine
 3. Normetanephrine
 4. Dopamine
 5. Metanephrine
 6. 3-Hydroxytyramine

Improve Polar Separations with 100% Aqueous Stable Phases

Dipeptides in 100% Aqueous Conditions



100% Aqueous Conditions Increase Retention and Selectivity of Polar Analytes



Column: Luna Omega 1.6 μm Polar C18
Dimensions: 50 x 2.1 mm
Part No.: 00B-4748-AN
Mobile Phase: A: Water with 0.1 % TFA
 B: Acetonitrile with 0.1 % TFA
Gradient:

Time (min)	% B
0	0
3	75



Flow Rate: 0.6 mL/min
Injection Volume: 1 μL
Temperature: 40 °C
Detection: UV @ 210 nm (ambient)
Sample: 1. Arg-Glu
 2. Gly-Tyr
 3. Trp-Gly
 4. Gly-Trp
 5. Pro-Trp

Column: Competitor 1.7 μm C18
Dimensions: 50 x 2.1 mm
Mobile Phase: A: Water with 0.1 % TFA
 B: Acetonitrile with 0.1 % TFA
Gradient:

Time (min)	% B
0	3
3	75

Flow Rate: 0.6 mL/min
Injection Volume: 1 μL
Temperature: 40 °C
Detection: UV @ 210 nm (ambient)
Sample: 1. Arg-Glu
 2. Gly-Tyr
 3. Trp-Gly
 4. Gly-Trp
 5. Pro-Trp

Water Soluble Vitamins under 100% Aqueous Conditions



Conditions same for both separations:

Column: Luna Omega 1.6 μm Polar C18
Dimensions: 50 x 2.1 mm
Part No.: 00B-4748-AN
Mobile Phase: A: 10 mM Ammonium Formate with 0.1 % Formic Acid
 B: Acetonitrile with 0.1 % Formic Acid
Gradient:

Time (min)	% B
0	0
4	90
4.1	0
7	0

Flow Rate: 0.4 mL/min
Temperature: 40 °C
Detection: MS/MS (SCIEX API 4000™) @ 450 °C
Sample: 1. Pyridoxamine
 2. Thiamine
 3. Nicotinic acid
 4. Pyridoxine
 5. Pantothenic acid
 6. Folic acid
 7. Riboflavin

Comparative separations may not be representative of all applications.

Phenomenex | WEB: www.phenomenex.com

Step
4

Applying Selectivity




17


Utilizing Unique Selectivities

The elution order of analytes can change depending upon the predominate selectivity characteristics of a UHPLC column and the utilized mobile phase. Therefore, by altering the stationary phase or mobile phase conditions, we can observe a unique elution order and extend retention of more polar analytes.

Luna® Omega 1.6µm C18



Kinetex® 1.7µm Biphenyl



Conditions for both columns:

Columns: Luna Omega 1.6µm C18
Kinetex 1.7µm Biphenyl

Dimensions: 50 x 2.1 mm

Part No.: 00B-4742-AN

00B-4628-AN

Mobile Phase: A: Water with 0.1 % Formic Acid
B: Acetonitrile with 0.1 % Formic Acid

Gradient:	Time (min)	% B
	0	5
	4	95
	5	95
	5.1	5

Flow Rate: 0.4 mL/min


Temperature: 40°C

Detection: MS/MS (SCIEX API 4000™) (ambient)

Sample: Drugs of Abuse


Aromatic Based Selectivity

Compounds with aromatic ring structures offer a specific type of selectivity associated with pi-pi bond interaction. The compound's aromaticity provides pi electrons that have the potential to interact with pi bonds, which can be found on phenyl-based stationary phases. This provides a unique, orthogonal selectivity compared to a traditional C18 phase.




Increased Separation of Aromatic Compounds Using a Biphenyl Phase

Kinetex 1.7 µm C18



Kinetex 1.7 µm Biphenyl



Conditions for both columns:

Columns: Kinetex® 1.7 µm C18
 Kinetex® 1.7 µm Biphenyl
Dimensions: 50 x 2.1 mm
Part No.: 00B-4628-AN
 00B-4475-AN
Mobile Phase: A: Water with 0.1% Formic Acid
 B: Methanol with 0.1% Formic Acid
Gradient:

Time (min)	% B
0	10
4	100
4.1	10

Flow Rate: 400 µL/min
Injection Volume: 10 µL
Temperature: 50 °C
Backpressure: 450 Bar
Detection: MS/MS (SCIEX API 4000™) (ambient)
Sample: 1. Morphine
 2. Hydromorphone



Protect Your Column's Selectivity



Save Time and Money

It's a fact! Chemical contaminants and particulates are a natural part of any chromatographic analysis. The easiest way to extend column performance is to remove these contaminants and particulates with SecurityGuard ULTRA before they reach your UHPLC column and degrade your chromatography.

With SecurityGuard ULTRA, you will experience:



- Increased UHPLC column lifetime
- Better column performance
- More reproducible chromatography
- Fewer wasted columns

SecurityGuard ULTRA

For all core-shell and/or < 3 μm particle columns (< 20,000 psi / 1,378 bar)





With SecurityGuard ULTRA



(24000 times magnification)

Without SecurityGuard ULTRA



(24000 times magnification)



“ We used to have to change out our columns every 2 to 3 months and ever since we started using the SecurityGuard cartridges we can do at least 6 months before changing a column out. ”

T. Serviss

The opinions stated herein are solely those of the speaker and not necessarily those of any company or organization.



▶ **Chat Now!** ◀
www.phenomenex.com/chat

A Phenomenex Technical Specialist is here to help nearly 24 hours a day!

As Voted by You

Kinetex® Core-Shell LC Columns Earned the SelectScience Gold Seal of Quality Award



Amgen®

“ The Kinetex column has worked great for our validated assays. We easily converted our HPLC methods to UPLC methods using the Kinetex column and have enjoyed being able to run fast UPLC chromatography... ”

GlycosBIO® Food Sciences

“ I really love the Kinetex columns. I now have shorter HPLC method runs and nice peak resolution. Methods that used to take 30 to 40 minutes on other columns now take 15–20 minutes. I’m enjoying the fact that my samples are analyzed more quickly without compromising the quality of the peaks in the chromatograms. ”

University of Texas MD Anderson Cancer Center

“ We have drastically improved sensitivity, reproducibility, and lifetime on our column after switching to the Kinetex technology. ”



KINETEX
Core-Shell Technology

The opinions stated herein are solely those of the speaker and not necessarily those of any company or organization.

Kinetex Core-Shell UHPLC Column Ordering Information

1.7 µm Minibore Columns (mm)					SecurityGuard™
	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	ULTRA Cartridges [‡]
Phases					3/pk
EVO C18	—	00B-4726-AN	00D-4726-AN	00F-4726-AN	AJ0-9298
F5	—	00B-4722-AN	00D-4722-AN	00F-4722-AN	AJ0-9322
Biphenyl	—	00B-4628-AN	00D-4628-AN	00F-4628-AN	AJ0-9209
XB-C18	00A-4498-AN	00B-4498-AN	00D-4498-AN	00F-4498-AN	AJ0-8782
C18	00A-4475-AN	00B-4475-AN	00D-4475-AN	00F-4475-AN	AJ0-8782
C8	00A-4499-AN	00B-4499-AN	00D-4499-AN	00F-4499-AN	AJ0-8784
HILIC	00A-4474-AN	00B-4474-AN	00D-4474-AN	—	AJ0-8786
Phenyl-Hexyl	—	00B-4500-AN	00D-4500-AN	00F-4500-AN	AJ0-8788

for 2.1 mm ID

1.7 µm MidBore™ Columns (mm)				SecurityGuard
	30 x 3.0	50 x 3.0	100 x 3.0	ULTRA Cartridges [‡]
Phases				3/pk
XB-C18	00A-4498-YO	00B-4498-YO	00D-4498-YO	AJ0-8775
C18	—	00B-4475-YO	00D-4475-YO	AJ0-8775
C8	00A-4499-YO	00B-4499-YO	00D-4499-YO	AJ0-8777
HILIC	—	00B-4474-YO	—	AJ0-8779

for 3.0 mm ID



1.7 µm Microbore Columns (mm)			
	50 x 1.0	100 x 1.0	150 x 1.0
Phases			
EVO C18	00B-4726-AO	00D-4726-AO	00F-4726-AO
Biphenyl	00B-4628-AO	00D-4628-AO	00F-4628-AO

1.3 µm Minibore Columns (mm)		
	30 x 2.1	50 x 2.1
Phases		
C18	00A-4515-AN	00B-4515-AN


[‡]SecurityGuard ULTRA cartridges require holder, Part No.: AJ0-9000

Luna® Omega Offers Better UHPLC Column Lifetimes

Phenomenex columns are engineered for durability and are able to withstand high system pressure. For example, lifetime studies are conducted to ensure superior performance and long column lifetimes.

Accelerated Lifetime Study

Less Increase in Pressure Over Time



explore

LUNA®
OMEGA

Conditions for both columns:

Columns: Luna Omega 1.6 µm C18
ACQUITY BEH 1.7 µm C18

Dimensions: 50 x 2.1 mm

Mobile Phase: A: Water with 0.1 % Formic Acid
B: Acetonitrile with 0.1 % Formic Acid

Gradient	Time (min)	% B
	0	5
	4	95
	4.1	5

Flow Rate: 0.4 mL/min

Temperature: 25 °C

Detection: UV @ 210 nm

Sample: Protein Matrix

App ID 23610

Luna Omega UHPLC Column Ordering Information

1.6 µm Microbore Columns (mm)			
Phases	50 x 1.0	100 x 1.0	150 x 1.0
Polar C18	00B-4748-A0	00D-4748-A0	00F-4748-A0
C18	00B-4742-A0	00D-4742-A0	00F-4742-A0
PS C18	00B-4752-A0	00D-4752-A0	-

1.6 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges†
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
Polar C18	00A-4748-AN	00B-4748-AN	00D-4748-AN	00F-4748-AN	AJ0-9505
PS C18	00A-4752-AN	00B-4752-AN	00D-4752-AN	00F-4752-AN	AJ0-9508
C18	00A-4742-AN	00B-4742-AN	00D-4742-AN	00F-4742-AN	AJ0-9502

for 2.1 mm ID



†SecurityGuard ULTRA cartridges require holder, Part No.: AJ0-9000



If Phenomenex analytical columns do not provide at least an equivalent separation as compared to a competing column of the same particle size, similar phase and dimensions, return the Phenomenex column with comparative data within 45 days for a FULL REFUND.

The Chromatographer's Guide To Improving UHPLC Column Selectivity

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auserinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 400-606-8099
f: +86 (0)22 2532-1033
phen@agela.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain

t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

Switzerland

t: +41 61 692 20 20
f: +41 61 692 20 22
swissinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Luna and Kinetex are registered trademarks and MidBore and SecurityGuard are trademarks of Phenomenex. Waters and ACQUITY are registered trademarks of Waters Technologies Corporation. HALO is a registered trademark of Advanced Materials Technology, Inc. Hypersil GOLD is a registered trademark of Thermo Hypersil-Keystone LLC. API 4000 is a trademark of AB SCIEX Ltd. AB SCIEX™ is being used under license. SelectScience is a registered trademark of SelectScience LLC.

Disclaimer

Phenomenex is not affiliated with Advanced Materials Technology, Thermo, or Waters. Comparative separations may not be representative of all applications.

Kinetex EVO is patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

The opinions stated herein are solely those of the speaker and not necessarily those of any company or organization.

© 2018 Phenomenex, Inc. All rights reserved.