

APPLICATIONS

USP Dexamethasone Assay and Organic Impurities by LC-UV using the Kinetex® 1.7 µm C18 Core-Shell UHPLC Column

Zeshan Aqeel and Phil Koerner
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

Dexamethasone is a corticosteroid widely used for its anti-inflammatory and immunosuppressant effects. This application note for the LC-UV assay and organic impurities for Dexamethasone is based on the official USP method. We compare the results obtained utilizing a Kinetex 1.7 µm C18 100 x 2.1 mm core-shell LC column with the original Waters® ACQUITY® BEH™ 1.7 µm C18 column indicated in the USP Chromatographic Columns¹. The standard and system suitability solutions for the assay and organic impurities were prepared in accordance with the USP monograph.

Both columns met the system suitability requirements for the assay: peak symmetry for Dexamethasone NMT 2.0, RSD NMT 0.73 % for Dexamethasone, and resolution NLT 1.5 between Betamethasone and Dexamethasone. Both columns also met the system suitability requirements for the organic impurities: RSD NMT 5.0 % for Dexamethasone and each of the three impurities, and resolution NLT 1.5 between Betamethasone and Dexamethasone.

1. USP Chromatographic Columns
(<https://www.uspchromcolumns.com/chrom/display>).

LC-UV Conditions

Column: Kinetex 1.7 µm C18 ([00D-4475-AN](#))
Waters ACQUITY UPLC BEH 1.7 µm C18

Dimension: 100 x 2.1 mm

Mobile Phase: A = 3.4 g/L of monobasic potassium phosphate in water. Adjust with phosphoric acid to a pH of 3.0
B = Acetonitrile

Gradient: Time (min)	% B
0	24
10	24
15	55
16	90
16.1	24
20	24

Flow Rate: 0.4 mL/min

Temperature: 35 °C


Detector: UV @ 240 nm

Injection: 2 µL


System: Agilent 1290 Infinity (binary)

Sample: As indicated in each chromatogram

Figure 1: System Suitability Solution (same for Assay and Organic Impurities) on Kinetex® 1.7 µm C18 (top) and Waters® ACQUITY® BEH™ 1.7 µm C18 (bottom)




Sample: 1. Betamethasone
2. Dexamethasone



Results – System Suitability Solution (Assay and Organic Impurities)


	Kinetex 1.7 µm C18		ACQUITY BEH 1.7 µm C18	
Analyte	Retention Time, min	Resolution, NLT 1.5	Retention Time, min	Resolution, NLT 1.5
Betamethasone	6.63		8.16	
Dexamethasone	7.06	2.18	8.71	2.40

Figure 2: Standard Solution (Assay) on Kinetex® 1.7 µm C18 (top) and Waters® ACQUITY® BEH™ 1.7 µm C18 (bottom)



Sample: 1. Dexamethasone

App ID 26286




App ID 26283


Results – Standard Solution, Assay

	Kinetex 1.7 µm C18		ACQUITY BEH 1.7 µm C18	
Analyte	Retention Time, min	RSD, NMT 0.73 %	Retention Time, min	RSD, MNT 0.73 %
Dexamethasone	7.04	0.53 %	8.69	0.43 %

Figure 3: Standard Solution (Organic Impurities) on Kinetex® 1.7 µm C18 (top) and Waters® ACQUITY® BEH™ 1.7 µm C18 (bottom)



- Sample:** 1. Betamethasone
2. Dexamethasone
3. Desoximetasone
4. Dexamethasone Acetate



Results – Standard Solution, Organic Impurities

Analyte	Kinetex 1.7 µm C18		ACQUITY BEH 1.7 µm C18	
	Retention Time, min	RSD, NMT 5.0%	Retention Time, min	RSD, MNT 5.0%
Betamethasone	6.62	1.35%	8.15	1.80%
Dexamethasone	7.05	3.98%	8.70	4.15%
Desoximetasone	12.30	1.09%	12.87	0.62%
Dexamethasone Acetate	13.81	1.15%	14.16	0.33%

APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

☎ **All other countries/regions**
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.
www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex is a registered trademark and BE-HAPPY is a trademark of Phenomenex. Waters and ACQUITY are registered trademarks, and BEH is a trademark of Waters Technologies Corporation.

Disclaimers

Comparative separations may not be representative of all applications.

Phenomenex is in no way affiliated with Waters.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.