

APPLICATIONS

Effect of Perfluorinated Alcohols in Impurity Analysis for Synthetic Oligonucleotides

Dr. Ivan Lebedev¹, Brian Rivera¹, and Dr. Helen Whitby²

¹Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

²Phenomenex, UK, Macclesfield

Overview

Therapeutic oligonucleotides represent a recent breakthrough in the pharmaceutical industry. Improved reliability to synthetic routes coupled with new drug delivery technologies has given a heightened level of interest in these clinically interesting targets. However, characterization of oligos, specifically by ion-pair reversed phase liquid chromatography (IP-RPLC), can be quite challenging. Oligos are manufactured by solid phase synthesis, where nucleotides are added in a step-wise phase. As such, n-1 and n+1 must be characterized, and this may require extensive method development to optimize. Further, it may be necessary for characterization and quantitation of other closely related impurities like oxidized phosphorothioates.

Typically, mobile phase composition for the analysis of oligos consists of a perfluorinated alcohol such as HFIP, as well as an alkylamine. To demonstrate the effect of changing the concentration of perfluorinated alcohols used in the impurity analysis of synthetic oligonucleotides we evaluated each, a single- and double-stranded oligonucleotide (BNA and siRNA).

HFIP is used as an acidic modifier for oligonucleotide analysis and facilitates ionization leading to a better signal in the mass spec. Previously it was reported that increasing concentrations of HFIP enhances the MS signal for oligonucleotides¹, however the effect to HFIP has not been fully explored for chromatographic performance.

In **Figure 1**, we can see improved separation of sense/anti-sense of an siRNA when decreasing HFIP concentration from 100 mM to 25 mM, with some loss of resolution from 25 to 12.5 mM. Nonetheless, the decreasing of HFIP does not seem to adversely effect chromatography and further can expand the experimental design in method optimization. A nominal effect is observed with BNA, as seen in **Figure 2**, with separation and impurity profiles being relatively similar.

Although the impurity profile of BNA by LC-UV might not show a marked improvement, the differences become much more apparent when the oligos are analyzed by high-resolution mass spectrometry. Decreasing the HFIP improves ionization efficiency 3x, as shown in differences in peak height in **Figure 3**. Interestingly, TICs do show the later eluting peaks with the 12.5 mM HFIP concentration. Spectra reveal to have similar spectral quality as the main peak (data not shown), but they differ enough in physicochemical properties to show sample heterogeneity that the MS could not determine if separation was not obtained.

In summary, chromatographic performance is often improved by decreasing HFIP concentration, with acceptable or even optimal chromatography being obtained at 12.5 mM. Even more compelling is the improved ionization with the decrease in HFIP yielding significantly higher peak heights, facilitating characterization and quantitation.

LC Conditions

Columns:	bioZen™ 2.6 µm Oligo
Dimension:	100 x 2.1 mm
Part No.:	OOD-4790-AN
Mobile Phase:	A: 4 mM TEA in Water; HFIP as noted B: 4 mM TEA in Methanol; HFIP as noted
Gradient:	5-30 % B in 14 minutes
Flow Rate:	0.3 mL/min
Injection:	200 ng (Figure 1 and 2) 12.5 ng (Figure 3)
Temperature:	65 °C (Figure 1 and 2) 55 °C (Figure 3)
Detection:	UV @ 260 nm (Figure 1 and 2) TOF-MS (Figure 3)
Sample:	As Noted

1. J Am Soc Mass Spectrom. 2017 Jan; 28(1): 190–199

Figure 1. Effect of HFIP Concentration on siRNA




Figure 2. Effect of HFIP Concentration on BNA





Figure 3. Effect of HFIP on Ionization Efficiency



APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/ChatNow.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 (12) 881 0121
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 All other countries/regions Corporate Office USA

t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

bioZen is a trademark of Phenomenex.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.