

Supported Liquid Extraction (SLE)

For Food Applications

Novum™ Synthetic SLE

Unique, Synthetic Sorbent that
Provides Consistent, Reliable Results

Strata® DE Diatomaceous Earth SLE

Cost-Effective, Guaranteed Alternative to
Traditional Diatomaceous Earth SLE Sorbents



 **phenomenex**[®]
...breaking with traditionSM

www.phenomenex.com/FoodSLE



Choose the Best SLE Solution for Your Analysis

Phenomenex provides two options for your SLE analysis, a unique synthetic sorbent and a traditional diatomaceous earth sorbent, each of which provides its own benefits. Select the best option for your work by using the chart below.



Synthetic	Sorbent	Diatomaceous Earth
Lot-to-lot consistency and reproducibility	Advantages	Cost effective and large volume capabilities
Ethyl Acetate, Methyl Tert-Butyl Ether (MTBE)	Extraction Solvents	Dichloromethane (DCM), Hexane, MTBE, Ethyl Acetate
MINI 96-Well Plates, MAX 96-Well Plates	Plate Formats	200 μ L 96-Well Plates, 400 μ L 96-Well Plates
1 cc, 3 cc, 6 cc, 12 cc	Tube Formats	12 cc and 60 cc

Need Help Selecting the Right SLE Option?



Call us

or



Live Chat

www.phenomenex.com/LiveChat

Explore Your SLE Options and Applications

Novum™ SLE

A Unique, Synthetic Sorbentp. 5

Strata® Diatomaceous Earth (DE) SLE

A Cost-Effective Alternative to Traditional SLE Sorbentsp. 5

Applications

Acrylamide from Coffee.....pp. 6–7

Sterols from Olive Oil.....pp. 8–10

Ordering Informationp. 11


Unique SLE Sorbents

SLE produces more reproducible results, increased accuracy, and higher throughput possibilities than Liquid-Liquid Extraction by utilizing a solid support to mimic LLE. With very little method development, both SLE options remove unwanted matrix interferences to provide cleaner samples than LLE.

Two Simple Steps for a Cleaner Extraction:

STEP
01


Load your sample in aqueous solvent.



- ▲ Interferences (i.e. phospholipids, proteins, salts, etc.)
- ◆ Target Analytes

STEP
02

Collect target analytes in water immiscible solvent for analysis.




Novum Synthetic SLE

As a unique, synthetic SLE sorbent, expect Novum to provide reliable, more consistent lot-to-lot cleanup that excels under extraction conditions where traditional diatomaceous earth SLE falls short, such as ethyl acetate.

Recommended Uses:

- Ethyl acetate extractions
- Smaller volume samples
- Low variation lot-to-lot




Strata DE Diatomaceous Earth SLE

Packed with diatomaceous earth, Strata DE is a cost-effective alternative to traditional SLE products such as Biotage® ISOLUTE® SLE+, Thermo Fisher® Hypersep™ SLE, and Agilent® Chem Elut™ SLE that won't require you to sacrifice your results.

Recommended Uses:

- Direct alternative to traditional diatomaceous earth SLE sorbents
- Non-oxygenated extraction solvents such as DCM and Hexane
- Large volume samples, up to 20 mL capacity with 60 cc tubes!



Acrylamide from Coffee


Pre-Treatment

Prepare ground coffee as directed on the coffee product using 1,500 mL of boiling water and 60 g of ground coffee. Prepare instant coffee using 250 μ L of boiling water and 2 g of instant coffee. Allow both samples to reach room temperature before adding 20 μ L Acrylamide-13C3 (4 μ g/mL in water) to 800 μ L of the prepared coffee. (100 ng/mL for ground coffee and 200 ng/mL for instant coffee). Then add 150 μ L 2% Ammonium hydroxide in water to the spiked samples and vortex for 30 seconds.

SLE Protocol

Cartridge:	Novum™ SLE 6 cc Tube
Part No.:	8B-S138-JCH
Load:	Pre-treated sample and pulse vacuum (~5-10" Hg) for 5-10 seconds or until sample has completely entered the sorbent
Wait:	5 minutes
Elute:	2x 2.5 mL Ethyl acetate/Tetrahydrofuran (1:1) and allow to elute by gravity into a tube containing 10 μ L Ethylene glycol*
Apply:	Vacuum or apply positive pressure at 5" Hg for 20-30 seconds to complete the extraction
Dry:	Sample under slow stream of Nitrogen at 45 °C
Reconstitute:	300 μ L Water

Achieve Ultra-Clean Samples with Novum SLE



*Ethylene glycol was added to prevent the sample from drying completely during the dry down step.

Recovery of Acrylamide


Sample ID	Ground coffee (100 ng/mL)	Instant coffee (200 ng/mL)
Mean of area ratio	1.89	3.75
STDV	0.01	0.06
CV (%)	0.78	1.61
Absolute Recovery (%)	94.9	92.8
n=	6	6

Excellent Recoveries!



Acrylamide from Coffee

Acrylamide from Coffee (100ng/mL)




LC-MS/MS Conditions

Column: Synergi™ 4 μm Hydro-RP
Dimensions: 50 x 2.0 mm
Part No.: 00B-4375-B0
Guard Cartridge: SecurityGuard™ Hydro-RP
Part No.: AJ0-7510
Mobile Phase: A: 0.1 % Formic acid in Water
B: 0.1 % Formic acid in Methanol
Gradient:

Time (min)	B (%)
0	0
0.6	0
0.85	100
3	100
3.01	0
5	0

Flow Rate: 0.3 mL/min
Injection Volume: 5 μL
Temperature: Ambient
Detection: MS/MS (SCIEX API 5000™), ESI+

App ID 22851



Sterols from Olive Oil

Pre-treatment

Add 40 μL of 1 mg/mL cholestanol in chloroform to a clean, dry 20 mL screw-top test tube and evaporate to dryness under a nitrogen flow.

Saponification

1. Add 200 mg of olive oil sample to the test tube containing the internal standard.
2. Add 1.5 mL of 2M Potassium hydroxide in 95 % Ethanol.
3. Cap the tube and heat in an 80 °C oven for 25 minutes.
4. Mix sample gently to ensure homogeneity (sample should appear as a clear solution) and continue heating for an additional 25 minutes.
5. After heating, add 13.5 mL of deionized water and mix. The entire diluted volume is now ready to load onto the SLE cartridge.

SLE Protocol

Cartridge:	Strata® DE SLE cartridge, 20mL loading capacity, 60 cc Tube
Part No.:	8B-S325-VFF
Load:	Diluted sample (from saponification step 5) plus 2x 1 mL DI water rinse (17 mL total volume, gravity flow)
Wait:	15 minutes
Extract:	3x 15 mL Diethyl ether (gravity flow)
Evaporate:	Dry under N_2 at 40 °C (greenish-yellow, oily residue)
Reconstitute:	5 mL of Hexane

SPE (Solid Phase Extraction) Protocol and Derivatization


Cartridge:	Strata Si-1 (1 g/6 mL) tube
Part No.:	8B-S012-JCH
Condition:	1. 2x 6 mL Hexane 2. 1 mL 0.2M Potassium hydroxide in 95 % Ethanol
Equilibrate:	5 mL Hexane immediately after potassium hydroxide elution
Load:	Reconstituted SLE extract (5 mL) followed by 2x 1 mL Hexane rinses
Wash:	85 mL Hexane/Diethyl ether (98:2) under 3" Hg vacuum, flow rate of 2 mL/min.*
Elute:	10 mL Hexane/Diethyl ether (60:40)
Dry:	Dry under N_2 at 50 °C. After evaporating to dryness, add 3-4 drops of acetone and then re-evaporate under N_2 to remove any occluded water. Place in 100 °C oven for 10 minutes.
Derivatization:	250 μL Pyridine/BSTFA (3:1) at 80 °C for 30 minutes

*To handle the large volume of eluant, a 60 mL empty reservoir tube was attached to the 6 mL SPE tube.




Sterols from Olive Oil

Extra Virgin Olive Oil Sterols



Adulterated Olive Oil:Extra Virgin Olive Oil/Canola Oil (50:50)



Sterols from Olive Oil

IOC Sterol Criteria for Virgin Olive Oil Classification

Standard Name	IOC Standard Criteria for Virgin Olive Oil	Extra Virgin Olive Oil		Adulterated Olive Oil	
		% Recovery	% RSD (n=3)	% Recovery	% RSD (n=2)
Apparent β-Sitosterol	$\geq 93.0\%$ of total sterols	94.6%	0.3	60.3%	1.2
Cholesterol	$\leq 0.5\%$ of total sterols	not detected	-	0.3%	13.3
Brassicasterol	$\leq 0.1\%$ of total sterols	not detected	-	8.7%	0.7
Campesterol	$\leq 4.0\%$ of total sterols	3.8%	6.8	29.1%	1.3
Stigmasterol	\leq Campesterol ($\leq 4.0\%$ of total sterols)	1.0%	9.0	0.6%	45.0
Δ-7-Stigmastenol	$\leq 0.5\%$ of total sterols	not detected	-	0.7%	2.9
Uvaol + Erythrodiol	$\leq 4.5\%$ of total sterols	1.8%	31	0.3%	1.0
Total Sterols	≥ 1000 mg/kg	1324 mg/kg	6	4221 mg/kg	1.0

Apparent β -sitosterol = β -sitosterol + Δ -5-avenasterol + Δ -5,23-stigmastadienol + clerosterol + sitostanol + Δ -5,24-stigmastadienol. Total sterols = cholesterol + 24-methylene cholesterol + brassicasterol + campesterol + campestanol + stigmasterol + Δ -7-campesterol + Δ -5,23-stigmastadienol + apparent β -sitosterol + Δ -7-avenasterol.

To read full technical note, go to www.phenomenex.com and search TN-0114



Start Your New SLE Method

Synthetic SLE

Novum™ SLE 96-Well Plates

Novum Simplified Liquid Extraction SLE Well Plates		
Part No.	Description	Unit
8E-S138-FGA	Novum SLE MINI 96-Well Plate	1/pk
8E-S138-5GA	Novum SLE MAX 96-Well Plate	1/pk

Novum Tubes

Novum Simplified Liquid Extraction SLE Tubes		
Part No.	Description	Unit
8B-S138-FAK	Novum SLE 1 cc Tubes	100/pk
8B-S138-5BJ	Novum SLE 3 cc Tubes	50/pk
8B-S138-JCH	Novum SLE 6 cc Tubes	30/pk
8B-S138-KDG	Novum SLE 12 cc Tubes	20/pk



Diatomaceous Earth SLE

Strata® DE 96-Well Plates

Strata DE Diatomaceous Earth SLE Well Plates		
Part No.	Description	Unit
8E-S325-FGB	Strata DE SLE 200 µL 96-Well Plate	2/pk
8E-S325-5GB	Strata DE SLE 400 µL 96-Well Plate	2/pk



Strata DE Tubes

Strata DE Diatomaceous Earth SLE Tubes		
Part No.	Description	Unit
8B-S325-KDG	Strata DE SLE 12 cc Tubes	20/pk
8B-S325-VFF	Strata DE SLE 60 cc Tubes	16/pk



If Phenomenex SLE products do not perform as well or better than your current SLE product, return the product with comparative data within 45 days for a FULL REFUND.



 **presston™ 100**

Presston 100 Positive Pressure Manifold

Presston 100		
Part No.	Description	Unit
AHO-9342	Presston 100 Positive Pressure Manifold, 1 mL Tube Complete Assembly	ea
AHO-9347	Presston 100 Positive Pressure Manifold, 3 mL Tube Complete Assembly	ea
AHO-9343	Presston 100 Positive Pressure Manifold, 6 mL Tube Complete Assembly	ea

The Presston 100 96-Well Positive Pressure Manifold can also process 1, 3, and 6 mL tubes using the following adapter kits:

Presston 100 Tube Adapter Kits (for AHO-9334)		
Part No.	Description	Unit
AHO-9344	1 mL Tube Adapter Kit	ea
AHO-9345	3 mL Tube Adapter Kit	ea
AHO-9346	6 mL Tube Adapter Kit	ea



Vacuum Manifolds

Vacuum Manifolds		
Part No.	Description	Unit
AHO-6023	12-Position Tube Vacuum Manifold Set	ea
AHO-6024	24-Position Tube Vacuum Manifold Set	ea



Phenomenex warrants that for a period of 12 months following delivery, the Presston 100 Positive Pressure Manifold you have purchased will perform in accordance with the published specifications and will be free from defects in materials or workmanship. In the event that the Presston 100 Positive Pressure Manifold does not meet this warranty, Phenomenex will repair or replace defective parts. Please visit www.phenomenex.com/Presston for complete warranty information.

Supported Liquid Extraction (SLE) For Food Applications

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auserinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 400-606-8099
f: +86 (0)22 2532-1033
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Portugal

t: +351 221 450 488
f: 34 91-413-2290
ptinfo@phenomenex.com

Spain

t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

Switzerland

t: +41 61 692 20 20
f: +41 61 692 20 22
swissinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com.

Trademarks

Strata is a registered trademark and Novum, Synergi, Presston, Zebtron, 5PLUS, Z-Liner, and SecurityGuard are trademarks of Phenomenex. Biotage and ISOLUTE are registered trademarks of Biotage AB Corp. Thermo Fisher and HyperSep are trademarks of Thermo Fisher Scientific. Agilent is a registered trademark and ChemElut is a trademark of Agilent Technologies, Inc. API 5000 is a trademark of AB SCIEX Pte. Ltd. AB SCIEX is being used under license.

Disclaimer

Phenomenex is not affiliated with Biotage AB Corp., Thermo Fisher Scientific, or Agilent Technologies, Inc.

Novum is patent pending.

© 2018 Phenomenex, Inc. All rights reserved.