



KINETEX
Core-Shell Technology



Tabletop HPLC/UHPLC Reference Tool Inside



Utilize your *Origami Skills* to maximize your chromatographic power!

Exclusive Offers
Just For You!



phenomenex
...breaking with traditionSM



www.phenomenex.com/Kinetex




KINETEX[®]
Core-Shell Technology



Year Anniversary

Guess what? I just turned 10 years old!

Over the last decade, I've been written into tens of thousands of methods, appeared in over 4,500 journals, and received several awards thanks to your confidence in me! I've also been cloned over 3.2 tresvigintillion times, packed into hundreds of thousands of columns, and sent all over the world. It's been a pretty cool ride so far.



But my proudest moment is being brought into YOUR lab to help you do all the great things you do.

I hope you'll take on this little "origami" project I've sent to you and construct your new tabletop core-shell reference tool. I want to ensure you get the absolute best performance out of me day after day—for 10 more years.

Yours in the Lab,

Kinetex Ball-Guy



THE RESOLUTION EQUATION

As You Know, It's All About Resolution



By modifying any one or all three of these important chromatographic factors (retention, selectivity, and efficiency), you can change or improve your resolution.

Utilize this tabletop tool to help make the modifications you need to supercharge your separations.

Questions?
www.phenomenex.com/Chat




The Resolution Equation:



COLUMN RETENTION

$$\left(\frac{\text{Retention } k}{k+1} \right)$$


Estimate of Column Void Volume



Core-Shell Column

Column Internal Diameter (mm)	Length (mm)	Void Volume (mL)
1.0	30	0.01
1.0	50	0.02
1.0	100	0.04
1.0	150	0.06
1.0	250	0.11
2.1	30	0.06
2.1	50	0.09
2.1	100	0.19
2.1	150	0.28
2.1	250	0.47

Column Internal Diameter (mm)	Length (mm)	Void Volume (mL)
3.0	30	0.11
3.0	50	0.19
3.0	100	0.38
3.0	150	0.57
3.0	250	0.95
4.6	30	0.27
4.6	50	0.45
4.6	100	0.90
4.6	150	1.35
4.6	250	2.24



Fully Porous Column

Column Internal Diameter (mm)	Length (mm)	Void Volume (mL)
1.0	30	0.02
1.0	50	0.03
1.0	100	0.05
1.0	150	0.08
1.0	250	0.13
2.1	30	0.07
2.1	50	0.12
2.1	100	0.24
2.1	150	0.35
2.1	250	0.59

Column Internal Diameter (mm)	Length (mm)	Void Volume (mL)
3.0	30	0.14
3.0	50	0.24
3.0	100	0.48
3.0	150	0.72
3.0	250	1.20
4.6	30	0.34
4.6	50	0.57
4.6	100	1.13
4.6	150	1.70
4.6	250	2.83

Estimates of core-shell and fully porous column void volumes are based on a particle volume exclusion factor of 68 % for fully porous and 54 % for core-shell particles.

To determine the empirical column void volume of a particular column installed on a particular system, we suggest an evaluation by injecting a non-retained compound as a void marker and utilizing the following formula to determine void volume:

Void Volume = Retention Time of Void Peak (Vo) * Flow Rate.

Material Characteristics


Packing Material Kinetex Phases	Particle Sizes (μm)	Pore Size (\AA)	Effective Surface Area (m^2/g)	Effective Carbon Load (%)	pH Range	Pressure Stability (bar)	Shipping Solvent
Kinetex Polar C18	2.6	100	200	9	1.5 – 8.5*	1,000/600**	Acetonitrile/Water (50:50)
Kinetex PS C18	2.6	100	200	9	1.5 – 8.5*	1,000/600**	Acetonitrile/Water (50:50)
Kinetex C18	1.3, 1.7, 2.6, 5	100	200	12	1.5 – 8.5*	1,000/600**	Acetonitrile/Water (50:50)
Kinetex EVO C18	1.7, 2.6, 5	100	200	11	1.5 – 12	1,000/600**	Acetonitrile/Water (45:55)
Kinetex XB-C18	1.7, 2.6, 3.5, 5	100	200	10	1.5 – 8.5*	1,000/600**	Acetonitrile/Water (50:50)
Kinetex C8	1.7, 2.6, 5	100	200	8	1.5 – 8.5*	1,000/600**	Acetonitrile/Water (45:55)
Kinetex Biphenyl	1.7, 2.6, 5	100	200	11	1.5 – 8.5*	1,000/600**	Acetonitrile/Water w/ 0.1 % Formic Acid (50:50)
Kinetex Phenyl-Hexyl	1.7, 2.6, 5	100	200	11	1.5 – 8.5*	1,000/600**	Acetonitrile/Water (45:55)
Kinetex F5	1.7, 2.6, 5	100	200	9	1.5 – 8.5*	1,000/600**	Acetonitrile/Water (40:60)
Kinetex HILIC	1.7, 2.6, 5	100	200	0	2.0 - 7.5*	1,000/600**	Acetonitrile/ 100 mM Ammonium Formate (93:7)
Kinetex PAH	3.5	—	—	12	1.5 – 8.5*	1,000/600**	Acetonitrile/Water (65:35)

* pH stability under gradient conditions. pH stability is 1.5 - 10 under isocratic conditions.

** 2.1 mm ID Kinetex columns are pressure stable up to 1,000 bar.

When using Kinetex 1.3 μm or 1.7 μm , increased performance can be achieved, however high pressure-capable instrumentation is required.

Visit www.phenomenex.com/Chat
to connect with our Technical Team
nearly 24/7!



MOBILE PHASE SELECTIVITY

$$\left(\frac{\alpha - 1}{\alpha} \right)$$

Table of Common HPLC/UHPLC Buffers and Additives

Name	Chemical Formula	Molecular Weight	pKa	UV (nm) Cutoff @ 1 AU	Buffer Range
Hydrochloric Acid	HCl	36.46	-7.00		
Sodium Hydroxide	NaOH	40.00	13.80		
Potassium Hydroxide	KOH	56.11	15.70		
Trifluoroacetic Acid**	C ₂ HF ₃ O ₂	114.03	0.23	210 nm (0.1 %)	< 1.5
Acetic Acid*	C ₂ H ₄ O ₂	60.05	4.76	210 nm (10 mM)	3.76 ~ 5.76
Sodium Acetate	C ₂ H ₃ NaO ₂	82.03	4.54		3.76 ~ 5.76
Potassium Acetate	C ₂ H ₃ KO ₂	98.13	4.54	210 nm (10 mM)	3.76 ~ 5.76
Phosphoric Acid	H ₃ O ₄ P	98.00	2.15	< 200 nm (0.1 %)	1.15 ~ 3.15
			7.20		6.20 ~ 8.20
			12.32		11.15 ~ 13.15
Monosodium Dihydrogen Phosphate	H ₂ NaO ₄ P	119.98	2.15		1.15 ~ 3.15
			7.20		6.20 ~ 8.20
			12.32		11.15 ~ 13.15
Monosodium Dihydrogen Phosphate Monohydrate	H ₄ NaO ₅ P	137.99	2.15		1.15 ~ 3.15
			7.20		6.20 ~ 8.20
			12.32		11.15 ~ 13.15
Monosodium Dihydrogen Phosphate Dehydrate	H ₃ NaO ₄ P	156.01	2.15		1.15 ~ 3.15
			7.20		6.20 ~ 8.20
			12.32		11.15 ~ 13.15
Disodium Monohydrogen Phosphate	HNa ₂ O ₄ P	141.96	2.15		1.15 ~ 3.15
			7.20		6.20 ~ 8.20
			12.32		11.15 ~ 13.15
Monopotassium Dihydrogen Phosphate	H ₂ KO ₄ P	136.08	2.15		1.15 ~ 3.15
			7.20		6.20 ~ 8.20
			12.32		11.15 ~ 13.15
Dipotassium Monohydrogen Phosphate	HK ₂ O ₄ P	174.18	2.15		1.15 ~ 3.15
			7.20		6.20 ~ 8.20
			12.32		11.15 ~ 13.15
Citric Acid	C ₆ H ₈ O ₇	192.12	3.13	230 nm (10 mM)	2.13 ~ 4.13
			4.76		3.76 ~ 5.76
			6.40		5.40 ~ 7.40
Formic Acid*	CH ₂ O ₂	46.02	3.74	210 nm (10 mM)	2.74 ~ 4.74
Sodium Formate	CHNaO ₂	68.01	3.74		2.74 ~ 4.74
Potassium Formate	CHKO ₂	84.11	3.74	210 nm (10 mM)	2.74 ~ 4.74
Ammonium Hydroxide	H ₅ NO	35.05	9.24	200 nm (10 mM)	8.24 ~ 10.24
Ammonium Chloride	ClH ₄ N	53.49	9.24		8.24 ~ 10.24
Ammonium Acetate*	C ₂ H ₇ NO ₂	77.08	4.76	205 nm (10 mM)	3.76 ~ 5.76
			9.24		8.24 ~ 10.24
Ammonium Formate*	CH ₃ NO ₂	63.06	3.74		2.74 ~ 4.74
			9.24		8.24 ~ 10.24
Tris	C ₄ H ₁₁ NO ₃	121.14	8.08		7.08 ~ 9.08
Tris Hydrochloride	C ₄ H ₁₂ ClNO ₃	157.60	8.08		7.08 ~ 9.08
Triethylamine*	C ₆ H ₁₅ N	101.19	10.72	< 200 nm (10 mM)	9.72 ~ 11.72
Triethylamine Hydrochloride	C ₆ H ₁₆ ClN	137.65	10.72		9.72 ~ 11.72
Pyrrolidine	C ₄ H ₉ N	71.12	11.30		10.30 ~ 12.30

* Volatile buffer, compatible with Mass Spectroscopy analysis.

** Compatible with Mass Spectroscopy but at low concentrations due to ion suppression effects.

Solvent Table

Solvent	Polarity Index	Refractive Index @ 20 °C	UV (nm) Cutoff @ 1 AU	Solubility in Water (% w/w)	Solvent	Polarity Index	Refractive Index @ 20 °C	UV (nm) Cutoff @ 1 AU	Solubility in Water (% w/w)
Acetic Acid	6.2	1.372	230	100	Di-Ethyl Ether	2.8	1.353	220	6.89
Acetone	5.1	1.359	330	100	Heptane	0.0	1.387	200	0.0003
Acetonitrile	5.8	1.344	190	100	Hexane	0.0	1.375	200	0.001
Benzene	2.7	1.501	280	0.18	Methanol	5.1	1.329	205	100
n-Butanol	4.0	1.394	254	0.43	Methyl-t-Butyl Ether ⁴	2.5	1.369	210	4.8
Butyl Acetate	3.9	1.399	215	7.81	Methyl Ethyl Ketone ⁵	4.7	1.379	329	24
Carbon Tetrachloride	1.6	1.466	263	0.08	Pentane	0.0	1.358	200	0.004
Chloroform	4.1	1.446	245	0.815	n-Propanol	4.0	1.384	210	100
Cyclohexane	0.2	1.426	200	0.01	Iso-Propanol ⁶	3.9	1.377	210	100
1,2-Dichloroethane ¹	3.5	1.444	225	0.81	Di-Iso-Propyl Ether	2.2	1.368	220	
Dichloromethane ²	3.1	1.424	235	1.6	Tetrahydrofuran	4.0	1.407	215	100
Dimethylformamide	6.4	1.431	268	100	Toluene	2.4	1.496	285	0.051
Dimethyl Sulfoxide ³	7.2	1.478	268	100	Trichloroethylene	1.0	1.477	273	0.11
Dioxane	4.8	1.422	215	100	Water	9.0	1.333	200	100
Ethanol	5.2	1.360	210	100	Xylene	2.5	1.500	290	0.018
Ethyl Acetate	4.4	1.372	260	8.7					

Synonym Footer

1. Ethylene Chloride
2. Methylene Chloride

3. Methyl Sulfoxide
4. tert-Butyl Methyl Ether

5. 2-Butanone
6. 2-Propanol

COLUMN EFFICIENCY

$$\left(\frac{\text{Efficiency}}{\sqrt{N}} \right) \frac{1}{4}$$

Typical HPLC/UHPLC Backpressure by Particle Size

Particle Size (µm)	Internal Diameter (mm)	Typical Flow (mL/min)	Typical Pressure (PSI)		
			50 mm**	150 mm**	250 mm**
	2.1	0.5	10,500	—	—
	1.0	0.1	2,600	—	—
	2.1	0.5	5,300	7,300	—
	3.0	0.3	7,000	—	—
	1.0	0.1	1,700	4,300	—
	2.1	0.5	3,000	6,800	—
	3.0	0.8	2,300	5,900	—
	4.6	1.8	2,800	5,300	6,900
	4.6	1.5	—	2,600	—
	2.1	0.3	580	1,200	—
	3.0	0.5	620	1,050	—
	4.6	1.3	1,100	1,600	1,700
	10	5	—	2,100	3,000
	21.2	25	630	1,300	2,300
	30	50	670	1,600	2,300
	50	80	—	1,000	—

** Column length in mm.

Typical Flow Rates by Column Internal Diameters

Column	Column Internal Diameter	Typical Flow Rates	Sample Load (Typical Estimate)	Application
Nano	50 - 75 µm	~200 - 500 nL/min	100 - 200 ng	Higher Sensitivity, Very Limited Sample, LC-MS/MS
Capillary	0.3 - 0.5 mm	0.001 - 0.02 mL/min	1 - 10 µg	High Sensitivity, Limited Sample, LC-MS/MS
Microbore	1.0 mm	0.02 - 0.1 mL/min	10 - 50 µg	Increased Sensitivity, LC-MS/MS
Analytical	2.0 - 4.6 mm	0.1 - 3.0 mL/min	0.1 - 1.5 mg	General Analytical LC Applications
Semi-Prep	7.8 w- 10.0 mm	5 - 10 mL/min	1 - 10 mg	Analytical and Limited Purification Applications
Preparative	21.2 - 100.0 mm	20 - 250 mL/min	20 - 250 mg	Higher Purification Yield

Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex is a registered trademark of Phenomenex.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2019 Phenomenex, Inc. All rights reserved.