

LC-MS/MS Separation of Doxepin and its Metabolite N-Desmethyldoxepin Using the Kinetex® 2.6 µm Biphenyl and Luna® Omega 1.6 µm C18 Columns

Roxana Eggleston-Rangel and Dr. Bryan Tackett
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

Doxepin hydrochloride belongs to a class of drugs referred to as tricyclic antidepressants (TCA) which are widely used to treat a broad range of psychotic depressive disorders. N-Desmethyldoxepin is the primary metabolite of Doxepin (**Figure 1**). This drives the interest in developing a reversed phase HPLC method able to separate and analyze Doxepin and its metabolites successfully and in a short amount of time. In this application note, we present an LC-MS/MS method to separate Doxepin and its metabolite by utilizing Kinetex 2.6 µm Biphenyl and Luna Omega 1.6 µm C18 columns.

The Kinetex Biphenyl column provided separation for Doxepin and its primary metabolite through a mixture of pi-pi and polar interactions; and the higher efficiency provided by the core-shell particle morphology also resulted in narrow peaks and increased MS sensitivity (**Figure 2a**).

The Luna Omega C18 column is a conventional C18 column with TMS endcapping, making this column a good choice for retention of analytes through hydrophobic interactions. Depending on the mobile phase used, there was a shift in elution order observed. Using Acetonitrile instead of Methanol in the mobile phase on the Luna Omega 1.6 µm C18 column gave narrower peaks and the same elution order observed with the Kinetex Biphenyl column (**Figure 2b, 2c**). There was a slightly shorter run time when using Methanol in the mobile phase on the Luna Omega 1.6 µm C18 column, which is attributed to the weaker elution strength of methanol vs. acetonitrile. Regardless of mobile phase, the run time was much shorter on the Luna Omega 1.6 µm C18 column compared to the Kinetex 2.6 µm Biphenyl column, which is likely due to the multiple interactions between the analytes and the Biphenyl stationary phase.

LC-MS/MS Conditions

Column: Kinetex 2.6 µm Biphenyl ([00B-4622-AN](#))
Luna Omega 1.6 µm C18 ([00B-4742-AN](#))

Dimension: 50 x 2.1 mm

Mobile Phase: A: 0.1 % Formic Acid in Water
B: See Chromatogram

Gradient: Time (min)	%B
0	40
6	70
6.1	40
8	40

Flow Rate: 500 µL/min

Injection Volume: 1 µL

Temperature: 40 °C

LC System: Agilent® 1200 Series


Detection: MS/MS

Detector: SCIEX® 4500

MRM Transitions

Analyte	Q1 (m/z)	Q3 (m/z)
Doxepin hydrochloride-1	280	107
Doxepin hydrochloride-2	280	115
N-Desmethyldoxepin-1	266.1	107
N-Desmethyldoxepin-2	266.1	235.1

Figure 1. Doxepin and N-Desmethyldoxepin



Interestingly, when using methanol both the Luna Omega C18 and Kinetex Biphenyl columns provided separation of the (E) and (Z) isomers of Doxepin, which are present at an approximate 85:15 ratio, with the Biphenyl providing complete separation. This was confirmed as the mass transitions for the (E) and (Z) isomers are the same. This separation could be due to the increased retention observed when using methanol, which is a weaker elution solvent than acetonitrile, with the increased interaction allowing for separation of the isomers.

Figure 2a. Kinetex® 2.6 µm Biphenyl




Figure 2b. Luna Omega 1.6 µm C18





Figure 2c. Luna® Omega 1.6 µm C18



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions
Corporate Office USA**

t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Luna and Kinetex are registered trademarks and BE-HAPPY is a trademark of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc. SCIEX and QTRAP are registered trademarks of AB SCIEX Pte. Ltd.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Agilent Technologies, Inc.
FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.
© 2021 Phenomenex, Inc. All rights reserved.

