

APPLICATIONS

Expanded Mycotoxins Analysis in Cannabis Matrices by LC-MS/MS

W. Maguire¹, C. Cerbu¹, K. Jambor¹, C. Reichmann¹, C. Call¹ and Scott Krepich²

¹ Columbia Food Labs, 12423 NE Whitaker Way, Portland, OR 97230 USA

² Phenomenex, Inc., 411 Madrid Avenue, Torrance, CA 90501, USA

Background

Mycotoxins are secondary fungal metabolites that can be harmful for human consumption, and the primary five (Ochratoxin A, Aflatoxin B1, B2, G1, and G2) are required testing in several states. Currently, Oregon does not mandate mycotoxin testing, however, given water activity requirements in Oregon along with the heat-stable nature of these compounds, they might persist in cured material. Processing of contaminated cannabis flower or trim with even trace levels of residues may concentrate these down to unsafe levels. Columbia Food Labs has taken the initiative to test for 15 common fungal toxins, above and beyond any current regulatory requirements, to ensure safety and scientific validity within the industry.

Introduction

An expanded list of 15 mycotoxins LC-MS/MS analysis is shown here with a simple solvent and modified QuEChERS salt extraction. The entire workflow is simple and accessible to most cannabis and food residue testing laboratories.

Mass Spec Parameters

Mass Spec Detector: Sciex[®] Triple Quad™ 5500

Samples were ionized using electrospray with positive/negative ion-mode polarity switching

MRM Transitions

| Compound Name | Q1 Mass | Q3 Mass | Retention Time (min) | Declustering Potential | Entrance Potential | Collision Energy | Collision Cell Exit Potential | Quantifier |
|-----------------------|---------|---------|----------------------|------------------------|--------------------|------------------|-------------------------------|------------|
| Patulin | 152.9 | 109 | 2.37 | -90 | -10 | -14 | -13 | |
| Patulin | 152.9 | 81 | 2.37 | -90 | -10 | -18 | -11 | X |
| Patulin | 152.9 | 53.1 | 2.37 | -90 | -10 | -22 | -9 | |
| Nivalenol | 311.02 | 281 | 2.51 | -105 | -10 | -14 | -17 | |
| Nivalenol (M-H+HCOOH) | 357.007 | 281.1 | 2.51 | -55 | -10 | -18 | -15 | X |
| Nivalenol (M-H+HCOOH) | 357.007 | 191 | 2.51 | -55 | -10 | -40 | -16 | |
| Deoxynivalenol (DON) | 296.9 | 248.9 | 3.08 | 71 | 10 | 17 | 18 | X |
| Deoxynivalenol (DON) | 296.9 | 203.1 | 3.08 | 71 | 10 | 21 | 14 | |
| Deoxynivalenol (DON) | 296.9 | 175 | 3.08 | 71 | 10 | 27 | 12 | |
| Deoxynivalenol (DON) | 296.9 | 91 | 3.08 | 71 | 10 | 65 | 4 | X |
| Aflatoxin G2 | 331.014 | 216.85 | 4.47 | 106 | 10 | 49 | 14 | |
| Aflatoxin G2 | 331.014 | 188.833 | 4.47 | 106 | 10 | 57 | 12 | X |
| Aflatoxin G1 | 329.059 | 243.35 | 4.71 | 81 | 10 | 39 | 16 | |
| Aflatoxin G1 | 329.059 | 311.2 | 4.71 | 81 | 10 | 31 | 10 | X |
| Aflatoxin B2 | 315.1 | 287.1 | 5.09 | 91 | 10 | 37 | 20 | |
| Aflatoxin B2 | 315.1 | 259.1 | 5.09 | 91 | 10 | 41 | 16 | X |
| Aflatoxin B1 | 313.039 | 285 | 5.41 | 86 | 10 | 31 | 20 | |

(Cont'd on next page)

Experimental Conditions

Sample Preparation:

Cannabis Flower Samples were ground in a Retch GM 200 knife mill, and 0.5 g sample soaked in 5 mL of 2 % ascorbic acid in water in a 50 mL falcon tube. 10 mL acetonitrile were added, followed by a modified EN 15662 QuEChERS salt extraction (4 g MgSO₄, 1 g NaCl, 1.5 g sodium citrate), shaken for 5 minutes, then centrifuged at 2500 RPM. The supernatant was diluted 5 x with aqueous ammonium formate buffer and filtered through a 0.45 µm syringe filter prior to injection to HPLC.

LC Method Parameters

Column: Luna[®] Omega 3 µm Polar C18
Dimensions: 100 x 2.1 mm
Part No.: [00D-4760-AN](#)
Flow Rate: 0.4 mL/min
Mobile Phase: A: 1 mM Ammonium formate + 0.1 % Formic acid in Water
 B: Methanol
Gradient:

| Time (min) | % B |
|------------|-----|
| 0 | 5 |
| 3 | 35 |
| 10 | 90 |
| 12 | 90 |


Temperature: 40 °C
Injection Volume: 3 µL
Detection: MS/MS – Sciex Triple Quad 5500

APPLICATIONS

MRM Transitions (Cont'd.)

| Compound Name | Q1 Mass | Q3 Mass | Retention Time (min) | Declustering Potential | Entrance Potential | Collision Energy | Collision Cell Exit Potential | Quantifier |
|-------------------|---------|---------|----------------------|------------------------|--------------------|------------------|-------------------------------|------------|
| Aflatoxin B1 | 313.039 | 241 | 5.41 | 86 | 10 | 49 | 18 | X |
| HT-2 Toxin | 442.4 | 263 | 6.66 | 46 | 10 | 19 | 18 | X |
| HT-2-Toxin | 442.4 | 214.9 | 6.66 | 46 | 10 | 17 | 14 | |
| HT-2-Toxin | 442.4 | 169 | 6.66 | 46 | 10 | 35 | 10 | |
| Fumonisin B1 | 722.2 | 334.2 | 7.02 | 146 | 10 | 53 | 24 | X |
| Fumonisin B1 | 722.2 | 352.2 | 7.02 | 146 | 10 | 49 | 30 | |
| Fumonisin B1 | 722.2 | 704.2 | 7.02 | 146 | 10 | 39 | 20 | |
| Fumonisin B1 | 722.2 | 67.1 | 7.02 | 146 | 10 | 129 | 10 | |
| Fumonisin B1 | 722.2 | 91 | 7.02 | 146 | 10 | 129 | 10 | |
| Ochratoxin B | 370.1 | 205 | 7.68 | 45 | 10 | 27 | 15 | X |
| Ochratoxin B | 370.1 | 187 | 7.68 | 45 | 10 | 35 | 12 | |
| T2-Toxin | 484.1 | 305.1 | 8.23 | 51 | 10 | 19 | 20 | X |
| T2-Toxin | 484.1 | 214.9 | 8.23 | 51 | 10 | 29 | 14 | |
| T2-Toxin | 484.1 | 185 | 8.23 | 51 | 10 | 31 | 12 | |
| Zearalenone (ZON) | 319 | 283.2 | 9.2 | 66 | 10 | 17 | 20 | X |
| Zearalenone (ZON) | 319 | 186.9 | 9.2 | 66 | 10 | 27 | 12 | |
| Zearalenone (ZON) | 319.4 | 185 | 9.2 | 66 | 10 | 35 | 8 | |
| Zearalenone (ZON) | 319 | 69.3 | 9.2 | 66 | 10 | 37 | 8 | |
| Ochratoxin A | 404.2 | 238.9 | 9.65 | 41 | 10 | 33 | 16 | X |
| Ochratoxin A | 404.2 | 358 | 9.65 | 41 | 10 | 21 | 10 | |
| Ochratoxin A | 404.2 | 220.8 | 9.65 | 41 | 10 | 49 | 14 | |
| Fumonisin B2 | 707.2 | 319.2 | 9.8 | 176 | 10 | 53 | 24 | X |
| Fumonisin B2 | 707.2 | 355.3 | 9.8 | 176 | 10 | 45 | 32 | |
| Fumonisin B2 | 707.2 | 55.1 | 9.8 | 176 | 10 | 129 | 26 | |
| Fumonisin B2 | 707.2 | 67 | 9.8 | 176 | 10 | 127 | 8 | |
| Fumonisin B2 | 707.2 | 69.1 | 9.8 | 176 | 10 | 107 | 10 | |

Chromatography of expanded Mycotoxins overlaid MRMs.



APPLICATIONS

Conclusion

Analysis of aflatoxins and mycotoxins can be performed using simple, well-characterized techniques and equipment readily available to most testing labs. As new regulations are developed across North America, mycotoxin testing analysis will likely become standard. Expanding the targeted scope to include a wider range of potential fungal contaminants within cannabis and hemp extracts will improve public safety as well as scientific credibility in the cannabis testing community.

Ordering Information



roQ Extraction Kits

Extraction kits contain fifty easy-pour salt packets and fifty 50 mL stand-alone centrifuge tubes

| Description | Unit | Part No. |
|--|-------|--------------------------|
| AOAC 2007.01 Method Extraction Kits | | |
| 6.0 g MgSO ₄ , 1.5 g NaOAc | 50/pk | KS0-8911 |
| EN 15662 Method Extraction Kits | | |
| 4.0 g MgSO ₄ , 1.0 g NaCl, 1.0 g SCTD, 0.5 g SCDS | 50/pk | KS0-8909 |
| Original Non-buffered Method Extraction Kits | | |
| 4.0 g MgSO ₄ , 1.0 g NaCl | 50/pk | KS0-8910 |
| 6.0 g MgSO ₄ , 1.5 g NaCl | 50/pk | KS0-8912 |

roQ Extraction Salt Packets

Salt packets only. Centrifuge tubes not included.

| Description | Unit | Part No. |
|--|-------|--------------------------|
| AOAC 2007.01 Method Extraction Packets | | |
| 6.0 g MgSO ₄ , 1.5 g NaOAc | 50/pk | KS0-9043 |
| EN 15662 Method Extraction Packets | | |
| 4.0 g MgSO ₄ , 1.0 g NaCl, 1.0 g SCTD, 0.5 g SCDS | 50/pk | KS0-9041 |
| Original Non-Buffered Method Extraction Packets | | |
| 4.0 g MgSO ₄ , 1.0 g NaCl | 50/pk | KS0-9042 |
| 6.0 g MgSO ₄ , 1.5 g NaCl | 50/pk | KS0-9044 |

APPLICATIONS

Luna[®] Omega Polar C18 LC Column Ordering Information

| 1.6 µm Minibore Columns (mm) | | | | | SecurityGuard [™] ULTRA Cartridges [‡] |
|------------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|---|
| Phases | 30 x 2.1 | 50 x 2.1 | 100 x 2.1 | 150 x 2.1 | 3/pk |
| Polar C18 | 00A-4748-AN | 00B-4748-AN | 00D-4748-AN | 00F-4748-AN | AJ0-9505 |

for 2.1 mm ID

| 3 µm Minibore and MidBore [™] Columns (mm) | | | | | | | | SecurityGuard [™] Cartridges (mm) |
|---|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|---|
| Phases | 30 x 2.1 | 50 x 2.1 | 100 x 2.1 | 150 x 2.1 | 50 x 3.0 | 100 x 3.0 | 150 x 3.0 | 4 x 2.0* |
| Polar C18 | 00A-4760-AN | 00B-4760-AN | 00D-4760-AN | 00F-4760-AN | 00B-4760-YO | 00D-4760-YO | 00F-4760-YO | AJ0-7600 |

for ID: 2.0 - 3.0 mm

| 3 µm Analytical Columns (mm) | | | | | SecurityGuard [™] Cartridges (mm) |
|------------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|---|
| Phases | 50 x 4.6 | 100 x 4.6 | 150 x 4.6 | 250 x 4.6 | 4 x 3.0* |
| Polar C18 | 00B-4760-E0 | 00D-4760-E0 | 00F-4760-E0 | 00G-4760-E0 | AJ0-7601 |

for ID: 3.2-8.0 mm

Kinetex[®] Core-Shell Column Ordering Information

| 1.7 µm Minibore Columns (mm) | | | | SecurityGuard [™] ULTRA Cartridges [‡] |
|------------------------------|-----------------------------|-----------------------------|-----------------------------|---|
| Phases | 50 x 2.1 | 100 x 2.1 | 150 x 2.1 | 3/pk |
| Biphenyl | 00B-4628-AN | 00D-4628-AN | 00F-4628-AN | AJ0-9209 |

for 2.1 mm ID

*SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#)

| 2.6 µm Minibore Columns (mm) | | | | SecurityGuard [™] ULTRA Cartridges [‡] | |
|------------------------------|-----------------------------|-----------------------------|-----------------------------|---|--------------------------|
| Phases | 30 x 2.1 | 50 x 2.1 | 100 x 2.1 | 150 x 2.1 | 3/pk |
| Biphenyl | 00A-4622-AN | 00B-4622-AN | 00D-4622-AN | 00F-4622-AN | AJ0-9209 |

for 2.1 mm ID

| 2.6 µm MidBore Columns (mm) | | | | SecurityGuard [™] ULTRA Cartridges [‡] |
|-----------------------------|-----------------------------|-----------------------------|-----------------------------|---|
| Phases | 50 x 3.0 | 100 x 3.0 | 150 x 3.0 | 3/pk |
| Biphenyl | 00B-4622-Y0 | 00D-4622-Y0 | 00F-4622-Y0 | AJ0-9208 |

for 3.0 mm ID

| 2.6 µm Analytical Columns (mm) | | | | SecurityGuard [™] ULTRA Cartridges [‡] |
|--------------------------------|-----------------------------|-----------------------------|-----------------------------|---|
| Phases | 50 x 4.6 | 100 x 4.6 | 150 x 4.6 | 3/pk |
| Biphenyl | 00B-4622-E0 | 00D-4622-E0 | 00F-4622-E0 | AJ0-9207 |

for 4.6 mm ID

* SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#)

* SecurityGuard Analytical Cartridges require holder, Part No.: [KJ0-4282](#)



Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions
Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks
Luna and Kinetex are registered trademarks, roQ, MidBore, and SecurityGuard are trademarks of Phenomenex. Sciex is a registered trademark and Triple Quad is a trademark of AB SCIEX Pte. Ltd. AB SCIEX[™] is being used under license.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 5, 863, 428
CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP, or ULTRA holders, or to any cartridges.
FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2019 Phenomenex, Inc. All rights reserved.

- Australia**
t: +61 (0)2-9428-6444
auinfo@phenomenex.com
- Austria**
t: +43 (0)1-319-1301
anfrage@phenomenex.com
- Belgium**
t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com
- Canada**
t: +1 (800) 543-3681
info@phenomenex.com
- China**
t: +86 400-606-8099
cninfo@phenomenex.com
- Denmark**
t: +45 4824 8048
nordicinfo@phenomenex.com
- Finland**
t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com
- France**
t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com
- Germany**
t: +49 (0)6021-58830-0
anfrage@phenomenex.com
- India**
t: +91 (0)40-3012 2400
indiainfo@phenomenex.com
- Ireland**
t: +353 (0)1 247 5405
eireinfo@phenomenex.com
- Italy**
t: +39 051 6327511
italiainfo@phenomenex.com
- Luxembourg**
t: +31 (0)30-2418700
nlinfo@phenomenex.com
- Mexico**
t: 01-800-844-5226
tecnicomx@phenomenex.com
- The Netherlands**
t: +31 (0)30-2418700
nlinfo@phenomenex.com
- New Zealand**
t: +64 (0)9-4780951
nzinfo@phenomenex.com
- Norway**
t: +47 810 02 005
nordicinfo@phenomenex.com
- Portugal**
t: +351 221 450 488
ptinfo@phenomenex.com
- Singapore**
t: +65 800-852-3944
sginfo@phenomenex.com
- Spain**
t: +34 91-413-8613
espinfo@phenomenex.com
- Sweden**
t: +46 (0)8 611 6950
nordicinfo@phenomenex.com
- Switzerland**
t: +41 (0)61 692 20 20
swissinfo@phenomenex.com
- Taiwan**
t: +886 (0) 0801-49-1246
twinfo@phenomenex.com
- United Kingdom**
t: +44 (0)1625-501367
ukinfo@phenomenex.com
- USA**
t: +1 (310) 212-0555
info@phenomenex.com
- All other countries Corporate Office USA** 
t: +1 (310) 212-0555
info@phenomenex.com