

explore

LUNA<sup>®</sup>

OMEGA

# Luna Omega C18

## Application eBook

### Your Go-To C18:

- Greater separation muscle
- Better peak shape through an inert foundation
- Extreme ruggedness and dependability



[www.phenomenex.com/LunaOmega](http://www.phenomenex.com/LunaOmega)

 **phenomenex<sup>®</sup>**

# Luna Omega

## Cutting Edge Fully Porous Silica Particle

Luna® is one of the most recognized HPLC brands on the market, delivering high efficiency, ruggedness, reproducibility, and dependability for a wide range of analyses. The new Luna Omega builds upon this legacy with an innovative yet rugged silica particle architecture, designed and manufactured by Phenomenex based on more than 20 years of applied knowledge, invention, and customer experience.


### Novel Design and Manufacturing Process

Within the novel manufacturing process of Luna Omega silica, we implement a proprietary processing technique to gain greater particle inertness, a stronger particle morphology, and more consistent porosity.





### Thermal Modified Pore Structure

Most importantly, through our proprietary process, we eliminate micropores, further improving column efficiency, inertness, and reproducibility.



Consistent Porosity

Absence of Micropores



# Enhanced with 20 Years of Technology, Innovation, and Experience

One of the world's leading HPLC brands, now enhanced for incredible LC performance! Luna Omega LC columns culminate 20 years of technological prowess, advancements, and innovation from Phenomenex!

With **astounding efficiency levels**, highly **versatile selectivities**, and **trusted accuracy**, Luna Omega columns will take your LC experience to a new level.



Luna® Omega UHPLC columns  
will boost your LC instrumentation!

## Table of Contents

<b>Omega</b>	Astounding Performance .....	p. 4
	Separation Muscle .....	p. 5
	Inert Foundation .....	p. 6
	Reproducible and Scalable .....	p. 7
<b>C18</b>	Selectivity Highlight Luna Omega C18 .....	p. 8
	Small and Large Compound Mixtures .....	p. 9
	Small Molecule Organic Compounds .....	p. 10
	C18 Application Highlights .....	pp. 11 - 13
	Application Highlights Overview .....	p. 14
	Ordering Information .....	p. 15






# Astounding Performance

The undeniably high efficiency levels found in each Luna Omega UHPLC/HPLC column provides you with the potential of huge gains in method performance. While traditional silica and hybrid fully porous particles claim high performance, when compared to Luna Omega, they drastically fall short and prevent UHPLC scientists from reaching their UHPLC potential.

## UHPLC Efficiency Comparison




**Conditions for all columns:**  
**Dimension:** 50 x 2.1 mm  
**Mobile Phase:** Acetonitrile/Water (65:35)  
**Flow Rate:** 0.5 mL/min  
**Temperature:** Ambient  
**Detection:** UV @ 254 nm  
**System:** ACQUITY UPLC®  
**Sample:** Naphthalene\*

Comparative separations may not be representative of all applications.


\*Efficiency (plates/meter) comparison is based on peak performance associated with the compound naphthalene for all columns.

Our industry leading bonding technologies in conjunction with high efficiency levels ensure excellent stationary phase coverage and improved separation power. Now, with Luna Omega 1.6, 3, or 5  $\mu\text{m}$ , you can turn difficult separations into resolution that is reliable and repeatable.


## Luna® Omega 1.6 $\mu\text{m}$ C18



## ACE® Excel™ 2 $\mu\text{m}$ C18-AR



## Agilent® ZORBAX® 1.8 $\mu\text{m}$ XDB-C18



**Conditions for all columns:**

**Columns:** Luna Omega 1.6 $\mu\text{m}$  C18  
ZORBAX 1.8 $\mu\text{m}$  XDB-C18  
ACE Excel/2 $\mu\text{m}$  C18-AR

**Dimension:** 50 x 2.1 mm

**Mobile Phase:** A: 0.1% Formic Acid in Water  
B: 0.1% Formic Acid in Methanol

Gradient	Time (min)	% B
	0	3
	7	90

**Flow Rate:** 0.3 mL/min

**Temperature:** 30 °C


**Detection:** MS/MS

- Sample:**
1. Succinic acid
  2. MMA
  3. Glutaric acid
  4. Methylsuccinate
  5. Ethylmalonic acid
  6. Hippuric acid
  7. Homovanillic acid
  8. Suberic acid


Comparative separations may not be representative of all applications.

Luna Omega HPLC and UHPLC columns contain a unique silica that is modified by using a proprietary, post-synthetic thermal treatment process to provide extraordinary mechanical strength and significantly greater inertness than traditional fully porous and hybrid materials. With this process, and our commitment to continuous improvement, the Luna Omega column is manufactured reproducibly from column-to-column, to batch-to-batch.


## Luna® Omega 1.6µm C18



## Waters® ACQUITY® BEH 1.7µm C18



## Thermo Scientific® Hypersil GOLD® 1.9µm C18



### Conditions for all columns:

**Columns:** Luna Omega 1.6µm C18  
ACQUITY BEH 1.7µm C18  
Hypersil GOLD 1.9µm C18

**Dimension:** 50 x 2.1 mm

**Mobile Phase:** A: 0.1% Formic Acid in Water  
B: 0.1% Formic Acid in Acetonitrile

Gradient:	Time (min)	% B
	0	5
	5	95
	6	95
	6.1	5
	8	5

**Flow Rate:** 0.4 mL/min  
**Temperature:** Ambient  
**Detection:** UV @ 254 nm

**Sample:** 1. Pindolol  
2. Chlorpheniramine  
3. Nortriptyline  
4. 3-Methyl-4-nitrobenzoic acid  
5. 5-Methyl salicylaldehyde  
6. Hexanophenone

Comparative separations may not be representative of all applications.




By setting a new standard for reliability, the Luna® Omega C18 spans UHPLC and HPLC with a scalable range of high-performance particle sizes that will ensure that your developed methods are easily transferred. From single compound identification to complex impurity profiles, the Luna Omega C18 will serve as a pillar for your lab to count on day in and day out.

## Batch-to-Batch Reproducibility Study


In this example, we compared three batches of Luna Omega C18 using all three different particle sizes on a complex QC Pharmaceutical representative sample.

### Conditions for all columns:

- Mobile Phase:** A: Water with 0.1 % Formic Acid  
B: Acetonitrile with 0.1 % Formic Acid
- Temperature:** 30 °C
- Detection:** UV @ 254 nm
- Injection Volume:** 5 µL
- Sample:** 5 mg/mL of Chlorhexidine and Related Substances



**Luna Omega 5 µm C18**  
Impurity Profile 3 Batch Comparison




**Column:** Luna Omega 5 µm C18  
**Dimension:** 250 x 4.6 mm  
**Part No.:** 00G-4785-E0

Gradient: Time (min)	% B
0	2
2.5	2
52.5	35
55	35
57.5	2
62.5	2



**Luna Omega 3 µm C18**  
Impurity Profile 3 Batch Comparison




**Column:** Luna Omega 3 µm C18  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4784-E0

Gradient: Time (min)	% B
0	2
1.5	2
31.5	35
34.5	35
36	2
42	2



**Luna Omega 1.6 µm C18**  
Impurity Profile 3 Batch Comparison



**Column:** Luna Omega 1.6 µm C18  
**Dimension:** 50 x 2.1 mm  
**Part No.:** 00B-4742-AN

Gradient: Time (min)	% B
0	2
0.5	2
10.5	35
11.5	35
12	2
14	2


# Ω Selectivity Highlight Luna Omega C18

## Selectivity Highlight Luna Omega C18

Luna® Omega C18 is an excellent first choice for chromatographers who are just starting method development or attempting to improve upon existing chromatographic results with other C18s. With its higher performance potential, excellent retention profile, and greater inertness, the Luna Omega C18 was designed to be the new reproducible, all-purpose HPLC/UHPLC solution for all industries.

<b>Phase</b>	C18
<b>Particle Size</b>	1.6, 3, 5 μm
<b>Pore Size</b>	100 Å
<b>pH Range</b>	1.5 - 8.5*
<b>Surface Area</b>	260 m <sup>2</sup> /g
<b>Carbon Load</b>	11 %
<b>Pressure Limit</b>	1034/600 bar**
<b>USP Listing</b>	L1


\*pH stability under gradient conditions. pH stability is 1.5-10 under isocratic conditions.  
\*\*1.6 μm Luna Omega columns are pressure stable up to 1,034 bar and 3 or 5 μm are stable up to 600 bar.



## Greater Retention and Better Results


Higher efficiency levels in combination with excellent stationary phase coverage and greater particle inertness, translates to improved separation power for you. Now you can utilize the greater retention of Luna Omega C18 to tackle both easy and difficult separations.

### Luna Omega 1.6 μm C18



Greater Peak Capacity and Retention

### Waters® ACQUITY® BEH 1.7 μm C18



#### Conditions for all columns:

**Columns:** Luna Omega 1.6 μm C18  
ACQUITY BEH 1.7 μm C18

**Dimension:** 50 x 2.1 mm

**Mobile Phase:** A: 0.1 % Formic Acid in Water  
B: 0.1 % Formic Acid in Acetonitrile

Gradient:	Time (min)	% B
	0	10
	3	55
	3.5	55
	3.51	10
	5	10

**Flow Rate:** 0.4 mL/min

**Temperature:** Ambient

**Detection:** UV @ 205 nm


**Sample:** 1. Acetaminophen  
2. 4-Aminobenzoic Acid  
3. 4-Hydroxybenzoic Acid  
4. 2-Acetaminophenol  
5. 3-Hydroxybenzoic Acid  
6. Salicylicamide  
7. Phenol  
8. Benzoic Acid  
9. Salicylic Acid

Comparative separations may not be representative of all applications.

# Small and Large Compound Mixtures

Strong and focused hydrophobic retention, incredible efficiency and valuable inertness of Luna® Omega C18 columns make them an excellent choice for small mixtures of compounds differing in hydrophobicity as well as large mixtures of compounds like impurity/degradation profiles and peptide maps.

## Phenols




**Columns:** Luna Omega 1.6µm C18  
**Dimension:** 50 x 2.1 mm  
**Part No.:** 00B-4742-AN  
**Mobile Phase:** A: 0.1% Formic Acid in Water  
 B: 0.1% Formic Acid in Acetonitrile  
**Gradient:**

Time (min)	% B
0	5
6	50
7	50
7.1	5
9	5

**Flow Rate:** 0.4 mL/min  
**Temperature:** 22°C  
**Detection:** UV @ 270 nm  
**Sample:** 1. 3-Hydroxyphenol  
 2. Phenol  
 3. 4-Nitrophenol  
 4. 4-Methylphenol  
 5. 2-Methylphenol  
 6. 2,4-Dimethylphenol  
 7. 1-Naphthol

## Synthetic Cannabinoids




**Columns:** Luna Omega 1.6µm C18  
**Dimension:** 50 x 2.1 mm  
**Part No.:** 00B-4742-AN  
**Mobile Phase:** A: 0.1% Formic Acid in Water  
 B: 0.1% Formic Acid in Acetonitrile  
**Gradient:**

Time (min)	% B
0	50
6	100

**Flow Rate:** 0.3 mL/min  
**Temperature:** 30°C  
**Detection:** MS/MS (SCIEX® API 4000™)  
**Sample:** 1. JWH-073 Butanoic acid metabolite  
 2. JWH-073 3-hydroxybutyl metabolite  
 3. JWH-018 Pentanoic acid metabolite  
 4. JWH-073 4-hydroxybutyl metabolite  
 5. JWH-018 4-hydroxypentyl metabolite  
 6. AM-2201 4-hydroxypentyl metabolite  
 7. JWH-018 5-hydroxypentyl metabolite  
 8. AM-694  
 9. AM-2201  
 10. JWH-073  
 11. JWH-018

## OTC Drug - Pill Formulation Profile




**Columns:** Luna Omega 1.6µm C18  
**Dimension:** 50 x 2.1 mm  
**Part No.:** 00B-4742-AN  
**Mobile Phase:** A: 20 mM Potassium Phosphate pH 7.2  
 B: Acetonitrile  
**Gradient:**

Time (min)	% B
0	5
5	70
6	70
6.1	5
8	5

**Flow Rate:** 0.4 mL/min  
**Temperature:** Ambient  
**Detection:** UV @ 254 nm  
**Sample:** OTC Drug Pill

## Peptide Map – Digested BSA




**Columns:** Luna Omega 1.6µm C18  
**Dimension:** 100 x 2.1 mm  
**Part No.:** 00D-4742-AN  
**Mobile Phase:** A: 0.1% Formic Acid in Water  
 B: 0.1% Formic Acid in Acetonitrile  
**Gradient:**

Time (min)	% B
0	3
50	50
50.1	3

**Flow Rate:** 0.4 mL/min  
**Temperature:** 40°C  
**Detection:** MS/MS (SCIEX API 4000)  
**Sample:** Tryptic digest of BSA

Comparative separations may not be representative of all applications.

## Sildenafil Citrate RS and Sildenafil N-Oxide USP System Suitability Solution on Luna® Omega 5 µm C18




**Columns:** Luna Omega 5 µm C18  
**Dimension:** 150 x 4.6 mm  
**Part No.:** 00F-4785-E0  
**Elution Type:** Isocratic  
 Find the full elution online at [www.phenomenex.com/Application/Detail/26200](http://www.phenomenex.com/Application/Detail/26200)  
**Gradient:**

Time (min)	% B
20	0

  
**Flow Rate:** 1.4 mL/min  
**Temperature:** 30 °C  
**Detection:** UV @ 290 nm  
**Sample:** 1. Sildenafil citrate  
 2. Sildenafil N-Oxide

## Lisinopril Ph. Eur Reference Solution A+B+C+D on Luna Omega 5 µm C18



**Columns:** Luna Omega 5 µm C18  
**Dimension:** 250 x 4.6 mm  
**Part No.:** 00G-4785-E0  
**Elution Type:** Gradient  
 Find the full elution online at [www.phenomenex.com/Application/Detail/26294](http://www.phenomenex.com/Application/Detail/26294)  
**Gradient:**


Time (min)	% B
0	0
2	0
37	100
62	100
62.5	0
66	0

  
**Flow Rate:** 1.4 mL/min  
**Temperature:** 30 °C  
**Detection:** UV @ 290 nm  
**Sample:** 1. Impurity A  
 2. Lisinopril  
 3. Impurity E  
 4. Impurity F  
 5. Impurity G


## Ph. Eur. Monograph 401: Levothyroxine Sodium

The suitability of this Luna Omega reversed phase HPLC column for the analysis of related substances according to Ph. Eur. monograph 401 and effect of adapting the flow rate according to the allowable adjustments of chapter 2.246 of the European Pharmacopoeia is shown here. System suitability requires a resolution of greater than 5 for the separation of impurity A and Levothyroxine.


**Figure 2a: Reference **a** on Luna® Omega 3 µm C18 with flow rate 1.0 mL/min**




**Figure 2b: Reference **d** on Luna Omega 3 µm C18 with flow rate 1.0 mL/min**



**Figure 3a: Reference **a** on Luna Omega 3 µm C18 with flow rate 1.32 mL/min**



**Figure 3b: Reference **d** on Luna Omega 3 µm C18 with flow rate 1.32 mL/min**



### HPLC Conditions

**Columns:** Luna Omega 3 µm C18

**Dimension:** 150 x 4.6 mm

**Part No.:** 00F-4874-E0

**Elution Type:** A: 1.97 g phosphoric acid in 2 L water  
B: 1.97 g phosphoric acid in 2 L acetonitrile

Gradient	Time (min)	% B
	0	30
	10	30
	40	80
	50	80

**Flow Rate:** as indicated on the chromatograms

**Temperature:** 25 °C

**Injection:** 25 µL

**System:** Shimadzu® Nexera® XR

**Detector:** UV @ 225 nm


**See the full application note online at**

[www.phenomenex.com/TN77471120](http://www.phenomenex.com/TN77471120)

# Ω C18 Application Highlight

## Rapid Analysis of 125 Pesticides from Groundwater by UHPLC-MS/MS

With this application, 125 pesticides and degradants can be quantified at approximately 10 ng/L in groundwater using UHPLC-MS/MS with a Luna® Omega 1.6 µm C18 column. The high resolution afforded by the Luna Omega C18 gives great separation, allowing the mass spectrometer to maximize dwell time and achieve very low detection limits in a short run of 8 minutes.



App ID 23803



**Column:** Luna Omega 1.6 µm C18  
**Dimensions:** 50 x 2.1 mm  
**Part No.:** 00B-4742-AN  
**Mobile Phase:** A: 2.5 mM Ammonium acetate + 0.05% Acetic acid in water  
B: 0.05% Acetic acid in acetonitrile  
**Gradient:**

Time (min)	% B
0	15
0.5	15
2.5	40
5.75	75
6.25	100
7	100
7.1	15
8	15

**Injection:** 75 µL- Direct Inject  
**Flow Rate:** 0.55 mL/min  
**Temperature:** 45 °C  
**Detection:** Agilent® 6495 Triple Quadrupole LC-MS  
**Sample:** 125 Pesticides and Degradants (see link below)

See the full application note online at


[www.phenomenex.com/TN45280816](http://www.phenomenex.com/TN45280816)



## Analysis of Synthetic Cathinones (Bath Salts) from Urine and Whole Blood by LC-MS/MS

This application illustrates the separation of 5 common synthetic cathinones using a Luna® Omega 1.6 µm C18 UHPLC column. Luna Omega columns contain a unique silica which is modified using a proprietary, post-synthetic thermal treatment process to provide extraordinary mechanical strength and significantly greater inertness than traditional fully porous and hybrid materials. This greatly minimizes any secondary interactions that negatively affect peak shape, allowing for greater method accuracy.

### Extracted Ion Chromatogram of Bath Salts in Whole Blood



APP ID 23948

### Sample Preparation

#### Solid Phase Extraction (SPE)

<b>Cartridge:</b>	Strata®-X-Drug B
<b>Part No.:</b>	8B-S128-UCH
<b>Condition:</b>	1 mL each of methanol, DI water, and 100 mM sodium acetate
<b>Load:</b>	Load pretreated sample
<b>Weak Wash:</b>	2 mL of 100 mM sodium acetate (pH 5.0)
<b>Strong Wash:</b>	1 mL of methanol
<b>Dry Down:</b>	2 minutes at >10" hg
<b>Elute:</b>	3 mL of ethyl acetate: IPA: ammonium hydroxide (70:20:10)
<b>Evaporate:</b>	to 500 µL and add 100 µL of HCl:methanol (20:80), evaporate to dryness under nitrogen
<b>Reconstitute:</b>	100 µL of methanol

#### Pretreatment

<b>Urine</b>	Add 2 mL of 100 mM sodium acetate buffer (pH 5.0) and 50 µL of internal standards (@10 ppm) to 2 mL of urine
<b>Whole blood</b>	1. Add 2 mL of ice cold methanol:acetonitrile (50:50) and 20 µL of internal standards (@10 ppm) and 2 mL of 100 mM sodium acetate to 1 mL of blood. 2. Centrifuge at 4700 rpm 10 °C for 5 min 3. Transfer supernatant for SPE

Q1	Q2	Analyte	Retention Time (min)
208.1	160.2	Methylone	1.56
222.3	174.1	Ethylone	1.69
222.3	174.1	Butylone	1.80
178.1	160.2	Mephedrone	1.84
232.4	91	α-PVP	2.15

See the full application note online at

[www.phenomenex.com/TN48401216](http://www.phenomenex.com/TN48401216)



# Application Highlights Overview

## Why chose Luna Omega columns?

- Greater separation muscle
- Better peak shape through an inert foundation
- Extreme ruggedness and dependability



## Application Highlights

### Environmental



**Batch-to-Batch Comparison: Steroid Mixture on Luna Omega 3 µm C18**



**Scalability, Reproducibility, and Increased Separation Power of Luna Omega C18**

### Forensic/Toxicology



**Ph. Eur. Monograph 401: Levothyroxine Sodium on Luna C18(2) and Luna Omega C18**



**Rapid Analysis of 125 Pesticides from Groundwater by UHPLC-MS/MS using Luna Omega C18**

### Pharmaceuticals



**Analysis of Synthetic Cathinones (Bath Salts) from Urine and Whole Blood by LC-MS/MS using Luna Omega C18**



**Luna Omega C18 Column-to-Column and Batch-to-Batch Reproducibility from HPLC to UHPLC**



**Demonstrating the Luna Omega C18's Reproducibility, Scalability – HPLC to UHPLC – 1.6 µm, 3 µm, and 5 µm**



**Comparison of Two High-Performance Particle Morphologies in the Separation of Hydrochlorothiazide and Chlorothiazide**

## Column Characteristics

Luna Omega Phases	Description	Particle Sizes (µm)	Pore Size (Å)	Surface Area (m <sup>2</sup> /g)	Carbon Load (%)	pH Stability	Reversed Phase	HILIC
C18	C18 ligand optimized for improved peak shape	1.6, 3, 5	100	260	11	1.5 - 8.5*	☾	
Polar C18	Enhanced selectivity/retention for polar analytes without diminishing useful non-polar retention	1.6, 3, 5	100	260	9	1.5 - 8.5*	☾	
PS C18	Mixed-mode functionality offering enhanced retention of polar acids along with improved peak shape for strong bases	1.6, 3, 5	100	260	9	1.5 - 8.5*	☾	☾
Sugar	Combined amide/amino stationary phase and polar endcapping offers enhanced HILIC retention of sugars through multiple interaction mechanisms.	3	100	260	<2	2.0-7.0	☾	

\* pH stability under gradient conditions. pH stability is 1.5 - 10 under isocratic conditions.

# Ordering Information



1.6 µm Microbore Columns (mm)			
Phases	50 x 1.0	100 x 1.0	150 x 1.0
<b>Polar C18</b>	00B-4748-A0	00D-4748-A0	00F-4748-A0
<b>PS C18</b>	—	00D-4752-A0	—
<b>C18</b>	00B-4742-A0	00D-4742-A0	00F-4742-A0

1.6 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges <sup>†</sup>
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
<b>Polar C18</b>	00A-4748-AN	00B-4748-AN	00D-4748-AN	00F-4748-AN	AJO-9505
<b>PS C18</b>	00A-4752-AN	00B-4752-AN	00D-4752-AN	00F-4752-AN	AJO-9508
<b>C18</b>	00A-4742-AN	00B-4742-AN	00D-4742-AN	00F-4742-AN	AJO-9502

for 2.1 mm ID

3 µm Micro LC Columns (mm)							Trap Column
Phases	50 x 0.30	100 x 0.30	150 x 0.30	50 x 0.50	100 x 0.50	150 x 0.50	20 x 0.30
<b>Polar C18</b>	00B-4760-AC	00D-4760-AC	00F-4760-AC	00B-4760-AF	00D-4760-AF	00F-4760-AF	—
<b>PS C18</b>	00B-4758-AC	00D-4758-AC	00F-4758-AC	00B-4758-AF	00D-4758-AF	00F-4758-AF	05M-4758-AC

3 µm Minibore Columns (mm)					SecurityGuard Cartridges (mm)
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	4 x 2.0* /10 pk
<b>Polar C18</b>	00A-4760-AN	00B-4760-AN	00D-4760-AN	00F-4760-AN	AJO-7600
<b>PS C18</b>	00A-4758-AN	00B-4758-AN	00D-4758-AN	00F-4758-AN	AJO-7605
<b>C18</b>	—	00B-4784-AN	00D-4784-AN	00F-4784-AN	AJO-7611
<b>SUGAR</b>	—	00B-4775-AN	00D-4775-AN	00F-4775-AN	AJO-4496

for ID: 2.0-3.0 mm

3 µm MidBore™ Columns (mm)				SecurityGuard Cartridges (mm)
Phases	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0* /10 pk
<b>Polar C18</b>	00B-4760-Y0	00D-4760-Y0	00F-4760-Y0	AJO-7600
<b>PS C18</b>	00B-4758-Y0	00D-4758-Y0	00F-4758-Y0	AJO-7605
<b>C18</b>	00B-4784-Y0	00D-4784-Y0	00F-4784-Y0	AJO-7611
<b>SUGAR</b>	—	—	00F-4775-Y0	AJO-4496

for ID: 2.0-3.0 mm

3 µm Analytical Columns (mm)					SecurityGuard Cartridges (mm)
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0* /10 pk
<b>Polar C18</b>	00B-4760-E0	00D-4760-E0	00F-4760-E0	00G-4760-E0	AJO-7601
<b>PS C18</b>	00B-4758-E0	00D-4758-E0	00F-4758-E0	00G-4758-E0	AJO-7606
<b>C18</b>	00B-4784-E0	00D-4784-E0	00F-4784-E0	00G-4784-E0	AJO-7612
<b>SUGAR</b>	—	00D-4775-E0	00F-4775-E0	00G-4775-E0	AJO-4495

for ID: 3.2-8.0 mm

5 µm Minibore and MidBore™ Columns (mm)							SecurityGuard Cartridges (mm)
Phases	50 x 2.1	100 x 2.1	150 x 2.1	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0* /10 pk
<b>Polar C18</b>	00B-4754-AN	00D-4754-AN	00F-4754-AN	00B-4754-Y0	00D-4754-Y0	00F-4754-Y0	AJO-7600
<b>PS C18</b>	00B-4753-AN	00D-4753-AN	00F-4753-AN	00B-4753-Y0	00D-4753-Y0	00F-4753-Y0	AJO-7605

for ID: 2.0 - 3.0 mm

5 µm Analytical Columns (mm)					SecurityGuard Cartridges (mm)
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0* /10 pk
<b>Polar C18</b>	00B-4754-E0	00D-4754-E0	00F-4754-E0	00G-4754-E0	AJO-7601
<b>PS C18</b>	00B-4753-E0	00D-4753-E0	00F-4753-E0	00G-4753-E0	AJO-7606
<b>C18</b>	00B-4785-E0	00D-4785-E0	00F-4785-E0	00G-4785-E0	AJO-7612

for ID: 3.2-8.0 mm

5 µm Semi-Preparative Columns (mm)			SecurityGuard Cartridges (mm)
Phases	250 x 10	10 x 10** /3 pk	
<b>Polar C18</b>	00G-4754-N0		AJO-9519
<b>PS C18</b>	00G-4753-N0		AJO-9520

for ID: 9-16 mm

5 µm Axia™ Packed Preparative Columns (mm)					SecurityGuard Cartridges (mm)
Phases	50 x 21.2	100 x 21.2	150 x 21.2	250 x 21.2	15 x 21.2** /ea
<b>Polar C18</b>	00B-4754-P0-AX	00D-4754-P0-AX	00F-4754-P0-AX	00G-4754-P0-AX	AJO-7603
<b>PS C18</b>	00B-4753-P0-AX	00D-4753-P0-AX	00F-4753-P0-AX	00G-4753-P0-AX	AJO-7608
<b>C18</b>	—	—	—	00G-4785-P0-AX	—

for ID: 18-29 mm

5 µm Axia™ Packed Preparative Columns (mm) (cont'd)					SecurityGuard Cartridges (mm)
Phases	100 x 30	150 x 30	250 x 30	250 x 50	15 x 30.0* /ea
<b>Polar C18</b>	00D-4754-U0-AX	00F-4754-U0-AX	00G-4754-U0-AX	00G-4754-V0-AX	AJO-7604
<b>PS C18</b>	00D-4753-U0-AX	00F-4753-U0-AX	00G-4753-U0-AX	00G-4753-V0-AX	AJO-7609

for ID: 30-49 mm

<sup>†</sup> SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000

<sup>\*\*</sup> SecurityGuard Analytical Cartridges require holder, Part No.: KJO-4282

<sup>\*\*\*</sup> SemiPREP SecurityGuard Cartridges require holder, Part No.: AJO-9281

<sup>††</sup> PREP SecurityGuard Cartridges require holder, Part No.: AJO-8223

<sup>\*\*††</sup> PREP SecurityGuard Cartridges require holder, Part No.: AJO-8277

# Luna Omega C18

## Application eBook

explore

LUNA®  
OMEGA

**Australia**  
t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

**Austria**  
t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**  
t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**  
t: +1 (800) 543-3681  
info@phenomenex.com

**China**  
t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**  
t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**  
t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**  
t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**  
t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**  
t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**  
t: +852 6012 8162  
hkinfo@phenomenex.com

**India**  
t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**  
t: +62 21 5010 9707  
indoinfo@phenomenex.com

**Ireland**  
t: +353 (0)1 247 5405  
eirinfo@phenomenex.com

**Italy**  
t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**  
t: +81 (0) 120-149-262  
jinfo@phenomenex.com

**Luxembourg**  
t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**  
t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**  
t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**  
t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**  
t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**  
t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**  
t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**  
t: +65 800-852-3944  
sginfo@phenomenex.com

**Slovakia**  
t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**  
t: +34 91-413-8613  
esinfo@phenomenex.com

**Sweden**  
t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**  
t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com


**Taiwan**  
t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**  
t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**  
t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**  
t: +1 (310) 212-0555  
info@phenomenex.com

☎ **All other countries/regions**  
**Corporate Office USA**  
t: +1 (310) 212-0555  
info@phenomenex.com



[www.phenomenex.com](http://www.phenomenex.com)

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

#### Terms & Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

#### Trademarks

Luna and Strata are registered trademarks and Axia, SecurityGuard, and MidBore are trademarks of Phenomenex. ACE is a registered trademark and Excel is a trademark of Advanced Chromatography Technologies Limited. Waters, ACQUITY, and UPLC are registered trademarks of Waters Technologies Corporation. Agilent and ZORBAX are registered trademarks of Agilent Technologies, Inc. Thermo Scientific is a registered trademark and Syncronis is a trademark of Thermo Fisher Scientific. Hypersil GOLD is a registered trademark of Thermo Hypersil-Keystone LLC. AkzoNobel and Kromasil are registered trademarks of AkzoNobel Pulp and Performance Chemicals AB. YMC is a registered trademark of YMC Co., Ltd. SCIEX is a registered trademark and API 4000 is a trademark of AB SCIEX Pte. Ltd. Shimadzu and Nexera are registered trademarks of Shimadzu Corporation.

#### Disclaimer

Phenomenex is not affiliated with Agilent Technologies, AkzoNobel, Thermo Fisher Scientific, YMC, or Waters Corporation.

Comparative separations may not be representative of all applicants.

Axia column and packing technology is patented by Phenomenex. U.S. Patent No. 7, 674, 383

Strata-X is patented by Phenomenex. U.S. Patent No. 7,119,145

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362

**CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridges.**

**FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures**

© 2021 Phenomenex Inc. All rights reserved.