

# APPLICATIONS

## Comparison of Three Unique and Complementary Micro LC Columns and Three Trap Selectivities Under Reversed Phase LC-MS/MS Conditions

Roxana Eggleston-Rangel, Ryan Splitstone, and Jason Anspach  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

### Overview

Micro LC uses columns with internal diameters (I.D.) that are typically within the range of 0.3 - 0.5-mm and packed with the same chromatographic materials used in traditional analytical column dimensions. As is the case in analytical scale LC separations, the chromatographic media impacts the chromatographic performance and selectivity of the separation and can be optimized. A common practice in micro LC is to use a trap-and-elute injection mode. The use of a trap-and-elute injection mode allows for both the significant decrease in injection time and the protection of the column from sample contamination. Micro LC traps can influence the recovery and analytes from the sample based upon the stationary phase chemistry chosen.

In this application, we investigate the influence of different micro LC column and trap stationary phase selectivity combinations on chromatographic performance, sample recovery, and selectivity. We compared micro LC columns in a 50 x 0.3 mm column dimension packed with three different micro LC stationary phase chemistries: Kinetex® XB-C18, Luna® Omega PS C18, and Luna Omega Polar C18.

The comparison was generated using a sample of 20 stable-isotope-labeled (SIL) peptides under general reversed-phase mobile phase conditions and MS/MS detection using a SCIEX® 5500 QTRAP®. This application highlights the differences in chromatographic performance and recovery when different combinations of micro LC column and trap selectivity are combined in trap-and-elute mode. Adjusting the column and trap selectivity combination is a useful method development tool for the optimization of your micro LC separation.

**Kinetex 2.6 µm XB-C18:** A C18 that is modified with protective iso-butyl side chains for improved analysis of polar compounds.

**Luna Omega PS C18:** The surface contains a positive charged ligand which improves peak shape for basic compounds, while the C18 ligand promotes general reversed phase hydrophobic retention.

**Luna Omega Polar C18:** The C18 ligand provides general hydrophobic interactions while a polar modified particle surface provides enhanced polar compound retention

### Micro LC Trap Phases & Dimension

**MicroTrap C18:** 10 x 0.3 mm

**MicroTrap Polar:** 10 x 0.3 mm

**MicroTrap PS:** 10 x 0.3 mm



# APPLICATIONS

## Micro LC Conditions


**Column:** Luna® Omega 3 µm Polar C18  
**Dimension:** 50 x 0.3 mm  
**Part No.:** [00B-4760-AC](#)  
**Mobile Phase:** A: Water with 0.1 % Formic Acid  
 B: Acetonitrile with 0.1 % Formic Acid  
**Gradient:**

Time (min)	% B
0	3
10	40
12	80
14	80
15	3
20	3

**Flow Rate:** 10 µL/min  
**Temperature:** Ambient (25 °C)  
**Detector:** MS/MS SCIEX® QTRAP® 5500  
**Injector Temp.:** 4 °C  
**Column Temp.:** 25 °C  
**Injection Volume:** 1 µL  
**Sample:** 20 stable-isotope-labeled (SIL) Peptide mix

Precursor Mz	Product Mz	Peptide Sequence
485.25	856.41	IGNEQGVSR
491.27	769.37	LVGTPAEER
408.55	593.35	AETSELHTSLK
473.26	555.32	GAYVEVTAK
583.31	753.41	AVGANPEQLTR
519.8	422.29	LDSTSIPVAK
593.8	729.38	SAEGLDASASLR
739.36	999.51	YDSINNTEVSGIR
533.32	711.41	AGLIVAEGVTK
657.34	724.37	YIELAPGVDNSK
768.9	725.39	ALENDIGVPSDATVK
636.35	759.44	VGNEIQYVALR
677.86	649.41	DGTFAVDGGPGVIAK
758.91	957.51	SPYVITGPGVVEYK
549.29	721.39	GFTAYYIPR
964.98	596.34	TVESLFPPEAETPGSAVR
613.35	878.51	VFTPLEVDVAK
540.27	796.35	LGLDFDSFR
883.47	713.41	AVYFYAPQIPLYANK
569.83	711.42	SGLLWQLVR


## Direct Injection on Luna Omega 3 µm Polar C18 Representative Chromatogram without Any Micro LC Trap



App ID 25918

# APPLICATIONS

## Luna® Omega 3 µm Polar C18 Representative Chromatogram with MicroTrap C18



**Column:** Luna Omega 3 µm Polar C18  
**Trap:** MicroTrap C18  
**Dimensions:** 50 x 0.3 mm + 10 x 0.3 mm  
**Part No.:** [00B-4760-AC](#)  
[05N-4252-AC](#)

**Loading Flow Rate:** 16 µL/min  
**Loading Conditions:**


Time (min)	A%
0	100
2	100

App ID 25918

### Micro LC Column and Trap Selectivity Configuration Comparison Summary

**Figure 1** shows the average peak height (blue), width (orange) and tailing (gray) of a Luna Omega Polar C18 column in line with a C18, Polar, or PS MicroTrap. LC conditions are the same as described above.


**Figure 1.** Luna Omega Polar Column with MicroTrap C18, Polar or PS Comparison




App ID 25920

# APPLICATIONS

**Figure 2.** Area Recovery for Kinetex<sup>®</sup> XB-C18




**Figure 3.** Area Recovery for Luna<sup>®</sup> Omega PS C18



App ID 25920

**Figure 4.** Area Recovery for Luna Omega Polar C18



## Micro LC Conditions

**Column:** Kinetex 2.6  $\mu$ m XB-C18 ([00B-4496-AC](#))  
Luna Omega 3  $\mu$ m PS C18 ([00B-4758-AC](#))  
Luna Omega 3  $\mu$ m Polar C18 ([00B-4760-AC](#))

**Trap:** MicroTrap C18 ([05N-4252-AC](#))  
MicroTrap Polar ([05N-4754-AC](#))  
MicroTrap PS ([05N-4753-AC](#))

**Dimension:** 50 x 0.3 mm + 10 x 0.3 mm

**Mobile Phase:** A: Water with 0.1 % Formic Acid  
B: Acetonitrile with 0.1 % Formic Acid

Gradient:	Time (min)	% B
	0	3
	10	40
	12	80
	14	80
	15	3
	20	3

**Flow Rate:** 10  $\mu$ L/min

**Temperature:** Ambient (25  $^{\circ}$ C)

**Detector:** MS/MS SCIEX<sup>®</sup> QTRAP<sup>®</sup> 5500

**Injector Temp.:** 4  $^{\circ}$ C

**Column Temp.:** 25  $^{\circ}$ C

**Injection Volume:** 1  $\mu$ L

**Sample:** 20 stable-isotope-labeled (SIL) Peptide mix

**Figures 2 – 4** show the average peptide peak area recovery when using a Kinetex XB-C18, a Luna Omega PS C18 or a Luna Omega Polar C18 column in line with a C18, Polar or a PS MicroTrap.

In **Figure 2**, the MicroTrap Polar in combination with the Kinetex XB-C18 column delivers the highest peak area recovery for this example. The specific combination and differences between the column and trap's media can result in differences in relative recovery, as seen in **Figures 3 – 4**. The relative difference in recovery illustrates the importance of optimizing micro LC column and trap selectivity combinations and demonstrates how using a different selectivity combinations can yield improved recoveries depending upon the chemical characteristics of your compounds of interest. **Figure 3** illustrates recovery improvements that can be obtained by varying the selectivity between the column and trap. Generally, unlike guard columns, a media difference between the two separation devices can be advantageous.

# APPLICATIONS

Figures 5 – 7 compare the trap performance of MicroTrap C18, MicroTrap Polar, and MicroTrap PS in combination with the Kinetex® 2.6 µm XB-C18 and under the same chromatographic conditions.

Figure 5. Area Recovery for Kinetex XB-C18




Figure 6. Height Recovery for Kinetex XB-C18





Figure 7. Peak Width & Tailing for Kinetex XB-C18



## Conclusion

In conclusion, similarly to analytical scale LC, the performance and optimization of your separation is affected by the chosen stationary phase chemistry. The specific combination of column and trap media can result in differences in relative recovery and separation performance as shown in Figures 5 – 6. Figure 5, the combination of the MicroTrap Polar and Kinetex XB-C18 yielded the highest overall average peak area. However, the MicroTrap PS and Kinetex XB-C18 combination provided the best peak shape and recovery performance for this given separation as seen in Figure 7. This application demonstrates the performance differences between different micro LC column and trap selectivity configurations and the impact on optimization of your micro LC separation performance.

# APPLICATIONS

## Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/ChatNow](http://www.phenomenex.com/ChatNow).

### Australia

t: +61 (0)2-9428-6444  
auiinfo@phenomenex.com

### Austria

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

### Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

### Canada

t: +1 (800) 543-3681  
info@phenomenex.com

### China

t: +86 400-606-8099  
cninfo@phenomenex.com

### Denmark

t: +45 4824 8048  
nordicinfo@phenomenex.com

### Finland

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

### France

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

### Germany

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

### India

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

### Ireland

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

### Italy

t: +39 051 6327511  
italiainfo@phenomenex.com

### Luxembourg

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

### Mexico

t: 01-800-844-5226  
tecnicomx@phenomenex.com

### The Netherlands

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

### New Zealand

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

### Norway

t: +47 810 02 005  
nordicinfo@phenomenex.com

### Poland

t: +48 (12) 881 0121  
pl-info@phenomenex.com

### Portugal

t: +351 221 450 488  
ptinfo@phenomenex.com

### Singapore

t: +65 800-852-3944  
sginfo@phenomenex.com

### Spain

t: +34 91-413-8613  
espinfo@phenomenex.com

### Sweden

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

### Switzerland

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

### Taiwan

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

### United Kingdom

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

### USA

t: +1 (310) 212-0555  
info@phenomenex.com

🌐 **All other countries/regions**  
**Corporate Office USA**  
t: +1 (310) 212-0555  
info@phenomenex.com

## [www.phenomenex.com](http://www.phenomenex.com)

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

### Trademarks

Luna and Kinetex are registered trademark of Phenomenex.  
SCIEX and QTRAP are registered trademarks and Scheduled MRM is a trademark of AB SCIEX Pte. Ltd. AB SCIEX™ is being used under license.  
*FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.*

© 2020 Phenomenex, Inc. All rights reserved.