

APPLICATIONS

A Fast and Effective Quantitation Method for Uracil, 5,6-Dihydrouracil, and 5-Fluorouracil from Human Serum by LC-MS/MS

Shahana Huq
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Introduction

Uracil, 5-Fluorouracil (5FU), and 5,6-Dihydrouracil (UH₂) offer unique markers for enzymatic activity. In order to determine the 5FU catabolic rate, uracil and its homologues are tested in correlation with each other to ultimately determine if sufficient metabolism is occurring. The primary goal of this study was to develop a sample preparation and LC-MS/MS method for quantitative analysis of uracil and its two homologues from human serum. This application focuses on developing a sample preparation method and LC-MS/MS analysis using a Strata[®]-X PRO solid phase extraction (SPE) and a Kinetex[®] PS C18 HPLC column respectively. The Kinetex PS C18 column is a positively charged, surface modified C18 phase that caters its unique selectivity to separating the very polar uracil and its homologues in this analysis, while the novel Strata-X PRO greatly reduces the phospholipids in the sample and provides cleaner extracts.


Materials and Methods

Reagents and Chemicals

Analytical reference standards and internal standards were purchased from Cerilliant[®] Corporation Corporation (Round Rock, TX, USA). Doubly stripped MS grade human serum and pooled human plasma K₂EDTA were obtained from Golden West Diagnostics and BioreclamationIVT[®] (Westbury, NY). All other reagents and chemicals were purchased from the Sigma-Aldrich[®] Company (St. Louis, MO). Ultrapure D.I. water was obtained from Sartorius arium[®] comfort I, courtesy of Sartorius Corporation (Bohemia, NY).


5,6-Dihydrouracil

Molar mass: 114.104 g/mol
Acidic pK_a: 11.73
LogP: -1.21




Uracil

Molar mass: 112.088 g/mol
Acidic pK_a: 8.8
LogP: -0.5



5-Fluorouracil

Molar mass: 130.078 g/mol
Acidic pK_a: 7.18
LogP: -0.66



SPE Procedure

Sample Pre-treatment: Add 100 µL Human serum (doubly stripped) to a tube. Add 900 µL of chilled (~0°C) Acetonitrile to the serum while vortexing. Centrifuge samples at 6000 rpm for 5 minutes.

96-Well Plate: Strata-X PRO, 30 mg/well

Part No.: 8E-S536-TGA

Condition: 1 mL Acetonitrile

Load: Pass the supernatant (From Pre-treatment) and apply vacuum to collect eluted extract

Dry Down: Under Nitrogen and heat at ~45 °C

Reconstitute: With 100 µL initial mobile phase

Note: For phospholipid analysis

(1) Human plasma K₂EDTA was utilized

(2) Direct injection (bypass dry down and reconstitute) of the eluted sample was made.

Quantitative Analysis for Uracil Homologues

Column: Kinetex 2.6 µm PS C18

Dimensions: 150 x 3.0 mm

Part No.: 00F-4780-YO

Mobile Phase: A: 0.1% Formic acid in Water
B: Methanol

Gradient:	Time (min)	% B
	0	7
	12	7

Flow Rate: 0.2 mL/min

Injection Volume: 10 µL

Column Temperature: 25 °C

Instrument: Agilent[®] 1260

Detection: MS/MS (SCIEX[®] Triple Quad[™] 4500, ESI Dual Polarity)
(Positive for U and UH₂, Negative for 5FU)

Qualitative Analysis for Phospholipids

Column: Kinetex 2.6 µm C18

Dimensions: 50 x 2.0 mm

Part No.: 00B-4462-AN

Mobile Phase: A: 0.1% Formic acid in Water
B: 0.1% Formic acid in Methanol

Gradient:	Time (min)	% B
	0	40
	0.5	95
	11.5	95
	11.51	40
	13.5	40

Flow Rate: 0.4 mL/min

Injection Volume: 2 µL

Column Temperature: 40 °C


Instrument: Agilent 1260

Detection: MS/MS (SCIEX[®] Triple Quad[™] 4500, ESI+)


Table 1.
Retention Time (RT), MRM Transition, and % Recovery for Analytes

Analyte	RT	Q1	Q3	Spiked conc. (ng/mL)	% Recovery	% CV
UH ₂	4.31	114.9	55.04	100	90	5.4
U	4.65	112.9	69.8	100	84	3.2
5FU	7.07	128.8	41.9	100	89	3.2
			86.1			
			58.9			
Uracil 1,3-15N ₂ (+Ve IS)	4.65	114.8	96.9	200	N/A	N/A
5 CI Uracil (-Ve IS)	10.82	145.1	42.1	200	N/A	N/A


Figure 1.
Representative Chromatogram of Extracted Human Serum Analyzed by a Kinetex® 2.6 µm PS C18 LC Column (Under ESI Positive Polarity)



Uracil 1,3-15N₂



Uracil



UH₂

App ID 25777

Figure 2.
Representative Chromatogram of Extracted Human Serum Analyzed by a Kinetex[®] 2.6µm PS C18 LC Column (Under ESI Negative Polarity)




Figure 3.
Representative Chromatogram for Qualitative Evaluation of Phospholipid in Extracted Human Plasma Samples With Two Different Clean-up Techniques Employed (A) Protein Precipitation (B) Strata[®]-X PRO




Figure 4.
Linearity Curve for Uracil Extracted Human Serum Over 200-fold (5 to 1000ng/mL) Dynamic Concentration Range




Figure 5.
Linearity Curve for UH₂ Extracted Human Serum Over 200-fold (5 to 1000ng/mL) Dynamic Concentration Range




Figure 6.
Linearity Curve for 5FU Extracted Human Serum Over 200-fold (5 to 1000ng/mL) Dynamic Concentration Range




Table 2.
Precision and Accuracy Data for QC Samples of Uracil

Expected Conc. (ng/mL)	Sample	Replicates (N = 4)	%CV	Accuracy
30	QC1-1	4	4.6	97.4
150	QC1-2	4	5.3	101.6
750	QC1-3	4	6.9	102

Table 3.
Precision and Accuracy Data for QC Samples of UH₂

Expected Conc. (ng/mL)	Sample	Replicates (N = 4)	%CV	Accuracy
30	QC1-1	4	11.8	89.3
150	QC1-2	4	8.1	96.2
750	QC1-3	4	5.3	94.8

Table 4.
Precision and Accuracy Data for QC Samples of 5FU

Expected Conc. (ng/mL)	Sample	Replicates (N = 4)	%CV	Accuracy
30	QC1-1	4	2.4	106.7
150	QC1-2	4	10.1	105.0
750	QC1-3	4	7.8	108.9

Discussion

The unique chemistry offered by the positively charged Kinetex[®] PS C18 LC column retains the uracil and its homologues (UH₂ and 5FU) adequately to provide great resolution between these 3 compounds (**Figure 1 and 2**). The column selectivity allows the mass spec to utilize a polarity switch from positive to negative and simultaneously analyze all compounds in one run. The qualitative analysis demonstrates the superiority in cleanliness while sample obtained from protein precipitation undergoes extraction, utilizing Strata[®]-X PRO SPE (**Figure 3**). The selectivity of Strata-X PRO chemistry promotes elimination of majority of the phospholipids present in plasma sample, resulting in cleaner extract with minimum background noise. The absolute recoveries of the above three analytes in extracted samples yield more than 80% (**Table 1**). The QC samples for replicate extraction at 3 different levels showed precision and accuracy data between 2.4-11.8% and 89.3-108.9% respectively, that are within acceptable industry standard (**Table 2, 3, & 4**). The dynamic range of this method was tested with eight calibrators from 5 to 1000 ng/mL concentration range with an acceptable linearity ($r \geq 0.9995$) value (**Figure 4-6**).

Conclusion

The prescribed sample preparation method utilizing a protein precipitation followed by a 2 step Strata-X PRO extraction results in a simple and rapid quantitation of uracil and its analogues from human serum. The unique selectivity of the Kinetex PS C18 column demonstrated good chromatographic performance and overcame the challenges of analysis of these extremely polar compounds that the reversed phase columns encounter.

References

1. Barbara Buchel, Peter Rhyn, Stefan Schurch, *Biomed. Chromatogr.* **2013**; 27 : 7-16.
2. Ruta Svobaite, Isabella Solassol, Frederic Pinguet, *Clinical Chemistry* 54:9,1463-1472 (**2008**).
3. Jenny P. Dai, Amra Tabakovic, Welley Loc, *Current Trends in Mass Spectrometry* **October 2013**.

Kinetex[®] Core-Shell LC Column Ordering Information


2.6 µm Minibore Columns (mm)					SecurityGuard [™] ULTRA Cartridges [†]
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
PS C18	00A-4780-AN	00B-4780-AN	00D-4780-AN	00F-4780-AN	AJ0-8951 for 2.1 mm ID


2.6 µm MidBore [™] Columns (mm)				SecurityGuard [™] ULTRA Cartridges [†]
Phases	50 x 3.0	100 x 3.0	150 x 3.0	3/pk
PS C18	00B-4780-Y0	00D-4780-Y0	00F-4780-Y0	AJ0-8950 for 3.0 mm ID


2.6 µm Analytical Columns (mm)					SecurityGuard [™] ULTRA Cartridges [†]
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
PS C18	00B-4780-E0	00D-4780-E0	00F-4780-E0	00G-4780-E0	AJ0-8949 for 4.6 mm ID

[†] SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#).

**Ordering Information
Strata[®]-X PRO SPE**

Format	Sorbent Mass	Part Number	Unit
Tube			
	10 mg	8B-S536-AAK	1 mL (100/box)
	30 mg	8B-S536-TAK	1 mL (100/box)
	30 mg	8B-S536-TBJ	3 mL (50/box)
	60 mg	8B-S536-UBJ	3 mL (50/box)
	200 mg	8B-S536-FBJ	3 mL (50/box)
	100 mg	8B-S536-ECH	6 mL (30/box)
	200 mg	8B-S536-FCH	6 mL (30/box)
	500 mg	8B-S536-HCH	6 mL (30/box)

96-Well Plate			
	10 mg/well	8E-S536-AGA	ea
	30 mg/well	8E-S536-TGA	ea
	60 mg/well	8E-S536-UGA	ea

96-Well Microelution Plate			
	2 mg/well	8M-S536-4GA	ea

Presston[™] 1000 Positive Pressure Manifold

Part No.	Description
AH1-7033	Presston 1000 Positive Pressure Manifold, 96-Well Plate



Phenomenex warrants the Presston 1000 will be free of defects in materials and workmanship under normal installation, use, and maintenance for a period of 12 months following delivery. Please visit www.phenomenex.com/Presstonwarranty for complete warranty information.



APPLICATIONS

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

**☎ All other countries/regions
Corporate Office USA**

t: +1 (310) 212-0555
info@phenomenex.com



Your happiness is our mission. Take
45 days to try our products. If you
are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex and Strata are registered trademarks, BE-HAPPY, MidBore, Presston and SecurityGuard are trademarks of Phenomenex. Cerilliant is a registered trademark of Cerilliant Corporation. Sigma-Aldrich is a registered trademark of Sigma-Aldrich Co., BioreclamationIVT a registered trademark of Bioreclamation-IVT Holdings, LLC. Sartorius and arium are registered trademarks of Sartorius AG. Agilent is a registered trademark of Agilent Technologies, Inc. SCIEX is a registered trademark and Triple Quad are trademarks of AB SCIEX Pte. Ltd. AB SCIEXTM is being used.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is not affiliated with Sigma-Aldrich, Inc. or Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.


www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com