

Ph. Eur. Monograph 0049: Paracetamol Related Substances on Core-Shell Kinetex® 5 µm C18 column

Zeshan Aqeel¹, Dirk Hansen² and Heiko Behr³

¹Phenomenex, Inc., 411 Madrid Avenue, Torrance, CA 90501 USA

²Phenomenex Ltd. Deutschland, Zeppelinstraße 5, 63741 Aschaffenburg, Germany

³Phenomenex Helvetia GmbH, Schäferweg 16, 4057 Basel, Switzerland

Overview

N-(4-hydroxyphenyl) acetamide, commonly referred to as paracetamol, is one of the most familiar analgesics and antipyretic therapeutics in today's drug market.

For this report, we focused on the complex related substances profile of paracetamol as identified in the Ph. Eur. monograph 0049 published in supplement 10.7 in October 2021. This monograph is based on Pharmeuropa draft monograph (PA/PH/Exp. 10A/T (19) 136 ANP – 32.1)¹, which was already investigated in our TN-1274² technical note.

Experimental

The experiments were performed on an Agilent® 1290 binary UHPLC system equipped with a UV-VIS detection set at 254 nm (no reference wavelength was utilized). Analytical reference standards for paracetamol, paracetamol impurity K (4-Aminophenol) and impurity J ((4-Chloroacetanilide) were obtained from Sigma-Aldrich® (St. Louis, Missouri, USA) and evaluated with the Kinetex 5 µm C18 column (Phenomenex, Torrance, California, USA).

LC Conditions

Column: Kinetex 5 µm C18

Dimension: 150 x 4.6 mm

Part No.: [00F-4601-E0](#)

Mobile Phase: A: 1.7 g KH₂PO₄ + 1.8 g K₂HPO₄ in 1000 mL water
B: Methanol

Gradient:	Time(min)	%B
	0.0	5
	2.3	5
	15.2	10
	29.6	10
	58.4	34
	60.8	34

Flow Rate: 1.5 mL/min

Injection: 50 µL

Column Temperature: 30 °C

System: Agilent 1290 binary UHPLC system

Detection: UV @ 254 nm

Solution	Step 1	Step 2	Step 3	Step 4	Final Conc.
Test Solution	50.0 mg paracetamol	dissolve with 0.75 mL methanol	dilute to 5.0 mL with water		10 mg/mL paracetamol
Ref a	dilute 1.0 mL Test Solution to 100.0 mL with solvent mixture	dilute 1.0 mL of that solution to 20.0 mL with solvent mixture			5 µg/mL paracetamol
Ref b	5.0 mg of imp J	dissolve with 25 mL of methanol	dilute to 250.0 mL with solvent mixture	dilute 1.0 mL of that solution to 200.0 mL with solvent mixture	0.1 µg/mL imp J
Ref c	5.0 mg of imp K	dissolve with solvent mixture	dilute to 100.0 mL with solvent mixture	dilute 1.0 mL of that solution to 10.0 mL with solvent mixture	5 µg/mL imp K
Ref d	dilute 1.0 mL Ref c to 10.0 mL with solvent mixture				0.5 µg/mL imp K
Ref e (system suitability)	1 mL Ref a	+ 1 mL Ref c	dilute to 10 mL with solvent mixture		0.5 µg/mL paracetamol 0.5 µg/mL imp K




LC Conditions (continued)

The dwell volume of the HPLC system used for the development was 1.13 mL. Therefore, we adjusted the isocratic hold according to Ph. Eur. Chapter 2.2.46³ to reflect the dwell volume of 150 µL of the Agilent® 1290 binary UHPLC system used for this study (for details see [TN-1274²](#)).


Chromatograms and Data

Figure 1. Reference solution (a) on Kinetex® 5 µm C18



#	Analyte	t _R (min)	Area	Height	Width	Area %	Symmetry
1	Paracetamol	6.148	303.9	44.2	0.1057	100.0	1.023

Figure 2. Reference solution (b) on Kinetex 5 µm C18




#	Analyte	t _R (min)	Area	Height	Width	Area %	Symmetry
1	Impurity J	47.755	9.6	0.42	0.3834	100.0	0.864




Chromatograms and data (continued)

Figure 3. Reference solution (d) on Kinetex® 5 µm C18



#	Analyte	t _R (min)	Area	Height	Width	Area %	Symmetry
1	Impurity K	2.026	5.4	2	0.0409	100.0	0.73

Figure 4. Reference solution (e) on Kinetex 5 µm C18 (system suitability requirement R ≥ 5.0)




#	Analyte	t _R (min)	Resolution
1	Impurity K	2.01	
2	Paracetamol	6.24	33.55



Chromatograms and data (continued)

Figure 5. Test solution on Kinetex® 5 µm C18



#	Analyte	t _r (min)	Area	Height	Width	Area %	Symmetry
1	Paracetamol	5.890	156222.4	3892.5	0.6689	99.966	1.683
2	Impurity B	11.972	10.7	9.6E-1	0.1852	0.007	0.947
3	Impurity D	16.872	33.2	2.2	0.1808	0.021	0.957
4	Impurity H	34.05	8.8	3.9E-1	0.3783	0.006	1.158

The details on the peak identification are described in our technical note [TN-1274](#)².

Figure 6. Batch-to-batch reproducibility analysis of reference solution (e) on Kinetex 5 µm C18




Table 1. Reproducibility data reference solution (e) on 4 different batches of Kinetex® 5 µm C18

Batch	t _R Impurity K (min)	t _R Paracetamol (min)	Resolution
5701-0070	2.05	6.33	33.23
5701-0075	2.01	6.24	32.20
5701-0077	1.99	6.20	32.27
5701-0079	1.99	6.24	33.55

Conclusion

The above experiments show Kinetex 5 µm C18 is suitable under the conditions outlined in the monograph for paracetamol and even gives increased resolution for the system suitability solution, reference solution (e) (**Figure 4**). With the Kinetex 5 µm C18 column we also demonstrated batch-to-batch reproducibility (**Table 1**, retention times and resolution) across multiple (4) batches. Therefore, Kinetex 5 µm C18 is a reliable solution for the analysis of paracetamol in routine laboratories following the Ph. Eur. regulations. Please also refer to our technical note [TN-1274](#)² discussing the draft of the monograph for paracetamol. The changes in elution order for 3 impurities (I, J, and L) observed during the work has been communicated to the EDQM. As a result the elution order for these impurities has been corrected with the release of the new monograph.

References

1. Paracetamol monograph draft (PA/PH/Exp. 10A/T (19) 136ANP published in Pharmeuropa 32.1 (01/2020)
2. Zeshan Aqeel, Dirk Hansen, and Heiko Behr (2020) [TN-1274 European Pharmacopoeia Paracetamol Monograph Draft Method: Achieving Improved Sensitivity, Resolution, and Separation for Paracetamol and All 14 Related Impurities using Kinetex® 5 µm C18 Core-Shell Columns](#)
3. European Pharmacopoeia; Supplement 10 – Chapter 2.2.46 Chromatographic Separation Techniques.



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions
Corporate Office USA**

t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

**BE-HAPPY™
GUARANTEE**

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of Phenomenex, Inc. or their respective owners in the United States and/or certain other countries.

Disclaimer

Phenomenex is not affiliated with Agilent Technologies, Inc. or Sigma-Aldrich Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.

