




The Ultimate Pre-Packed Preparative Column for HPLC and SFC

GUARANTEED!

Axia PREP LC columns offer:

- Increased Performance
- Groundbreaking Lifetimes
- Optimized Loadability
- Increased Reproducibility



The Axia™ Advantage

Available in over 40 unique achiral and chiral selectivities, Axia advanced preparative column packing and column hardware design offer several advantages. Unlike traditional column packing methods, the Axia packing method offers increased sorbent bed density for increased performance and eliminates media bed collapse as a source of premature column failure in preparative HPLC/SFC columns.



guarantee

If Axia packed columns do not provide at least an equivalent separation as compared to a competing preparative column of the same particle size, same phase, and dimensions, return the column with comparative data within 45 days for a FULL REFUND. Only applies to 21.2 mm ID columns.

“

“I find Axia Columns to be very robust and durable. I often use the prep column for much longer than predicted with reproducible peaks. This saves us a significant amount of money.”

*David Wisnoski
GlaxoSmithKline, USA*

”

Axia Technology pp. 4-7

Award winning column packing technology

Media Selection Chart pp. 8-9

Quickly find the media for your purification needs

Media Selectivity Details pp. 10-25

Achiral and Chiral Chemistries

Kinetex® pp. 10-13

Gemini® pp. 14-15

Synergi™ pp. 16-17

Luna® pp. 18-19

Jupiter® p. 20

Clarity® p. 21

Lux® (chiral) pp. 22-25

Preparative Column Protection pp. 26-27

SecurityGuard™

PhenoLogix pp. 28-29

Screening, Optimization and Scale-up Services

Ordering Information pp. 30-31

We routinely use Axia packed columns from Phenomenex for peptide purifications. Among various preparative HPLC columns we have used, the Axia packed Luna columns (5µm) stand out. We have been very satisfied with the increased loading capacity and excellent performance.



Guangcheng Jiang
Ferring Research Institute, Inc., USA

Axia™ Technology vs. Traditional "OBD" Prep Column Packing

Traditional Slurry Packing

Traditional slurry packing processes, like the Waters® OBD™ (Optimum Bed Density) column packing approach, involve the column being removed from the column packing station once it is packed.

Several potential problems with this packing method are:

- Variability in column performance due to increased number of manual operations required for assembly
- Potential silica media damage during recompression
- Level of process control is based on traditional slurry packing technology



Conventional packing process involves:
 Compression → Decompression → Recompression → Final Column




Diagram from Waters Corporation U.S. Patent No. 7,399,410


Axia Packing Technology

Axia packed preparative columns involve a single axial compression step unlike conventional packed preparative columns. The ideal column bed density is custom calculated and automated for each specific media and column size. Computer control of the entire process ensures both proper bed density and column uniformity every time.

During the Axia packing process, the packing piston is locked in place, eliminating any decompression and then recompression of the media sorbent, thus maintaining media and column bed integrity. This solves common lifetime and performance problems associated with conventional packing processes for preparative columns.




Axia Packing Process Involves:
 Compression → Final Column



Traditional packed preparative columns produce non-uniform media beds with sheared and crushed particles


Decompression and then recompression during packing can damage the media and lead to increased column-to-column variability, flow disturbances, and decreased column lifetimes.



Axia packed columns produce uniform media bed with intact particles

The highly tuned patented process and hardware eliminates potential decompression ensuring bed stability and optimal packing density.

The media found on the inlet frit of the Axia packed column shows no signs of damage unlike the media found on inlet frit of traditionally packed prep columns.




“

We are using chromatography media from Phenomenex for GPL/GMP purposes, therefore we audited Phenomenex USA as a manufacturer. From the beginning, we were impressed with Phenomenex and the attitude of their employees. Phenomenex is a unique company in many aspects. Their degree of dedication to customer service, to the organization of the QMS system and last but not least the positive atmosphere in the company is impressive. The outcome of the audit was to our fullest satisfaction.

Major Generic Pharma Company, Europe

”



View an animated packing process comparison at www.AxiaPrep.com

Axia™ Technology Outperforms Traditional Packing Processes!

Because of the constant pressure placed on the integrated packing piston, Axia packed columns possess the dynamic capability of maintaining a consistent, homogeneous media bed. This results in superior column performance no matter which media selectivity you choose.


To better understand how much Axia technology improves column performance over traditionally slurry packed preparative columns we scaled-up a 5 µm Lux® Cellulose-1 chiral media analytical column and packed the same media into two different

150 x 21.2 mm I.D. columns. One column was packed using Axia technology and the other prep column was packed using the traditional slurry packing process.


The Axia packing technology had a substantial increase in column efficiency resulting in increased resolution over traditionally packed preparative columns. With increased resolution you are able to increase your sample load enabling you to purify more target compound(s) per purification run. This equates to better throughput and economics.

Warfarin Chiral Purification in Normal Phase Mode


Analytical




Column: Lux 5 µm Cellulose-1
Dimensions: 150 x 4.6 mm
Mobile Phase: Hexane/Ethanol (75:25)
Flow Rate: 1 mL/min
Temperature: Ambient
Inj. Volume: 100 µL



Standard Packing and Hardware



Axia Technology and Hardware



Conditions for both PREP columns:

Media: Lux 5 µm Cellulose-1
Dimensions: 150 x 21.2 mm
Mobile Phase: Hexane / Ethanol (75:25)

Flow Rate: 20 mL/min
Temperature: Ambient
Inj. Volume: 2 mL

Column (mm)	Analytical 150 x 4.6	Standard 150 x 21.2	Axia 150 x 21.2
Mass Loaded (mg)	2	40	40
Resolution*	1.5	2.85	3.72
Plates (N)	117	535	760

42% Increase in Efficiency

* Resolution calculated with peak width at baseline and center retention time due to the overloaded peaks being off-scale

Tip:


For more detailed information on this warfarin application, request technote:

TN-9002: Scaling from Analytical to Preparative Chiral Chromatography While Balancing Purity, Yield, and Throughput under HPLC and SFC Conditions


Unmatched Column Reproducibility

The completely automated Axia™ packing system provides feed-back control and infinite tuning of packing density for specific media characteristics such as mechanical strength and porosity. An optimum bed density can be consistently reproduced column-to-column.


This directly translates into consistent efficiency and peak asymmetry measurements and decreases the column variability seen in traditionally packed preparative columns.



Average Efficiency (N) with Synergi™
4 µm Hydro-RP 100 x 21.2 mm



Average Peak Asymmetry with Gemini®
5 µm C18 50 x 21.2 mm



Axia columns provide me with first rate quality and engineering. Reliability, reproducibility, and durability are provided with all Axia columns that I use. I can literally purify 2500 samples per column. The time and cost savings are tremendous.

Derrick Miyao
—Large Biotech Manufacturer, USA




We have used Phenomenex Axia prep-HPLC columns for several years and they consistently provide excellent separation and reproducibility for a variety of different compounds.


Jeremy R. Wolf
ABC Laboratories, USA



Phase Selection Chart

Quickly find the media for your purification needs.





Kinetex

1st Core-Shell Preparative Column Ever!
Pages 10-13



Gemini

High pH Separations
Pages 14-15



Synergi

Unique Chemistries for Complex Mixtures
Pages 16-17



Luna

Proven Purification Performance
Pages 18-19



Jupiter

Increase Loadability for Biomolecule Separations
Page 20



Clarity

Purification of Synthetic Oligonucleotides
Page 21



Lux

Polysaccharide Supports with Excellent Enantioselectivity
Pages 22-25



First Core-Shell Preparative HPLC/SFC Column Ever!

Kinetex® Core-Shell Technology produces increased efficiencies over traditional, fully porous columns, yielding remarkable chromatographic resolution, higher peak capacities, and greater sensitivity, so labs can get even more out of their HPLC analyses!

The benefits of Kinetex Core-Shell Technology include:


- Increased efficiencies over traditional fully porous columns
- Seamless scalability from HPLC/UHPLC to Preparative LC
- Kinetex 5 μm provides better performance than traditional fully porous 5 and 3 μm materials



High Column Efficiency

Combining 5 μm Kinetex core-shell and Axia™ technology can provide the highest separation efficiency of any pre-packed preparative HPLC column.


Waters® XBridge® 5 μm C18 Prep OBD™



Conditions for both columns:

Columns: Kinetex 5 μm XB-C18 Axia Packed
Waters XBridge 5 μm C18 Prep OBD
Dimensions: 150 x 21.2 mm (Kinetex)
150 x 19 mm (XBridge)
Mobile Phase: Water/ Acetonitrile (50:50)
Injection Volume: 10 μL

Kinetex 5 μm XB-C18 Axia Packed



Flow Rate: 25 mL/min
Temperature: Ambient
Detection: UV @ 254 nm
Sample: 1. Uracil
2. Acetophenone
3. Toluene
4. Naphthalene

Key: ● Best Suited ○ Very Good


Key: <div><div></div> Best Suited</div> <div><div></div> Very Good</div>						Applications					Type of Compounds				Loading				
	Small Molecules		Peptides		Proteins		Chiral		Oligonucleotides		Acids		Polar		Hydrophobic		Bases		Available Surface Area
Packing Material	Particle Size (μm)	Pore Size (Å)	Surface Area (m²/g)	Carbon Load (%)	pH Range														
Kinetex C18	1.3, 1.7, 2.6, 5	100	200	12	1.5-8.5*	<div></div>	<div></div>					<div></div>		<div></div>		<div></div>		<div></div>	
Kinetex XB-C18	1.7, 2.6, 5	100	200	10	1.5-8.5*	<div></div>	<div></div>					<div></div>		<div></div>		<div></div>		<div></div>	
Kinetex C8	1.7, 2.6, 5	100	200	8	1.5-8.5*	<div></div>	<div></div>					<div></div>		<div></div>		<div></div>		<div></div>	
Kinetex Phenyl-Hexyl	1.7, 2.6, 5	100	200	11	1.5-8.5*	<div></div>	<div></div>					<div></div>		<div></div>		<div></div>		<div></div>	

*Columns are pH stable from 1.5-10 under isocratic conditions. Columns are pH stable 1.5-8.5 under gradient conditions. Comparative separations may not be representative of all applications.

Excellent Loadability!

With narrower peak widths than fully porous columns across every sample load, Axia™ packed Kinetex 5 µm columns give you the capability of increased sample load and higher throughput for vastly improved purification performance and economics.

Waters® XBridge® 5 µm C18 Prep OBD™




Conditions for both columns:

Columns: Kinetex 5 µm C18 Axia Packed
 XBridge 5 µm C18 Prep OBD
Dimensions: 50 x 21.2 mm (Kinetex)
 50 x 19 mm (XBridge)
Mobile Phase: A: Water with 0.5 % Formic acid
 B: Acetonitrile with 0.5 % Formic acid
Gradient:

Time (min)	% B
0	20
8	50
11	50

Kinetex 5 µm C18 Axia Packed



Flow Rate: 30 mL/min
Temperature: Ambient
Detection: UV @ 254 nm
Sample: 200 mg/mL in DMSO
 1. Doxepin (From 1 - 500 mg on-column)
 2. Amitriptyline (From 1 - 500 mg on-column)

Kinetex Axia Preparative columns are fantastic! I currently use two Kinetex 5 µm C18 150 x 21.2 mm columns in parallel for high throughput purifications (<100 mg scale), and Kinetex core shell media delivers significantly improved peak shape and lower back pressure compared to many of the industry. I can also analyze quickly my purified fractions with the same core shell phase on my analytical UPLC® system.

Chris DeVore
 Neurocrine Biosciences, USA

Axia packed column has a great efficiency for the separation of several classes of natural compounds. Due to its low back pressure and therefore high flow work conditions, time for conditioning the columns is sped up greatly!


Sylvian Cretton
 -Europe

Seamless Scalability from HPLC/UHPLC to PREP

The recent addition of the Kinetex® 5 µm in the Axia™ packed format (21.2 mm ID) makes it the first core-shell sorbent commercially available for small-scale preparative applications. Combine this with the fact that the entire Kinetex core-shell line


is fully scalable from 1.3 µm to 5 µm, means that transferring high performance HPLC/UHPLC methods to preparative HPLC and SFC formats is fast and simple.

Analytical method — Kinetex 2.6 µm XB-C18




Column: Kinetex 2.6 µm XB-C18
Dimensions: 150 x 4.6 mm
Part No.: 00F-4496-E0
Mobile Phase: A: 0.1 % TFA in Water
 B: 0.1 % TFA in Acetonitrile
Gradient: Linear 85:15 (A/B) to 5:95 (A/B) over 10 minutes
Injection Volume: 10 µL
Flow Rate: 1 mL/min
Temperature: Ambient
Detection: UV @ 210 nm
Sample: Crude peptide mix

Preparative scale-up and fraction collection — Kinetex 5 µm XB-C18



Column: Kinetex 5 µm XB-C18 Axia Packed
Dimensions: 150 x 21.2 mm
Part No.: 00F-4605-P0-AX
Mobile Phase: A: 0.1 % TFA in Water
 B: 0.1 % TFA in Acetonitrile
Gradient: Linear 85:15 (A/B) to 5:95 (A/B) over 10 minutes
Injection Volume: 1 mL
Flow Rate: 20 mL/min
Temperature: Ambient
Detection: UV @ 210 nm
Sample: Crude peptide mix

Analytical fraction analysis — Kinetex 2.6 µm XB-C18



Column: Kinetex 2.6 µm XB-C18
Dimensions: 150 x 4.6 mm
Part No.: 00F-4496-E0
Mobile Phase: A: 0.1 % TFA in Water
 B: 0.1 % TFA in Acetonitrile
Gradient: Linear 85:15 (A/B) to 5:95 (A/B) over 10 minutes
Injection Volume: 10 µL
Flow Rate: 1 mL/min
Temperature: Ambient
Detection: UV @ 210 nm
Sample: Purified Fractions


Tip:

For more information on the power of Kinetex core-shell scalability, request technical note TN-1135 at:
www.phenomenex.com/Kinetex/AxiaRequest

A Broad Spectrum of Column Selectivities

Kinetex® core-shell columns are available in a wide range of stationary phases, allowing you to optimize your separation for maximum resolution and loadability across HPLC, UHPLC, and Preparative HPLC and SFC applications.

C18



Conditions for all columns:

Column: Kinetex 2.6 µm C18

Kinetex 2.6 µm XB-C18

Kinetex 2.6 µm Phenyl-Hexyl

Kinetex 2.6 µm C8

Dimensions: 50 x 2.1 mm

Mobile Phase: A: 0.1 % Formic acid in Water

B: 0.1 % Formic acid in Acetonitrile

Gradient	Time (min)	% B
	0.0	5
	0.2	5
	4.2	95
	4.21	5
	5.5	5


Flow Rate: 0.8 mL/min

Temperature: 30 °C


Detection: UV @ 254 nm (ambient)

- Sample:**
1. Pyridine
 2. Acetaminophen
 3. Sulfathiazole
 4. Pindolol
 5. Quinidine
 6. Benzyl alcohol
 7. Phenol
 8. Acebutolol
 9. Chlorpheniramine
 10. Triprolidine
 11. Prednisolone
 12. 3-Methyl-4-nitrobenzoic acid
 13. Nortriptyline
 14. 2-Hydroxy-5-methylbenzaldehyde
 15. Diflunisal
 16. Hexanophenone


XB-C18



Phenyl-Hexyl



C8



Setting the Standard for pH Method Development



Gemini features a pH stability from 1-12, making it optimal for high alkaline washes and high pH purifications of basic drugs.


Optimized parameters include:

- Innovative surface layer for increased pH stability
- High-surface area for increased loading
- Silica smoothness for stable packing beds
- Bonding density for excellent reproducibility

TWIN (Two-In-One) Technology™

Gemini C18 and C6-Phenyl


During the final stage of silica manufacturing a unique silica-organic layer is grafted to create a completely new composite particle. Since the internal base silica is unaltered by this manufacturing process, the particle retains its mechanical strength and rigidity of the silica. This provides excellent efficiency, while the silica-organic shell protects the particle from chemical attack at extreme pH conditions.



Second-Generation TWIN-NX™ Technology

Gemini NX-C18

TWIN-NX technology uses an improved patented organo-silica grafting process which incorporates highly stabilizing ethane cross-linking. These organic groups are evenly incorporated into the grafted layers on the silica surface while maintaining a pure silica core. This not only provides resistance to high pH attack, but also maintains the high efficiency and mechanical strength of a silica particle.




Key: <div><div></div> Best Suited</div> <div><div></div> Very Good</div>						Applications					Type of Compounds				Loading
						Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic	Bases	Available Surface Area
Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m²/g)	Carbon Load (%)	pH Range										
Gemini C18	3, 5, 10	110	375	14	1.0-12.0	<div></div>	<div></div>				<div></div>	<div></div>	<div></div>	<div></div>	<div></div>
Gemini C6-Phenyl	3, 5	110	375	12	1.0-12.0	<div></div>	<div></div>				<div></div>		<div></div>	<div></div>	<div></div>
Gemini NX-C18	3, 5, 10	110	375	14	1.0-12.0	<div></div>	<div></div>				<div></div>		<div></div>	<div></div>	<div></div>

Key: ● Best Suited ○ Very Good

Flexibility in pH adjustments allows for increased purification performance.


Separating basic compounds at higher pH levels produces dramatic changes when compared to low pH conditions. At pH 10.5, the basic compounds become neutralized and are more hydrophobic. The retention for the uncharged basic compounds increases providing an increase in separation along with superior peak shapes.

Gemini® NX-C18 with 0.5% TFA



Column: Gemini NX-C18 5 µm
Dimensions: 50 x 21.2 mm
Mobile Phase: A: 0.5% TFA in Water
 B: Acetonitrile
Gradient: 5% to 95% in 5 min
Flow Rate: 30 mL/min
Detection: UV @ 254 nm
Sample: 1. Diphenhydramine
 2. Propranolol


Gemini NX-C18 with 0.2% NH₄OH



Column: Gemini NX-C18 5 µm
Dimensions: 50 x 21.2 mm
Mobile Phase: A: 0.2% NH₄OH in Water
 B: Acetonitrile
Gradient: 5% to 95% in 5 min
Flow Rate: 30 mL/min
Detection: UV @ 254 nm
Sample: 1. Diphenhydramine
 2. Propranolol


Excellent resolution at high pH

Gemini NX-C18 Low pH 2 with 0.5% TFA



Column: Gemini NX-C18 5 µm
Dimensions: 50 x 21.2 mm
Mobile Phase: A: 0.5% TFA in Water
 B: Acetonitrile
Gradient: 5% to 95% in 5 min
Flow Rate: 30 mL/min
Detection: UV @ 254 nm
Sample: 1. Diphenhydramine
 2. Oxybutynin
 3. Terfenadine

Gemini NX-C18 with 0.2% NH₄OH



Column: Gemini NX-C18 5 µm
Dimensions: 50 x 21.2 mm
Mobile Phase: A: 0.2% NH₄OH in Water
 B: Acetonitrile
Gradient: 5% to 95% in 5 min
Flow Rate: 30 mL/min
Detection: UV @ 254 nm
Sample: 1. Diphenhydramine
 2. Oxybutynin
 3. Terfenadine

Separation shape improvement provides opportunity for increased loading

“Our Phenomenex Gemini and Luna Axia packed columns are the work-horses in our lab. These columns exhibit outstanding performance for challenging separations while also handling a high workload for standard separations. Longevity has also been excellent with some columns lasting 2 years or more. Dependability is so important in my line of work and these columns never disappoint!!

—Major Pharmaceutical Company, USA”

Tip:

If you want longer Gemini NX-C18 Axia packed column lifetimes, request:

TN-1138 Increase Column Performance and Lifetime in Peptide and Protein Purification using Aggressive Wash Conditions

Increased Loading with Unique Selectivities



Synergi is available in four unique phases, each offering dramatic differences in:


- Selectivity
- Retention time
- Resolution

The unique selectivity profiles found within the Synergi product line offer complementary selectivity to the standard C18, C8, or silica phases traditionally employed in preparative HPLC.

Synergi Polar-RP

For Polar and Aromatic Mixtures


(100 % Aqueous Stable)



Synergi Fusion-RP

Balanced Non-polar and Polar Performance


(100 % Aqueous Stable)



Synergi Hydro-RP


Strong Non-polar and Polar Retention

(100 % Aqueous Stable)



Synergi Max-RP

Excellent for Basic Compounds at Neutral pH




Packing Material	Particle Size (μm)	Pore Size (Å)	Surface Area (m ² /g)	Carbon Load (%)	pH Range	Applications					Type of Compounds				Loading Available Surface Area
						Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic	Bases	
Synergi Fusion-RP	4, 10	80	475	12	1.5-10.0*	●	●				●	●	●	●	●
Synergi Max-RP	4, 10	80	475	17	1.5-10.0*	●	●				●		●	●	●
Synergi Hydro-RP	4, 10	80	475	19	1.5-7.5	●	●				●	●	●	●	●
Synergi Polar-RP	4, 10	80	475	11	1.5-7.0	●	●				●	●		●	●

*pH range is 1.5-10 under isocratic conditions and 1.5-9.0 under gradient conditions.

Selectivity Like No Other

Offering a balanced combination of hydrophobic and polar selectivity, Synergi™ Fusion-RP separates compounds exhibiting moderately polar and hydrophobic characteristics.

Hydrophobic basic compounds




Conditions for all columns:

Columns: Synergi 4 µm Fusion-RP
Typical C18
Dimensions: 150 x 4.6 mm
Mobile Phase: 20 mM Potassium Phosphate,
pH 2.5 / Acetonitrile (75:25)
Flow Rate: 1.0 mL/min
Detection: UV @ 210 nm
Sample: 1. Maleic acid
2. Chlorpheniramine
3. Triprolidine
4. Diphenhydramine

The slightest variations in compound polarity and aromaticity are exploited by Synergi Polar-RP to achieve the greatest separation between polar and/or aromatic compounds.

Increased resolution of polar compounds with Synergi Polar-RP compared to traditional C18 phases



Conditions for all columns:

Columns: Synergi 4 µm Polar-RP
Waters 5 µm SymmetryShield RPC18
Waters 5 µm Symmetry C18
Waters 5 µm XTerra RP18
Dimensions: 150 x 4.6 mm
Mobile Phase: 20 mM Potassium phosphate pH 3 / Methanol (50:50)
Flow Rate: 1.0 mL/min
Detection: UV @ 230 nm
Temperature: Ambient
Injection: 2 µL
Sample: 1. Metoprolol (0.15 µg) 4. Alprenolol (0.3 µg)
2. Pindolol (0.6 µg) 5. Propranolol (0.04 µg)
3. Metoprolol (0.15 µg) 6. Ethylparaben (0.4 µg)



“We regularly use RP stationary phases from Phenomenex for our separation problems. Especially Synergi Polar-RP was found to often show the desired selectivity, distinguishing this phase from other RP phases.”

CARBOGEN AMCIS, Switzerland




Comparative separations may not be representative of all applications.

Media for One of the World's Leading PREP HPLC Columns

Luna® high surface area (400 m²/g) silica packing materials provide optimized parameters specifically designed for the purification of small molecules and peptides. This media allows high loading with excellent lifetimes.

Optimized loading parameters include:

- Silica smoothness for stable packed beds
- Optimum pore size/distribution provide outstanding performance
- High pore volume offers increased surface area
- Fine tuned bonding density for excellent reproducibility
- Greater loading capacity with an extended pH range of 1.5 to 10.0*



We use the Phenomenex Luna HPLC as our standard purification media to purify our customer's peptides. In addition to the excellent loadability and selectivity of the media itself, the Phenomenex PREP Team supports their entire line of products very effectively.

Major Biotech Manufacturer, USA







Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m ² /g)	Carbon Load (%)	pH Range	Applications					Type of Compounds				Loading
						Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic	Bases	
Luna C18(2)	3, 5, 10, 10-PREP, 15	100	400	17.5	1.5-10.0*	●	●				●		●	●	●
Luna C8(2)	3, 5, 10, 10-PREP, 15	100	400	13.5	1.5-10.0*	●	●				●		●	●	●
Luna C5	5, 10	100	440	12.5	1.5-10.0*	●	●				●		●	●	●
Luna Phenyl-Hexyl	3, 5, 10, 10-PREP, 15	100	400	17.5	1.5-10.0*	●	●				●	●	●	●	●
Luna Silica(2)	3, 5, 10, 10-PREP, 15	100	400	-	-	●					●	●		●	●
Luna CN	3, 5, 10	100	400	7	1.5-7.0	●					●	●		●	●
Luna NH ₂	3, 5, 10	100	400	9.5	1.5-11.0	●					●	●		●	●
Luna SCX	5, 10	100	400		2-7.0	●								●	●
Luna HILIC	3, 5	200	200	5.7	1.5-8.0	●	●				●	●		●	●
Luna PFP(2)	5, 10	100	400	11.5	1.5-9.0	●						●		●	●

*pH range is 1.5-10 under isocratic conditions and 1.5-8.5 under gradient conditions.

Simple Scale-Up

Axia™ column technology provides the same high efficiency chromatographic performance for preparative scale columns (21.2, 30, and 50 mm ID) as obtained in 4.6 mm ID analytical columns. This improvement in preparative column performance across

all lengths and internal diameters makes it easier to select the appropriate column size to achieve the desired purity and yield without having to compromise on performance.



Columns: Luna 5 µm C18(2)
Dimensions: As Noted
Mobile Phase: A: 0.5 % TFA in Water
B: 0.5 % TFA in Acetonitrile
Gradient: A/B (95:5) to A/B (5:95) in 5 minutes


Flow Rate: As Noted
Detection: UV @ 254 nm
Temperature: Ambient
Injection: See chromatograms
Sample: Suzuki Reaction Mixture

Proven Media for Peptide Purifications

Optimal compromise between throughput, recovery, yield. Ability to perform high loading (0.74 g on column) and achieve high purity (>98 %) in a single purification run.


Preparative Purification of Bivalirudin (20 amino acid peptide also known as Angiomax)

Purification Elution Profile at 1.5 % Specific Load



Column: Luna 10 µm-PREP C8(2)
Dimensions: 250 x 21.2 mm
Part No.: 00G-4324-P0-AX
Mobile Phase: A: 100 mM Ammonium acetate pH 4.7 in Water
B: Acetonitrile
Gradient: 10 to 50 % B in 40 min; hold at 80 % B and 20 % C for 5 min; re-equilibration at 90 % A and 10 % for 10 min
Flow Rate: 21 mL/min
Temperature: Ambient
Detection: UV @ 280 nm
Injection Volume: 105 mL
Sample Concentration: 7 mg/mL in water

Purity Confirmation of Combined Fractions



**11 Combined fractions 27.8 – 29.8 min;
Recovery 80.5 % with purity ≥ 98.5 %**

Peak No.	Time (min)	Area	Area %
1	12.74	73.7	0.35
2	13.83	40.6	0.19
3	14.37	21118.7	98.53
4	15.858	200.5	0.93

Column: Luna 5 µm C18(2)
Dimensions: 250 x 4.6 mm
Mobile Phase: A: 0.1 % TFA in Water
B: 0.1 % TFA in Acetonitrile
Gradient: 20 % to 50 % B in 30 min
Flow Rate: 1 mL/min
Temperature: 25 °C
Detection: UV @ 220 nm
Injection Volume: 2 µL
Sample: Combined Fractions

Media for Biomolecules

The Jupiter HPLC column portfolio, including Jupiter 300 and Jupiter Proteo, offers optimized reversed phase solutions for peptide and protein purification. Identify, purify, and analyze almost any protein with Jupiter columns.




Jupiter Proteo 90Å

- For separation of proteins and peptides < 10,000 MW
- C12 bonded onto an ultra-high surface area (475 m²/g) silica for increased peak capacity and resolution of peptide separations
- Direct scale-up from analytical to preparative and bulk materials

Resolve Peptides with Similar Hydrophobicity

Jupiter Proteo is able to fully resolve peptides that differ in hydrophobicity by one methyl group.



Columns: Phenomenex Jupiter 4µm Proteo 90 Å
Agilent Technologies ZORBAX® 5µm SB-C18 300 Å
Alltech Associates, Inc. Vydac® 5µm MS54 300 Å
Alltech Associates, Inc. Vydac 5µm TP54 300 Å

Dimensions: 250 x 4.6 mm

Mobile Phase: A: 0.1 % TFA in Water
B: 0.085 % TFA in Acetonitrile

Gradient: A/B (95:5) to A/B (55:45) in 20 minutes

Flow Rate: 1 mL/min

Temperature: 22 °C

Detection: UV @ 214 nm


Sample: 1. NH₂-Arg-Gly-Gly-Ala-Gly-Gly-Leu-Gly-Leu-Gly-Lys-Amide
2. Ac-Arg-Gly-Gly-Gly-Gly-Gly-Leu-Gly-Leu-Gly-Lys-Amide
3. Ac-Arg-Gly-Ala-Gly-Gly-Gly-Leu-Gly-Leu-Gly-Lys-Amide
4. Ac-Arg-Gly-Val-Gly-Gly-Gly-Leu-Gly-Leu-Gly-Lys-Amide
5. Ac-Arg-Gly-Val-Gly-Gly-Gly-Leu-Gly-Leu-Gly-Lys-Amide

Jupiter 300Å

- For separation of proteins > 10,000 MW
- Available with C18, and C4 bonded phases
- 1.5 – 10 pH stability for method ruggedness and easy protein removal
- Direct scale up to preparative and bulk materials

Compare PEGylated vs. Native Forms of Proteins

Reversed phase separation of PEGylated and native proteins on a Jupiter 300 C4 column. Note the good resolution of multiple PEGylated forms for all proteins tested.



Columns: Jupiter 300 5µm C4 300 Å
Dimensions: 150 x 4.6 mm
Part No.: 00F-4167-E0
Mobile Phase: A: 2 % Acetonitrile / 0.1 % TFA in Water
B: 70 % Acetonitrile / 20 % IPA / 0.08 % TFA in Water
Gradient: A/B (85:15) to A/B (30:70) in 25 min
Flow Rate: 1 mL/min
Temperature: 45 °C
Detection: UV @ 214 nm
Sample: PEGylated and Native Proteins

“ We purchased the Jupiter 300 C18 300 Å column a few months ago and have been quite impressed with its performance. The Jupiter 300 column provides better separation of the proteins. As for reproducibility, the control profiles have not changed since day one of its use.

Major Biotech Company, Europe ”

Key: ● Best Suited ○ Very Good

Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m ² /g)	Carbon Load (%)	pH Range	Applications					Type of Compounds				Loading
						Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic	Bases	
Jupiter C18	5, 10, 15	300	170	13.3	1.5-10.0			●			○		●	○	●
Jupiter C4	5, 10, 15	300	170	5	1.5-10.0			●			○		●	○	●
Jupiter Proteo	4, 10	90	475	15	1.5-10.0	●	●				○		●	○	●

Comparative separations may not be representative of all applications.

Purification of Synthetic Oligonucleotides




Clarity Oligo-RP™

Unique media specifically designed for reversed phase purification of oligonucleotides with balanced hydrophobicity and polar selectivity. The media is based on composite particle TWIN™ technology that provides improved selectivity and efficiency for oligonucleotides when compared to competing hybrid, polymer, and silica media.

RP-HPLC Preparative Purification

- Easily separate N-1 failure sequences from target oligo with >90 % purities
- Purify oligos up to 60 nt in length
- Trityl-off purification of DNA, RNA, thioates, and modified/ labeled oligonucleotides
- 3 µm, 5 µm, 10 µm particles for seamless scaling

Preparative 20nt DNA Oligo-RP Purification



Clarity Oligo-WAX™


Clarity Oligo-WAX is a crosslinked weak anion exchanger media designed for successful ion-exchange purification of synthetic DNA/RNA. Oligo-WAX is an advantageous combination of purity, capacity, mechanical strength, cost, and efficiency.

- Excellent efficiency column results in > 90 % purities due to good fractionation of closely eluting compounds
- High loading capacity due to very high density ligand
- Increase productivity by running at higher flow rates and pressures

Purify Failure Sequences and Contaminants from Target Sequence

Ion-exchange is an excellent separation mode for purifying contaminants and failure sequences from target sequences. Clarity Oligo-WAX, due to its increased efficiency compared to other ion-exchange columns, has the ability to recognize minute charge differences in nucleotide sequences such as failure sequences or base substitutions.

DNA Purification of N-1 Sequence from Target N Sequence



We have used the Axia prep columns and not had problems with them. I have never had to adjust for retention gaps. This speaks directly to the quality of Phenomenex's phases and the quality of their PREP columns.

-Major Biotech Company, USA



Key: ● Best Suited ○ Very Good

	Applications					Type of Compounds				Loading
	Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic	Bases	
Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m ² /g)	Carbon Load (%)	pH Range					Available Surface Area
Oligo-RP	3, 5, 10	110	375	14	1.0-12.0	●	○	●	○	●
Oligo-Wax	10	360	150	-	1.0-11.0	●	●			○


Complete Chiral Solutions



Achieving optimal chiral separation is easier than ever with five unique Lux polysaccharide stationary phases to screen. Choose a phase, then transfer the method to lab scale, process, pilot, and commercial scale.

Lux chiral preparative columns simplify the separation process:

- Unique and traditional phases that increase the success rate of the chiral screen
- Consistent particle size distribution so performance is maintained
- Mechanically strong media for increased stability
- Available in multiple particle sizes for direct scale up (3 μ m and 5 μ m packed columns for screening and small scale purifications; 10 μ m and 20 μ m bulk media for process scale purifications)



Cellulose-O-CONH

Lux Cellulose-1

Cellulose tris(3,5-dimethylphenylcarbamate)
Guaranteed Alternative to
 CHIRALCEL® OD®, OD-H®, OD-3, OD-RH, and OD-3R



Cellulose-O

Lux Cellulose-3

Cellulose tris(4-methylbenzoate)
Guaranteed Alternative to
 CHIRALCEL OJ®, OJ-H®, OJ-3, OJ-RH, and OJ-3R



Amylose-O-CONH

Lux Amylose-2


Amylose tris(5-chloro-2-methylphenylcarbamate)
Guaranteed Alternative to
 CHIRALPAK® AY®, AY-H®, AY-3, AY-RH, and AY-3R

Resolve Your Enantiomers with Five Unique Phases

The Lux family of bulk cellulose chiral selectors provides a variety of complementary selectivities.

Screen for the most effective chiral separation under the following conditions:


- Reversed Phase
- Polar Organic
- Normal Phase
- Supercritical Fluid Chromatography (SFC)



Cellulose-O-CONH

Lux Cellulose-2

Cellulose tris(3-chloro-4-methylphenylcarbamate)
Guaranteed Alternative to
 CHIRALCEL OZ, OZ-H®, OZ-3, OZ-RH, and OZ-3R



Cellulose-O-CONH

Lux Cellulose-4

Cellulose tris(4-chloro-3-methylphenylcarbamate)
Guaranteed Alternative to
 CHIRALCEL OX-H, OX-3, OX-RH, and OX-3R



Free Screening Services

see page 28

Key: ● Best Suited ○ Very Good

Packing Material	Particle Size (μm)	Pore Size (Å)	Surface Area (m²/g)	Carbon Load (%)	pH Range	Chiral Applications									
						Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic	Bases	Surface Area
Lux Cellulose-1	5, 10, 20	1,000	-	-	2-9.0	●			●		○	○	○	○	●
Lux Cellulose-2	5, 10, 20	1,000	-	-	2-9.0	●			●		○	○	○	○	●
Lux Cellulose-3	5, 20	1,000	-	-	2-9.0	●			●		○	○	○	○	●
Lux Cellulose-4	5, 20	1,000	-	-	2-9.0	●			●		○	○	○	○	●
Lux Amylose-2	5	1,000	-	-	2-9.0	●			●		○	○	○	○	●

Column Screening for Optimal Chiral Resolution


Being able to utilize differences in selectivity in each of the five Lux® columns can help develop methods more efficiently by offering broad and complementary chiral recognition abilities.

In the example below, a simple screen determined which column gave the best separation.


Etozolin

Based on a five phase screen under reversed phase conditions, the optimal chiral stationary phase for resolving Etozolin is Lux Cellulose-3.


Lux 5 µm Cellulose-3




Conditions for all columns:
Column: As noted
Dimension: 250 x 4.6 mm
Mobile Phase: Acetonitrile /
 20 mM Ammonium
 bicarbonate with 0.1 % Diethylamine (60:40)
Flow Rate: 1 mL/min
Temperature: Ambient
Detection: UV @ 220 nm




Lux 5 µm Amylose-2




Lux 5 µm Cellulose-2



Lux 5 µm Cellulose-1




Lux 3 µm Cellulose-4




Innovative chiral selector will succeed where others fail.


Lux 5 µm Cellulose-4




Conditions for all columns:
Dimensions: 250 x 4.6 mm
Mobile Phase: 0.1 % Diethylamine in Hexane /
 0.1 % Diethylamine in Isopropanol (90:10)
Flow Rate: 1 mL/min
Detection: UV @ 220 nm
Temperature: Ambient



CHIRALCEL® 5 µm OD-H®



CHIRALPAK® 5 µm AD®-H




Comparative separations may not be representative of all applications.

Load More with an Increase in Column Length

Axia™ column technology allows separation to scale up directly based on column length. With the 100 mm length column a 32 mg/load separation was achieved and an increased sample


load of 80 mg/load was achieved on the longer 250 mm length column. As expected when increasing the load, the peak width and tailing increased but there was no loss of resolution.




Conditions for all columns:

Columns: Lux® 5 µm Cellulose-1
Dimensions: as noted
Mobile Phase: Methanol / Isopropanol (90:10)
Flow Rate: as noted
Detection: as noted
Sample: Dissolved in mobile phase as noted

Dimensions: 100 x 4.6 mm
Flow Rate: 1 mL/min
Detection: UV @ 220 nm
Sample: 5 µg in 2 µL




No resolution loss with increased sample load



Dimensions: 100 x 21.2 mm
Flow Rate: 20 mL/min
Detection: UV @ 220 nm and 254 nm
Sample: 32 mg in 640 µL

2.5x Load Increase



Dimensions: 250 x 21.2 mm
Flow Rate: 20 mL/min
Detection: UV @ 220 nm and 254 nm
Sample: 80 mg in 1600 µL



Lux Axia preparative column are wonderful! I regularly use Lux chiral stationary phase Cellulose-2 and Cellulose-4 and less frequently, the Lux Amylose-2. In our community of chiral analysis/purification scientists, there are some who use the CC4 column instead of the *equivalent* Lux Cellulose-4. On several occasions we've seen separation and good peak shape on the Lux Cellulose-4 that was completely missing from the CC4. Customer support and delivery times are always within a few days.


Julia G. Christie
 GlaxoSmithKline, USA






Analytical and Axia™ packed columns have been extensively tested on various SFC systems and all column ID's and lengths are SFC compatible.


Baseline Separation of Enantiomers



Dimensions: 250 x 4.6 mm
Flow Rate: 2.5 mL/min
Detection: UV @ 220 nm
Load: 300 µg 10 µL



5x Load Increase



Dimensions: 250 x 4.6 mm
Flow Rate: 2.5 mL/min
Detection: UV @ 254 nm
Load: 1.5 mg in 50 µL

Conditions for all columns:

Columns: Lux® 5 µm Cellulose-1


Mobile Phase: Methanol with 0.1 % DEA/
Carbon Dioxide (25:75)

Column Temperature: 35 °C

Polarimeter: ALP-PDR-Chiral

Sample: Terfenadine with ethanol
dissolution solvent

70x Load Increase



High loading capacity media along with stacking injections allow for increased yields and productivity

Closer stacked injections can not be used due to the impurities eluting after the major enantiomers

Dimensions: 250 x 21.2 mm
Flow Rate: 50 mL/min
Detection: UV @ 220 nm
Load: 105 mg in 3.5 mL

Tip:

For more information on this warfarin application, request technote:
TN-9002: Scaling from Analytical to Preparative Chiral Chromatography While Balancing Purity, Yield, and Throughput under HPLC and SFC Conditions

Want to Extend the Lifetime of Your Axia™ Column?



Use the SecurityGuard PREP column protection system


SecurityGuard:

- Extends prep column lifetime by as much as 5x
- Protects column from samples that precipitate out of solution
- Protects column from contaminants
- Stable and leak-free up to 60 mL/min

The SecurityGuard PREP system was designed to effectively (and inexpensively) protect your valuable prep columns from the damaging effects of mobile phase and sample chemical contaminants and particulates, without altering your chromatographic results.


Forced Degradation Study

Injection 1: Axia Packed column with SecurityGuard PREP cartridge





Conditions

Column: Luna® 10 µm C18(2) Axia Packed
Dimension: 50 x 21.2 mm
Part No.: 00B-4253-P0-AX
Mobile Phase: A: 0.1 % TFA in Water
 B: 0.1 % TFA in Water/Acetonitrile (25:75)
Gradient: Linear 93:7 (A/B) to 100% B over 5 minutes
Injection Volume: 420 µL
Flow Rate: 60 mL/min
Temperature: Ambient
Detection: UV @ 270 nm
Sample: 1. Nadolol
 2. Metoprolol
 3. Propranolol




Injection 240: Axia Packed column with SecurityGuard PREP cartridge



Guard Cartridge

Injection 241: Axia Packed column after removing SecurityGuard column protection system



guarantee

If the SecurityGuard PREP cartridge protection system does not perform as well or better than your current guard cartridge system of similar phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.

SecurityGuard™ PREP System

Simply match the appropriate SecurityGuard PREP cartridge and cartridge holder to your Axia column's chemistry and ID to greatly extend your column's lifetime

Cartridges and Holders

Step 1: Choose column ID

Step 2: Match column phase

Ordering Information

Material	Description	pH Stability	/ea	/ea
Cartridges for General Purpose/Pharmaceutical				
Kinetex C18	(C18, XB-C18 Core-Shell Technology)	1.5 - 10	AJ0-9145	AJ0-9204
Kinetex C8	(C8 Core-Shell Technology)	1.5 - 10	AJ0-9205	AJ0-9217
Kinetex Phenyl-Hexyl	(Phenyl-Hexyl Core-Shell Technology)	1.5 - 10	AJ0-9147	AJ0-9216
C18	(ODS, Octadecyl)	1.5 - 10	AJ0-7839	AJ0-8301
C12	(Dodecyl)	1.5 - 10	AJ0-7842	AJ0-8304
C8	(MOS, Octyl)	1.5 - 10	AJ0-7840	AJ0-8302
Silica	—	—	AJ0-7229	AJ0-8312
NH ₂	(Amino, Aminopropyl)	1.5 - 11	AJ0-8162	AJ0-8309
CN	(Cyano, Cyanopropyl)	2 - 7.5	AJ0-8220	AJ0-8311
Phenyl	(Phenylhexyl)	1.5 - 10	AJ0-7841	AJ0-8303
PFP(2)	Pentafluorophenyl	1.5 - 8	AJ0-8377	AJ0-8378
SCX	(SA, Strong Cation Exchanger)	2.5 - 7.5	AJ0-8595	AJ0-8596
RP-1	(Reversed Phase - Polymer)	0 - 14	AJ0-8358	—
Polar-RP	(Ether-linked Phenyl)	1.5 - 7	AJ0-7845	AJ0-8307
Fusion-RP	(C18 Polar Embedded)	1.5 - 10	AJ0-7844	AJ0-8306
AQ C18	(Polar Endcapped C18)	1.5 - 7.5	AJ0-7843	AJ0-8305
Gemini® NX-C18	(C18 TWIN NX™ Technology)	1 - 12	AJ0-8370	AJ0-8371
Gemini C18	(C18 TWIN™ Technology)	1 - 12	AJ0-7846	AJ0-8308
Oligo-RP™	(C18 TWIN Technology)	1 - 12	AJ0-8210	AJ0-8310
Oligo-WAX™	(WA, Weak Anion Exchanger)	1.5 - 11	AJ0-8639	AJ0-8420
Cartridges for Protein and Polypeptide Reversed Phase			/ea	/ea
<i>For use with silica columns for separation of proteins & peptides, such as Jupiter® (Phenomenex); Vydac® 218TP, 214TP (Alltech Associates, Inc.); Nucleosil® 300 Å C18, C4 (Macherey-Nagel); HYPERSIL® 300 Å (Thermo-Hypersil-Keystone), and other widepore or 300 Å brands.</i>				
Widepore C18	(ODS, Octadecyl)	1.5 - 10	AJ0-7230	AJ0-8313
Widepore C4	(Butyl)	1.5 - 10	AJ0-7231	AJ0-8314
Cartridges for Chiral			/ea	/ea
<i>For use with chiral columns, such as Lux® Cellulose-1, -2, -3, -4, & Amylose-2 (Phenomenex); CHIRALCEL® OD-H®, CHIRALCEL® OJ-H® & CHIRALPAK® AD®-H (DAICEL Corporation)</i>				
Lux Cellulose-1	Cellulose tris (3, 5-dimethylphenylcarbamate)	2 - 9	AJ0-8405	AJ0-8406
Lux Cellulose-2	Cellulose tris (3-chloro-4-methylphenylcarbamate)	2 - 9	AJ0-8400	AJ0-8401
Lux Cellulose-3	Cellulose tris (4-methylbenzoate)	2 - 9	AJ0-8624	AJ0-8625
Lux Cellulose-4	Cellulose tris (4-chloro-3-methylphenylcarbamate)	2 - 9	AJ0-8629	AJ0-8630
Lux Amylose-2	Amylose tris (5-chloro-2-methylphenylcarbamate)	2 - 9	AJ0-8473	AJ0-8474
Prep Guard Cartridge Holders (one-time purchase only)			/kit	/kit
			AJ0-8223	AJ0-8277

HPLC Holder Kit; includes column coupler

SFC Holder Kit; includes column coupler



Choose your HPLC or SFC Column Internal Diameter, ID (mm):

18.0 - 29.0

30.0 - 49.0

Use Cartridges (mm):

15 x 21.2

15 x 30.0



I used about 10-12 SecurityGuard Cartridges and received extremely long Axia column lifetimes (8000 injections). During this time the columns showed extremely good efficiency with no significant changes concerning the backpressure.

Claudia Oschwald
Bayer, Germany



A New Era of Technical Support Services

Let Us Do the Work for You

PhenoLogix, our in-house application support lab, saves you time and money by screening multiple scout columns and solvent strategies for new purification methods or revalidating your current methods. We work together to make you successful by minimizing your process purification development time and optimizing your purification method.

Chiral Screening

- Normal Phase
- Reversed Phase
- Polar Organic
- SFC

Method Optimization Services

- Fast Turnaround
- Easy Method Transfer
- Continued Support

Preparative and Process Scale-Up

- Media Screening
- Small Scale Purification
- DAC Packing Assistance



Your Method, Our Scientists

Quality Products, Advanced Performance, Complete Support

For more information or to begin a project today, please contact your local Phenomenex representative or email us at phenologix@phenomenex.com

You can also visit us online:

www.phenomenex.com/phenologix



“

Our scientists at American Peptide have taken advantage of Phenomenex's column packing services, application development, and project-specific consultation services for some of our most challenging separations.

”

American Peptide Company, USA



Ordering Information

Achiral Phases

Gemini®

Phase	Length	ID	Part No.	Price
5 µm				
NX-C18	50	21.2	00B-4454-P0-AX	
	50	30	00B-4454-U0-AX	
	75	30	00C-4454-U0-AX	
	100	21.2	00D-4454-P0-AX	
	100	30	00D-4454-U0-AX	
	150	21.2	00F-4454-P0-AX	
	150	30	00F-4454-U0-AX	
	250	21.2	00G-4454-P0-AX	
	250	30	00G-4454-U0-AX	
C18	50	21.2	00B-4435-P0-AX	
	50	30	00B-4435-U0-AX	
	75	30	00C-4435-U0-AX	
	100	21.2	00D-4435-P0-AX	
	100	30	00D-4435-U0-AX	
	150	21.2	00F-4435-P0-AX	
	150	30	00F-4435-U0-AX	
	250	21.2	00G-4435-P0-AX	
	250	30	00G-4435-U0-AX	
C6-Phenyl	75	30	00C-4444-U0-AX	
	100	21.2	00D-4444-P0-AX	
	100	30	00D-4444-U0-AX	
	150	21.2	00F-4444-P0-AX	
	250	21.2	00G-4444-P0-AX	
10 µm				
NX-C18	50	21.2	00B-4455-P0-AX	
	50	50	00B-4455-V0-AX	
	100	21.2	00D-4455-P0-AX	
	100	30	00D-4455-U0-AX	
	100	50	00D-4455-V0-AX	
	150	21.2	00F-4455-P0-AX	
	150	30	00F-4455-U0-AX	
	150	50	00F-4455-V0-AX	
	250	21.2	00G-4455-P0-AX	
	250	30	00G-4455-U0-AX	
	250	50	00G-4455-V0-AX	
C18	50	21.2	00B-4436-P0-AX	
	100	21.2	00D-4436-P0-AX	
	100	30	00D-4436-U0-AX	
	100	50	00D-4436-V0-AX	
	150	21.2	00F-4436-P0-AX	
	150	30	00F-4436-U0-AX	
	150	50	00F-4436-V0-AX	
	250	21.2	00G-4436-P0-AX	
	250	30	00G-4436-U0-AX	
	250	50	00G-4436-V0-AX	

Jupiter®

Phase	Length	ID	Part No.	Price
10 µm				
Proteo	250	21.2	00G-4397-P0-AX	
	250	30	00G-4397-U0-AX	
C18 300 Å	100	21.2	00D-4055-P0-AX	
	250	30	00G-4055-U0-AX	
C4 300 Å	100	21.2	00D-4168-P0-AX	
	150	21.2	00F-4168-P0-AX	
	250	21.2	00G-4168-P0-AX	

Luna®

Phase	Length	ID	Part No.	Price
5 µm				
C18(2)	50	21.2	00B-4252-P0-AX	
	50	30	00B-4252-U0-AX	
	100	21.2	00D-4252-P0-AX	
	100	30	00D-4252-U0-AX	
	150	21.2	00F-4252-P0-AX	
	250	21.2	00G-4252-P0-AX	
	250	30	00G-4252-U0-AX	
C8(2)	50	21.2	00B-4249-P0-AX	
	50	30	00B-4249-U0-AX	
	100	21.2	00D-4249-P0-AX	
	100	30	00D-4249-U0-AX	
C5	150	21.2	00F-4043-P0-AX	
	250	21.2	00G-4043-P0-AX	
CN	50	21.2	00B-4255-P0-AX	
	100	30	00D-4255-P0-AX	
	150	21.2	00F-4255-P0-AX	
	250	21.2	00G-4255-P0-AX	
Phenyl-Hexyl	50	21.2	00B-4257-P0-AX	
	100	21.2	00D-4257-P0-AX	
	100	30	00D-4257-U0-AX	
	150	21.2	00F-4257-P0-AX	
NH₂	100	21.2	00D-4378-P0-AX	
	150	21.2	00F-4378-P0-AX	
	250	21.2	00G-4378-P0-AX	
HILIC	50	21.2	00B-4450-P0-AX	
	100	21.2	00D-4450-P0-AX	
	150	21.2	00F-4450-P0-AX	
	250	21.2	00G-4450-P0-AX	
	250	30	00G-4450-U0-AX	
PFP(2)	50	21.2	00B-4448-P0-AX	
	100	21.2	00D-4448-P0-AX	
	100	30	00D-4448-U0-AX	
	150	21.2	00F-4448-P0-AX	
	250	21.2	00G-4448-P0-AX	
	250	30	00G-4448-U0-AX	
	250	50	00G-4448-V0-AX	
Silica (2)	50	21.2	00B-4274-P0-AX	
	100	21.2	00D-4274-P0-AX	
	150	21.2	00F-4274-P0-AX	
	250	21.2	00G-4274-P0-AX	
	250	30	00G-4274-U0-AX	
10 µm				
C18(2)	50	21.2	00B-4253-P0-AX	
	50	30	00B-4253-U0-AX	
	100	21.2	00D-4253-P0-AX	
	250	21.2	00G-4253-P0-AX	
	250	30	00G-4253-U0-AX	
	250	50	00G-4253-V0-AX	
	250	50	00G-4253-V0-AX	
C8(2)	50	21.2	00B-4250-P0-AX	
	250	21.2	00G-4250-P0-AX	
	250	50	00G-4250-V0-AX	
C5	100	21.2	00D-4092-P0-AX	
	250	21.2	00G-4092-P0-AX	
	250	50	00G-4092-V0-AX	
CN	250	21.2	00G-4300-P0-AX	
Phenyl-Hexyl	250	21.2	00G-4285-P0-AX	
	250	30	00G-4285-U0-AX	
NH₂	250	21.2	00G-4379-P0-AX	
Silica (2)	250	21.2	00G-4091-P0-AX	
	250	30	00G-4091-U0-AX	
	250	50	00G-4091-V0-AX	
15 µm				
C18(2)	250	21.2	00G-4273-P0-AX	
	250	30	00G-4273-U0-AX	
	250	50	00G-4273-V0-AX	
C8(2)	250	21.2	00G-4272-P0-AX	
Silica (2)	250	21.2	00G-4271-P0-AX	
	250	30	00G-4271-U0-AX	

Achiral Phases (continued)

Kinetex®

Phase	Length	ID	Part No.	Price
5 µm				
XB-C18	50	21.2	00B-4605-P0-AX	
	50	30	00B-4605-U0-AX	
	100	21.2	00D-4605-P0-AX	
	100	30	00D-4605-U0-AX	
	150	21.2	00F-4605-P0-AX	
	150	30	00F-4605-U0-AX	
	250	21.2	00G-4605-P0-AX	
	250	30	00G-4605-U0-AX	
C18	50	21.2	00B-4601-P0-AX	
	50	30	00B-4601-U0-AX	
	100	21.2	00D-4601-P0-AX	
	100	30	00D-4601-U0-AX	
	150	21.2	00F-4601-P0-AX	
	150	30	00F-4601-U0-AX	
	250	21.2	00G-4601-P0-AX	
	250	30	00G-4601-U0-AX	
C8	50	21.2	00B-4608-P0-AX	
	50	30	00B-4608-U0-AX	
	100	21.2	00D-4608-P0-AX	
	100	30	00D-4608-U0-AX	
	150	21.2	00F-4608-P0-AX	
	150	30	00F-4608-U0-AX	
	250	21.2	00G-4608-P0-AX	
	250	30	00G-4608-U0-AX	
Phenyl-Hexyl	50	21.2	00B-4603-P0-AX	
	50	30	00B-4603-U0-AX	
	100	21.2	00D-4603-P0-AX	
	100	30	00D-4603-U0-AX	
	150	21.2	00F-4603-P0-AX	
	150	30	00F-4603-U0-AX	
	250	21.2	00G-4603-P0-AX	
	250	30	00G-4603-U0-AX	

Synergi™

Phase	Length	ID	Part No.	Price
4 µm				
Fusion-RP	50	21.2	00B-4424-P0-AX	
	100	21.2	00D-4424-P0-AX	
	100	30	00D-4424-U0-AX	
	150	21.2	00F-4424-P0-AX	
	250	21.2	00G-4424-P0-AX	
Hydro-RP	50	21.2	00B-4375-P0-AX	
	75	30	00C-4375-U0-AX	
	100	21.2	00D-4375-P0-AX	
	100	30	00D-4375-U0-AX	
	150	21.2	00F-4375-P0-AX	
	250	21.2	00G-4375-P0-AX	
	250	30	00G-4375-U0-AX	
	250	30	00G-4375-U0-AX	
Max-RP	50	21.2	00B-4337-P0-AX	
	100	21.2	00D-4337-P0-AX	
	100	30	00D-4337-U0-AX	
	150	21.2	00F-4337-P0-AX	
	250	21.2	00G-4337-P0-AX	
Polar-RP	50	21.2	00B-4336-P0-AX	
	50	30	00B-4336-U0-AX	
	75	30	00C-4336-U0-AX	
	100	21.2	00D-4336-P0-AX	
	100	30	00D-4336-U0-AX	
	150	21.2	00F-4336-P0-AX	
	250	21.2	00G-4336-P0-AX	
	250	30	00G-4336-U0-AX	
10 µm				
Fusion-RP	250	21.2	00G-4425-P0-AX	
Hydro-RP	250	21.2	00G-4376-P0-AX	
	250	30	00G-4376-U0-AX	
Max-RP	100	21.2	00D-4350-P0-AX	
	250	21.2	00G-4350-P0-AX	
Polar-RP	100	30	00G-4350-U0-AX	
	250	21.2	00G-4351-P0-AX	

For additional sizes not displayed, please contact your Phenomenex technical consultant or local distributor.

Clarity®

Phase	Length	ID	Part No.	Price
5 µm				
Oligo-RP	100	21.2	00D-4442-P0-AX	
	250	21.2	00G-4442-P0-AX	
10 µm				
Oligo-RP	150	21.2	00F-4445-P0-AX	
	150	30	00F-4445-U0-AX	
	250	21.2	00G-4445-P0-AX	
Oligo-WAX	150	21.2	00F-4451-P0-AX	
	150	30	00F-4451-U0-AX	
	250	21.2	00G-4451-P0-AX	

Chiral Phases

Lux®

Phase	Length	ID	Part No.	Price
5 µm				
Amylose-2	150	21.2	00F-4472-P0-AX	
	250	21.2	00G-4472-P0-AX	
	250	30	00G-4472-U0-AX	
Cellulose-1	150	21.2	00F-4459-P0-AX	
	250	21.2	00G-4459-P0-AX	
	250	30	00G-4459-U0-AX	
	250	50	00G-4459-V0-AX	
Cellulose-2	150	21.2	00F-4457-P0-AX	
	250	21.2	00G-4457-P0-AX	
	250	30	00G-4457-U0-AX	
	250	50	00G-4457-V0-AX	
Cellulose-3	150	21.2	00F-4493-P0-AX	
	250	21.2	00G-4493-P0-AX	
	250	30	00G-4493-U0-AX	
	250	50	00G-4493-V0-AX	
Cellulose-4	150	21.2	00F-4491-P0-AX	
	250	21.2	00G-4491-P0-AX	
	250	30	00G-4491-U0-AX	
	250	50	00G-4491-V0-AX	

Tip:

Want to stay on top of the latest information and product releases in the Prep world?

Sign up for our Prep Newsletter.


It's Free and Informative!

www.Phenomenex.com/goPREP

guarantee

If Axia™ packed columns do not provide at least an equivalent separation as compared to a competing preparative column of the same particle size, same phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND. Only applies to 21.2mm ID columns.

The Ultimate Pre-packed Preparative HPLC/SFC Column



Australia

t: 02-9428-6444
f: 02-9428-6445
auinfo@phenomenex.com

Austria

t: 01-319-1301
f: 01-319-1300
anfrage@phenomenex.com

Belgium

t: 02 503 4015 (French)
t: 02 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: (800) 543-3681
f: (310) 328-7768
info@phenomenex.com

Denmark

t: 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: 09 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: 01 30 09 21 10
f: 01 30 09 21 11
franceinfo@phenomenex.com

Germany

t: 06021-58830-0
f: 06021-58830-11
anfrage@phenomenex.com

India

t: 040-3012 2400
f: 040-3012 2411
indiainfo@phenomenex.com

Ireland

t: 01 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: 051 6327511
f: 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 001-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: 030-2418700
f: 030-2383749
nlinfo@phenomenex.com

New Zealand

t: 09-4780951
f: 09-4780952
nzinfo@phenomenex.com

Norway

t: 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: (800) 541-HPLC
f: (310) 328-7768
info@phenomenex.com

Sweden

t: 08 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: 01625-501367
f: 01625-501796
ukinfo@phenomenex.com

United States

t: (310) 212-0555
f: (310) 328-7768
info@phenomenex.com

All other countries: Corporate Office USA



t: (310) 212-0555
f: (310) 328-7768
info@phenomenex.com



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at www.Phenomenex.com/TermsAndConditions

Trademarks

Clarity, Gemini, Kinetex, and Luna are registered trademarks, and Oligo-RP, Oligo-WAX, SecurityGuard, Synergi, TWIN-NX, TWIN Technology, and Two-In-One Technology are trademarks of Phenomenex. UPLC, XBridge, Symmetry, XTerra, and Waters are registered trademarks and OBD and SymmetryShield are trademarks of Waters Corporation. ZORBAX is a registered trademark of Agilent Technologies, Inc. Vydac is a registered trademark of Alltech Associates, Inc. CHIRALCEL and CHIRALPAK are registered trademarks of DAICEL Corporation.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is not affiliated with Waters Corporation, Agilent Technologies, Inc., Alltech Associates, Inc. or DAICEL Corporation. Columns used for comparison studies were manufactured by and purchased from Waters Corp., Agilent Technologies, Inc., Alltech Associates, Inc. and DAICEL Corporation.

The opinions stated herein are solely those of the speaker and not necessarily those of any company or organization.

Axia is patented by Phenomenex. U.S. Patent No. 7,674,383

Gemini is patented by Phenomenex. U.S. Patent No. 7,563,367

© 2014 Phenomenex, Inc. All rights reserved.