

Ph. Eur. Monograph 2918: Test for Everolimus Impurity A on Kinetex® 2.7 µm C18

Zeshan Aqeel¹ and Heiko Behr²

¹Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 1430, USA

²Phenomenex Ltd. Deutschland, Zeppelinstr. 5, 63741 Aschaffenburg, Germany

Overview

Everolimus is an immunosuppressant used to prevent rejection of organ transplants and in the treatment of renal cell cancer and other tumors.

In this application note we demonstrate the separation of Everolimus from Sirolimus (Impurity A) according to Ph. Eur. Monograph 2918 for the determination of Impurity A.

The chromatogram obtained for reference solution (a) was used to confirm the retention time for Impurity A. Three batches of the Kinetex 2.7 µm C18 column were used to demonstrate excellent batch-to-batch reproducibility, with each column meeting the system suitability requirement for the separation between Impurity A and Everolimus, as determined by the peak-to-valley ratio (minimum 2.0); where H_p = height above the baseline of the peak due to Impurity A and H_v = height above the baseline of the lowest point of the curve separating this peak from the peak due to Everolimus.

All reference solutions were prepared as indicated in Ph. Eur. monograph 2918 for Everolimus. (NOTE: Plastic labware was not used for solution preparation, as recommended, to minimize the formation of the everolimus tautomer; therefore, the observed peak for the tautomer was significantly larger than expected). The Everolimus (catalog no. Y0002048) certified reference standard (CRS) was purchased from the European Directorate for the Quality of Medicines & HealthCare (EDQM) Council of Europe; and the Sirolimus (Impurity A) reference standard (catalog no. 1612765) was obtained from US Pharmacopeia (USP).

LC-UV Conditions

Column: Kinetex 2.7 µm C18

Dimension: 150 x 2.1 mm (00F-4783-AN)

Mobile Phase: A = Formic acid/Methanol/0.7 g/L Ammonium Formate/Acetonitrile (0.05:100:450:450 v/v/v/v)

B = Formic acid/Methanol/0.96 g/L Ammonium Formate/Acetonitrile (0.05:100:330:570 v/v/v/v)

Gradient: Time (min)	% B
0	0
14	0
23	100
28	100
30	0
34	0

Flow Rate: 0.9 mL/min

Temperature: 60 °C

Detector: UV @ 278 nm

Injection: 3.5 µL

System: Waters® ACQUITY® I-Class


Sample: 1. Sirolimus (Impurity A)
2. Everolimus
3. Everolimus Tautomer

Table 1: Preparation of Test and Reference Solutions


Solution	Composition
Test Solution	Dissolve 50.0 mg of Everolimus CRS (or the substance to be examined) in acetonitrile and dilute to 25.0 mL with acetonitrile.
Reference solution (a)	Dissolve 1 mg of USP Sirolimus RS (impurity A) in 5 mL of acetonitrile. Dilute 0.8 mL of the solution to 10 mL with the test solution.




Original chromatogram published in EDQM for reference solution (a)




Kinetex® 2.7 µm C18 150 x 2.1 mm (00F-4783-AN); reference solution (a) on 3 different batches



$H_{p,r}$ (Impurity A)	$H_{v,r}$ (Everolimus)	$p/v = H_v/H_p$
17458	5668	3.08



$H_{p,r}$ (Impurity A)	$H_{v,r}$ (Everolimus)	$p/v = H_v/H_p$
17299	5751	3.01



$H_{p,r}$ (Impurity A)	$H_{v,r}$ (Everolimus)	$p/v = H_v/H_p$
17141	5470	3.13



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia t: +61 (0)2-9428-6444 auinfo@phenomenex.com	India t: +91 (0)40-3012 2400 indiainfo@phenomenex.com	Slovakia t: +420 272 017 077 sk-info@phenomenex.com
Austria t: +43 (0)1-319-1301 anfrage@phenomenex.com	Ireland t: +353 (0)1 247 5405 eireinfo@phenomenex.com	Spain t: +34 91-413-8613 espinfo@phenomenex.com
Belgium t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com	Italy t: +39 051 6327511 italiainfo@phenomenex.com	Sweden t: +46 (0)8 611 6950 nordicinfo@phenomenex.com
Canada t: +1 (800) 543-3681 info@phenomenex.com	Luxembourg t: +31 (0)30-2418700 nlinfo@phenomenex.com	Switzerland t: +41 (0)61 692 20 20 swissinfo@phenomenex.com
China t: +86 400-606-8099 cninfo@phenomenex.com	Mexico t: 01-800-844-5226 tecnicomx@phenomenex.com	Taiwan t: +886 (0) 0801-49-1246 twinfo@phenomenex.com
Czech Republic t: +420 272 017 077 cz-info@phenomenex.com	The Netherlands t: +31 (0)30-2418700 nlinfo@phenomenex.com	Thailand t: +66 (0) 2 566 0287 thaiinfo@phenomenex.com
Denmark t: +45 4824 8048 nordicinfo@phenomenex.com	New Zealand t: +64 (0)9-4780951 nzinfo@phenomenex.com	United Kingdom t: +44 (0)1625-501367 ukinfo@phenomenex.com
Finland t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com	Norway t: +47 810 02 005 nordicinfo@phenomenex.com	USA t: +1 (310) 212-0555 info@phenomenex.com
France t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com	Poland t: +48 22 104 21 72 pl-info@phenomenex.com	Ⓜ All other countries/regions Corporate Office USA t: +1 (310) 212-0555 info@phenomenex.com
Germany t: +49 (0)6021-58830-0 anfrage@phenomenex.com	Portugal t: +351 221 450 488 ptinfo@phenomenex.com	
Hong Kong t: +852 6012 8162 hkinfo@phenomenex.com	Singapore t: +65 800-852-3944 sginfo@phenomenex.com	

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex is a registered trademark and BE-HAPPY is a trademark of Phenomenex. Waters and ACQUITY are registered trademarks of Waters Technologies Corporation..

Disclaimer

Comparative separations may not be representative of all applications.

Phenomenex is in no way affiliated with Waters.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.

