

Ph. Eur. Monograph 2244: Tacrolimus Monohydrate Related Substances on Luna® 3 µm C18(2) and Luna Omega 3 µm C18

Janja Vasilčić¹, Amra Perva Uzunalić, PhD¹, Dirk Hansen, PhD², Heiko Behr, PhD², and Bryan Tackett, PhD³

¹University of Maribor, Faculty of Chemistry and Chemical Engineering Group for Separation Analysis, Smetanova 17, 2000 Maribor, Slovenia

²Phenomenex Ltd. Deutschland, Zeppelinstr. 5, 63741 Aschaffenburg, Germany

³Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

Tacrolimus is an immunosuppressive drug. The molecule is an inhibitor of calcineurin, a phosphatase regulating the interleukin-2 production.

In this application note we show the separation of Tacrolimus from its related substances following Ph. Eur. monograph 2244. The Luna 3 µm C18(2) and Luna Omega 3 µm C18 used for this study met the system suitability criteria of a $R \geq 3$ between Tacrolimus ($t_R \cong 23$ min) and impurity A ($t_R \cong 20$ min).

All reference solutions were prepared as indicated in Ph. Eur. monograph 2244 for Tacrolimus monohydrate. All solutions were prepared in glassware protected from light and allowed the test and reference solutions to stand for 3 hours at room temperature before use. The following certified reference standards (CRS) were purchased from the European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: Allee Kastner CS 30026 F - 67081 Strasbourg (France):

- Y0001926, batch 1 Tacrolimus monohydrate CRS
- Y0001933, batch 1 Tacrolimus for system suitability CRS

Table 1. Preparation of Buffer Solution, Mobile Phase Solution, and Solvent

Solution	Composition
Buffer Solution	Add 0.4 mL phosphoric acid to 500 mL water, then dilute to 1000 mL
Mobile Phase 1	1,1-Dimethylether/acetonitrile (19:81, v/v)
Solvent	Water/acetonitrile (30:70, v/v)

LC-UV Conditions

Columns: Luna 3 µm C18(2) ([00F-4251-E0](#))
Luna Omega 3 µm C18 ([00F-4784-E0](#))

Dimension: 150 x 4.6 mm

Mobile Phase: A = **Mobile Phase 1/Buffer Solution** (20:80, v/v)
B = **Mobile Phase 1/Buffer Solution** (80:20, v/v)

Gradient: Time (min)	%B
0	28
30	28
53	85
55	28 System equilibration
60	28 System equilibration

Injection: 20 µL

Temperature: 60 °C

Detector: UV @ 220 nm

System: Agilent® 1260 Infinity I

Table 2. Preparation of Test and Reference Solutions

Solution	Composition
Reference Solution (a) = Sample/Test Solution	Dissolve 30.0 mg reference substance Tacrolimus monohydrate CRS in Solvent and dilute to 10.0 mL with the same Solvent . Mix well (SaS). Allow the sample solution to stand for 3 h at room temperature before use (RSa).
Reference Solution (b)	Dissolve approx. 3.0 mg reference substance Tacrolimus for system suitability CRS (containing impurity A) in Solvent and dilute to 1.0 mL with the same Solvent . Mix well.
Reference Solution (c)	Transfer 1.0 mL of Reference Solution (a) into a 100 mL volumetric flask, fill up to volume with Solvent and mix well (RSa-1). Transfer 1.0 mL of RSa-1 into a 10 mL volumetric flask, fill up to volume with Solvent and mix well.



Figure 1. System Suitability Test using Reference Solution (b) on Luna® 3 µm C18(2)




Table 3. System Suitability Test using Reference Solution (b) on Luna 3 µm C18(2)

Inj. No.	Impurity B		Tacrolimus Compound II		Impurity A		Tacrolimus			R Tacrolimus/ Impurity A
	t _R	Area	t _R	Area	t _R	Area	t _R	Area	Tailing Factor	
1	12.592	142433	14.887	1136649	20.496	440305	23.394	12843556	1.24	3.73
2	12.596	142624	14.894	1136618	20.503	439897	23.402	12842176	1.24	3.73
3	12.607	141228	14.907	1136898	20.523	440548	23.426	12836480	1.24	3.73
4	12.612	139698	14.916	1136346	20.535	439356	23.441	12830299	1.24	3.73
5	12.618	139487	14.926	1137175	20.549	439833	23.454	12833823	1.24	3.73
6	12.631	138135	14.942	1136265	20.572	439412	23.483	12822469	1.24	3.73
Average	12.609	140601	14.912	1136659	20.530	439892	23.433	12834801	1.24	3.73
% RSD	0.11	1.27	0.14	0.03	0.14	0.11	0.14	0.06	0.08	0.06



Figure 2. Reference Solution (a) on Luna® 3 µm C18(2)




Table 4. Reference Solution (a) on Luna 3 µm C18(2)

Inj. No.	Impurity B		Tacrolimus Compound II		Tacrolimus		Tailing Factor
	t _R	Area	t _R	Area	t _R	Area	
1	12.631	121958	14.952	1238528	23.485	13837607	1.26
2	12.640	120184	14.964	1197193	23.506	13829646	1.25
3	12.638	119503	14.959	1199238	23.490	13822786	1.26
4	12.606	119101	14.922	1197651	23.430	13802206	1.26
5	12.606	118229	14.932	1196083	23.452	13799546	1.26
6	12.652	116425	14.990	1193951	23.547	13796236	1.26
Average	12.629	119233	14.953	1203774	23.485	13814671	1.26
% RSD	0.15	1.56	0.16	1.42	0.18	0.13	0.12



Figure 3. System Suitability Test using Reference Solution (b) on Luna® Omega 3 µm C18




Table 5. System Suitability Test using Reference Solution (b) on Luna Omega 3 µm C18

Inj. No.	Impurity B		Tacrolimus Compound II		Impurity A		Tacrolimus			R Tacrolimus/ Impurity A
	t _R	Area	t _R	Area	t _R	Area	t _R	Area	Tailing Factor	
1	11.953	55691	14.223	1306845	19.634	466575	22.376	13764864	1.32	3.85
2	11.952	55444	14.223	1306289	19.632	467499	22.374	13759471	1.32	3.85
3	11.957	55624	14.228	1307104	19.642	466395	22.385	13746057	1.32	3.85
4	11.967	55414	14.243	1303581	19.660	464397	22.407	13720535	1.32	3.85
5	11.973	55491	14.253	1303009	19.677	464309	22.428	13717911	1.32	3.86
6	11.983	55342	14.267	1305552	19.695	464679	22.448	13728662	1.32	3.85
Average	11.964	55501	14.240	1305397	19.657	465642	22.403	13739583	1.32	3.85
% RSD	0.10	0.24	0.13	0.13	0.13	0.29	0.13	0.15	0.13	0.09



Figure 4. Reference Solution (a) on Luna® Omega 3 µm C18




Table 6. Reference Solution (a) on Luna Omega 3 µm C18

Inj. No.	Impurity B		Tacrolimus Compound II		Tacrolimus		Tailing Factor
	t _R	Area	t _R	Area	t _R	Area	
1	11.994	102163	14.281	1380416	22.455	14784437	1.35
2	12.005	101120	14.295	1380206	22.475	14777072	1.35
3	12.030	99050	14.328	1379529	22.536	14776594	1.35
4	12.061	96605	14.368	1378994	22.602	14778160	1.35
5	12.088	94767	14.400	1378935	22.659	14768606	1.35
6	12.131	92101	14.456	1378421	22.747	14770702	1.35
Average	12.052	97634	14.355	1379417	22.579	14775928	1.35
% RSD	0.44	3.96	0.46	0.06	0.50	0.04	0.07

Conclusion

Both columns used for this study showed good separation of all Tacrolimus related compounds according to Ph. Eur. monograph 2244. The Luna 3 µm C18(2) column met the system suitability criteria with a resolution of 3.73 between Impurity A and Tacrolimus (Table 3). The Luna Omega 3 µm C18 column achieved a resolution of 3.85 between Impurity A and Tacrolimus (Table 5). Therefore, both columns are suitable for the analysis of Tacrolimus and related substances following the Ph. Eur. monograph 2244.



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

**🌐 All other countries/regions
Corporate Office USA**

t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Luna is a registered trademark and BE-HAPPY is a trademark of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.

