

Separation of Diazoxide and its Organic Impurities per USP Monograph

Lauren Nakasone, Zeshan Aqeel, and Bryan Tackett, PhD
Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501, USA

Overview


Diazoxide is a member of the thiazide family of drugs and prevents insulin release from the pancreas, helping to return blood sugar to normal levels. The development of a quick and efficient analysis of Diazoxide and its related organic impurities is of interest for generic drug manufacturers. In this application note, we report the separation of Diazoxide and its related organic impurities using a Kinetex™ 2.6 µm Biphenyl column according to the USP monograph for Diazoxide.

System suitability per USP Monograph for the Diazoxide Assay is a symmetry factor no more than (NMT) 2.0 and a percent relative standard deviation (%RSD) of NMT 0.73 % for Diazoxide. System suitability requirements for Diazoxide Assay were met by the Kinetex 2.6 µm Biphenyl column (**Figure 2**).

System suitability per USP Monograph for the Diazoxide Organic Impurities is a symmetry factor NMT 2.0, a %RSD of NMT 5.0 %, and a signal-to-noise (S/N) ratio no less than 10 for Diazoxide. All requirements for System Suitability for Diazoxide Organic Impurities were met by the Kinetex 2.6 µm Biphenyl column (**Figures 3 and 4**).

All solutions were prepared as indicated in the USP Monograph for Diazoxide. USP Diazoxide RS (Catalog No. 1186000) was purchased from USP.

Figure 1. Diazoxide



LC-UV Conditions - Assay

Column: Kinetex 2.6 µm Biphenyl
Dimension: 50 x 4.6 mm
Part No.: [00B-4622-E0](#)
Mobile Phase: Acetonitrile/Methanol/Buffer (10:10:80, v/v/v)
 Buffer: Dissolve 1.56 g of Sodium Phosphate Monobasic in 1 L of water. Adjust with Phosphoric Acid to a pH of 2.5.
Flow Rate: 1 mL/min (Isocratic)
Injection Volume: 20 µL
Temperature: 30 °C
Detector: UV @ 254 nm
System: Agilent® 1260 Binary UHPLC

LC-UV Conditions - Organic Impurities

Column: Kinetex 2.6 µm Biphenyl
Dimension: 50 x 4.6 mm
Part No.: [00B-4622-E0](#)
Mobile Phase: A: Methanol/Buffer (15:85, v/v)
 B: Acetonitrile/Methanol (70:30, v/v)

Buffer: Dissolve 1.56 g of Sodium Phosphate Monobasic in 1 L of water. Adjust with Phosphoric Acid to a pH of 2.5.

Gradient:	Time (min)	%B
	0	5
	10	15
	20	65
	25	65
	25.1	5
	30	5


Flow Rate: 0.8 mL/min
Injection Volume: 20 µL
Temperature: 30 °C
Detector: UV @ 220 nm
System: Agilent 1260 Binary UHPLC



Table 1. Preparation of Solutions

Solution	Composition
Diluent	Methanol/Buffer (50:50, v/v)
Standard Solution – Assay	0.05 mg/mL of USP Diazoxide RS in Diluent Sonicate to dissolve if needed
Standard Solution – Organic Impurities	0.5 µg/mL of USP Diazoxide RS in Diluent
Sensitivity Solution – Organic Impurities	0.25 µg/mL of USP Diazoxide RS in Diluent


Figure 2. Standard Solution – Assay



Peak No.	Analyte	Retention Time (min)	Symmetry Factor	Area %RSD
1	Diazoxide	5.43	0.96	0.11
N = 6 Injections				




Figure 3. Standard Solution – Organic Impurities



Peak No.	Analyte	Retention Time (min)	Symmetry Factor	Area %RSD
1	Diazoxide	8.74	0.93	1.59

N = 6 Injections

Figure 4. Sensitivity Solution – Organic Impurities



Peak No.	Analyte	Retention Time (min)	S/N Ratio
1	Diazoxide	8.62	25.87

N = 3 Injections



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

☉ All other countries/regions

Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex and BE-HAPPY are trademarks of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc.

Disclaimer

Phenomenex is in no way affiliated with Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.

